COAMATIC® Antithrombin - 82 1991 63

Intended use

For the quantitative determination of the heparin cofactor activity of Antithrombin (AT) in human citrated plasma.

Summary and principle

Antithrombin is the most important natural inhibitor of the coagulation cascade. By inhibiting the coagulation proteases, especially thrombin, factor Xa, and factor IXa, AT prevents uncontrolled coagulation and thrombosis. Plasma is incubated with an excess of factor Xa (FXa) in the presence of heparin. The residual quantity of FXa is determined by the rate of hydrolysis of the chromogenic substrate S-2765. The pNA release measured at 405 nm is inversely proportional to the AT level in the range from 0-120% of normal plasma.¹

1. AT + Heparin → [AT · Heparin]
2. [AT · Heparin] + FXa (excess) → [AT · Heparin · FXa] + FXa (residual)

FXa (residual) → Peptide + pNA

Composition

1. **S-2765, 8 mg** Lyophilized chromogenic substrate N-α-Cbo-D-Arg-Gly-Arg-pNA·2HCl. **1 vial**
2. **Buffer with heparin, 25 mL** Tris buffer, pH 8.2, ionic strength 0.25. **1 vial**
3. **Factor Xa, 29 nkat** Lyophilized purified bovine Factor Xa and bovine albumin. **1 vial**

PRECAUTION AND WARNINGS

Hazard class: **none**
 Hazard statements: **none**
 Precautionary statements: **none**
 This product is for *in vitro* diagnostic use.

Preparation

The reagents are reconstituted according to the specific instrument application. For microplate and test tube techniques:

1. **Substrate S-2765:** reconstitute with 10 mL of CLSI Type II water or equivalent. ²
2. **Factor Xa:** reconstitute with 10 mL of buffer with heparin. Replace the stoppers and swirl gently. Make sure of the complete reconstitution of the product. Keep reagent at 15-25°C for 10-30 min and invert before use.

NOTE: Other reagent reconstitution volumes may apply for automated methods. The reagents are not interchangeable between lots.

Storage conditions and stability

Unopened reagents are stable until the expiration date shown on the vial when stored at 2-8°C.

1. **S-2765:** Stability after reconstitution: 6 months at 2-8°C in the original vial.
2. **Buffer:** with heparin. Stability when opened: 3 months at 2-8°C in the original vial.
3. **Factor Xa:** Stability after reconstitution: 3 months at 2-8°C in the original vial.

WARNING: Do not use reagents beyond the expiry date printed on the package label. Discard if the substrate solution appears yellow. Avoid contamination by microorganisms.

Specimen collection and preparation

Nine parts of freshly drawn venous blood are collected into one part trisodium citrate. Centrifugation: 2000 x g for 10-20 minutes at 20-25°C. Refer to CLSI document H21-A5 for further instructions on specimen collection, handling and storage.⁵

Additional reagents and control plasmas

1. Deionized water, filtered through 0.22 mm or CLSI type II water.²
2. Calibration plasma
3. Control Plasma Abnormal and Normal
4. Saline (0.9% NaCl).
5. Acetic acid 20% or citric acid 2% (end-point methods).

Materials required but not provided

- Spectrophotometer, 405 nm (and 490 nm for microplate procedure)
- Microplate* or semi-micro cuvettes (1 cm)
- Centrifuge, 2000 x g
- Incubator 37°C ±0.2°C
- Vortex mixer
- Stopwatch
- Calibrated pipettes
- Plastic test tubes

*NOTE: Do not use microplates intended for coating!

Quality controls

Normal and abnormal controls are recommended for a complete quality control program.³ Assigned values of Controls should be traceable to the International Standard. Each laboratory should establish its own mean and standard deviation and should establish a quality control program to monitor laboratory testing. Controls should be analyzed at least every 8 hours in accordance with good laboratory practice. Refer to Westgard et al⁴ for identification and resolution for out of control situations.

Results

Antithrombin results are reported in activity (%).

Expected values

106±18% (2 SD, n=109; M=53 F=56) in a normal healthy population evaluated with Coamatic Antithrombin. Due to many variables which may affect results, each laboratory should establish its own normal range.

Procedures

All conditions and performance characteristics included in this package insert are referred to microplate method. Detailed instrument settings including instructions for preparation of the reagents for a variety of automated instruments are available on request from Chromogenix.

*NOTE: Not all instrument applications are available in all countries.

Calibration

A standard curve is obtained by analyzing different dilutions in saline of Calibration Plasma, to obtain antithrombin levels of 0, 25, 50, 75 and 100%. A standard about of 120% can be prepared by diluting 25 µL Calibration Plasma with 2500 µL of saline. The assigned value of calibrator plasma should be traceable to the International Standard.

Assay conditions for microplate and test tube techniques:

Samples/controls/standards	25 µL
Saline	3000 µL
Mix well	

Microplate method

Diluted samples/controls/standards to wells
 Factor Xa
 Incubate at 37°C, 90 sec
 S-2765 (pre-heat at 37°C)
 50 µL
 50 µL
 A. Kinetic method: read DA/min at 405 nm for 30-90 sec.
 B. End-point method: proceed as described below.
 Incubate at 37°C, 90 sec
 Acetic acid 20% or 2% citric acid
 50 µL

Mix
 Read the absorbance against water at 405 nm. If possible, read and subtract the absorbance at 490 nm in order to compensate for differences in the material of the microplate wells.

Test tube method


Use 200 µL instead of 50 µL for all pipetting steps.

ENGLISH - Insert revision 03/2015

Calculation

Plot the change in absorbance per minute (ΔA/min) or absorbance (A) for the standards against their antithrombin activity on linear graph paper. Plot ΔA/min or A on the Y axis and % antithrombin on the X axis. Connect the standard points with the best fitting straight line. Samples are evaluated based on this standard curve. An example of a typical standard curve (microplate method) is shown below.

Performance Characteristics

Limitations/interfering substances

No drug interference reported. Since the method is based on inhibition of factor Xa there is no influence from heparin cofactor II, α₂-macroglobulin or α₁-antitrypsin.^{6,7,8} In plasmas where contact activation has occurred, a contribution to the substrate activity might be produced. Thus an underestimation of the antithrombin level might follow. To improve the validity of the assay the value obtained in absence of factor Xa can be subtracted from the sample value. Bilirubin, haemoglobin and plasma from hyperlipaemic patients interfere in absorbance readings. In these cases individual plasma blanks are necessary when acid stopped methods are used.

Precision:

Microplate method	CV% (Within day)	n	CV% (Between days)	n
Mean (%AT)				
100	3.1	20	2.6	5
50	4.8	20	4.3	5

Correlation:

The assay has been compared with Coatest Antithrombin and Coacut Antithrombin R using normal plasmas, patient plasmas and mixtures of normal and deficiency plasmas.:

System	Slope	Intercept	r	Reference method	n
Microplate	1.03	-2.6	0.97	Coatest	45
ACL	1.01	-1.6	0.96	Coatest	43
Microplate	1.13	-7.5	0.99	Coacut	43
ACL	0.95	2.3	0.99	Coacut	43

The precision and correlation results were obtained using specific lots of reagents and controls.

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