COATEST® SP4 FVIII - 82 4094 63

Intended use
For the photometric determination of factor VIII activity in citrated plasma, such as when identifying factor VIII deficiency or monitoring patients on replacement therapy, as well as for potency estimation of FVIII concentrates.

Measurement principle
1. Factor X $\xrightarrow{\text{Factor IXa, Ca}^{2+}, \text{phospholipid}}$ factor Xa
factor VIII
2. S-2765 $\xrightarrow{\text{factor Xa}}$ Peptide + pNA

In the presence of calcium and phospholipids, factor X is activated to factor Xa by factor IXa. This generation is greatly stimulated by factor VIII, which may be considered as a cofactor in this reaction. By using optimal amounts of Ca²⁺ and phospholipids and an excess of factors IXa and X, the rate of activation of factor X is solely dependent on the amount of factor VIII. Factor Xa hydrolyses the chromogenic substrate S-2765 thus liberating the chromophoric group, pNA. The color is then read photometrically at 405 nm. The generated factor Xa and thus the intensity of color is proportional to the factor VIII activity in the sample. Hydrolysis of S-2765 by thrombin formed is prevented by the addition of the synthetic thrombin inhibitor, I-2581, together with the substrate.

- Composition**
- S-2765 15.4 mg + I-2581** 1 vial
Chromogenic substrate (N-a-Z-D-Arg-Gly-Arg-pNA), 15.4 mg, synthetic thrombin inhibitor, 0.4 mg, and mannitol added as a bulking agent.
 - Factor IXa + factor X 2.7 IU** 4 vials
Lyophilized bovine factors IXa and X with bovine albumin added as a stabilizing agent.
 - CaCl₂ 6 ml** 1 vial
Calcium chloride solution, 0,025 mol/L.
 - Buffer, stock solution 20 mL** 1 vial
20 mL concentrated Tris buffer containing NaCl and BSA. Characteristics of tenfold diluted buffer: Tris 0,05 mol/L, pH 7,3, 10 mg/L Ciprofloxacin and 1,0% BSA.
 - Phospholipid 2 mL** 1 vial
Mixture of highly purified phospholipids and 10 mg/L Ciprofloxacin.

PRECAUTION AND WARNINGS
Hazard class: **none**
Hazard statements: **none**
Precautionary statements: **none**
This product is for *in vitro* diagnostic use.

- Preparation**
The reagents are reconstituted according to the specific instrument application. For microplate and test tube techniques:
- S-2765 + I-2581:** Reconstitute with 12.0 mL of sterile water or NCCLS type II water¹¹, to obtain a concentration of 1.8 mmol/L.
 - Factor IXa + Factor X:** Reconstitute with 3.0 mL of sterile water or NCCLS type II water¹¹
 - CaCl₂:** Ready to use.
 - Buffer, stock solution:** Dilute 1:10 (1+9) with sterile water or NCCLS type II water¹¹. Prepare a new buffer working solution each day.
 - Phospholipid:** Ready to use.

- Reagent storage and stability**
When kept at 2-8°C the sealed reagents are stable until the expiry date printed on the label. Contamination by microorganisms should be avoided once the vials are opened.
- S-2765 + I-2581:** Stability after reconstitution: 3 months at 2-8°C.
 - Factor IXa + Factor X:** Stability after reconstitution: 12 hours at 2-8°C. The solution can be stored frozen in aliquots at -20°C (or at lower temperature) for 3 months. Do not refreeze.
 - CaCl₂ (0,025 mol/L):** Opened vial is stable 3 months at 2-8°C.
 - Buffer, stock solution (Tris 0,05 mol/L, pH 7,3, 10 mg/L Ciprofloxacin and 1,0% BSA):** Once opened the buffer is stable 3 months at 2-8°C. Prepare a new buffer working solution each day.
 - Phospholipid:** Opened vial is stable for 3 months at 2-8°C. Shake gently before use.

- Reagents and materials required but not provided**
- Deionized water, filtered through 0.22 µm or NCCLS type II water.¹¹
 - Acetic acid 20% or citric acid 2%.
 - Control Plasma Abnormal and Normal calibrated against an International Standard for Factor VIII
 - Calibration plasma calibrated against an International Standard
 - Photometer, 405 nm (and 490 nm for microplate procedure)
 - Heat incubator 37°C ± 0.2°C
 - Semi-micro cuvettes
 - Centrifuge, 2000xg
 - Plastic test tubes
 - Stopwatch
 - Vortex mixer
 - Calibrated pipettes

Specimen collection
Nine parts of freshly drawn venous blood are collected into one part trisodium citrate. Centrifugation: 2000 x g for 10-20 minutes at 20-25°C. Refer to NCCLS document H21-A4 for further instructions on specimen collection, handling and storage.¹²

Quality Control
Normal and abnormal controls for plasma or concentrates are recommended for reliable quality control.¹³ Assigned values of Controls should be traceable to the International Standard. Periodically within each run a control should be analyzed. The control material should be treated in the same way as a test sample. A range of allowable variation should be established for controls in each laboratory. If a value outside the established control range is obtained, a complete check of calibration, reagents and instrument performance should be made.

Results
Factor VIII results are reported in % activity (100% factor VIII activity is equivalent to 1.0 IU/mL).

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Expected values
Range: 49 - 126 % (2 SD, n=121) in a normal healthy population evaluated with Coatest SP Factor VIII (test tube method).

Due to the many variables that may affect results, each laboratory should establish its own normal range, avoiding inadvertent losses of factor VIII activity.

Procedure
NOTE: All conditions included in this package insert refer to the manual method. Detailed settings for the ACL 8000/9000/10000 including instructions for preparation of the reagents are available on request from Chromogenix.

Two ranges of factor VIII are defined (20-150% and 1-20%).
Range 20-150%:
Prepare a solution of phospholipid+factor IXa and factor X reagent by mixing:
• 1 volume of phospholipid
• 5 volumes of factor IXa+factor X reagent
Keep on ice or at 2-8°C
Shake gently just before use.
The mixture is stable for 4 hours at 2-8°C or 12 hours on ice.

Calibration
A standard curve is required for each Coatest SP Factor VIII kit. Normal human plasma, calibrated against an International Standard, is used for preparation of standard dilutions in plastic tubes using pre-cooled buffer working solution according to the table below:

Standard	Predilution		Buffer working-solution		Final dilution	
	%	Plasma µL	µL	Predilution µL	µL	Buffer working-solution µL
150	-	undiluted	25	25	2000	
120	200	50	20	25	2000	
100	100	50	25	25	2000	
75	100	100	25	2000		
50	100	200	25	2000		
21	100	600	25	2000		

The assigned percentage values of the standard dilutions are those obtained from a normal plasma containing 1.0 IU factor VIII/mL. In case the factor VIII content of the normal plasma differs from this value, an appropriate correction factor should be used.

Preparation of plasma sample
Use plastic test tubes.
Test plasma or concentrate 25 µL
Buffer working solution (2-8°C) 3000 µL
Mix well. Keep at 2-8°C.
The assay must be performed within 30 minutes after dilution because of the lability of factor VIII.

Assay
NOTE: The assay should be performed in plastic material.

Standard	Predilution		Buffer working-solution		Final dilution	
	%	Plasma µL	µL	Predilution µL	µL	Buffer working-solution µL
Phospholipid+ FIXa + FX (2-8°C)			200 µL		200 µL	
Test plasma or standard dilution (2-8°C)			100 µL		100 µL	
CaCl ₂ (37°C)			100 µL		100 µL	
Mix and incubate at 37°C exactly 5 min						
S-2765+I-2581 (37°C)			200 µL		200 µL	
Mix and incubate at 37°C exactly 5 min						
Acetic acid 20% or citric acid 2% (20-25°C)			100 µL			

Acid-stopped method: Read the absorbance of the sample against a reagent blank (buffer working solution instead of sample) within 4 hours.
Acid-stopped method: Because of the large dilution of the plasma, no sample blanks have to be included.
Initial rate method: Transfer immediately to a 1 cm semi-micro cuvette (pre-heated to 37°C) and measure the absorbance change at 405 nm.
Range 1-20%:
Prepare a solution of phospholipid+factor IXa and factor X reagent by mixing:
• 1 volume of phospholipid
• 5 volumes of factor IXa+factor X reagent
Keep on ice or at 2-8°C.
Shake gently just before use.
The mixture is stable for 4 hours at 2-8°C or 12 hours on ice.

Calibration
A standard curve is required for each Coatest SP Factor VIII kit. Normal human plasma, calibrated against the International Standard, is used for preparation of standard dilutions in plastic tubes using pre-cooled buffer working solution according to the table below:

Standard	Predilution		Buffer working-solution		Final dilution	
	%	Plasma µL	µL	Predilution µL	µL	Buffer working-solution µL
20	50	200	25	25	2000	
14,3	50	300	25	25	2000	
9,1	50	500	25	25	2000	
4,8	25	500	25	25	2000	
1,2	25	2000	25	25	2000	

The assigned percentage values of the standard dilutions are those obtained from a normal plasma containing 1.0 IU FVIII/mL. In case the FVIII content of the normal plasma differs from this value, an appropriate correction factor should be used.

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ENGLISH - Insert revision 02/2015


Preparation of plasma sample
Use plastic test tubes.
Test plasma or concentrate 25 µL
Buffer working solution (2-8°C) 2000 µL
Mix well. Keep at 2-8°C.
The assay must be performed within 30 minutes after dilution because of the lability of factor VIII.

Assay
Because of the fairly small generation of FXa in samples with <5% factor VIII, the acid-stopped method is preferred for this range of factor VIII.
NOTE: The assay should be performed in plastic material.
Phospholipid+ FIXa + FX (2-8°C) 200 µL
Test plasma or standard dilution (2-8°C) 100 µL
Mix and incubate at 37°C 4-5 min
CaCl₂ (37°C) 100 µL
Mix and incubate at 37°C exactly 10 min
S-2765+I-2581 (37°C) 200 µL
Mix and incubate at 37°C exactly 10 min
Acetic acid 20% or citric acid 2% (20-25°C) 100 µL

Read the absorbance of the sample against a reagent blank (buffer working solution instead of sample) within 4 hours. Because of the large dilution of the plasma, no sample blanks have to be included.
NOTE: The above described assay can also be conveniently performed in microplates by a four-fold reduction of all volumes, and keeping all other conditions such as incubation and hydrolysis times identical. In this case read the absorbance at 405 nm and 490 nm. Subtract the A₉₀ from the A₀₅ to correct for differences in microplate wells.

LIMITATIONS OF PROCEDURE
The activation reaction should be performed in plastic material since glass surfaces may interfere with the generation of factor Xa. Factor VIII is a labile coagulation factor and in order to obtain the accuracy which the method offers it is important to work in a carefully standardized manner throughout the assay procedure.

Calculation
Plot the change in absorbance per minute (ΔA/min) or absorbance (A) for the standards against their concentrations of factor VIII on linear graph paper. Read the % FVIII value for the corresponding absorbance for the unknown sample from the standard curve.

Specificity and Interfering Factors

FVIII results are not affected by Triglycerides at concentrations of 700 mg/dL. Bilirubin at concentrations of 20 mg/dL, Hemoglobin at concentrations of 100 mg/dL and unfractionated (UF) Heparin at concentrations of 1.0 IU/mL.
NOTE: Hemolyzed samples in the low range should not be analyzed.
Due to the high dilutions used, there is no underestimation of factor VIII activity in samples containing Lupus anticoagulant.

Precision
Within run and total precision was assessed over multiple runs.

System	%CV (Within run)	n	%CV (Total)	N
Test Tube method Mean FVIII				
83 %	3,4	2	5,3	80
14 %	4,3	2	5,6	80

Correlation:

System	Slope	Intercept	r	Reference method	n
Test Tube method	1,0851	-4,80	0,9873	Coatest Factor VIII (natural porcine phospholipid)	181

This study (n=181) was performed using samples from healthy individuals, as well as samples from patients with various levels of FVIII deficiency, von Willebrand's disease and other disorders.

Linearity System
Test Tube method: 0 -150% factor VIII

Detection Limit System
Test Tube method: The assay allows detection of 1% factor VIII activity.

Sensitivity: System
Test Tube method ΔA₄₀₅ per 1% of FVIII activity: Low range 0.034 Normal range 0.009

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