

S-2251™

For Laboratory Use Only

For General Laboratory Use

S-2251™

S-2251 is a chromogenic substrate for plasmin and streptokinase-activated plasminogen.

COMPOSITION

Each vial contains chromogenic substrate S-2251 25 mg and mannitol 60 mg as a bulking agent.

CHEMISTRY

Chemical name: H-D-Valyl-L-leucyl-L-lysine-p-Nitroaniline dihydrochloride

Formula: H-D-Val-Leu-Lys-pNA · 2HCl


Mol. wt: 551.6

$\epsilon_{316 \text{ nm}}$: $1.27 \cdot 10^4 \text{ mol}^{-1} \cdot \text{L} \cdot \text{cm}^{-1}$

Solubility: > 40 mmol/L in H₂O

Stability: Substance: Stable until expiry date if stored at 2-8°C. Avoid exposure to light. The substance is hygroscopic and should be stored dry.
Solution: 3 mmol/L in H₂O is stable for at least 6 months at 2-8°C. Contamination by microorganisms may cause hydrolysis.

Suitable stock solution: 3-4 mmol/L in H₂O.

PRINCIPLE

The method for the determination of activity is based on the difference in absorbance (optical density) between the pNA formed and the original substrate. The rate of pNA formation, i.e. the increase in absorbance per second at 405 nm, is proportional to the enzymatic activity and is conveniently determined with a photometer.

KINETIC DATA

Plasmin (human): $K_m = 3 \cdot 10^{-4} \text{ mol/L}$,
 $V = 0.5 \cdot 10^{-6} \text{ mol/min} \cdot \text{CU}$

Plasminogen · SK: $K_m = 2 \cdot 10^{-4} \text{ mol/L}$,
 $V = 1 \cdot 10^{-6} \text{ mol/min} \cdot \text{mL plasma}$.
Determined at 37°C in 2.5 mL
0.05 mol/L Tris buffer pH 7.4, I 0.5.

CHROMOGENIX

STANDARDIZATION

An activity of $\Delta A/\text{min}=0.05$ (37°C) is obtained by using a substrate concentration of $2 \cdot k_m$ and:

1. 0.010 CU/mL of human plasmin from Chromogenix AB.
2. 0.0011 U/mL of the plasmin standard from NIBSC, Potters Bar, London.
3. 0.0078 CU/mL of SK-activated human plasminogen from Chromogenix AB.

The substrate is insensitive to kallikrein (glandular and plasma) and urokinase.

APPLICATIONS

The substrate has been used for the determination of:

1. Antiplasmin in plasma (1,2,3,5)
2. Plasminogen in plasma (4,5,6,7)



1. EDY J et al.: Inhibition of plasmin by normal and antiplasmindepleted human plasma. *Thromb Res* 8, 513-518 (1976)
2. TEGER-NILSSON A-C et al.: Determination of a new rapid plasmin inhibitor in human blood by means of a plasmin specific tripeptide substrate. *Scand J Clin Lab Invest* 37, 403-409, (1977).
3. Chromogenix AB. Determination of antiplasmin in plasma with S-2251. Laboratory Instruction.
4. SORIA J et al.: Dosage du plasminogène a l'aide d'un substrat chromogène tripeptidique. *Pathologie Biologie* 24, 725-729 (1976).
5. FRIBERGER P.: Methods for the determination of plasmin, antiplasmin and plasminogen by means of the substrate S-2251. *Haemostasis* 7, 138-145, (1975)
6. PHILO R D and GAFFNEY P J.: Some observation on the assay of plasminogen using streptokinase and S-2251. *Haemostasis* 11 suppl. 1, 66 (1982).
7. Chromogenix AB. Determination of plasminogen in plasma with S-2251. Laboratory Instruction.

CHROMOGENIX

S-2251

Printed Insert Sheet: 301933
Revision: R2
Issued: December 2008
C.O.: 406499

LANGUAGES

ENGLISH

TECHNICAL SPEC'S

PAPER: White paper,
50-60 g/m² weight.
SIZE: 4.1 x 5.9" (104 x 150 mm.).
PRINT: Front/Back.
PRINT COLOR: All type in black.