
# TECHNOZYM<sup>®</sup> D-Dimer ELISA

For research  
use only

REF 2599006 TECHNOZYM<sup>®</sup> D-Dimer ELISA

symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Klíčova slova / Značenje simbola	
	manufacturer / Hersteller / fabbricante / fabricante / fabricante / fabricant / Tillverkaren / Fabrikanten / Produzent / Κατασκευαστής / Производитель / Производител / výrobce / Proizvođač
	storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage / lagringstemperatur / opbevaringstemperatur / Oppbevaringstemperatur / θερμοκρασία αποθήκευσης / съхранение на / teplota skladování / температура хранения / teplota skladování / Temperatura lagerovanja
	 determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations / bestämmingar / bestemmelser / Bestemmelser / προσδιορισμοί / брой тестове / stanovení / определенний / роčet stanovení / Definicija
<b>AQUA</b>	distilled water / destilliertes Wasser / acqua distillata / agua destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / дестилирана вода / destilovaná voda / дистиλλированная вода / destilovaná voda / Destilovaná Voda
<b>LOT</b>	lot / Charge / lotto / lote / lote / lot / sats / serie / Parti / ηαρία / партия номер / šarže / lot / šarže / Serija
<b>BUF</b>	Reaction buffer / Reaktionspuffer / tampone di reazione / tampón de reacció / Tampão de reação / tampon de réaction / Reaktionsbuffert / Reaktionsbuffer / Reaktionsbuffer / διάλυμα αντίδρασης / Реакционен буфер / Рабочий буферный раствор / Reakční pufr / Reakcioni pufer
<b>MTP</b>	microtiter plate / Mikrotiterplatte / placa microtiter / microplaca / microplaca / microplaques sensibilisées / Mikrotiterplatta / Mikrotiterplade / mikrotiterplate / πλάκα μικροπιλοδότησης / Микротитерная плака / Микропланшет / Mikrotitrační destička / Mikrotitracione ploče
<b>CAL</b>	Calibrator / Kalibrator / Calibratore / calibrador / calibrador / calibreur / Kalibrator / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrator
<b>REF</b>	catalogue number / Katalognummer / numero di catalogo / número de catálogo / número de referência / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / καταложен номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj
<b>CONJ</b>	Conjugate / Konjugat / Coniugato / conjugado / conjugado / conjugaté / Konjugerad / Konjugat / Konjugat / συνδεδετικό / Конюгат / Конъюгат / Konjugát / Konjugat
<b>RTU</b>	ready to use / gebrauchsfertig / pronto all'uso / listo para usar / pronto a usar / prêt à l'emploi / færdig att användas / færdig til brug / klar til bruk/ έτοιμο προς χρήση / Готов за употреба / готов к использованию / k přímému použití / Razrediti ili rastvoriti
<b>CONT</b>	Control / Kontrolle / controllo / control / control / contrôle / Kontroll / Kontroll / Kontroll / διάλυμα ελέγχου / Контрол / Контрольный образец / Kontrola / Kontrola
<b>STOP</b>	stop solution / Stoppløsning / Soluzione di arresto / solución de parada / solução de paragem / solution d'arrêt / Stoppløsning / Stop-opløsning / Stoppløsning / διάλυμα παύσης / Стоп раствор / Стоп-раствор / Zastavovací roztok / Stop solucija
<b>DIL</b>	dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o dissolver / diluir ou dissolver em / diluer ou dissoudre dans / spädd eller upplöst i / fortyndes eller opløses i / Fortyndes eller opløses i / αρωση ή διάλυση σε / растворите или разрежете с / zředit anebo rozpuštit v / разбавить или растворить в / nafedte nebo rozpustit v / razrediti ili rastvoriti u
<b>SUB</b>	substrate / Substrat / substrato / substrato / substrato / substrat / Substrat / Substrat / Substrat / υπόστρωμα / Субстрат / Субстрат / Substrát / Substrat
<b>INC</b>	incubation buffer / Inkubationspuffer / tampone di incubazione / tampón de incubación / tampão de incubação / tampon d'incubation / Inkubationsbuffert/ Inkubationsbuffer/ Vaskebufferkonsentrat / διάλυμα επώασης / Инкубационен буфер / Буфер для инкубации / Inkubační pufr / Inkubacioni pufer
<b>WASH</b>	washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado / solução de lavado concentrada / tampão de lavagem concentrado / Tampon de lavage concentré / Vattenlösningkoncentrat / Vaskeopløsningskoncentrat / vaskeløsningskonsentrat / συμπυκνωμένο διάλυμα πλύσης / Концентриран миещ разтвор / Концентрат промывочного раствора / Koncentrát promývacího roztoku / Koncetrat solucije za ispiranje
<b>RUO</b>	For research use only

**PRODUCT DESCRIPTION**

**INTENDED USE**

TECHNOZYM® D-Dimer ELISA can be used to determine D-Dimer concentrations in plasma. D-Dimer is generated as a specific soluble degradation product during fibrinolysis. Elevated D-Dimer levels are found in cases of disseminated intravascular coagulation (DIC), deep vein thrombosis (DVT) and pulmonary embolism (PE) but other circumstances may also lead to high D-Dimer levels such as old age, pregnancy, cancer, liver disease and infection. The ELISA system has a higher sensitivity when compared with latex agglutination tests thus giving an improved aid in the diagnosis for the exclusion of venous thromboembolism (VTE).

**COMPOSITION**

1. ELISA test strips (12), with 8 wells each, coated with anti D-Dimer monoclonal antibody; the drying agent is supplied in an aluminium bag.
2. Washing buffer concentrate (PBS; pH 7.3); containing detergent; 0.01 merthiolate; 1 bottle, 80 mL.
3. Incubation buffer (PBS; pH 7.3); contains stabiliser protein; 0.05% proclin; and blue dye, 1 bottle, 90 mL, ready for use.
4. Calibrators (Standards) numbered 1-5; lyophilised, 1 vial each. **Concentrations are lot-dependent; consult label on the vial.**
5. Control plasmas "low level" and "high level", lyophilised; 1 vial each. **Concentrations are lot-dependent; consult the label on the vial.**
6. Conjugate monoclonal Anti-D-Dimer-POX; dyed blue; 1 vial, 0.3 mL.
7. Chromogen TMB (tetramethylbenzidine); 1 bottle, 12 mL; ready for use.
8. Stopping solution: sulphuric acid 0.45 mol/L; 1 bottle 12 mL; ready for use.
9. Adhesive film for ELISA test strips (2).

**MATERIAL REQUIRED** (but not supplied with the kit)

1. Distilled water
2. Test tubes for diluting standard and samples
3. Measuring cylinder (1000 mL)
4. Precision pipettes (10, 100 and 1000 µL)
5. Variable pipette (1000 µL)
6. Multichannel and/or dispensing pipettes (100 and 200 µL)
7. ELISA washer or multichannel pipette
8. ELISA reader with 450 nm filter, with a 620 nm reference filter if available.

**WARNING AND PRECAUTIONS**

- All human blood or plasma products as well as samples must be considered as potentially infectious. They have to be handled with appropriate care and in strict observance of safety regulations. The rules pertaining to disposal are the same as applied to disposing hospital waste.
- Calibrators and control plasmas made from human blood and any individual plasma involved in the procedure is HBs Ag, HIV 1/2 Ab and HCV-Ab-negative (see labels on kit and/or bottles).
- Stopping solution (sulphuric acid) may irritate the skin. Should acid get into your eyes, wash out immediately with water and consult a doctor.
- The reagents sometimes contain preserving agents (merthiolate). Beware of swallowing! Avoid contact with skin or mucous membranes

**STABILITY AND STORAGE**

The expiry date printed on the labels applies to storage of the unopened bottles at + 2...8 °C. Stability after reconstitution/opening:

Material/Reagent	State	Storage	Stability
Calibrators, control plasmas	after reconstitution	-20 °C	6 months
ELISA test strip	after opening	2 ... 8 °C with adhesive film in plastic bag with drying agent	expiry date
Washing buffer conc.	after opening	2 ... 8 °C	6 months
Washing buffer	1+11,5 dilution of concentrate	2 ... 8 °C	3 weeks
Incubation buffer (= sample dilution buffer)	after opening	2 ... 8 °C	2 months
Conjugate	after opening	2 ... 8 °C	6 months
	working solution	room temperature	60 minutes
Chromogen TMB	after opening	2 ... 8 °C	expiry date

**TEST PROCEDURE**

**PREPARATION OF SAMPLES**

Sample material: Plasma.  
Citrate or EDTA plasma should be used. Centrifuge for 15 minutes at a minimum of 2500 g (DIN 58905). The plasma sample may be stored for 3 hours at room temperature; otherwise the sample ought to be frozen immediately after centrifugation. Stable at -20°C for 6 months. Avoid repeated freeze-thaw cycles. Use of serum samples is not recommended.

**PREPARATION OF REAGENT**

1. Before starting the test, all the required components are to be brought to room temperature.
2. Preparing the washing buffer: Dilute 1 part by volume washing buffer concentrate with 11.5 parts by volume distilled water (1+11.5). Mix well! (Diluted washing buffer concentrate = washing buffer). There may be crystalline precipitations which will dissolve at 37°C within 10 minutes.
3. Reconstituting calibrators and control plasmas: Calibrators and control plasmas are reconstituted with 0.5 mL distilled water and mixed for 10 seconds after a reconstitution time of 15 minutes (vortex mixer). Reconstituted components are clear to slightly turbid. Control plasmas, normal and low level samples are used undiluted. For high level samples (up to 15-20 µg/mL) a dilution of 1:20 (20 µL plasma + 380 µL incubation buffer) is recommended, for abnormally high samples dilute 1:200 (40 µL plasma +360 µL incubation buffer).
4. Preparing the conjugate working solution (1+50): Dilute 1 part by volume conjugate with 50 parts by volume incubation buffer

**For 8 test wells: Mix 20 µL conjugate with 1000 µL incubation buffer**

**PERFORMANCE OF THE TEST**

<b>SAMPLE INCUBATION</b> (reference 1,2)	Pipette calibrators, control plasmas, diluted samples into test wells; cover test strips with film	100 µL
	Incubate at 37°C	60 minutes
<b>WASHING</b> (reference 1,3,4)	Washing buffer	3 x 200 µL
<b>CONJUGATE REACTION</b> (reference 1,2)	Pipette conjugate working solution into wells, cover test strip with film	100 µL
	Incubate at 37°C	60 minutes
<b>WASHING</b> (reference 1,3,4)	Washing buffer	3 x 200 µL
<b>SUBSTRATE REACTION</b> (reference 1,2)	Pipette Substrate solution into test wells, cover test strip with film	100 µL
	Incubate at room temperature	10 minutes
<b>STOPPING</b> (reference 1,2)	Pipette stopping solution into wells	100 µL
<b>MEASURING</b> (reference 5)	ELISA-Reader, 450 nm	shake 10 sec., measure within 10 min.

**References**

1. Reagents of different lots must not be combined
2. Precision and performance, among others, essentially depend on the following factors:
  - Thorough mixing of all substances used for dilution, 10 sec. with Vortex Mixer
  - Test calibrators, controls and samples in duplicates
  - Incubate at indicated temperature (RT: room temperature, 20...25°C)
  - Strict observance of the order of pipetting and of the time element as indicated
  - The time for sample incubation, conjugate and substrate reaction as indicated starts after pipetting the last sample. Incubation times should not vary by more than 5%.
  - During sample incubation and conjugate reaction, the time for pipetting calibrators/control plasmas/samples and/or conjugate solutions must not exceed 60 seconds per ELISA test strip (8 wells).
  - During substrate reaction and at stopping, the time needed for pipetting the substrate and/or the stopping solution must not exceed 10 seconds per ELISA test strip. Short pipetting times may be secured by using multichannel- and dispensing pipettes.
3. Label/number strips with a water resistant pen in case the strips accidentally fall out of the frame during testing.
4. After the last washing, wells must be aspirated thoroughly, turned upside down and positioned on a blotting paper; by gentle tapping, the last remnants must be removed.
5. By measuring the difference in wave lengths at 450 and 620 nm the precision of the test is increased.

**LIMITATION OF THE TEST**

The test cannot be used with serum samples. Interference with other Fibrin degradation products (FDP) can be mostly excluded as the antibody used is specific for a neo-antigen on the D-Dimer structure.

**ANALYSIS RESULTS**

**CALCULATION OF THE RESULTS**


Setting up a reference curve: X axis: Concentration D-Dimer ng/mL  
Y axis: Extinction

Graph plot is linear-linear with a linear, point to point fit or cubic spline.

**Assessment of reference curve**

- The extinction coefficient of the highest calibrator should be between 1.0 and 2.5.
- The validity of the test may be checked on the basis of the calculated control values.

**Example of standard curve**

**Measuring concentration of samples**

- Read off the concentration from the reference curve.
- If there are samples with extinction coefficients higher than that of the highest point on the curve, they have to be prediluted with incubation buffer. The measured concentration then has to be multiplied with the dilution factor.

**REFERENCE RANGE**

Normal range for D-Dimer was determined between 0 – 250 ng/mL (n = 148). It is recommended that individual laboratories establish their own normal range.

**STANDARDISATION**

There is no international standard for D-Dimer. Technoclone's standard was calibrated against the standard of a commercial D-Dimer Latex assay.

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