

TECOmedical Group

TECO®

Human Alpha GST ELISA

Serum & Urine

Instructions for use
English

Catalogue No. TE1056
For **Research Use Only**

Symbol Description



Kit Instructions



Lot Number



Expiry Date



Storage Temperature



Manufacturer



TE 1056



Caution: caustic



Intended use



96

Tests



diapharma[®]

 **TECO**medical AG
Headquarters **TECO**medical Group

Gewerbestrasse 10
4450 Sissach
Switzerland

phone + 41(0) 985 81 00
fax + 41(0) 61 985 81 09

info@tecomedical.com
www.tecomedical.com

DIAPHARMAgroup INC
Distributor: DIAPHARMA Group

8948 Beckett Road
West Chester, OH 45069
USA

phone 800 526 5224
fax 513 860 9635

info@diapharma.com
www.diapharma.com

Technical Services

Germany phone 0800 985 99 99
France phone 0800 100 437
Benelux phone +31(0)33 4951 473


Technical Services

United States & Canada

phone 800 447 3846
email techsupport@diapharma.com

TECO® Human Alpha GST Urine ELISA Kit

CONT Reagents and Materials Supplied:

SYMBOL	DESCRIPTION	FORMAT
1	αGST Antibody Coated Microtiter Plate 12 strips of 8 wells (96 break apart wells in total), in a frame with cover plate. Ready to use.	1 plate
S	Standard Stock 2 mg/L	1 x 0.2 mL
C1	Control C1 Ready to use.	1 x 1.0 mL
C2	Control C2 Ready to use.	1 x 1.0 mL
2	Wash Buffer 50x	1 x 30 mL
3	Sample Diluent Ready to use.	1 x 30 mL
6	Enzyme Conjugate Ready to use.	1 x 12 mL
7	TMB Substrate Ready to use.	1 x 12 mL
8	Stop Solution – 1 M HCl 1 M hydrochloric acid. Ready to use.	1 x 12 mL
	Kit instructions	1 x

Storage

Store kit at 2–8 °C. Do not freeze. Store unused reagents at 2–8 °C.

Intended Use

The Alpha GST ELISA provides a method for the quantitative determination of alpha glutathione S-transferase (α GST) in human urine, serum and plasma. To assay α GST in other media or assay other GST subclasses, contact us for further information. The Alpha GST ELISA is for research use only and not for use in diagnostic procedures.

Background

URINE

In kidney, alpha glutathione S-transferase (α GST) is found in the proximal tubule region whereas pi glutathione S-transferase (π GST) is confined mainly to the distal tubules¹. Low levels of α GST are released into the urine in normal individuals, as confirmed by immunoassay and Western blot analysis². Any event which precipitates proximal tubular damage may cause increased release of α GST into urine and elevations of urinary α GST levels have been shown to be indicative of proximal tubule damage in nephrotoxicity³⁻⁵, environmental toxicity⁶, surgery⁷, acute renal failure⁸ and transplantation⁹⁻¹². The release of α GST has been shown to be associated with distal tubular damage⁶, thus simultaneous measurement of α GST and π GST may allow discrimination between proximal and tubular damage^{5, 9-11}.

SERUM

In liver, alpha glutathione S-transferase is located in the hepatocytes whereas pi GST (π GST) is confined to the intrahepatic bile duct cells^{1, 13-14}. This heterogeneous GST subclass distribution suggests that the isoenzymes have unique in vivo functions in different hepatic regions and that the detection of GST subclass levels in biological fluids would be of significant use in monitoring the integrity of specific hepatic regions. Currently, liver injury is studied by the measurement of liver enzymes such as alanine aminotransferase (ALT) and aspartate aminotransferase (AST). A disadvantage of these markers is that they are not distributed uniformly throughout the liver, the periportal concentration being greater than the centrilobular¹⁵. In contrast, α GST has been found to be equally distributed in both the centrilobular and periportal regions¹³⁻¹⁴. Since the centrilobular hepatocytes are very susceptible to damage in a variety of conditions including Allograft Rejection¹⁶⁻¹⁸, Viral Hepatitis¹⁹, and Hepatotoxicity²⁰, α GST is a more sensitive indicator of hepatic status.

Alpha GST ELISA is a specific, precise immunoassay for α GST^{21,22} and, being a quantitative test, is unaffected by modulators of enzyme activity (e.g. bile salts and bilirubin)²¹. Thus, it is now possible to use α GST quantitation to study the hepatocellular status of individuals at risk of hepatic damage.

Assay Principle

Alpha GST ELISA is a quantitative enzyme immunoassay. The test procedure is based on the sequential addition of sample, antibody-enzyme conjugate and substrate to microassay wells coated with anti- α GST IgG. The resultant color intensity is proportional to the amount of α GST present in the sample.

Materials Required and not Supplied

- Pipettes 10 μ L – 1000 μ L
- Multichannel pipettes for 100 μ L
- Graduated cylinders for reconstituting or diluting reagents
- Manual Aspiration System or Automatic washer for ELISA plates
- Aqua dest
- Vortex mixer
- ELISA plate reader suitable for 96 well formats and capable of measuring at 450 and 405 nm (Reference: 590-650 nm)
- ELISA plate shaker (500 rpm) (orbital shaker)
- Software package for data generation and analysis

Warnings and Precautions

This kit is intended for research use by professional persons only.

Follow the instructions carefully.

Observe expiration dates stated on the labels and the specified stability for reconstituted reagents. Refer to "Materials Safety Data Sheet" for more detailed safety information.

Material of animal origin used in the preparation of this kit has been obtained from animals certified as healthy but these materials should be handled as potentially infectious.

Material of human origin used in the preparation of this kit has been tested and found non reactive for HIV-1 and HIV-2 as well as for HCV antibodies and HbsAg but should, nonetheless, be handled as potentially infectious.

TECOmedical AG is not liable for loss or harm caused by non-observance of the Kit instructions.

- 1 For research use only.
- 2 Treat all specimen samples as potentially biohazardous material.
Follow General Precautions when handling contents of this kit and any patient samples.
- 3 Disposal of containers and unused contents should be done in accordance with federal and local regulatory requirements.
- 4 Use the supplied reagents as an integral unit prior to the expiration date indicated on the package label.
- 5 Store assay reagents as indicated.
- 6 Do not use coated strips if pouch is punctured.
- 7 Test each sample in duplicate.
- 8 Use of multichannel pipettes is recommended to ensure the timely delivery of liquids.
- 9
 - a. 1 M hydrochloric acid is caustic and can be harmful for skin, eyes and mucosae.
 - b. Handle TMB with care. Do not ingest. Avoid contact with skin, eyes, or clothing. Should there be any contact, wash with water. If ingested, call a physician.
- 10 A mercury-free preservative is used. Incidental contact with or ingestion of buffer solutions may cause irritation of skin, eyes or mouth. Should there be any contact, wash with water. If ingested, call a physician.

Reagents Stability and Storage

1 αGST Antibody Coated Microtiter Plate

12 break apart strips of 8 wells (96 in total) in a frame and sealed in a foil bag. Fit strip wells firmly into the frame. After opening, return any unused wells to the original foil package and seal.

Store at 2–8 °C until expiration date.

Cover for microtiter plates.

S Standard Stock

1 vial of standard stock containing αGST (2 mg/L).

Store at 2–8 °C until expiration date.

C1 Control 1

1 vial of low control. Concentration see Certificate of Analysis.

Ready to use. Store at 2–8 °C until expiration date.

C2 Control 2

1 vial of high control. Concentration see Certificate of Analysis.

Ready to use. Store at 2–8 °C until expiration date.

2 Wash Buffer 50x

1 vial of 30 ml Wash Buffer concentrate. Dilute the 1:50 concentrate with deionized or distilled water up to 1500 ml. Prepare only the volume of Wash Buffer required for the assay. Store undiluted at 2–8 °C until expiration date. The diluted washing solution is stable for 4 weeks at 2–8 °C.

3 Sample Diluent

1 vial of 30 ml.

Ready to use. Store at 2–8 °C until expiration date.

6 Enzyme Conjugate

1 vial of 12 ml.

Ready to use. Store at 2–8 °C until expiration date.

7 TMB Substrate

1 vial of 12 ml of H₂O₂ stabilized tetramethylbenzidine.

Ready to use. Store at 2–8 °C until expiration date.

8 Stop Solution – 1 M HCl

1 vial of 12 ml of 1 M hydrochloric acid.

Ready to use. Store at 2–8°C until expiration date.

Sample Collection and Storage

SERUM / PLASMA

Alpha GST ELISA can be used to measure α GST in serum, EDTA or sodium-heparin plasma samples. Collect all blood samples in an appropriate tube and observe routine precautions for venipuncture. Mix the tube immediately after collection by inverting several times. Centrifuge within 3 hours from time of collection and transfer the sample from the original tube for storage at 2-8°C. If not tested within 24 hours, aliquot the sample and store at -20°C or -80°C. Inspect samples for turbidity. Turbid samples should be centrifuged and aspirated again to remove remaining insoluble matter.

Serum and plasma samples can be stored at 20-25°C for up to 48 hours, at 2-8°C for up to one week or at -20°C for >1 year. Repeated freeze thawing of samples should be avoided to prevent loss of α GST (up to 20% drop in α GST concentration observed after 3 freeze-thaw cycles as measured by ELISA).

GST (up to

URINE

Alpha GST ELISA can be used to measure α GST in any urine sample but, due to the diurnal variation in proteinuria²³, it is important for optimal results that timed, quantitative, urine samples are collected and the collection period and volume recorded. This will enable α GST excretion to be expressed as rate (ng/min), refer to Appendix 1. Overnight or 24 hour urine samples are recommended. For the use of other collection methods and periods, contact for advice.

Do not store urine samples without the addition of Urine Stabilizing Buffer (USB). USB must be added within 12 hours of sample collection.

As soon as possible after sample collection, add 100 μ L of Urine Stabilizing Buffer (TE1050 or TE1055) to 400 μ L urine (4/5 dilution of sample), even if the samples are not to be stored. The presence of blood will not affect α GST measurements.

It is recommended that samples are assayed as soon as possible after collection. After the addition of USB, samples can be stored at 20-25°C for up to 48 hours, at 2-8°C for up to one week or at -20°C for >1 year. Repeated freeze thawing of samples should be avoided to prevent loss of α GST (up to 20% drop in α GST concentration observed after 3 freeze-thaw cycles as measured by ELISA).

Sample Preparation

SERUM / PLASMA

Immediately prior to the assay, dilute samples dilute 1/5 by adding 50µL sample to 200µL Sample Diluent.

URINE

Immediately prior to the assay, dilute samples 1/2 by adding 125µL stabilized urine sample to 125µL Sample Diluent.

Kit Controls (C1 and C2)

The control sample are ready to use and does not require dilution

Preparation Standard Curve

Mix Standard (A) by vortexing for 5 - 10 seconds. Using labelled tubes prepare further standards as follows:

αGST Standard Concentration (µg/L)	Standard Volume (µL)	Sample Diluent Volume (µL)
64 (A)	32 (S)	968
32 (B)	300 (A)	300
16 (C)	300 (B)	300
8 (D)	300 (C)	300
4 (E)	300 (D)	300
2 (F)	300 (E)	300
1 (G)	300 (F)	300
0 (H)	0	300

Note:

Please only use Standard A in the assay, if high values (>32 µg/L) are expected in the samples and if the additional measurement at 405 nm is performed (Range extension).

Alpha GST standards must be used within 30 minutes of preparation.

Assay Procedure

All determinations (standards, controls and samples) should be assayed in duplicate. When performing the assay, the standards, controls and samples should be pipetted as fast as possible (<15 minutes).

To avoid distortions due to differences in incubation times, Enzyme Conjugate, Substrate Solution and Stop Solution should be added to the plate in the same order and with the same time interval as the samples. A multichannel pipette is essential.

Allow all reagents to stand at 20–25°C for at least 30 minutes. During all incubation steps, plates should be sealed with the adhesive foil or a plastic cover. For light protection, incubate in a dark chamber or cover plate with aluminium foil.

- 1 Allocate the wells of the Microtiter plate **1** for standards, controls and samples.
- 2 Pipette 100 µl of each standards (A until H), controls (C1 and C2) and diluted samples into the corresponding wells.
- 3 Cover the wells with a plastic cover and incubate the plate for 1h ± 5 min at (20–25°C) on a shaker (500 rpm).
- 4 After incubation, aspirate the wells by using a plate washer or manually decant by inverting the plate. Wash the wells 4 times with 350 µl diluted Wash Buffer per well. After the last wash cycle tap the inverted wells on a dry absorbent surface to remove excess Wash Solution. The use of an automatic plate washer is recommended.
- 5 Following the last washing step, pipette 100 µl of the Enzyme Conjugate **6** in each well (multichannel pipette).
- 6 Cover the wells with a plastic cover and incubate the plate for 1h ± 5 min at 20–25°C on a shaker (500 rpm).
- 7 After incubation wash the wells 4 times with Wash Buffer as described in step 4.
- 8 Pipette 100 µl of the TMB Substrate **7** in each well (multichannel pipette).
- 9 Incubate the plate for 15-30 min, in the dark, at 20–25°C on a shaker (500 rpm).
- 10 Stop the reaction by adding 100 µl of Stop Solution **8** (multichannel pipette).
- 11 Read the absorbance of the wells (450 and/or 405 nm).
Reference filter at 590–650 nm.

Result Analysis

A standard curve can be established by plotting standard concentration on the x-axis (linear scale) against the absorbance of the standards on the y-axis (linear scale). The α -GST concentrations can then be read off the standard curve. A 4-parameter curve fit should be used for automatic data reduction.

If both, the normal and the extended range should be used, read both wavelengths (see below). Start with the calculation of the 450 nm results. Thereafter recalculate values $>32 \mu\text{g/L}$ by using the 405 nm results. Samples with values above $64 \mu\text{g/L}$ have to be measured again pre-diluted.

Please note:

Multiply the calculated α GST by the appropriate dilution factor in order to obtain the actual α GST. Results for stabilized urine samples should be multiplied by an additional factor of 1.25 to compensate for the dilution of sample with Urine Stabilizing Buffer.


Typical Results 450 nm (Normal range: 1 to 32 $\mu\text{g/L}$)

(Example only, not for use in calculation of actual results)

Standard	Absorbance at	
	450 nm	$\mu\text{g/L}$
B	2.975	32
C	1.347	16
D	0.548	8
E	0.229	4
F	0.113	2
G	0.070	1
H	0.042	0
C1	0,202	3.6 (2.7-4.5)
C2	1.721	19.5 (14.6-24.4)

Typical Results 405 nm (Extended range: 4 to 64 µg/L)
 (Example only, not for use in calculation of actual results)

Standard	Absorbance at 405 nm	
	Absorbance	µg/L
A	1.657	64
B	0.956	32
C	0.448	16
D	0.203	8
E	0.104	4
H	0.048	0
C1	---	---
C2	0.562	19.5 (14.6-24.4)



Reference Ranges

SERUM / PLASMA

Samples were obtained from apparently healthy donors without any clinical abnormal indications. **ELISA** in order to establish the α GST concentration in the normal population.

The reference interval (5th to 95th percentiles) for Alpha GST ELISA is 0-12.0 μ g/L in serum (n=120). The reference intervals reflect the donor population of this study group. It is recommended that each laboratory determine their own reference range appropriate for their study group.

URINE

Samples were obtained from apparently healthy donors without any clinical abnormal indications. α GST levels were determined using the Alpha GST ELISA in order to establish the α GST concentration in the normal population.

The reference interval (5th to 95th percentiles) for Alpha GST ELISA is 0-29.0 μ g/L in urine (n=120). The reference intervals reflect the donor population of this study group. It is recommended that each laboratory determine their own reference range appropriate for their study group.

Performance Characteristics

SPECIFICITY

Alpha GST ELISA is highly specific for α GST. No cross-reactivity was observed with μ GST at 500 μ g/L, or π GST at 500 μ g/L.

MEASURING RANGE AND SENSITIVITY

The limit of detection (LoD) of Alpha GST ELISA was found to be <0.5 μ g/L α GST, which corresponds to 1.25 μ g/L in a stabilized urine sample diluted 1/2. For serum/plasma, the LoD is 2.5 μ g/L for a sample diluted 1/5.

The measurement Range for serum/plasma is 5 – 320 μ g/L and for urine 2.5 to 160 μ g/L.

RECOVERY OF ADDED α GST

The mean recovery of α GST was 94% in different urine samples and 105% in serum samples of the expected value.

Sample	Added $\mu\text{g/L}$	Measured $\mu\text{g/L}$	Expected $\mu\text{g/L}$	Spike Recovery %
Urine #1	0	5,7		
	20	21,0	25,7	82
	5	9,8	10,7	92
Urine #2	0	0,0		
	20	17,1	20,0	85
	5	5,7	5,0	115
Urine #3	0	3,7		
	20	21,5	23,7	91
	5	8,7	8,7	100
Urine #4	0	4,0		
	20	21,0	24,0	87
	5	8,7	9,0	97

Sample	Added $\mu\text{g/L}$	Measured $\mu\text{g/L}$	Expected $\mu\text{g/L}$	Spike Recovery %
Serum #1	0	0,0		
	20	17,9	20,0	90
	5	5,8	5,0	116
Serum #2	0	0,0		
	20	18,5	20,0	92
	5	6,1	5,0	121
Serum #3	0	0,0		
	20	18,8	20,0	94
	5	5,9	5,0	117
Serum #4	0	0,0		
	20	19,2	20,0	96
	5	5,8	5,0	117

INTERFERENCE

Potentially interfering endogenous substances were evaluated to determine their effect on α GST recovery using Alpha GST ELISA. The endogenous substances listed below were spiked into urine pools containing endogenous α GST at a concentration of $\sim 300\mu\text{g/L}$ and assayed to determine the degree of interference. The degree of interference with each test substance is presented in the table below. The percentage bias for each interferent was calculated as:

$$\% \text{ Bias} = \left[\frac{[\alpha\text{GST}] \mu\text{g/L interferent-spiked urine}}{[\alpha\text{GST}] \mu\text{g/L non-spiked urine}} \right] \times 100 - 100$$

Interfering Substance	Interferent Concentration (mg/dL)	Interference in Urine (% Bias)
Bilirubin (conjugated)	20	0%
Bilirubin (unconjugated)	20	1%
Hemoglobin	2000	-7%
Albumin	6000	3%
Lipid*	1500	-5%
Human IgG	4	3%
Tamm-Horsfall Protein**	5	-22%

* Performed with 20% intralipid.

** The endogenous concentration of the urine sample pool used for testing was unknown. The average THP concentration in healthy subjects is estimated at 6.1 – 9.0 mg/dL²⁴ and thus, the final concentration is likely to be in excess of 11.1mg/dL.

No significant interference was observed in this assay with EDTA up to 3.4µmol/L or sodium heparin up to 3,000U/L. Studies also indicated that samples with rheumatoid factor do not cause interference.

APPENDIX 1 - Urine

EXPRESSING α GST RELEASE RATE

Excretion of α GST is constant with time, not urine volume. This means that it may be more relevant to express α GST release in terms of rate (ng/min) rather than concentration. This can be important in situations of unusual diuresis, such as oligo or polyuria. The rate of release is obtained as follows:

URINE COLLECTION

Collect urine samples as described in 'Sample Collection'. Note the time of urination (T2), time of the previous urination (T1) and the total urine volume (V).

CALCULATION OF α GST EXCRETION RATE

1. Determine urinary α GST levels (µg/L) using Alpha GST ELISA.
2. Calculate the period over which the urine was collected (T = T2 - T1) in minutes.
3. Note the urine volume in mL (V).
4. Calculate the rate of release as follows:

$$\text{ng } \alpha\text{GST/min} = \frac{[\alpha\text{GST}] \mu\text{g/L} \times V}{T}$$

REFERENCES

1. **Campbell, J.A.H. et al.** (1991). Immunohistologic localization of alpha, mu and pi class glutathione S-transferase in human tissues. *Cancer (Philadelphia)* **67 (6)**, 1608-1613.
2. **Hassett, B. and Doyle, S.** (1995). Biotrin International internal research.
2. **Goldberg, M.E. et al.** (1999). Dose of compound A, not sevoflurane, determines changes in the biochemical markers of renal injury in healthy volunteers. *Anesthesia and Analgesia* **88(2)**, 437-445.
3. **Eger II, E.I. et al.** (1997). Nephrotoxicity of sevoflurane versus desflurane anesthesia in volunteers. *Anesthesia and Analgesia* **84**, 160-168.
5. **Kirby K.B. et al.** (1997). Urinary glutathione transferase as an early marker of renal impairment in psoriasis patients treated with Cyclosporin A (CsA). Paper presented at the XIVth International Congress of Nephrology 25-29 May 1997, Sydney, Australia.
1. **Sundberg, A.G.M. et al.** (1994). Glutathione transferases in the urine: sensitive methods for detection of kidney damage induced by nephrotoxic agents in humans. *Environmental Health Perspectives* **102 (Suppl 3)**, 293-296.
2. **Cressey G et al.** (2002). Renal tubular injury after infrarenal aortic aneurysm repair. *Journal Cardiothoracic and Vascular Anesthesia* **16(3)**, 290-3.
8. **Cakalaroski, K. et al.** (1999). α -glutathione S transferases as markers of tubular cell dysfunction in acute renal failure patients. Abstract from the Third Congress of the Balkan Cities Association of Nephrology, Dialysis and Artificial Organs (BANTAO) Belgrade, Yugoslavia, 1998. *Nephrol. Dial. Transplant* **14**, 2978.
9. **Sundberg, A.G.M. et al.** (1994). Quantitation of glutathione transferase-pi in the urine by radioimmunoassay. *Nephron* **66(2)**, 162-169.
10. **Stegeman, C.A. et al.** (1996). Differential diagnosis of early graft dysfunction by urinary excretion of α GST **7(9)**, 1986.
11. **Kievit, J.K. et al.** (1997). Release of alpha-glutathione S-transferase (α GST) and pi-glutathione S-transferase (π GST) from ischemic damaged kidneys into the machine perfusate - relevance to viability assessment *Transplantation Proceedings* **29(8)**, 3591-3593.
12. **Daeman, J.W.H.C. et al.** (1997). Glutathione S-transferase as predictor of functional outcome in transplantation of machine-preserved non-heart beating donor kidneys. *Transplantation* **63(1)**, 89-93.
13. **Sundberg, A.G. et al.** (1993). Immunohistochemical localisation of alpha and pi class glutathione transferases in normal human tissues. *Pharmacology and Toxicology* **72(4-5)**, 321-331.
14. **Manning, F. et al.** (1995). Argutus Medical International Internal Research.
15. **Beckett, G. J. and Hayes, J.D.** (1993) Glutathione S-transferases: biomedical applications. *Advances in Clinical Chemistry* **30**, 281-380.
16. **Trull, A.K. et al.** (1994). Serum alpha-glutathione S-transferase: a sensitive marker of hepatocellular damage associated with acute liver allograft rejection. *Transplantation* **58(12)**, 1345-51.
17. **Platz K.P. et al.** (1997). Determination of alpha- and Pi-glutathione-S-transferase will improve monitoring after liver transplantation. *Transplant Proc.* **29 (7)**, 2827-2829.
18. **Hughes, V.F. et al.** (1997) Randomized trial to evaluate the clinical benefits of serum α -glutathione S-transferase concentration monitoring after liver transplantation. *Transplantation* **64(10)**, 1446-1452.
19. **Nelson D.R. et al.** (1995). alpha-glutathione S-transferase as a marker of hepatocellular damage in chronic hepatitis C virus infection. *Am. J. Clin. Pathol.* **104(2)**, 193-198.
20. **Murray, J. M. et al.** (1992). Indocyanine green clearance and hepatic function during and after prolonged anaesthesia: comparison of halothane with isoflurane *Br. J. Anaesth.* **68(2)**, 168-171.
21. **Rees, G.W. et al.** (1995). Evaluation of an enzyme-immunometric assay for serum alpha-gultathione S-transferase. *Ann. Clin. Biochem.* **32**, 575-583.
22. **Doyle, S. et al.** (1994). Detection of serum α glutathione S-transferase by enzyme immunoassay. *International Symposium on Liver and Drugs*, Bratislava, Slovakia.
23. **Jung, K.** (1994). Urinary enzymes and low molecular weight proteins as markers of tubular dysfunction. *Kidney International Suppl.* **47**, S29-33.
24. **Lau, W-H. et al.** (2008). Qualification and application of an ELISA for the determination of Tamm Horsfall protein (THP) in human urine and its use for screening of kidney stone disease. *Journal of Biological Sciences* **4(4)**, 215-22.

TECO® Human Alpha GST Urine ELISA

Assay Procedure – Quick Guide

- Bring samples and reagents to 20-25°C. Mix the samples well.
- Wash Buffer **2** : Dilute 1:50 with deionized or distilled water.
- Prepare standards and samples as described in 'Preparation of Reagents' and 'Sample Preparation'.

Prepare the required number of Assay Strips **1**

Pipette 100 µl standards, controls and (diluted) samples into assay wells

Incubate 1 h at 20-25°C on a shaker (500 rpm)

Aspirate and wash 4 times with 350 µl Wash Buffer, aspirate and tap the inverted wells on a clean dry absorbent surface

Pipette 100 µl Enzyme Conjugate into **6** each well

Incubate 1 h at 20-25°C on a shaker (500 rpm)

Aspirate and wash 4 times with 350 µl Wash Buffer, aspirate and tap the inverted wells on a clean dry absorbent surface

Pipette 100 µl TMB Substrate Solution **7**

Incubate 15-30 min at 20-25°C in the dark on a shaker (500 rpm)

Pipette 100 µl Stop Solution **8**

Read the Optical Density at 450 and 405 nm using a reference filter between 590-650 nm. Analyze the assay results using a 4-parameter curve fit: $y = (A-D)/(1+(x/C)^B) + D$.

 **Please read Kit instruction before using the Quick Guide**