# Coagulation Reference

GB

IT

PT

SE

NO

SRB

DE

ES

FR

DK

GR

HR

REF 5220110 Coagulation Reference 5 x 1 mL

REF 5220120 Coagulation Reference 50 x 1 mL

symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Ključova slova / Značenje simbola			
	manufacturer / Hersteller / fabbricante / fabricante / fabricante / fabricant / Tillverkaren / Fabrikanten / Produzent / Κατασκευαστής / Производител / Производитель / výrobc / Proizvođač		expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration / utgångsdatum / udløbsdato / Utløpsdato / Ημερομηνία λήξης / срок на годност / datum expirace / срок годности / datum expirace / Rok trajanja
	storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage / lagringstemperatur / oppbevaringstemperatur / Oppbevaringstemperatur / θερμοκρασία αποθήκευσης / съхранение на / teplota skladování / температура хранения / teplota skladování / Temperatura lagerovanija		consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation / se användarinstruktioner / følg bruksveiledning / Følg bruksanvisningen / συμβουλευθείτε τις οδηγίες για τη χρήση / прочетете инструкцията за работа / potfeba fídit se instrukcemi / перед использованием читайте инструкцию / sledujte návod k použití / Pročitaj uprustvo pre upotrebe
			determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations / bestämmingar / bestemmelser / Bestemmelser / προσδιορισμοί / брой тестове / stanovení / определний / роčet stanovení / Definicija
<b>AQUA</b>	distilled water / destilliertes Wasser / acqua distillata / agua destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / destilirana voda / destilovaná voda / дистиллированная вода / destilovaná voda / Destilisaná Voda	<b>LOT</b>	lot / Charge / lotto / lote / lote / lot / sats / serie / Parti / παρτία / партида номер / šarže / lot / šarže / Serija
<b>BUF</b>	Reaction buffer / Reaktionspuffer / tampone di reazione / tampón de reacció / Tampão de reação / tampon de réaction / Reaktionsbuffert / Reaktionsbuffer / Reaktionsbuffer / διάλυμα αντίδρασης / Реакционен буфер / Рабочий буферный раствор / Reakční pufr / Reakcioni pufer	<b>MTP</b>	microtiter plate / Mikrotiterplatte / placa microtiter / microplaca / microplaca / microplaques sensibilisées / Mikrotiterplatta / Mikrotiterplade / mikrotiterplate / πλάκα μικροπιλοδότησης / Микротитърна плака / Микропланшет / Mikrotitrační destička / Mikrotitracione ploče
<b>CAL</b>	Calibrator / Kalibrator / Calibratore / calibrador / calibrador / calibreur / Kalibrator / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrator	<b>REF</b>	catalogue number / Katalognummer / numero di catalogo / número de catálogo / número de referência / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / каталожен номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj
<b>CONJ</b>	Conjugate / Konjugat / Coniugato / conjugado / conjugado / conjugaté / Konjugerat / Konjugat / Konjugat / συνδεδετικό / Конюгат / Конъюгат / Konjugát / Konjugat	<b>RTU</b>	ready to use / gebrauchsfertig / pronto all'uso / listo para usar / pronto a usar / prêt à l'emploi / færdig att användas / færdig til brug / klar til bruk / έτοιμο προς χρήση / Готов за употреба / готов к использованию / k přímému použití / Razrediti ili rastvoriti
<b>CONT</b>	Control / Kontrolle / controllo / control / control / contrôle / Kontroll / Kontroll / Kontroll / διάλυμα ελέγχου / Контрол / Контрольный образец / Kontrola / Kontrola	<b>STOP</b>	stop solution / Stoppløsning / Soluzione di arresto / solución de parada / solução de paragem / solution d'arrêt / Stoppløsning / Stop-opløsning / Stoppløsning / διάλυμα πάυσης / Стоп разтвор / Стоп-разтвор / Zastavovací roztok / Stop solucija
<b>DIL</b>	dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o disolver / diluir ou dissolver em / diluer ou dissoudre dans / spädd eller upplöst i / fortyndes eller opløses i / Fortyndes eller opløses i / αραιωση ή διάλυση σε / разтворете или разреждете с / zředit anebo rozpusťt v / разбавить или растворить в / nafedte nebo rozpusťte v / razrediti ili rastvoriti u	<b>SUB</b>	substrate / Substrat / substrato / substrato / substrato / substrat / Substrat / Substrat / Substrat / υπόστρωμα / Субстрат / Субстрат / Substrát / Substrat
<b>INC</b>	incubation buffer / Inkubationspuffer / tampone di incubazione / tampón de incubación / tampão de incubação / tampon d'incubation / Inkubationsbuffert / Inkubationsbuffer / Vaskebufferkonsentrat / διάλυμα επώασης / Инкубационен буфер / Буфер для инкубации / Inkubacijski pufr / Inkubacioni pufer	<b>WASH</b>	washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado / solução de lavado concentrada / tampão de lavagem concentrado / Tampon de lavage concentré / Vattenlösingskonsentrat / Vaskeopløsningskonsentrat / vaskeløsningskonsentrat / συμπυκνωμένο διάλυμα πλύσης / Концентриран миеш разтвор / Концентрат промывочного раствора / Koncentrát promývacieho roztoku / Koncetrat solucije za ispiranje
<b>RUO</b>	For Research Use Only		

**PRODUCT DESCRIPTION****INTENDED USE**

Calibration plasma for all screening tests and coagulation factors. Coagulation Reference can be used for:

- Establishing calibration curves for screening tests.
- Establishing calibration curves for all coagulation factors given in the table.
- Accuracy control for all screening tests, i.e. prothrombin time (PT), partial thromboplastin time (aPTT), thrombin time (TT) and determination of fibrinogen (Clauss) within the normal range.
- Precision control of all parameters listed in the table.

**COMPOSITION**

Coagulation Reference is prepared from selected citrated plasma donations of healthy donors. The clotting activity is normally distributed. Thereby the "average" presence of all coagulation factors and inhibitors is secured. Coagulation Reference contains stabilizers but no bactericide additives.

**MATERIAL REQUIRED** (not supplied with the kit)

- Pipettes
- Distilled water


**WARNING AND PRECAUTIONS**

- For Research Use Only
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.
- Each single donor plasma and each lot of Coagulation Reference are tested and found negative for HbsAg, HIV 1/2 Ab and HCV Ab. However, universal precautions (treating all human source materials as if potentially infectious) should be exercised.

**STABILITY AND STORAGE**

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8 °C.

Stability after reconstitution, without loss of activity:

Freezing of the reconstituted reagent is not recommended due to possible loss of activity.

\*=room temperature

**TEST PROCEDURE****PREPARATION AND PERFORMANCE OF THE TEST**

The lyophilized Coagulation Reference is reconstituted in the volume of distilled water specified for the respective single factor determination in the table. The Coagulation Reference contains 1 Unit/mL of the respective single factor to be determined.

For the indications in g/L the Coagulation Reference is reconstituted with 1 mL of distilled water.

- Open the vial carefully and reconstitute the contents in the volume of distilled water specified in the table for the respective test by carefully rotating the vial until the product is completely reconstituted (avoid frothing).  
**Please consider that Coagulation Reference has to be reconstituted in a different volume of distilled water, when using certain prothrombin time methods such reagents from Axis Shield (see footnotes of the batch table).**
- Allow the reconstituted Coagulation Reference to stand for 10 minutes at room temperature before use<sup>1</sup>. Invert to mix before use (avoid foaming).
- Treat the reconstituted plasma as a citrated patient sample according to the instructions of the respective test.

**TABLE OF ASSIGNED VALUES**

The stated expected value, the confidence range and the coagulation time have been determined by the following institutions as part of proficiency tests (see table column "Determined by ☉").

**A** INSTAND (Institute for standardisation and documentation in medical laboratories, Duesseldorf)

**B** UK Quality Assurance Scheme

Or with reagents and methods of the indicated companies using one or more of the following standard preparations:

**C** Diagnostic Reagents Ltd

**D** Technoclone GmbH, Vienna

**E** Institute f. Clinical Chemistry and Haematology, St.Gallen

1. international standard (WHO)
2. national standard
3. fresh citrated plasma pool from 8 to 10 selected donors
4. frozen plasma pools of at least 100 donors
5. highly-purified protein preparations

**LIMITATION OF THE TEST**

The figures in the table are only applicable to the indicated lot number of Coagulation Reference and the method has to adhere strictly as prescribed for the respective reagent. The guidelines for evaluation given in this leaflet have been established assuming that a constant sensitivity of different batches is guaranteed by the manufactures of each reagent.

**ANALYSES RESULTS****EVALUATION OF THE RESULTS**

The figures in the table are only applicable to the indicated lot number of Coagulation Reference and the method has to adhere strictly as prescribed for the respective reagent.

**STANDARDIZATION**

See "Table of assigned values".

**LITERATURE**

Please contact Technoclone or your local distributor.

<sup>1</sup> For standardization a reconstitution time of 30 minutes is recommended.

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