# Technoclot Reference

REF 5220170 Technoclot Reference 10 x 1 mL

REF 5220175 Technoclot Reference 50 x 1 mL

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**PRODUCT DESCRIPTION****INTENDEND USE**

Calibration plasma for research tests and coagulation factors. Technoclot Reference can be used for:

- Establishing calibration curves for research tests.
- Establishing calibration curves for all coagulation factors and inhibitors given in the table.
- Accuracy control for research tests, prothrombin time (PT), partial thromboplastin time (aPTT), thrombin time (TT) and determination of fibrinogen (Clauss) within the normal range.
- Precision control of all parameters listed in the table.

**COMPOSITION**

Technoclot Reference is prepared from selected citrated plasmas of healthy donors. By establishing this pool the clotting activity of the single donor plasmas is normally distributed. Thereby the "average" presence of all coagulation factors and inhibitors is secured. Technoclot Reference contains stabilizers but no bactericide additives.

**MATERIAL REQUIRED** (not supplied with the kit)

- Pipettes
- Distilled water

**WARNING AND PRECAUTIONS**

- RUO for *research* use only
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.
- Each single donor plasma and each lot of Technoclot Reference are tested and found negative for Hb<sub>s</sub>Ag, HIV 1/2 Ab and HCV Ab. However, universal precautions (treating all human source materials as if potentially infectious) should be exercised.

**STABILITY AND STORAGE**

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8 °C.

Stability after reconstitution, without loss of activity:

----- RT* ----- 2 hours
----------------------------------

Freezing of the reconstituted reagent is not recommended due to possible loss of activity.

\*=room temperature

**TEST PROCEDURE****PREPARATION AND PERFORMANCE OF THE TEST**

The lyophilized Technoclot Reference is reconstituted in the volume of distilled water specified for the respective single factor determination in the table. The Technoclot Reference contains 1 Unit/mL of the respective single factor to be determined.

For the indications in g/L the Technoclot Reference is reconstituted with 1 mL of distilled water.

- Open the vial carefully and reconstitute the contents in the volume of distilled water specified in the table for the respective test by carefully rotating the vial (avoid frothing).
- Allow the reconstituted Technoclot Reference to stand for 10 minutes at room temperature before use <sup>1</sup>.
- Treat the reconstituted plasma as a citrated patient sample according to the instructions of the respective test.

**TABLE OF ASSIGNED VALUES**

The stated expected value, the confidence range and the coagulation time have been determined by Technoclone.

**LIMITATION OF THE TEST**

The figures in the table are only applicable to the indicated lot number of Technoclot Reference and the method has to adhere strictly as prescribed for the respective reagent. The guidelines for evaluation given in this leaflet have been established assuming that a constant sensitivity of different batches is guaranteed by the manufactures of each reagent.

**ANALYSES RESULTS****EVALUATION OF THE RESULTS**

The figures in the table are only applicable to the indicated lot number of Technoclot Reference and the method has to adhere strictly as prescribed for the respective reagent.

**STANDARDIZATION**

See "Table of assigned values".

**LITERATURE**

Please contact Technoclone or your local distributor.

<sup>1</sup> For standardization a reconstitution time of 30 minutes is recommended.

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