TECHNOCHROM[®] FXIII

GB

DE

FR

REF 5360010 TECHNOCHROM[®] FXIII

3 x 3 mL

Symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Ključova slova / Značenje simbola			
	Manufacturer / Hersteller / fabbricante / fabricante / fabricante / fabricant / Tillverkaren / Fabrikanten / Produzent / Κατασκευαστής / Производител / Производител / výrobce / Proizvođač		Expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration / utgångsdatum / udløbsdato / Utløpsdato / Ημερομηνία λήξης / срок на годност / datum expirace/ срок годности / datum expirace / Rok trajanja
	Storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage / lagringstemperatur / opbevaringstemperatur / Oppbevaringstemperatur / θερμοκρασία αποθήκευσης / съхранение на / teplota skladování / температура хранения / teplota skladování / Temperatura lagerovanja		Consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation / se användarinstruktioner / følg brugsvejledning / Følg bruksanvisningen / συμβουλευθείτε τις οδηγίες για τη χρήση / прочетете инструкцията за работа / potfeba řidit se instrukcemi / перед использованием читайте инструкцию / sledujte návod k použití / Pročitaj upustvo pre upotrebe
			Determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations / bestämmingar / bestemmelse / Bestemmelse / προσδιορισμοί / брой тестове / stanovení / определний / počet stanovení / Definicija
AQUA	Distilled water / destilliertes Wasser / acqua distillata / agua destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / destilipirana voda / destilovaná voda / дистиллированная вода / destilovaná voda / Destiliana Voda	LOT	Lot / Charge / lotto / lote / lote / lot / sats / serie / Parti / ηαρία / партида номер / šarže / лот / šarže / Serija
BUF	Reaction buffer / Reaktionspuffer / tampone di reazione / tampón de reacció / Tampão de reação / tampon de réaction / Reaktionsbuffert / Reaktionsbuffer / Reaksjonsbuffer / διάλυμα αντίδρασης / Реакционен буфер / Рабочий буферный раствор / Reakční pufr / Reakcioni pufer	MTP	Microtiter plate / Mikrotiterplatte / placa microtiter / microplaca / microplaca / microplaques sensibilisées / Mikrotiterplatta / Mikrotiterplade / mikrotiterplate / πλάκα μικροτιτλοδότησης / Микротитърна плака / Микропланшет / Mikrotitrační destička / Mikrotitracione ploče
CAL	Calibrator / Kalibrator / Calibratore / calibrador / calibrador / calibreur / Kalibrator / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrator	REF	Catalogue number / Katalognummer / numero di catalogo / número de catálogo / número de referencia / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / каталожен номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj
CONJ	Conjugate / Konjugat / Coniugato / conjugado / conjugado / conjugaté / Konjugerad / Konjugat / Konjugat / συνδεδετικό / Конюгат / Конъюгат / Konjugát / Konjugat	RTU	Ready to use / gebrauchsfertig / pronto all'uso / listo para usar / pronto a usar / prêt à l'emploi / færdig att användas / færdig til brug / klar til bruk / έτοιμο προς χρήση / Готов за употреба / готов к использованию / k přímému použití / Razrediti ili rastvoriti
CONT	Control / Kontrolle / controllo / control / control / contrôle / Kontroll / Kontroll / Kontroll / διάλυμα ελέγχου / Контрол / Контрольный образец / Kontrola / Kontrola	STOP	Stop solution / Stopplösung / Solución de arresto / solución de parada / solução de paragem / solution d'arrêt / Stoppløsning / Stop-opløsning / Stoppløsning / διάλυμα παύσης / Стоп разтвор / Стоп-разтвор / Zastavovací roztok / Stop solucija
DIL	Dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o dissolver / diluir ou dissolver em / diluer ou dissoudre dans / späd eller upplös i / fortyndes eller opløses i / Fortyndes eller opløses i / αραιωση ή διάλυση σε / разтворете или разрежете с / zředit anebo rozpustit v / разбавить или растворить в / naředite nebo rozpustte v / razrediti ili rastvoriti u	SUB	Substrate / Substrat / substrato / substrato / substrato / substrat / Substrat / Substrat / Substrat / υπόστρωμα / Субстрат / Субстрат / Substrát / Substrat
INC	Incubation buffer / Inkubationspuffer / tampone di incubazione / tampón de incubación / tampão de incubação / tampon d'incubation / Inkubationsbuffert/ Inkubationsbuffer/ Vaskebufferkonsentrat / διάλυμα επώασης / Инкубационен буфер / Буфер для инкубации / Inkubační pufr / Inkubacioni pufer	WASH	Washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado / solución de lavado concentrada / tampão de lavagem concentrado / Tampon de lavage concentré / Vattenlösningsskoncentrat / Vaskeopløsningskonsentrat / vaskeløsningskonsentrat / συμπυκνωμένο διάλυμα πλύσης / Концентриран миеш разтвор / Концентрат промывочного раствора / Koncentrát promývaciho roztoku / Koncetrat solucije za ispiranje
FI C	For Research Use Only		


PRODUCT DESCRIPTION

INTENDED USE

The TECHNOCHROM® FXIII kit is a reagent kit for the determination of blood coagulation Factor XIII (FXIII) activity to detect inherited or acquired FXIII deficiencies, abnormal FXIII with decreased activity and elevated FXIII level. Inherited FXIII deficiency is a rare, but severe bleeding diathesis with occasional wound healing impairment and in women with habitual abortion. Acquired FXIII deficiency due to an anti-FXIII autoantibody is also a very severe haemorrhagic diathesis. Consumption of FXIII in various diseases (malignant diseases, Chron disease Henoch-Schoenlein purpura, major surgery, etc...) usually results in moderate decrease of FXIII level. The assay can also be used for monitoring FXIII replacement therapy. Factor XIII (FXIII) is transformed into an active transglutaminase (active FXIII; FXIIIa) during its activation. The determination of FXIII activity is based on the measurement of ammonia released during the transglutaminase reaction. FXIII present in the plasma sample is activated by thrombin and Ca²⁺. Fibrin formed from fibrinogen by thrombin accelerates this reaction. The polymerisation of fibrin is prevented by a tetra peptide. The formed FXIIIa then cross-links the amine substrate glycine ethyl ester (GEE) to the glutamine residue of specific peptide substrate PI(1-12), and ammonia is released. In the indicator reaction the amount of released ammonia is monitored in a glutamate dehydrogenase catalysed NADPH-dependent reaction. The consumption of NADPH is measured spectrophotometrically by the decrease of absorbance at 340 nm. Within a time window the decrease of absorbance is directly proportional to the FXIII activity.

TEST PRINCIPLE

COMPOSITION

The TECHNOCHROM® FXIII Kit for photometric FXIII determinations contains:

mL	reagent	other data
3 x 3	FXIII Activator Reagent	Hepes, Thrombin, Polybrene, Fibrin pol. inh tetrapept.
3 x 3	FXIII Detection Reagent	Hepes, Glycine, α-ketoglutarate, PI, ADP
3 x 3	FXIII NADPH Solution	NADPH
3 x 1	FXIII Inhibitor Reagent	2-iodoacetamide
1 x 6	Stabilizer Solution	Na-acid

MATERIAL REQUIRED (not supplied with the kit)

- Pipettes
- Distilled water
- Reference Plasma**

REF	5220110	Coagulation Reference	5 x 1 mL
-----	---------	-----------------------	----------
- Control Plasma Normal and Abnormal**

REF	5020040	Coagulation Control N	5 x 1 mL
REF	5021055	Coagulation Control A	5 x 1 mL

** or any other package sizes, special Ceveron® alpha or TECHNOCLOT® Control and Calibration reagents of Technoclone.

WARNING AND PRECAUTIONS

- For Research Use Only
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.
- Each single donor plasma and each lot of Ref. Standards are tested and found negative for Hb_sAg, HIV 1/2 Ab and HCV Ab. However, universal precautions (treating all human source materials as if potentially infectious) should be exercised.
- The Stabilizer solution contains sodium azide. Do not ingest and avoid contact with the skin and mucous membranes. It can be removed with plenty of water. Sodium azide can form explosive azides when contacting heavy metals (e.g. copper or lead)!

STABILITY AND STORAGE

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8 °C.

Stability after reconstitution :

Stability	+2...8°C	-20°C
FXIII Activator Reagent	3 days	3 months
FXIII Detection Reagent	3 days	3 months
FXIII NADPH Solution	3 days	3 months
FXIII Inhibitor Reagent	3 days	3 months

Stability Working Solution	+2...8°C	-20°C
Working Solution	2 days	1 month

TEST PROCEDURE

PREPARATION OF PLASMA SAMPLES

Plasma separation:

Mix 9 parts of venous blood and 1 part sodium citrate solution (0.11 mol/L) and centrifuge for 15 minutes at a RCF of at least 2500 g (corresponding to DIN 58905). The plasma sample can be stored at room temperature up to 24 hours, or at 4°C up to 3 days; otherwise the sample has to be frozen immediately after centrifugation. Stability at -20°C or -70°C is 6 month.

PREPARATION OF REAGENTS

All reagents including distilled water should have reached room temperature before use. The lyophilized reagents are dissolved in the volume of distilled water, NADPH or Stabilizer solutions indicated on the vial and are ready for use after 10 minutes. For standardization test a reconstitution time of 30 min is recommended.

- 1 FXIII Activator Reagent: DIL 3 mL NADPH 3
- 2 FXIII Detection Reagent: DIL 3 mL AQUA
- 3 FXIII NADPH Solution: DIL 3 mL AQUA
- 4 FXIII Inhibitor Reagent: DIL 1 mL STAB 5
- 5 FXIII Stabilizer Solution: RTU

WORKING SOLUTION

Prepare the working solution by mixing equal parts of reagent 1 FXIII Activator Reagent and 2 FXIII Detection Reagent and use this mixture for the preparation of sample reagent solution and blank reagent solution.

SAMPLE REAGENT SOLUTION

Add 1/20 part of 5 Stabilizer solution to 1 part of the working solution. The reagent solutions prepared this way are used for the assays.

BLANK REAGENT SOLUTION

Add 1/20 part of 4 FXIII Inhibitor Reagent to 1 part of the working solution.

PERFORMANCE OF THE TEST

Wavelength: 340 nm Measurements are done at + 37°C
In temperature controlled spectrophotometer, microplate readers (ELISA) or clinical laboratory analyser.

Pipetting scheme manual method:

Reagent / Solution	Sample	Blank
Plasma	25 µL	25 mL
Sample reagent solution	250 µL	-
Blank reagent solution	-	250 µL

Read absorbance after 5 min.
Repeat reading after exactly 2, 4 and 6 min.

If the assay is adapted to a clinical laboratory analyzer, sample and reagent volumes are proportionally decreased, and the measuring interval is set to 5-10 min.

It is important to use plasma blank to correct for reactions leading to FXIIIa independent NH₃ production and NADPH consumption. Although in most cases such reactions cause only a few percent of alteration in the FXIII activity results, in patient with severe FXIII deficiency the omission of such correction might lead to clinically significant overestimation of FXIII activity.

For each new lot of TECHNOCHROM® FXIII, the FXIII Activity of the calibration plasma should be determined – the result will be used for the analysis.

APPLICATIONS FOR INSTRUMENTS

Application Sheets for auto analyzers are available from Technoclone or your local distributor upon request.

LIMITATIONS OF THE TEST

Bilirubin up to 200 µmol/L and triglycerides up to 7.5 mmol/L do not interfere with the assay. At higher bilirubin or triglyceride concentrations the assay may be repeated with diluted plasma sample.

High concentration of ammonia (>200 µmol/L) may result in significant NADPH consumption, i.e. the absorbance significantly decreases during the lag phase of the reaction. This may lead to an underestimation of FXIII activity. In these rare cases dilute the plasma sample with physiological saline solution and repeat the assay.

ANALYSES RESULTS

From the absorbance changes during the measuring periods mean ΔA/min values of the sample and the sample blank are calculated. Using automatic kinetic photometers or clinical laboratory analysers mean ΔA/min is directly obtained. FXIII activity can be evaluated by using a factor.

The factor (F) is calculated from the ΔA/min value measured for the calibration plasma and from the assigned FXIII activity of the calibration plasma according to the following formula:

$$F = \frac{\text{FXIII}_{\text{cal.}}}{\Delta A_{\text{cal.}} / \text{min} - \Delta A_{\text{cal. blank}} / \text{min}}$$

ΔA/min value of the sample is multiplied by the factor to obtain results as percentage of normal average.

$$\text{FXIII sample} = F \times (\Delta A_{\text{sample}} / \text{min} - \Delta A_{\text{sample blank}} / \text{min})$$

STANDARDIZATION

The Coagulation Reference is calibrated against the FXIII Reference Preparation of the WHO.

REFERENCE RANGE

69 – 143 % of normal

LINEARITY AND DETECTION LIMIT

The method is linear up to 300% FXIII activity.
The detection limit is 0.6%.

LITERATURE

Please contact Technoclone or your local distributor.

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