
# TECHNOCHROM® F VIII:C

REF 5344101 TECHNOCHROM® F VIII:C

symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Ključova slova / Značenje simbola			
	manufacturer / Hersteller / fabbricante / fabricante / fabricant / Tillverkaren / Fabrikanten / Produzent / Κατασκευαστής / Производитель / Производител / výrobc / Proizvođač		expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration / utgångsdatum / udløbsdato / Utløpsdato / Ημερομηνία λήξης / срок на годност / datum expirace/ срок годности / datum expirace / Rok trajanja
	storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage / lagringstemperatur / orbevaringstemperatur / Oppbevaringstemperatur / θερμοκρασία αποθήκευσης / съхранение на / teplota skladování / температура хранения / teplota skladování / Temperatura lagerovanja		consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation / se användarinstruktioner / følg brugsvejledning / Følg bruksanvisningen / συμβουλευθείτε τις οδηγίες για τη χρήση / прочетете инструкцията за работа / potřeba řídit se instrukcemi / перед использованием читайте инструкцию / sledujte návod k použití / Pročitaj uputstvo pre upotrebu
	CE-mark / CE-Kennzeichnung / marchio CE / marca de CE / Simbolo da CE / marquage CE / CE-märkning / CE-mærket / CE-merke / CE-σημάδι / CE марка / CE-označení / маркировка CE / značka CE / CE-marka		determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations / bestämmingar / bestemmelser / Bestemmelser / προσδιορισμοί / брой тестове / stanovení / определний / počet stanovení / Definicija
<b>AQUA</b>	distilled water / destilliertes Wasser / acqua distillata / agua destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / destilirana voda / destilovaná voda / дистиллированная вода / destilovaná voda / Destilovaná Voda	<b>LOT</b>	lot / Charge / lotto / lote / lote / lot / sats / serie / Parti / παρτίδα / партида номер / šarže / lot / šarže / Serija
<b>BUF</b>	Reaction buffer / Reaktionspuffer / tampone di reazione / tampón de reacció / Tampão de reação / tampon de réaction / Reaktionsbuffert / Reaktionsbuffer / Reaktionsbuffer / διάλυμα αντίδρασης / Реакционен буфер / Рабочий буферный раствор / Reakční pufr / Reakcioni pufer	<b>MTP</b>	microtiter plate / Mikrotiterplatte / placa microtiter / microplaca / microplaca / microplaques sensibilisées / Mikrotiterplatta / Mikrotiterplade / mikrotiterplate / πλάκα μικροπιλοδότησης / Микротитерная плака / Микропланшет / Mikrotitrační destička / Mikrotitracione ploče
<b>CAL</b>	Calibrator / Kalibrator / Calibratore / calibrador / calibrador / calibreur / Kalibrator / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrator	<b>REF</b>	catalogue number / Katalognummer / numero di catalogo / numéro de catálogo / número de referência / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / каталожен номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj
<b>CONJ</b>	Conjugate / Konjugat / Coniugato / conjugado / conjugado / conjugaté / Konjugerad / Konjugat / Konjugat / συνδεδετικό / Конюгат / Конъюгат / Konjugát / Konjugat	<b>RTU</b>	ready to use / gebrauchsfertig / pronto all'uso / listo para usar / pronto a usar / prêt à l'emploi / færdig att användas / færdig til brug / klar til bruk / έτοιμο προς χρήση / Готов за употреба / готов к использованию / k přímému použití / Razrediti ili rastvoriti
<b>CONT</b>	Control / Kontrolle / controllo / control / control / contrôle / Kontroll / Kontroll / Kontroll / διάλυμα ελέγχου / Контрол / Контрольный образец / Kontrola / Kontrola	<b>STOP</b>	stop solution / Stopplösung / Soluzione di arresto / solución de parada / solução de paragem / solution d'arrêt / Stopplösning / Stop-opløsning / Stoppløsning / διάλυμα παύσης / Стоп разтвор / Стоп-раствор / Zastavovací roztok / Stop solucija
<b>DIL</b>	dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o dissolver / diluir ou dissolver em / diluer ou dissoudre dans / späd eller upplös i / fortyndes eller opløses i / Fortyndes eller opløses i / αραίωση ή διάλυση σε / разтворете или разрежете с / zředit nebo rozpustit v / разбавить или растворить в / naředte nebo rozpustte v / razrediti ili rastvoriti u	<b>SUB</b>	substrate / Substrat / substrato / substrato / substrato / substrat / Substrat / Substrat / Substrat / υπόστρωμα / Субстрат / Субстрат / Substrát / Substrat
<b>INC</b>	incubation buffer / Inkubationspuffer / tampone di incubazione / tampón de incubación / tampão de incubação / tampon d'incubation / Inkubationsbuffert/ Inkubationsbuffer/ Vaskebufferkonsentrat / διάλυμα επώασης / Инкубационен буфер / Буфер для инкубации / Inkubací pufr / Inkubacioni pufer	<b>WASH</b>	washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado / solución de lavado concentrada / tampão de lavagem concentrado / Tampon de lavage concentré / Vattenlösningkoncentrat / Vaskeopløsningskoncentrat / vaskeløsningskoncentrat / συμπυκνωμένο διάλυμα πλύσης / Концентриран миеш разтвор / Концентрат промывочного раствора / Koncentrat promývavého roztoku / Koncentrat solucije za ispiranje
<b>RUO</b>	for research use only		


**PRODUCT DESCRIPTION**

**INTENDED USE**

TECHNOCHROM® F VIII:C contains reagents for the photometric determination of Factor VIII activity in plasma and for plasma derivatives. It can be used for assaying Factor VIII deficiencies as well as monitoring Factor VIII substitution therapies.

**TEST PRINCIPLE**

**COMPOSITION**

Reagent kit for 2x20 photometric FVIII:C determinations

mL	reagent	other data
2 x 2	Substrate FXa-1+ αNAPAP	10 μmol / vial F Xa-1, 0.012 μmol / vial αNAPAP
2 x 2	Reagent A	Phospholipid, albumin
2 x 2	Reagent B	F IX <sub>a8</sub> , F X, Ca <sup>++</sup> , albumin, thrombin
1 x 1	Ref. Stand. FVIII 1	~ 130 % or 1.30 IU FVIII / mL
1 x 1	Ref. Stand. FVIII 2	~ 70 % or 0.70 IU FVIII / mL
1 x 1	Ref. Stand. FVIII 3	~ 10 % or 0.1 IU FVIII / mL
1 x 1	Ref. Stand. FVIII 4	≤ 0.5 % or 0.005 IU FVIII / mL
3 x 30	F VIII dilution buffer	3.4 g/L Imidazole, 5.85 g/L NaCl, 0.2% Albumin, pH 7.4
2 x 8	F VIII reaction buffer	6.06 g/L Tris, 3.03 g/L Na <sub>2</sub> EDTA, 25 g/L NaCl, pH 8.3

**MATERIAL REQUIRED** (not supplied with the kit)

- Pipettes - Distilled water
- For the endpoint method: 20 % acetic acid.
- Control Plasma Normal and Abnormal

<b>REF</b>	5020040	Coagulation Control N	5 x 1 mL
<b>REF</b>	5021055	Coagulation Control A	5 x 1 mL

**WARNING AND PRECAUTIONS**

- RUO for research use only
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.
- Each single donor plasma and each lot of Reference Standards are tested and found negative for Hb<sub>s</sub>Ag, HIV 1/2 Ab and HCV Ab. However, universal precautions (treating all human source materials as if potentially infectious) should be exercised.

**STABILITY AND STORAGE**

The expiry date printed on the labels applies to storage of the unopened bottles at +2... 8 °C.

Stability after reconstitution:

Stability	+37°C	+20°C	+15°C	+4°C	-20°C
substrate	8 hours	1 month	1 month	6 months	1 year
Reagent A	-	17 hours	24 hours	2 days	14 days
Reagent B	-	8 hours	24 hours	2 days	14 days
Ref. Standards F VIII	-	8 hours	8 hours	2 days	1 month
substrate-buffer-mixture	8 hours	8 hours	-	-	-

**TEST PROCEDURE**

**PREPARATION OF PLASMA SAMPLES**

Plasma separation:

Mix 9 parts of venous blood and 1 part sodium citrate solution (0.11 mol/L) and centrifuge for 15 minutes at a RCF of at least 2500 g (corresponding to DIN 58905). The plasma sample may not be stored at room temperature for more than three hours; otherwise the sample has to be frozen immediately after centrifugation. Stability at -20°C is 1 month.

Sample preparation:

Before testing the plasma samples are diluted with FVIII dilution buffer at a ratio of 1:41 (0.05 ml sample + 2.00 ml FVIII dilution buffer). For plasma samples with F VIII activity of 0 – 20% the samples should be dissolved 1:10 (100 μL plasma and 900 μL FVIII dilution buffer).

**PREPARATION OF REAGENT**

All reagents including distilled water should have reached room temperature before use. The lyophilized reagents are dissolved in the volume of distilled water indicated on the vial and are ready for use after 10 minutes. For standardization test a reconstitution time of 30 min is recommended. The reconstituted substrate (2 mL) is diluted with 8 mL FVIII reaction buffer (1+4).

**PERFORMANCE OF THE TEST**

Wavelength: 405 nm

Reagent A and B and the diluted sample are kept at room temperature, the substrate-buffer-mixture at + 37°C. Measurements are done at + 37°C using air as blank.

Pipetting scheme: Pipette into plastic tubes or cuvettes.

Kinetic determination		Endpoint determination
100 μL	Reagent A	100 μL
100 μL	Reagent B	100 μL
100 μL	diluted sample	100 μL
5 minutes	incubation 37°C Substrate-buffer-mixture	5 minutes
500 μL		500 μL
Record the increase of absorption at 405 nm per minute at 37°C. The reaction is linear for 3min.	incubation 37°C Acetic acid 20%	5 minutes 200 μL
		The extinction is measured at 405 nm against blank.

Pipetting scheme: microtiter plates

Kinetic determination		Endpoint determination
25 μL 25 μL 25 μL	Reagent A Reagent B diluted sample	25 μL 25 μL 25 μL
5 Minuten 125 μL	incubation 37°C Substrate-buffer-mixture	5 minutes 125 μL
Record the increase of absorption at 405 nm per minute at 37°C. The reaction is linear for 3min.	incubation 37°C Acetic acid 20%	5 minutes 50 μL

**LIMITATION OF THE TEST**


- The values found when testing Coagulation Control N and Coagulation Control A should be compared to the value given on the data-key for the corresponding lot of Control plasma.
- If the results obtained are outside the recommended range, avoid measuring patient samples until the problem is solved.
- A new calibration is required for each batch of TECHNOCHROM® F VIII:C and for each instrument used. Also a new calibration is recommended, if software changes are introduced or following a major service of either instruments or equipment.

**ANALYSES RESULTS**

**CALCULATION OF THE RESULTS**

To establish a calibration curve each Reference Standard is dissolved in the volume of distilled water indicated on the vial, diluted 1:41 (0.05 mL plasma + 2.00 mL F VIII dilution buffer) and tested in the same way as a patient sample in the assay. The increasing absorption (ΔA/min) in the kinetic method or the absorptions (A) in the endpoint method are plotted on linear graph paper against the % F VIII values indicated on the label of the respective Reference Standard to give a calibration curve. The blank reading may be used in the reference curve as the 0 % F VIII value.

Example: manual method (kinetic and endpoint):

All samples diluted 1:41 can be directly read off from the calibration curve. For samples diluted other than 1:41, the % activity read off from the calibration curve has to be converted as follows:

$$\frac{\% F VIII (\text{calibration curve})}{41} \times \text{actual dilution ratio} = \% F VIII \text{ sample}$$

**NORMAL RANGE**

0.6 – 1.5 IU F VIII (60 - 150% of normal)<sup>1</sup>

**STANDARDIZATION**

The Reference Standards F VIII were calibrated against the WHO plasma standard. Concentrations are lot-dependent, consult the label on the vials.

**PERFORMANCE CHARACTERISTICS**

Performance data are given below. Results obtained in individual laboratories may differ.

**PRECISION**

Reproducibility was determined with different samples (in series and day to day). The following results were obtained:

Sample	Intra assay		Inter assay	
	Sample 1	Sample 2	Sample 1	Sample 2
n	29	18	8	5
MV	42.2 Ext	81.0 Ext	90.3 %	22.5%
SD (%)	1.87	1.69	2.786	1.46
CV (%)	2.5%	2.1%	3.08	6.5%

**COMPARISON OF METHODS OR CORRELATION**

Following correlation (%) was obtained in comparing TECHNOCHROM® FVIII:C (Technoclone) with the FVIII deficient Plasma Imidazole buffer method

$$n=42 \quad y = 1,1071x - 5,9113 \quad R^2 = 0,957$$

**LINEARITY**

0.0 – 144 (Activity %)

**DETECTION LIMIT**

0% (Activity %)

**LITERATURE**

Please contact Technoclone or your local distributor for literature or technical applications for the test.


(1) H.Lang, M. Oberreiter, B. Moritz: A new Chromogenic Factor VIII:C-Assay Based on a System of Human Proteins. Thromb. Haemostas. Abstracts XIII<sup>th</sup> Congress ISTH, Th 65 (6), 943 (1991).

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REF	5020040	Coagulation Control N	5 x 1 mL
REF	5021055	Coagulation Control A	5 x 1 mL

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