
TECHNOZYM® ADAMTS-13 INH

For Research
Use Only

REF	5450401	TECHNOZYM® ADAMTS-13 INH	
REF	5450451	TECHNOZYM® ADAMTS-13 INH	
REF	5450461	TECHNOZYM® ADAMTS-13 INH Calibrator Set	5 x 0.5 mL
REF	5450463	TECHNOZYM® ADAMTS-13 INH Control Set	2 x 0.5 mL

symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Klíčova slova / Značenje simbola

	manufacturer / Hersteller / fabbricante / fabricante / fabricant / Tillverkaren / Fabrikanten / Produzent / Κατασκευαστής / Производитель / Производител / Угобвец / Proizvođač		expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration / utgångsdatum / udløbsdato / Utløpsdato / Ημερομηνία λήξης / срок на годност / datum expirace / срок годности / datum expirace / Rok trajanja
	storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage / lagringstemperatur / oppbevaringstemperatur / Oppbevaringstemperatur / θερμοκρασία αποθήκευσης / съхранение на / teplota skladování / температура хранения / teplota skladování / Temperatura lagerovnja		consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation / se användarinstruktioner / følg brugsvejledning / Følg bruksanvisningen / συμβουλευθείτε τις οδηγίες για τη χρήση / прочетете инструкцията за работа / potfeba fidiť se instrukcemi / перед использованием читайте инструкцию / sledujte návod k použití / Pročitaj upustvo pre upotrebe
			determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations / bestämmingar / bestemmelse / Bestemmelse / προσδιορισμοί / брой тестове / stanovení / определний / роčet stanovení / Definicija
AQUA	distilled water / destilliertes Wasser / acqua distillata / agua destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / destilirana voda / destilovaná voda / destilированная вода / destilovaná voda / Destilissana Voda	LOT	lot / Charge / lotto / lote / lote / lot / sats / serie / Parti / партиа / партида номер / šarže / lot / šarže / Serija
BUF	Reaction buffer / Reaktionspuffer / tampone di reazione / tampón de reacció / Tampão de reação / tampon de réaction / Reaktionsbuffert / Reaktionsbuffer / Reaksjonsbuffer / διάλυμα αντίδρασης / Реакционен буфер / Рабочий буферный раствор / Reakčni pufr / Reakcioni pufer	MTP	microtiter plate / Mikrotiterplatte / placa microtiter / microplaca / microplaca / microplaques sensibilisées / Mikrotiterplatta / Mikrotiterplade / mikrotiterplate / πλάκα μικροτιτλοδότησης / Микротитерная плака / Микропланшет / Mikrotitrační destička / Mikrotitrationsplöche
CAL	Calibrator / Kalibrator / Calibratore / calibrador / calibrador / calibreur / Kalibrator / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrator	REF	catalogue number / Katalognummer / numero di catalogo / número de catálogo / número de referència / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / καταλογен номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj
CONJ	Conjugate / Konjugat / Coniugato / conjugado / conjugado / conjugaté / Konjugerad / Konjugat / Konjugat / συνδεδετικό / Конюгат / Конъюгат / Konjugát / Konjugat	RTU	ready to use / gebrauchsfertig / pronto all'uso / listo para usar / pronto a usar / prêt à l'emploi / færdig att användas / færdig til brug / klar til bruk / έτοιμο προς χρήση / Готов за употреба / готов к использованию / k přímému použití / Razrediti ili rastvoriti
CONT	Control / Kontrolle / controllo / control / control / contrôle / Kontroll / Kontroll / Kontroll / διάλυμα ελέγχου / Контрол / Контрольный образец / Kontrola / Kontrola	STOP	stop solution / Stopplösung / Soluzione di arresto / solución de parada / solução de paragem / solution d'arrêt / Stopplösning / Stop-opløsning / Stoppløsning / διάλυμα παύσης / Стоп раствор / Стоп-раствор / Zastavovací roztok / Stop solucija
DIL	dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o dissolver / diluir ou dissolver em / diluer ou dissoudre dans / späd eller upplös i / fortyndes eller opløses i / Fortyndes eller opløses i / αραιωση ή διάλυση σε / растворите или разредете с / zředit anebo rozpustit v / разбавить или растворить в / naředte nebo rozpustite v / razrediti ili rastvoriti u	SUB	substrate / Substrat / substrato / substrato / substrato / substrat / Substrat / Substrat / Substrat / υπόστρωμα / Субстрат / Субстрат / Substrát / Substrat
INC	incubation buffer / Inkubationspuffer / tampone di incubazione / tampón de incubación / tampão de incubação / tampon d'incubation / Inkubationsbuffert / Inkubationsbuffer / Vaskebufferkoncentrat / διάλυμα επώασης / Инкубационен буфер / Буфер для инкубации / Inkubační pufr / Inkubacioni pufer	WASH	washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado / solución de lavado concentrada / tampão de lavagem concentrado / Tampon de lavage concentré / Vattenlösningkoncentrat / Vaskeoppløsningskoncentrat / vaskeløsningskoncentrat / συμπυκνωμένο διάλυμα πλύσης / Концентриран мицк разтвор / Концентрат промывочного раствора / Koncentrat promyvachho roztoku / Koncetrat solucije za ispiranje

PRODUCT DESCRIPTION

INTENDED USE

The TECHNOZYM® ADAMTS-13 INH ELISA is a test for detection of human autoantibodies (IgG) in serum or plasma against ADAMTS-13, the protease which is responsible for cleaving unusually large vWF (ULvWF). Most of these autoantibodies inhibit ADAMTS-13 activity and thus non-cleaved ULvWF accumulates in plasma. This is believed to be the major cause for thrombotic thrombocytopenic purpura (TTP). TECHNOZYM® ADAMTS-13 INH test makes it possible to differentiate between congenital (gene polymorphisms) and acquired (autoantibodies) TTP when coupled to an activity assay and to control efficacy of plasma exchange therapy.

COMPOSITION

1. ELISA test strips (6 for 48 T. and 12 for 96 T.), with 8 wells each, coated with a recombinant form of ADAMTS-13 protease; the drying agent is supplied in an aluminium bag.
2. Washing buffer concentrate (PBS; pH 7.3); containing detergent; 0.01% merthiolate; 1 vial, 80 mL.
3. Incubation buffer (= sample dilution buffer) (PBS; pH 7.3); contains stabiliser protein; 0.05% proclin; and dye, 1 vial, 90 mL, ready for use.
4. Calibrators (Standards) numbered from 1 to 5; lyophilised; 1 vial each; 0.5 mL. **Concentrations are lot-specific; consult the label on the vial.**
5. Positive and negative control plasma; lyophilised; 1 vial each; 0.5 mL. **Concentrations are lot-specific; consult the label on the vial.**
6. Conjugate: anti-human IgG POX; dyed blue; 1 vial, 0.3 mL.
7. Chromogenic substrate TMB (tetramethylbenzidine); 1 vial; 12 mL; ready to use.
8. Stopping solution sulphuric acid 0.45 mol/L; 1 vial; 12 mL; ready to use.
9. Adhesive film: for ELISA test strips; 2 pieces.

MATERIAL REQUIRED (not supplied with the kit)

1. Distilled water
2. Test tubes for diluting samples
3. Measuring cylinder (1000 mL)
4. Precision pipettes (10, 100 and 1000 µL)
5. Variable pipette (1000 µL)
6. Multichannel and/or dispensing pipettes (100 and 200 µL)
7. ELISA washer or multichannel pipette
8. ELISA reader with 450 nm filter, with a 620 nm reference filter if available.

STABILITY AND STORAGE

The expiry date printed on the labels applies to storage of the unopened vial at +2...+8 °C. Stability after reconstitution/opening:

Material/Reagent	State	Storage	Stability
Calibrators, control plasmas	after reconstitution	-20 °C	6 months
ELISA test strip	after opening	2 ... 8 °C with adhesive film in plastic bag with drying agent	expiry date
Washing buffer conc.	after opening	2 ... 8 °C	6 months
Washing buffer	1+11.5 dilution of concentrate	2 ... 8 °C	3 weeks
Incubation buffer (= sample dilution buffer)	after opening	2 ... 8 °C	2 months
Conjugate	after opening	2 ... 8 °C	6 months
	working solution	room temperature	60 minutes
Chromogen TMB	after opening	2 ... 8 °C	expiry date

TEST PROCEDURE

PREPARATION OF THE SAMPLES

Sample material: Human serum or citrated plasma. Samples maybe stored for three hours at room temperature. At -20°C they can be stored for several months. Samples may not be frozen and thawed several times.

PREPARATION OF REAGENT

1. Before starting the test, all the required components are to be brought to room temperature.
2. Preparing the washing buffer: Dilute 1 part by volume washing buffer concentrate with 11.5 parts by volume distilled water (1+11.5). Mix well! (Diluted washing buffer concentrate = washing buffer). There may be crystalline precipitations which will dissolve at 37°C within 10 minutes.
3. Reconstituting calibrators and control plasmas: Calibrators and control plasmas are reconstituted with 500 µL distilled water and mixed for 10 seconds after a reconstitution time of 15 minutes (vortex mixer).

No dilution is necessary for calibrators and controls!

Reconstituted components are clear to slightly turbid.

4. Sample Dilution: Samples are diluted with incubation buffer (1+100) : Dilute 1 part by volume sample with 100 parts by volume incubation buffer

10 µL sample + 1000 µL incubation (=sample dilution) buffer

5. Preparing the conjugate working solution (1+50): Dilute 1 part by volume conjugate with 50 parts by volume incubation buffer.

For 8 test wells: Mix 20 µL conjugate with 1000 µL incubation (=sample dilution) buffer

PERFORMANCE OF THE TEST

SAMPLE INCUBATION (reference 1,2)	Pipette calibrators, control plasmas, diluted samples into test wells; cover test strips with film	100 µL
	Incubate at room temperature	60 minutes
WASHING (reference 1,3,4)	Washing buffer	3 x 200 µL
CONJUGATE REACTION (reference 1,2)	Pipette conjugate working solution into wells, cover test strip with film	100 µL
	Incubate at room temperature	60 minutes
WASHING (reference 1,3,4)	Washing buffer	3 x 200 µL
SUBSTRATE REACTION (reference 1,2)	Pipette Substrate solution into test wells, cover test strip with film	100 µL
	Incubate at room temperature	10 minutes
STOPPING (reference 1,2)	Pipette stopping solution into wells	100 µL
MEASURING (reference 5, 6)	ELISA-Reader, 450 nm	

References

1. Reagents of different lots must not be combined
2. Precision and performance, among others, essentially depend on the following factors:
 - Thorough mixing of all substances used for dilution, 10 sec. with Vortex Mixer
 - Test calibrators, controls and samples in duplicates
 - Incubate at indicated temperature (RT: room temperature, 20...25°C)
 - Strict observance of the order of pipetting and of the time element as indicated
 - The time for sample incubation, conjugate and substrate reaction as indicated starts after pipetting the last sample. Incubation times should not vary by more than ± 5%.
 - During sample incubation and conjugate reaction, the time for pipetting calibrators/control plasmas/samples and/or conjugate solutions must not exceed 60 seconds per ELISA test strip (8 wells).
 - During substrate reaction and at stopping, the time needed for pipetting the substrate and/or the stopping solution must not exceed 10 seconds per ELISA test strip. Short pipetting times may be secured by using multichannel- and dispensing pipettes.
 - **Use incubation buffer from actual kit box, do not use incubation buffer left from previous boxes. Keep incubation buffer free from contaminants**
3. Label/number strips with a water resistant pen in case the strips accidentally fall out of the frame during testing.
4. After the last washing, wells must be aspirated thoroughly, turned upside down and positioned on a blotting paper; by gentle tapping, the last remnants must be removed.
5. By measuring the difference in wave lengths at 450 and 620 nm the precision of the test is increased.
6. shake 10 sec., measure within 10 min.

WARNING AND PRECAUTIONS

- All human blood or plasma products as well as samples must be considered as potentially infectious. They have to be handled with appropriate care and in strict observance of safety regulations. The rules pertaining to disposal are the same as applied to disposing hospital waste.
- Calibrators and control plasmas are made from human blood and any individual plasma involved in the procedure is HBsAg, HIV 1/2 Ab and HCV-Ab-negative (see labels on the vials). However, all human blood products should be handled as potentially infectious material.
- Stopping solution (sulphuric acid) may irritate the skin. Should acid get into your eyes, wash out immediately with water and consult a doctor.
- The reagents sometimes contain preserving agents (merthiolate). Beware of swallowing! Avoid contact with skin or mucous membranes

LIMITATION OF THE TEST

Samples with high concentrations of other than anti ADAMTS-13 autoantibodies may result in weak positive or borderline results.

ANALYSES RESULTS

CALCULATION OF THE RESULTS


Setting up a reference curve: X axis: concentration ADAMTS-13 IgG (U/mL)
Y axis: Extinction at 450nm

Graph plot is linear-linear with a linear or point to point fit.

Assessment of reference curve

- The extinction coefficient of the highest calibrator should be between 1.0 and 2.5.
- The extinction coefficient of the lowest Calibrator should be <0.1.
- The extinction difference between these two values should be at least 1.0.
- The validity of the test may be checked on the basis of the calculated control values.

Example of standard curve

Measuring concentration of samples

- Read off the concentration from the reference curve.
- If there are samples with extinction coefficients higher than that of the highest point on the curve, they have to be prediluted with incubation buffer (1+1, or 1+3). The measured concentration then has to be multiplied with the dilution factor 2 or 4, respectively.

INTERPRETATION OF RESULTS / REFERENCE RANGE

negative samples: < 12 units/mL

borderline: 12 – 15 units/mL

positive samples: > 15 units/mL

(n=193)

It is recommended that individual laboratories establish their own normal range. When interpreting the serological results the history of the patient has to be taken into account.

STANDARDISATION

Standards are calibrated against a plasma with a very high titre of anti ADAMTS-13 IgG. A 1:200 dilution of this reference plasma is defined to contain an antibody concentration of 100 Units/mL. (arbitrary units).

PERFORMANCE CHARACTERISTICS

Performance data are given below. Results obtained in individual laboratories may differ.

PRECISION

Reproducibility was determined with different samples (in series and day to day). The following results were obtained.

Sample	Intra assay variation		Inter assay variation	
	Sample 1	Sample 2	Sample 3	Sample 4
N	24	24	20	20
Mean (U/mL)	62,9	31,2	75,4	12,8
SD (U/mL)	4,9	0,9	2,7	0,7
CV (%)	7,86%	2,78%	3,54%	5,52%

ASSAY RANGE

2 – 104 U/mL

DETECTION LIMIT

1,68 U/mL

LITERATURE

- Sadler JE. A new name in thrombosis, ADAMTS13. Proc Natl Acad Sci U S A. 2002;99: 11552-11554.
- Tsai HM, Li A, Rock G. Inhibitors of von Willebrand factor-cleaving protease in thrombotic thrombocytopenic purpura. Clin Lab. 2001;47:387-392.
- Fontana S, Hovinga JA, Studd JD et al. Plasma therapy in thrombotic thrombocytopenic purpura: review of the literature and the Bern experience in a subgroup of patients with severe acquired ADAMTS-13 deficiency. Semin Hematol. 2004;41:48-59.
- Klaus C, Plaimauer B, Studd JD et al. Epitope mapping of ADAMTS13 autoantibodies in acquired thrombotic thrombocytopenic purpura. Blood. 2004;103:4514-4519.

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1,68 U/mL

LITERATUR

- Sadler JE. A new name in thrombosis, ADAMTS13. Proc Natl Acad Sci U S A. 2002;99: 11552-11554.
- Tsai HM, Li A, Rock G. Inhibitors of von Willebrand factor-cleaving protease in thrombotic thrombocytopenic purpura. Clin Lab. 2001;47:387-392.
- Fontana S, Hovinga JA, Studt JD et al. Plasma therapy in thrombotic thrombocytopenic purpura: review of the literature and the Bern experience in a subgroup of patients with severe acquired ADAMTS-13 deficiency. Semin Hematol. 2004;41:48-59.
- Klaus C, Plaimauer B, Studt JD et al. Epitope mapping of ADAMTS13 autoantibodies in acquired thrombotic thrombocytopenic purpura. Blood. 2004;103:4514-4519.

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2 – 104 U/mL

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1,68 U/mL

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- Sadler JE. A new name in thrombosis, ADAMTS13. Proc Natl Acad Sci U S A. 2002;99: 11552-11554.
- Tsai HM, Li A, Rock G. Inhibitors of von Willebrand factor-cleaving protease in thrombotic thrombocytopenic purpura. Clin Lab. 2001;47:387-392.
- Fontana S, Howings JA, Studd JD et al. Plasma therapy in thrombotic thrombocytopenic purpura: review of the literature and the Bern experience in a subgroup of patients with severe acquired ADAMTS-13 deficiency. Semin Hematol. 2004;41:48-59.
- Klaus C, Plaimauer B, Studd JD et al. Epitope mapping of ADAMTS13 autoantibodies in acquired thrombotic thrombocytopenic purpura. Blood. 2004;103:4514-4519.

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1,68 U/mL

LITTÉRATURE

- Sadler JE. A new name in thrombosis, ADAMTS13. Proc Natl Acad Sci U S A. 2002;99: 11552-11554.
- Tsai HM, Li A, Rock G. Inhibitors of von Willebrand factor-cleaving protease in thrombotic thrombocytopenic purpura. Clin Lab. 2001;47:387-392.
- Fontana S, Hovinga JA, Studt JD et al. Plasma therapy in thrombotic thrombocytopenic purpura: a review of the literature and the Bern experience in a subgroup of patients with severe acquired ADAMTS-13 deficiency. Semin Hematol. 2004;41:48-59.
- Klaus C, Plaimauer B, Studt JD et al. Epitope mapping of ADAMTS13 autoantibodies in acquired thrombotic thrombocytopenic purpura. Blood. 2004;103:4514-4519.