
# Protein C Deficient Plasma immunads.

FOR RESEARCH  
USE ONLY

REF 5346806 Protein C Deficient Plasma, immunads.

2 x 1 mL

**symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Ključova slova / Značenje simbola**

	manufacturer / Hersteller / fabbricante / fabricante / fabricante / fabricant / Tillverkaren / Fabrikanten / Producent / Κατασκευαστής / Производител / Производител / výrobce / Proizvođač		expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration / utgångsdatum / udløbsdato / Utløpsdato / Ημερομηνία λήξης / срок на годност / datum expirace/ срок годности / datum expirace / Rok trajanja
	storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage / lagringstemperatur / opbevaringstemperatur / Oppbevaringstemperatur / θερμοκρασία αποθήκευσης / съхранение на / toplota skladowání / температура хранения / toplota skladowání / Temperatura lagerovanja		consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation / se användarinstruktioner / følg bruksvejledning / Følg bruksanvisningen / συμβουλευθείτε τις οδηγίες για τη χρήση / прочетете инструкцията за работа / potěba řídit se instrukcemi / перед использованием читайте инструкцию / sledujte návod k použití / Pročitaj upustvo pre upotrebe
	CE-mark / CE-Kennzeichnung / marchio CE / marca de CE / Simbolo da CE / marquage CE / CE-märkning / CE-mærket / CE-merke / CE-σημάδι / CE марка / CE-označení / маркировка CE / značka CE / CE-marka		determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations / bestämmingar / bestemmelser / Bestemmelser / προσδιορισμοί / брой тестове / stanovení / определний / ročet stanovení / Definicija
<b>AQUA</b>	distilled water / destilliertes Wasser / acqua distillata / agua destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / дестилирана вода / destilovaná voda / дистиллированная вода / destilovaná voda / Serija	<b>LOT</b>	lot / Charge / lotto / lote / lote / lot / sats / serie / Parti / партиа / партида номер / šarže / lot / šarže / in vitro diagnostika
<b>BUF</b>	Reaction buffer / Reaktionspuffer / tampone di reazione / tampón de reacció / Tampão de reação / tampon de réaction / Reaktionspuffer / Reaktionsbuffer / Reaktionsbuffer / διάλυμα αντίδρασης / Реакционен буфер / Рабочий буферный раствор / Reakčni pufr / Reakcioni pufer	<b>MTP</b>	microtiter plate / Mikrotiterplatte / placa microtiter / microplaca / microplaca / microplaques sensibilisées / Mikrotiterplatta / Mikrotiterplade / mikrotiterplate / πλάκα μικροτιτλοδότησης / Микротитрѐна плака / Микропланшет / Mikrotitrační destička / Mikrotitracione ploče
<b>CAL</b>	Calibrator / Kalibrator / Calibratore / calibrador / calibrador / calibreur / Kalibrator / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrator	<b>REF</b>	catalogue number / Katalognummer / numero di catalogo / número de catálogo / número de referência / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / καταλογен номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj
<b>CONJ</b>	Conjugate / Konjugat / Coniugato / conjugado / conjugado / conjugaté / Konjugerad / Konjugat / Konjugat / συνδεδετικό / Конюгат / Конъюгат / Konjugát / Konjugat	<b>RTU</b>	ready to use / gebrauchsfertig / pronto all'uso / listo para usar / pronto a usar / prêt à l'emploi / færdig att användas / færdig til brug / klar til bruk / έτοιμο προς χρήση / Готов за употреба / готов к использованию / k přímému použití / Razrediti ili rastvoriti
<b>CONT</b>	Control / Kontrolle / controllo / control / control / contrôle / Kontroll / Kontroll / Kontroll / διάλυμα ελέγχου / Контрол / Контрольный образец / Kontrola / Kontrola	<b>STOP</b>	stop solution / Stopplösung / Soluzione di arresto / solución de parada / solução de paragem / solution d'arrêt / Stoppløsning / Stop-opløsning / Stoppløsning / διάλυμα παύσης / Стоп разтвор / Стоп-раствор / Zastavovací roztok / Stop solucija
<b>DIL</b>	dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o dissolver / diluir ou dissolver em / diluer ou dissoudre dans / späd eller upplös i / fortyndes eller opløses i / Fortyndes eller opløses i / αρωση ή διάλυση σε / растворите или разредете с / zředit nebo rozpustit v / разбавить или растворить в / naředte nebo rozpustte v / razrediti ili rastvoriti u	<b>SUB</b>	substrate / Substrat / substrato / substrato / substrato / substrat / Substrat / Substrat / Substrat / υπόστρωμα / Субстрат / Субстрат / Substrát / Substrat
<b>INC</b>	incubation buffer / Inkubationspuffer / tampone di incubazione / tampón de incubación / tampão de incubação / tampon d'incubation / Inkubationspuffer/ Inkubationsbuffer/ Vaskebufferkoncentrat / διάλυμα επώασης / Инкубационен буфер / Буфер для инкубации / Inkubační pufr / Inkubacioni pufer	<b>WASH</b>	washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado / solución de lavado concentrada / tampão de lavagem concentrado / Tampon de lavage concentré / Vattenlösningskoncentrat / Vaskeopløsningskoncentrat / vaskeløsningskoncentrat / συμπυκνωμένο διάλυμα πλύσης / Концентриран миеш разтвор / Концентрат промывочного раствора / Koncentrát promývacího roztoku / Koncetrat solucije za ispiranje
<b>RUO</b>	for research use only		

## PRODUCT DESCRIPTION

### INTENDED USE

For functional Protein C determination. Protein C is activated by the Protein C Activator Protac® (a component of the snake venom from Agkistrodon contortrix). Using Protein C Deficient Plasma the anticoagulant effect of Activated Protein C (APC) on the activated Partial Thromboplastin Time (aPTT) is determined.

### COMPOSITION

Protein C Deficient Plasma is a lyophilised, stabilised and immunadsorbed human plasma with a Protein C content ≤ 1%.

### MATERIAL REQUIRED (not supplied with the kit)

- Pipettes
- Distilled water
- Solutions/buffers:

REF 5410012 Imidazolpuffer 90 mL

REF 5277015 CaCl<sub>2</sub> 25 mmol/L solution 100 mL

- Reagents

aPTT Reagent

REF 5346212 Protac® 5 x 1 mL

- Control Plasmas and Calibrators

REF 5020040 Coagulation Control N 5 x 1 mL

REF 5021055 Coagulation Control A 5 x 1 mL

REF 5220110 Coagulation Reference 5 x 1 mL

\*\*\* or any other package sizes, special Ceveron® alpha or Technoclote Control and Calibration reagents of Technoclone.

### WARNING AND PRECAUTIONS

- RUO for research use only
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.
- This lot of reagents prepared from human blood and each single plasma used for this lot are HB<sub>s</sub>Ag, HIV 1/2 Ab and HCV Ab negative (see package label and vial label).

### STABILITY AND STORAGE

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8 °C.

Stability after reconstitution:

RT* (Ceveron® alpha)	-20 °C
2 hours	1 month

Upon storage, caps should be screwed tightly.

\*=room temperature

- When storing, the tubes must be securely capped.
- Deep frozen reagent must be thawed for at least 10 minutes at 37°C and mixed thoroughly before use
- Repeated freezing is not recommended

### TEST PROCEDURE

#### PREPARATION OF PLASMA SAMPLES

Mix 9 parts of venous blood and 1 part of Sodium Citrate solution (0.11 mol/L) and centrifuge for 15 min at a RCF of at least 2500 (corresponding to DIN 58905). Store the plasma at room temperature (up to 1 day).

All samples are diluted with Imidazole buffer.

The patients plasmas are determined after diluting 1:10. Other sample dilutions (1:3 to 1:20) can be used (see evaluation!).

Within the scope of quality assurance, Coagulation Control A is to be used in order to check correct operation. Coagulation Control N confirms precision.

### PREPARATION OF REAGENT

Reconstitute the lyophilised reagents in the prescribed quantities of distilled water and allow them to stand for 10 min at room temperature<sup>1</sup>. Prewarm the CaCl<sub>2</sub> 25 mmol/L solution to 37°C.

### PERFORMANCE OF THE TEST

#### CEVERON

Technoclone provides application sheets for Ceveron® alpha. The application sheets contain analyser/assay specific handling and performance information which may differ from that provided in this instruction for use. In this case the information contained in the application sheet supersedes the information in this instruction for use. Please consult the instruction manual of the Ceveron® alpha.

### MANUAL

Pipetting scheme:

0.1 mL	diluted plasma
+ 0.1 mL	Protein C deficient plasma, immune-adsorbed
+ 0.1 mL	Protac®
+ 0.1 mL	aPTT reagent
(incubate for at 37 °C in accordance with the aPTT reagent used)	
+ 0.1 mL	CaCl <sub>2</sub> 25 mmol/L solution 37°C
determine the point of coagulation	

### ANALYSES RESULTS

#### REFERENCE RANGE

70-130 % of normal

#### CALCULATION OF THE RESULTS

The sample dilutions (generally 1:10) are chosen to ensure that the values in seconds obtained can be read from the calibration curve. Plasma levels in the lower range of the calibration curve (≤ 20% protein C) and levels outside of the calibration curve ought to be measured once again with a different factor of dilution.

Samples diluted 1:10 can be read off directly from the calibration curve. When dilutions other than 1:10 have been prepared, the Protein C content of the sample is calculated from the value read from the calibration curve by means of the following formula:

$$\frac{\% \text{ Protein C}}{10} \times \text{actual dilution} = \% \text{ Protein C of the sample}$$

#### CALIBRATION CURVE

Reconstitute the Coagulation Reference according to the table accompanying the plasma.

Prepare from this solution a geometric series of dilutions (1:5 to 1:80).

Determine the coagulation time of the different dilutions and plot them in double-log paper (x-axis: activity / y-axis: seconds)

Example of a calibration curve

Dilution	1:5	1:10	1:20	1:40	1:80
% Protein C	200	100	50	25	12.5

#### QUALITY CONTROL

In order to verify the accuracy of the results, Coagulation Control A and N should always be tested in the same way as the patient's plasma.

#### LIMITATION OF THE TEST

Incorrect sample handling can lead to partial activation of the coagulation factors and to falsely elevated single factor determinations.

A new calibration is required for each batch of reagents where a calibration curve is necessary and for each instrument used. Also a new calibration is recommended, if software changes are introduced or following a major service of either instruments or equipment.

#### LITERATURE

Please contact Technoclone or your local distributor.

<sup>1</sup> For standardisation a reconstitution time of 30 min is recommended.

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