Control Plasma LMW Heparin - 82 3492 63

Intended use
For the quality control of chromogenic heparin assays.
Summary and principle
Heparin is the most frequently used antithrombotic therapeutic drug. The biological activity of this sulfated glycosaminoglycan resides in its ability to accelerate (up to 2000-fold) the inhibitory effect of antithrombin (AT) on the coagulation proteases. The Low and High Levels of Control Plasma LMW Heparin are prepared by means of a dedicated process and contain Low and High concentration of LMW Heparin. They are intended for the assessment of precision and accuracy of the assay at Heparin low and high levels.
Composition
The Control Plasma LMW Heparin kit consists of: Control Plasma LMW Heparin, Low: 4 vials of a lyophilized powder of human plasma containing low molecular weight (LMW) heparin Control Plasma LMW Heparin, High: 4 vials of a lyophilized powder of human plasma containing low molecular weight (LMW) heparin
PRECAUTIONS AND WARNINGS:
Each donor unit used in the preparation of human source reagent has been tested by FDA approved methods for the presence of Hepatitis B surface antigen and antibodies to HIV 1 and 2 and Hepatitis C and found to be negative. However, since no test can completely rule

out the presence of these blood borne diseases, the handling and disposal of human source reagents from this product should be made with care. Handle as potentially infectious. ¹ Hazard class: None Risk phrases: None Safety phrases: None This product is for <i>in vitro</i> diagnostic use
Preparation
Dissolve the contents of each vial with 1.0 mL of NCCLS Type II water or equivalent. ² Replace the stopper and swirl gently. Make sure of the complete reconstitution of the product. Keep the control at 15-25°C for 20 minutes and invert to mix before use. Do not shake. Avoid foam formation.
Reagent storage and stability
Unopened controls are stable until the expiration date shown on the vial when stored at 2-8°C. Stability after reconstitution: 24 hours at 20-25°C, or 14 days at 2-8°C in the original vial.
Limitations
These products are designed as controls for monitoring Heparin assays. These controls are subjected to the limitations of the assay system. Deviations may indicate possible problems with one or more components in the test system.

Performance characteristics
Refer to the Chromogenix Heparin reagent's inserts for performance characteristics.
Traceability of calibrators and control materials
Each lot of the Control Plasma LMW Heparin is calibrated against the International WHO Standard for LMW heparin
Expected values
The enclosed Low and High LMW Heparin Controls concentration ranges were determined over multiple runs on an ACL Family instrument using a specific lot of COAMATIC® Heparin reagent: these values must be considered only as a general indication once used with other Heparin kits. Each laboratory should establish its own mean and standard deviation and should establish quality control program to monitor laboratory testing. Controls should be analyzed at least once every 8 hours shift in accordance with good laboratory practice. Refer to Westgard <i>et al</i> for identification and resolution for out-of-control situations. ³

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- RICHARDSON J H and BACKLEY W E. Eds.Biosafety in Microbiological and Biomedical Laboratories. US. Dept. of Health and Human Services, Public Health Service, HHS Publication No (CDC) 84-8395, Washington, D.C. 1984.
- National Committee for Clinical Laboratory Standards. Specifications for reagent water used in the clinical laboratory, NCCLS Approved Standard: ASC-3.
- WESTGARD J O, BARRY P L. Cost-effective Quality Control: Managing the Quality and Productivity of Analytical Process, AACC Press 1986.

CONTROL PLASMA LMW 82 3492 63

Printed Insert Sheet: 302063
Revision: R5
Issued: October 2013
C.O.: 436617

LANGUAGES

ENGLISH
GERMAN
SPANISH
FRENCH
ITALIAN
PortuguêsE
DANSK
SVENSK
GREEK

TECHNICAL SPEC'S

PAPER: White paper, 50-60 g/m² weight.
SIZE: 11 x 17" (280 x 432 mm.).
PRINT: Front/Back.
PRINT COLOR: Front - Top rule Orange Pantone 137, all remaining type in black.
Back - All type in black.