ENGLISH

REAADS® Anti-Cardiolipin IgG/IgM Semi-Quantitative Test Kit

For *In Vitro* Diagnostic Use

An enzyme-linked immunosorbent assay (ELISA) for the semi-quantitative determination of IgG and IgM anti-cardiolipin antibodies in human serum or plasma.

INTENDED USE

For the detection and semi-quantitation of anti-cardiolipin antibodies in individuals with systemic lupus erythematosus (SLE) and lupus-like disorders (anti-phospholipid syndrome).

SUMMARY AND EXPLANATION OF THE TEST

Anti-phospholipid antibodies are autoantibodies that react with most negatively charged phospholipids, including cardiolipin (CL).^{1,2} Additionally, anti-phospholipid antibodies are known to prolong *in vitro* phospholipid-dependent coagulation tests and have been historically referred to as the “lupus anticoagulant”.^{1,3,4} Paradoxically, patients with the lupus anticoagulant do not present with abnormal bleeding except in the presence of other hemostatic abnormalities.³

Anti-cardiolipin (aCL) antibodies are frequently found in patients with systemic lupus erythematosus (SLE). They are also found in patients with other autoimmune diseases, as well as in some individuals with no apparent previous underlying disease.^{1,5,6} Elevated levels of aCL antibodies have been reported to be significantly associated with the presence of both venous and arterial thrombosis, thrombocytopenia, and recurrent fetal loss. The term “anti-phospholipid syndrome” (APS) has been introduced to describe patients who present these clinical manifestations, in association with aCL antibodies or the lupus anticoagulant.^{6,7,8,9,10}

The REAADS Anti-Cardiolipin Test Kit uses a well known ELISA format to detect aCL antibodies in human serum or plasma. Solid-phase immunoassays are generally considered more sensitive¹⁰ and more specific¹¹ for detecting aCL antibodies than coagulation assays. The REAADS Anti-Cardiolipin Test Kit provides rapid, highly reproducible, accurate, and objective results in units that are traceable to an internationally recognized reference preparation. The values for IgG aCL antibodies and IgM aCL antibodies are reported separately in GPL and MPL units, respectively.

PRINCIPLE OF THE TEST

The test is performed as an indirect ELISA. Diluted serum samples, calibrator sera, and controls are incubated in cardiolipin coated microwells, allowing aCL antibodies present in the samples to react with the immobilized antigen. After the removal of unbound serum proteins by washing, antibodies specific for human IgG or IgM labeled with horseradish peroxidase (HRP) are added forming complexes with the cardiolipin bound antibodies. Two enzyme-conjugated antibody solutions are provided, one specific for human IgG antibodies and one specific for human IgM antibodies. The concentration of IgG aCL antibodies and IgM aCL antibodies must be determined separately. Following another washing step, the bound enzyme-antibody conjugate is assayed by the addition of tetramethylbenzidine (TMB) and hydrogen peroxide (H₂O₂) as the chromogenic substrate. Color develops in the wells at an intensity proportional to the serum concentration of aCL antibodies.

Results are obtained by reading the O.D. (optical density or absorbance) of each well with a spectrophotometer. Calibrator sera are provided, with the IgG and IgM aCL concentrations expressed in GPL or MPL units, respectively. These units are traceable to the internationally recognized preparations of the Phospholipid Standardization Laboratory, University of Louisville. The user has the option of running either a single point calibrator or a four-point calibration curve. For single point calibration, dividing the concentration value of the calibrator sera by the O.D. value of the calibrator provides a conversion factor (one for IgG and one for IgM aCL). The O.D. values of all the other samples are multiplied by the conversion factors to obtain IgG and IgM aCL antibody concentrations in standard units. One GPL unit is equivalent to 1 µg/mL of an affinity purified standard IgG sample, and one MPL unit is equivalent to 1 µg/mL of an affinity purified standard IgM sample.¹² For multipoint calibration, perform a linear regression analysis with calibrator values against calibrator O.D.s. Control and patient results are determined from the calibration curve.

REAGENTS

Store at 2 - 8°C. Do Not Freeze.

Each REAADS Anti-Cardiolipin 96-microwell Test Kit contains the following reagents (volumes may vary depending on kit size and configuration):

- 96 stabilized beef heart cardiolipin (diphosphatidyl glycerol) coated microwells (12 strips of 8 breakaway wells), with frame.
- 1 bottle (60 mL) Sample Diluent* (green solution); contains bovine calf serum.
- 3 vials (0.250 mL) aCL IgG Calibrator Sera* (1-high, 2-moderate, 3-low)(human); see vial label for antibody concentration in GPL units. Calibrator 2 should be used when performing single point calibration.
- 3 vials (0.250 mL) aCL IgM Calibrator Sera* (1-high, 2-moderate, 3-low)(human); see vial label for antibody concentration in MPL units. Calibrator 2 should be used when performing single point calibration.
- 1 vial (0.250 mL) aCL IgG Positive Control Serum* (human); see vial label for expected GPL range.
- 1 vial (0.250 mL) aCL IgM Positive Control Serum* (human); see vial label for expected MPL range.
- 1 vial (0.250 mL) aCL Normal Control Serum* (human); see vial label for expected GPL and MPL ranges.
- 1 bottle (15 mL) anti-human IgG (goat) HRP-Conjugated Antibody Solution (blue solution).
- 1 bottle (15 mL) anti-human IgM (goat) HRP-Conjugated Antibody Solution (red solution).
- 1 bottle (15 mL) One Component Substrate (TMB and H₂O₂); ready to use.
- 1 bottle (15 mL) Stopping Solution (0.36 N sulfuric acid).
- 2 bottles (30 mL) Wash Concentrate (33X PBS).

***CAUTION: Contains sodium azide**

WARNINGS AND PRECAUTIONS

For *In Vitro* Diagnostic Use

1. Human source material used to prepare the calibrators and controls included in this kit has been tested and shown to be negative for antibodies to HBsAg, HCV, and HIV I & 2 by FDA required tests. However, all human blood derivatives, including patient samples, should be handled as potentially infectious material.
2. Do not pipette by mouth.
3. Do not smoke, eat, or drink in areas where specimens or kit reagents are handled.
4. Wear disposable gloves while handling kit reagents and wash hands thoroughly afterwards.
5. Certain components of this product contain sodium azide as a preservative. Sodium azide has been reported to form lead and copper azides when left in contact with these metals. These metal azides are explosive. Any solutions containing azide must be thoroughly flushed with copious amounts of water to prevent the build-up of explosive metal azides in the plumbing system.
6. Certain components are labeled with the following:
Irritating to eyes (R 36). Avoid contact with skin (S 24). Avoid contact with eyes (S 25). In case of contact with eyes, rinse immediately with plenty of water and seek medical advice (S 26). If swallowed, seek medical advice immediately and show this container or label (S 46).

Warning  . Biological Risk .

SPECIMEN COLLECTION AND PREPARATION

Serum is the preferred sample matrix. Blood should be collected by venipuncture, and the serum separated from the cells by centrifugation after clot formation. If not tested immediately, specimens should be stored at 2 to 8°C. If specimens are to be stored for more than 72 hours, they should be frozen at -20°C or below. Avoid repeated freezing and thawing. Do not use hemolyzed, icteric, or lipemic serum or plasma as these conditions may cause aberrant results. Specimens containing visible particulate matter should be clarified by centrifugation before testing.

Plasma collected with most anticoagulants except heparin may be used. Blood should be collected by venipuncture and the plasma immediately separated from the cells by centrifugation at 1500g for 10 minutes. The supernatant must be carefully removed after centrifugation to avoid contamination with platelets. Repeating the centrifugation and separation steps may be advisable to minimize platelet contamination. Lysed or aged platelets can react with anti-phospholipid antibodies leading to aberrant results. If not tested immediately, plasma samples should be stored as described for serum.

INSTRUCTIONS FOR USE

Materials Provided:

REAADS Anti-Cardiolipin Test Kit; see "Reagents" for a complete listing.

Materials Required but not Provided:

- Reagent grade water to prepare PBS wash solution and to zero or blank the plate reader during the final assay step
- Graduated cylinders
- Precision pipettors capable of delivering between 5 μ L and 1000 μ L, with appropriate tips
- Miscellaneous glassware appropriate for small volume handling
- Flask or bottle, 1 liter
- Wash bottles, preferably with the tip partially cut back to provide a wide stream, or an automated or semi-automated washing system
- Disposable gloves
- Plate reading spectrophotometer capable of reading absorbance at 450 nm (with a 650 nm reference if available)
- Multichannel pipettors capable of delivering to 8 wells simultaneously

Procedural Notes

1. Bring serum samples and kit reagents to room temperature (18-26°C) and mix well before using; avoid foaming. Return all unused samples and reagents to refrigerated storage as soon as possible.
2. All dilutions of calibrators, controls, and test sera must be made just prior to use in the assay.
3. A single water blank well can be set up on each plate with each run. No sample or kit reagents are to be added to this well. Instead, add 200 μ L of reagent grade water to the well immediately prior to reading the plate in the spectrophotometer. The plate reader should be programmed to zero or blank against an air or a water well.
4. Good washing technique is critical for optimal performance of the assay. Adequate washing is best accomplished by directing a forceful stream of wash solution into the bottom of the microwells from a plastic squeeze bottle with a wide tip. Wash solution in the water blank well will not interfere with the procedure. An automated microtiter plate washing system can also be used.
5. **IMPORTANT:** Failure to adequately remove residual PBS can cause inconsistent color development of the Substrate Solution.
6. Use a multichannel pipettor capable of delivering to 8 wells simultaneously when possible. This speeds the process and provides more uniform incubation and reaction times for all wells.
7. Careful controlled timing of all steps is critical. All calibrators, controls, and samples must be added within a five minute period. Batch size of samples should not be larger than the amount that can be added within this time period.
8. For all incubations, the start of the incubation period begins with the completion of reagent or sample addition.
9. Addition of all samples and reagents should be performed at the same rate and in the same sequence.
10. Incubation temperatures above or below normal room temperature (18 - 26°C) may contribute to inaccurate results.
11. Avoid contamination of reagents when opening and removing aliquots from the primary vials.
12. Do not use Tween 20 or other detergents in this assay.
13. Do not use kit components beyond expiration date.
14. Do not use kit components from different kit lot numbers.

REAGENT PREPARATION

Wash Solution (PBS): Measure 30 mL of Wash Concentrate (33X PBS) and dilute to 1 liter with reagent grade water. The pH of the final solution should be 7.35 ± 0.1 . Store unused PBS solution in the refrigerator at 2 - 8°C. Discard if the solution shows signs of microbial or cross-contamination.

Assay Procedure

1. The assay can be performed with a single point calibration (Calibrator 2) or a four-point calibration curve (Calibrators 1, 2, and 3 plus sample diluent/reagent blank as Calibrator 4 equal to 0 GPL or 0 MPL units). A reagent blank control should be run for each conjugate, IgG and IgM, with the single point and multipoint calibration method; Sample Diluent without serum is added to the well. This well will be treated the same as sample wells in subsequent assay steps.
2. Remove any microwell strips that will not be used from the frame and store them in the bag provided.
3. Prepare a 1:50 dilution of the calibrators, controls, and patient samples in Sample Diluent (green solution); e.g., 10 μ L sample added to 490 μ L Sample Diluent equals a 1:50 sample dilution.
4. Add 100 μ L of diluted calibrators (including the reagent blank/Calibrator 4), controls, and patient samples to the appropriate microwells.
5. Incubate 15 minutes at room temperature. After the incubation is complete, carefully invert the microwells and empty the sample fluid. Do not allow samples to contaminate other microwells.

6. Wash 4 times with PBS. Each well should be filled with PBS per wash. Invert microwells between each wash to empty fluid. Use a snapping motion of the wrist to shake the liquid from the wells. The frame must be squeezed at the center on the top and bottom to retain microwell modules during washing. Blot on absorbent paper to remove residual wash fluid. Do not allow wells to dry out between steps.
7. Add 100 µL anti-human IgG HRP-Conjugated Antibody Solution (blue) to the wells corresponding to the IgG calibrator, controls, reagent blank, and patient samples. Add 100 µL anti-human IgM HRP Conjugated Antibody Solution (red) to the wells corresponding to the IgM calibrator, controls, reagent blank, and patient samples.
8. Incubate for 15 minutes at room temperature. After the incubation is complete, carefully invert the microwells and empty the conjugate solutions. Take care not to allow cross-contamination of anti-IgG and anti-IgM Antibody Conjugate Solutions.
9. Wash 4 times with PBS as in step 6. Use a snapping motion to drain the liquid and blot on absorbent paper after the final wash. Do not allow the wells to dry out.
10. Add 100 µL One Component Substrate to each well and incubate for 10 minutes at room temperature. Add the substrate to the wells at a steady rate. Blue color will develop in wells with positive samples.
11. Add 100 µL Stopping Solution (0.36 N sulfuric acid) to each well to stop the enzyme reaction. Be sure to add Stopping Solution to the wells in the same order and at the same rate as the Substrate was added. Blue substrate will turn yellow and colorless substrate will remain colorless. Blank or zero the plate reader against an air or a water blank well. Read the O.D. of each well at 450 nm. The O.D. values should be measured within 5 minutes after the addition of Stopping Solution.

RESULTS

Single Point Calibration


1. Calculate the mean O.D. values if duplicates of Calibrator 2, controls, and patient samples were performed.
2. Divide the concentration value of Calibrator 2 (printed on the vial label) by the O.D. or mean O.D. value of the calibrator serum to obtain the conversion factor.
3. Multiply the O.D. or mean O.D. value for each of the controls and patient samples by the appropriate conversion factor to obtain a concentration value in GPL or MPL units.

<p>Conversion Factor = <u>Anti-cardiolipin concentration of Calibrator 2 (GPL or MPL)</u> Absorbance value of the Calibrator 2 (O.D. – G or M)</p> <p>Anti-cardiolipin concentration of sample = Conversion Factor X Absorbance of the sample (O.D.)</p>
--

4. The conversion factor must be calculated for both calibrators for each assay run. Using a conversion factor from another assay or interchanging GPL and MPL conversion factors will invalidate the results.

Multi-Point Curve Calibration

1. Calculate the mean O.D. values if duplicates of the calibrators, controls and patient samples were performed.
2. Perform linear regression or quadratic/2nd order polynomial regression analysis with the four calibrator values against the mean O.D.s for each calibrator. (See vial labels for GPL or MPL units. Calibrator 4 [sample diluent] is equal to 0 GPL or 0 MPL units)
3. The calibrator curve can be plotted either automatically using a validated software program or manually with graph paper. It is recommended to use a zero intercept when generating the regression line to avoid negative values. If this option is not available, any negative values should be reported as zero units. When generating the curve manually, draw a best fit line through the plotted points using a zero intercept.
4. Determine the control and patient sample values from the calibrator curve.
5. Example of a multi-point curve calibration.

Using the example calibration curve provided, a specimen O.D. of 1.000 at 450 nm would correspond to a calculated value of 26.2 units. The calibration curve provided is an example only and should not be used to calculate patient results. A new calibration curve should be performed with every test run.

QUALITY CONTROL

1. The O.D. value of Calibrator 2 should be at least 0.600 to assure that the kit is functioning properly. Calibrator 2 O.D. readings of less than 0.600 may indicate that the kit is no longer suitable for use.
2. The O.D. of Calibrator 4 or reagent blank should be less than 0.050 when the spectrophotometer has been blanked against an air or water well. Readings greater than 0.050 may indicate possible reagent contamination or inadequate plate washing.
3. The anti-cardiolipin values obtained for the control sera should be within the ranges indicated on the container labels. Occasional small deviations outside these ranges are acceptable.
4. O.D. values for duplicates (if performed) of the controls or patient samples should be within 20% of the mean O.D. value for samples with absorbance readings greater than 0.200.
5. Each laboratory should periodically determine its own normal cut-off values for the appropriate population of patients. See Performance Characteristics, Clinical Specificity, as an example.
6. Samples with anti-cardiolipin values greater than 100 GPL or 60 MPL may be reported as “greater than 100 GPL or 60 MPL”.
7. Assure that all quality control parameters have been met (see Quality Control) before reporting test results.

EXPECTED VALUES

Normal Range:

Serum samples from 94 healthy blood donors were tested and the following normal range established (mean + 2 SD):

- Less than 23 GPL
- Less than 11 MPL

EXPECTED PREVALENCE

SLE:

Serum samples from 149 individuals with SLE were tested with the kit. Thirty-one of the samples (21%) were positive for IgG anti-cardiolipin antibodies. Fifteen of the samples (10%) were positive for IgM anti-cardiolipin antibodies. No correlation was found between anti-cardiolipin antibody levels and anti-dsDNA antibody levels or disease activity. Monoclonal antibodies specific for dsDNA and ssDNA have also been tested in the assay and were shown to be nonreactive with the cardiolipin coated microwells.

Other Disease States:

Twenty-four serum samples from patients with osteoarthritis (OA) were tested in the assay. Only one was positive for IgG aCL antibodies. None were positive for IgM aCL antibodies.

Ninety-three serum samples from patients with progressive systemic sclerosis (PSS) were tested in the assay. Twelve of the samples (13%) were positive for IgG aCL, and three of the samples (3%) were positive for IgM aCL.

The clinical significance of positive results in disease states other than SLE is still under investigation.

PERFORMANCE CHARACTERISTICS

Clinical Specificity

Normal Samples:

Serum samples from 94 healthy blood donors were assayed for the presence of IgG and IgM aCL antibodies. The normal cut-off was calculated as the mean GPL or MPL antibody level plus two standard deviations. The cut-off value was calculated separately for IgG and IgM antibodies. The calculated values are 23 GPL for IgG and 11 MPL for IgM. Using these cut-off values, the assay is 97% specific for IgG antibodies and 96% specific for IgM antibodies.

Disease Controls:

Serum samples from 14 patients with SLE or a lupus-like disorder with no history of thrombosis or other features of the antiphospholipid syndrome, were tested in the assay. None of the samples were positive for IgG aCL antibodies. One sample was positive for IgM aCL antibodies. The overall specificity of the assay for this sample population was 93%.

Clinical Sensitivity Serum samples from 18 patients with SLE or a lupus-like disorder, known to have had at least one thrombotic event, were tested in the assay. Eleven samples were positive for a sensitivity of 61% in this sample population. Seven of the samples had elevated levels of IgG aCL antibodies only, three samples had elevated levels of both IgG and IgM aCL antibodies, and one sample was positive for IgM aCL antibodies only.

Precision Two samples with known GPL values (one low and one high) and two samples with known MPL values (one low and one high) were assayed in 28 replicates on three different occasions. The intraassay and interassay coefficients of variation (CVs) are presented in the following table. The reported intraassay coefficient of variation is the mean of the three separate intraassay CVs.

Sample (aCL conc.)	Mean Intraassay CV	Interassay CV
Low (20.1 GPL)	7.7%	11.9%
Low (16.4 MPL)	10.5%	11.6%
High (48.8 GPL)	9.4%	13.2%
High (64.6 MPL)	7.4%	10.9%

Recovery Two samples with known GPL concentration values (one low, LS, and one high, HS) were diluted 1:1 with standards containing known amounts of IgG anti-cardiolipin antibodies. The calculated, observed, and recovery values are presented in the following table. Anti-cardiolipin concentration values are reported in GPL units.

Sample	<u>GPL Concentration</u>		Recovery Value (%)
	Calculated	Observed	
LS	--	11.0	--
LS + 6.3	8.7	8.8	101%
LS + 12.5	11.8	10.6	89.9%
LS + 25	18.0	17.0	94.2%
LS + 50	30.5	25.6	83.8%
HS	--	35.9	--
HS + 6.3	21.1	24.7	115%
HS + 12.5	24.2	27.7	115%
HS + 2.5	30.5	31.8	104%
HS + 50	43.0	39.8	92.6%

Two samples with known MPL concentration values (one low, LS, and one high, HS) were diluted 1:1 with standards containing known amounts of IgM anti-cardiolipin antibodies. The calculated, observed, and recovery values are presented in the table below. Anti-cardiolipin concentration values are reported in MPL units.

Sample	MPL Concentration		Recovery Value (%)
	Calculated	Observed	
LS	--	7.1	--
LS + 6.9	7.0	8.6	122%
LS + 13.8	10.5	12.4	118%
LS + 27.5	17.3	19.0	110%
LS + 55	31.1	30.6	98.4%
HS	--	47.1	--
HS + 6.9	27.0	35.2	130%
HS + 13.8	30.5	37.6	123%
HS + 27.5	37.3	39.4	106%
HS + 55	51.1	52.4	102%

LIMITATIONS OF THE TEST

The anti-cardiolipin antibody concentration values obtained from this assay are an aid to diagnosis only. Each physician must interpret these results in light of the patient's history, physical findings, and other diagnostic procedures. If clinical findings suggest the presence of anti-phospholipid antibodies and the patient is negative for anti-cardiolipin antibodies, some investigators recommend testing for the lupus anticoagulant to confirm the negative result. A patient is considered positive for anti-phospholipid antibodies if one or both of the tests give positive results.¹³

Patients with current or prior syphilis infections may have a positive result in the REAADS aCL assay without increased risk of thrombosis. In a study of 23 patients with confirmed syphilis (FTA-Abs positive), nine of the samples (39%) were positive in the REAADS aCL assay. If a patient's history is consistent with a diagnosis of syphilis infection, this should be confirmed or ruled out by an assay specific for anti-treponemal antibodies. Anti-cardiolipin antibodies can appear transiently at low levels during many infections. If a patient first tests positive while there are clinical signs of infection, the test should be repeated after an interval of six months.¹⁴

WARRANTY

This product is warranted to perform as described in this package insert. Corgenix, Inc. disclaims any implied warranty of merchantability or fitness for a particular use, and in no event shall Corgenix, Inc. be liable for consequential damage.

For Technical or Customer Service in the United States, phone 1-800-729-5661. Outside the United States, phone (303) 457-4345, fax (303) 457-4519, or contact a Corgenix authorized distributor.

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REFERENCES

1. McNiel HP, Chesterman CN, Krilis SA. Immunology and clinical importance of antiphospholipid antibodies. *Adv Immunol* 1991; 49:193-280.
2. Harris EN. Annotation. Antiphospholipid antibodies. *Br J Haematol* 1990; 74:1-9.
3. Love PE, Santoro SA. Antiphospholipid antibodies: Anticardiolipin and the lupus anticoagulant in systemic lupus erythematosus (SLE) and in non-SLE disorders. Prevalence and clinical significance. *Ann Intern Med* 1990; 112:682-698.
4. Sammaritano LR, Gharavi AE, Lockshin MD. Antiphospholipid antibody syndrome: immunologic and clinical aspects. *Seminars in Arthritis and Rheumatism* 1990; 20:81-96.
5. Hughes GRV, Harris EN, Gharavi AE. The anticardiolipin syndrome (Editorial). *J Rheumatol* 1986; 13:486-489.
6. Hughes GRV. The antiphospholipid syndrome: ten years on. *Lancet* 1993; 342:341-344.
7. Wilson, WA, Gharavi AE, Koike T, et al. International consensus statement on preliminary classification criteria for definite antiphospholipid syndrome. *Arthritis Rheum* 1999;42:1309-1311.
8. Loizou S, McCrean JD, Rudge AC, et al. Measurement of anti-cardiolipin antibodies by an enzyme linked immunosorbent assay (ELISA): standardization and quantitation of results. *Clin Exp Immunol* 1985; 62:738-745.
9. Triplett DA. Assays for detection of antiphospholipid antibodies. *Lupus* 1994; 3:281-287.
10. Schick PK, Kurica KB, Chacko GK. Location of phosphatidylethanolamine and phosphatidylserine in the human platelet plasma membrane. *J Clin Invest* 1976; 57:1221-1226.
11. Branch DW, Rote NS, Dostal DA, Scott JR. Association of lupus anticoagulant with antibody against phosphatidylserine. *Clin Immunol Immunopathol* 1987; 42:63-75.
12. Wilson WA, Faghiri Z, Taheri F, Gharavi AE. Significance of IgA antiphospholipid antibodies. *Lupus* 1998; 7:Suppl 2, S110-S113.
13. Bick RL, Baker WF. Antiphospholipid syndrome and thrombosis. *Seminars in Thrombosis and Hemostasis* 1999; 25: 333-350.
14. Barna LK, Triplett DA, Foster WM, Gaddis ML. A study of relationships among anti-phosphatidylserine and anticardiolipin antibodies, the lupus anticoagulant and clinical complications. *Lupus* 1994; 3:357.
15. Keedy K, Santos M, Lopez L. Prevalence of antibodies to anti-phosphatidylserine in autoimmune diseases. *Lupus* 1994; 3:344.
16. Keedy K, Santos M, Lopez L. Prevalence and clinical significance of anti-phosphatidylserine antibodies. *Arthritis Rheum* 1994; 37:S392.
17. Jansen NL, Snell JR, Moy JN. Myocardial infarction as the presenting manifestation of systemic lupus erythematosus with antiphosphatidylserine antibodies. *Annals of Allergy, Asthma, and Immunology* 1996; 76:266-268.
18. Blank M, Tincani A, Shoenfeld Y. Induction of experimental antiphospholipid syndrome in naive mice with purified IgG anti-phosphatidylserine antibodies. *J Rheumatol* 1994; 21:100-104.
19. McNiel H, Simpson R, Chesterman C, Krilis S. Antiphospholipid antibodies are directed against a complex antigen that includes a lipid-binding inhibitor of coagulation: β_2 -glycoprotein I (apolipoprotein H). *Proc Natl Acad Sci USA* 1990; 87:4120.
20. Kandiah D, Krilis S. Beta₂-glycoprotein I. *Lupus* 1994; 3:207-212.
21. Keedy K, Santos M, Lopez L. Anti-phosphatidylserine antibodies require β_2 -glycoprotein I as cofactor in ELISA. *Lupus* 1994; 3:327.
22. Vaarala O, Palosuo T, Kleemola M, Aho K. Anticardiolipin response in acute infections. *Clinical Immunology and Immunopathology* 1986; 41: 8-15.
23. Maclean C, Flegg PJ, Kilpatrick DC. Anti-cardiolipin antibodies and HIV infection. *Clinical Experimental Immunology*. 1990; 81: 263-266.

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Corgenix, Inc.
 11575 Main Street, Suite 400
 Broomfield, Colorado 80020, USA

READS® is a registered trademark of Corgenix, Inc.

© 2015, Corgenix, Inc.

THIS PAGE INTENTIONALLY LEFT BLANK

THIS PAGE INTENTIONALLY LEFT BLANK

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MT Promedt Consulting GmbH