



## KALLIKREIN INHIBITOR ACTIVITY

Determination of kallikrein inhibitor activity in plasma, with S-2302.

### Measurement Principle

Plasma is incubated with a purified plasma kallikrein preparation. The amount of kallikrein inhibited is proportional to the activity of the kallikrein inhibitor present in the plasma. The remaining amount of kallikrein activity is then determined by using the substrate H-D-Pro-Phe-Arg-pNA (S-2302). The rate at which p-nitroaniline (pNA) is released is measured photometrically at 405 nm. This can be followed on a recorder (initial rate method) or read after stopping the reaction with acetic acid (acid stopped method).



### Reagents

1. S-2302, 25 mg Art. No. S820340  
Reconstitute the substrate S-2302 (MW: 611.6) with 20 ml of distilled water.

### 2. Plasma Kallikrein

Use purified human plasma kallikrein (refer to Gallimore MJ et al., 1978). Prepare a solution of 1 nkat<sub>S-2302</sub>/ml human plasma kallikrein in Tris buffer pH 7.8.

1 nkat<sub>S-2302</sub> corresponds to 0.06 U or 0.017 PEU (refer to Friberger P et al. 1979).

### 3. Tris Buffer, pH 7.8 (25°C)

Tris	6.1 g	(50 mmol/l)
NaCl	21.1 g	(361 mmol/l)
Polybrene	20 mg	
Distilled water	800 ml	

Adjust the pH to 7.8 at 25°C by adding an appropriate amount (approx. 38 ml) of 1 mol/l HCl. Fill up to 1000 ml with distilled water. The buffer, if not contaminated, will remain stable for six months at 2 to 8°C.

### 4. Normal plasma

Blood samples are taken from at least 10 healthy donors. For the preparation of the samples, refer to the Specimen collection section.

### 5. Acetic acid, 20%

Acetic acid is used in the acid stopped method.

### Equipment

1. Spectro- or filter photometer, 405 nm
2. Siliconised semi-microcuvettes, 1 cm
3. Centrifuge
4. Thermostat, 37°C
5. Stop watch
6. Disposable plastic tubes
- Additional equipment for the initial rate method
7. Photometer with cuvette housing, thermostated at 37°C.

### Specimen collection

Blood (9 vol) is mixed with 0.1 mol/l sodium citrate (1 vol) and centrifuged at 2000 x g for 20 minutes at 15-25°C. In order to avoid low-temperature activation of plasma kallikrein inhibitor the plasma should be kept at 15-25°C for not more than a few hours or immediately frozen at -20°C or below. After thawing at 37°C the plasma should be kept at 15-25°C and used as soon as possible. Frozen plasma may lose some plasma kallikrein inhibitor activity on freezing or thawing, but is stable for several months at -20°C or below.

### Standard curve

Normal plasma has a kallikrein inhibitor concentration of 100% and is diluted according to the table below (see Note 1).

K-Inhibitor %	Normal plasma $\mu$ l	Buffer $\mu$ l
25	100	300
50	200	200
75	300	100
100	400	-

### Method

Sample dilution	Tube No. 1
Buffer	1900 $\mu$ l
Test plasma or standard (see note 1)	100 $\mu$ l
Mix	

Initial rate method	Tube No. 2
Sample from tube No. 1	200 $\mu$ l
Incubate at 37°C	3-4 min
Plasma kallikrein	200 $\mu$ l
Mix and incubate at 37°C	5 min
Substrate (37°C)	200 $\mu$ l
Mix	

Transfer sample immediately to a 1 cm siliconised semi-microcuvette (preheated to 37°C) for measurement of the absorbance change in a photometer at 405 nm and at 37°C. Calculate  $\Delta A/\text{min}$ .

Acid stopped method	Tube No. 2
Sample from tube No. 1	200 $\mu$ l
Incubate at 37°C	3-4 min
Plasma kallikrein	200 $\mu$ l
Mix and incubate at 37°C	5 min
Substrate (37°C)	200 $\mu$ l
Mix and incubate at 37°C	4 min
Acetic acid 20%	200 $\mu$ l
Mix	



Plasma blanks are prepared by adding the reagents in reverse order without incubation. Read the absorbance (A) of the sample against its blank in a photometer at 405 nm. The colour is stable for at least 4 hours.

### Calculation

#### Plasma Kallikrein inhibitor in percentage of normal plasma

Plot A or  $\Delta A/\text{min}$  for the standards against their concentration of kallikrein inhibitor on log-lin graph paper.

Read the kallikrein inhibitor value for the corresponding A or  $\Delta A/\text{min}$  of the unknown test sample from the standard curve.

#### Plasma Kallikrein inhibitor in enzyme activity units

In each test series a kallikrein activity determination with buffer instead of sample dilution must be performed.

The difference between this activity and the sample activity is then calculated.

Initial rate method:

$$\mu\text{kat/l} = (\Delta A/\text{min buffer} - \Delta A/\text{min sample}) \times 104$$

$$\text{U/l} = (\Delta A/\text{min buffer} - \Delta A/\text{min sample}) \times 6\,250$$

Acid stopped method:

$$\mu\text{kat/l} = (A_{\text{buffer}} - A_{\text{sample}}) \times 34.7$$

$$\text{U/l} = (A_{\text{buffer}} - A_{\text{sample}}) \times 2\,080$$

### Notes

1. A 150% standard is prepared by diluting 300  $\mu\text{l}$  normal plasma with 3700  $\mu\text{l}$  buffer. A 200% standard is prepared by diluting 100  $\mu\text{l}$  normal plasma with 900  $\mu\text{l}$  buffer. For 0% use the buffer only (note that the adsorbance to surfaces can result in lower readings when plasma is absent).

2. It is suggested, that the spontaneous kallikrein activity ( $\alpha$ 2-M complex) should be determined in patients in whom the kallikrein system is suspected to be activated.

See Determination of Kallikrein-like Activity in Plasma.

### Bibliography

Gallimore MJ et al. The purification of human plasma kallikrein with weak plasminogen activator activity. *Thromb Res* 2, 409-420 (1978).

Friberger P et al. Determination of prekallikrein in plasma by means of a chromogenic tripeptide substrate for plasma kallikrein. In: *KININS II. Biochemistry, Pathophysiology and Clinical aspects*, Eds Fujii S et al. Plenum Publishing Corp 67-82 (1979).

Gallimore MJ et al: Studies on plasma inhibitors of plasma kallikrein using chromogenic peptide substrate assay. *Thromb Res* 16, 695-703 (1979).

Aasen AO et al. Studies on components of the plasma kallikrein-kinin system in plasma samples from normal individuals and patients with septic shock. *Advances in Shock Research* 4. Alan R. Liss Inc., NY. 1-10 (1980).

Walker ID et al. The coagulation, fibrinolytic and plasma kallikrein systems in acute pancreatitis. In *Progress in Chemical Fibrinolysis and Thrombolysis*. Vol. V. Ed Davidson J F, Nilsson I M & Åstedt B. 291-293 (1980).

Gallimore MJ et al. Further studies on components of the plasma kallikrein system in plasma samples from cancer patients and normal individuals. In: *Progress in Chemical Fibrinolysis and Thrombolysis*. Vol. V. Ed Davidson J F, Nilsson I M & Åstedt B p 256-258 (1980).

Gallimore MJ et al. Activation of the coagulation, fibrinolytic and kallikrein-kinin systems during cardiopulmonary bypass applying extreme hemodilution. *Vth Int Conf on Synthetic Fibrinolytic and Thrombolytic Agents. Progress in Fibrinolysis*. June 17-20, Malmö 1980. Abstract 51.

Gallimore MJ & Friberger P: Simple chromogenic peptide substrate assays for determining prekallikrein, kallikrein inhibition and kallikrein "like" activity in human plasma. *Thromb Res* 25, 293-298, (1982).

Baskova IP et al. Inhibition of plasma kallikrein Kininase and kinin-like activities of preparation from the medicinal leeches. *Thromb Res* 67, 721-730 (1992).