


CHYMOTRYPSIN

Determination of chymotrypsin with S-2586.

Measurement Principle

The chymotrypsin activity is determined by its amidolytic effect on the substrate MeO-Suc-Arg-Pro-Tyr-pNA (S-2586). The rate at which p-nitroaniline (pNA) is released is measured photometrically at 405 nm. This can be followed on a recorder (initial rate method) or read after stopping the reaction with acetic acid (acid stopped method). The correlation between the change in absorbance per minute ($\Delta A/\text{min}$) or absorbance (A) and the chymotrypsin activity is linear in the 0.05-1.0 $\mu\text{kat/l}$ or 3-60 U/l range. The amidolytic activity of different chymotrypsin preparations does not necessarily parallel the protease activity.



Reagents

- S-2586, 25 mg** Art. No. S820894
Reconstitute the substrate S-2586 (MW: 705.3) with 60 ml of distilled water.
- Tris/Calcium Buffer, pH 8.3 (25°C)**

Tris	12.1 g	(100 mmol/l)
NaCl	56.2 g	(960 mmol/l)
Distilled water	800 ml	

Adjust the pH to 8.3 at 25°C by adding approximately 50 ml of 1 mol/l HCl. Add 10 ml of 1 mol/l CaCl_2 solution. Fill up to 1000 ml with distilled water. The buffer, if not contaminated, will remain stable for two months at 2-8°C.

- Acetic acid, 20%**
Acetic acid is used in the acid-stopped method.

Equipment

- Spectro- or filter photometer, 405 nm
- Semi-microcuvettes, 1 cm.
- Thermostat, 37°C
- Stop watch
- Disposable plastic tubes
- Photometer with cuvette housing, thermostated at 37°C (for the initial rate method)

Sample

The sample containing chymotrypsin is dissolved in or diluted with 1 mmol/l HCl to a concentration of 0.1 g/l. This stock solution is stable for more than two weeks at 2-8°C. Before assay, the solution is diluted 1:200 with 1 mmol/l HCl. If the sample is a pure protein, it is advisable to use 0.1% Carbowax 6000 (Union Carbide, NY) or 1% albumin (previously checked for amidolytic activity) to avoid adsorption to surfaces.

Method

Initial rate method	
Buffer	200 μl
Incubate at 37°C	3-4 min
Chymotrypsin sample	200 μl
Mix and incubate at 37°C	2-3 min
Substrate (37°C)	200 μl
Mix	

Transfer the sample immediately to a 1 cm semi-microcuvette (preheated to 37°C) for measurement of the absorbance change in a photometer at 405 nm and at 37°C. Calculate $\Delta A/\text{min}$.

Acid stopped method	Sample	Blank
Buffer	200 μl	200 μl
Incubate at 37°C	3-4 min	-
Chymotrypsin sample	200 μl	200 μl
Mix and incubate at 37°C	2-3 min	-
Substrate (37°C)	200 μl	-
Mix and incubate at 37°C	3 min	-
Acetic acid 20%	200 μl	200 μl
Mix	yes	-
Substrate (37°C)	-	200 μl
Mix	-	yes

Read the absorbance (A) of the sample against a water or sample blank in a photometer at 405 nm. The colour is stable for at least 4 hours.

Calculation

Calculate the chymotrypsin activity of the stock solution from the following formulas:

$$\begin{aligned} \text{Initial rate method: } \mu\text{kat/l} &= 5.19 \times \Delta A/\text{min} \times 200 \\ &\text{U/l} = 311 \times \Delta A/\text{min} \times 200 \\ \text{Acid stopped method: } \mu\text{kat/l} &= 2.31 \times A \times 200 \\ &\text{U/l} = 138 \times A \times 200 \end{aligned}$$

Bibliography

- Tözsér J et al. Active centre studies on bovine pancreatic chymotrypsin with tripeptidyl-p-nitroanilide substrates. Acta Biochim Biophys Hung 21, 335-348 (1986).
- Takayama TK et al. Characterization of the precursor of prostate-specific antigen. Activation by trypsin and by human glandular kallikrein. J Biol Chem 272, 21582-8 (1997).
- Pejler G, Sadler JE. Mechanism by which heparin proteoglycan modulates mast cell chyma activity. Biochemistry 38, 12187-95 (1999).