
Coagulation Control AK f. Ceveron® alpha

For Research Use Only

REF 5011030 Coagulation Control AK f. Ceveron® alpha 5 x 1 mL

REF 5011035 Coagulation Control AK f. Ceveron® alpha 50 x 1 mL

symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Ključova slova / Značenje simbola		
	manufacturer / Hersteller / fabbricante / fabricante / fabricante / fabricant / Tillverkaren / Fabrikanten / Produzent / Κατασκευαστής / Производитель / Производител / výrobce / Proizvođač	 expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration / utgångsdatum / udløbsdato / Utløpsdato / Ημερομηνία λήξης / срок на годност / datum expirace / срок годности / datum expirace / Rok trajanja
	storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage / lagringstemperatur / opbevaringstemperatur / Oppbevaringstemperatur / Θερμοκρασία αποθήκευσης / съхранение на / teplota skladování / температура хранения / teplota skladování / Temperatura lagerovanja	 consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation / se användarinstruktioner / følg brugsvejledning / Følg bruksanvisningen / συμβουλευθείτε τις οδηγίες για τη χρήση / прочетете инструкцията за работа / potfeba fídit se instrukcemi / перед использованием читайте инструкцию / sledujte návod k použití / Pročítaj upustvo pre upotrebe
		 determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations / bestämmingar / bestemmelse / Bestemmelse / προσδιορισμοί / брой тестове / stanovení / определний / роčet stanovení / Definicija
AQUA	distilled water / destilliertes Wasser / acqua distillata / agua destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / destilirana voda / destilovaná voda / дистиллированная вода / destilovaná voda / Serija	LOT lot / Charge / lotto / lote / lote / lot / sats / serie / Parti / παρτίδα / партида номер / šarže / лот / šarže / in vitro dijagnostika
BUF	Reaction buffer / Reaktionspuffer / tampone di reazione / tampón de reacció / Tampão de reação / tampon de réaction / Reaktionspuffer / Reaktionsbuffer / Reaktionsbuffer / διάλυμα αντίδρασης / Реакционен буфер / Рабочий буферный раствор / Reakční pufr / Reakcioni pufer	MTP microtiter plate / Mikrotiterplatte / placa microtiter / microplaca / microplaca / microplaques sensibilisées / Mikrotiterplatta / Mikrotiterplade / mikrotiterplate / πλάκα μικροτιτλοδότησης / Микротитрѐна плака / Микропланшет / Mikrotitrační destička / Mikrotitracione ploče
CAL	Calibrator / Kalibrator / Calibratore / calibrador / calibrador / calibreur / Kalibrator / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrator	REF catalogue number / Katalognummer / numero di catalogo / número de catálogo / número de referència / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / καταλογен номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj
CONJ	Conjugate / Konjugat / Coniugato / conjugado / conjugado / conjugaté / Konjugerad / Konjugat / Konjugat / συνδεδετικό / Конюгат / Конъюгат / Konjugát / Konjugat	RTU ready to use / gebrauchsfertig / pronto all'uso / listo para usar / pronto a usar / prêt à l'emploi / færdig att användas / færdig til brug / klar til bruk / έτοιμο προς χρήση / Готов за употреба / готов к использованию / k přímému použití / Razrediti ili rastvoriti
CONT	Control / Kontrolle / controllo / control / control / contrôle / Kontroll / Kontroll / Kontroll / διάλυμα ελέγχου / Контрол / Контрольный образец / Kontrola / Kontrola	STOP stop solution / Stoppløsning / Soluzione di arresto / solución de parada / solução de paragem / solution d'arrêt / Stoppløsning / Stop-opløsning / Stoppløsning / διάλυμα παύσης / Стоп разтвор / Стоп-раствор / Zastavovací roztok / Stop solucija
DIL	dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o dissolver / diluir ou dissolver em / diluer ou dissoudre dans / späd eller upplös i / fortyndes eller opløses i / Fortyndes eller opløses i / αραιωση ή διάλυση σε / разворете или разредете с / zředit anebo rozpustit v / разбавить или растворить в / nafedte nebo rozpustte v / razrediti ili rastvoriti u	SUB substrate / Substrat / substrato / substrato / substrato / substrat / Substrat / Substrat / Substrat / υπόστρωμα / Субстрат / Субстрат / Substrát / Substrat
INC	incubation buffer / Inkubationspuffer / tampone di incubazione / tampón de incubación / tampão de incubação / tampon d'incubation / Inkubationspuffer / Inkubationsbuffer / Vaskebufferkonsentrat / διάλυμα επώασης / Инкубационен буфер / Буфер для инкубации / Inkubacijski pufr / Inkubacioni pufer	WASH washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado / solución de lavado concentrada / tampão de lavagem concentrado / Tampon de lavage concentré / Vattenlösningsskoncentrat / Vaskeopløsningskoncentrat / vaskeløsningskoncentrat / συμπυκνωμένο διάλυμα πλύσης / Концентриран миеш разтвор / Концентрат промывочного раствора / Koncentrát promývacího roztoku / Koncetrat solucije za ispiranje
RUO	For Research Use Only	

PRODUCT DESCRIPTION

INTENDED USE

Coagulation Control AK can be used for

- Accuracy control of prothrombin time determination in the therapeutic range of oral anticoagulant therapy.
 - in INR
 - in % of normal
- Precision control of prothrombin time in the therapeutic range of oral anticoagulant therapy.
 - in INR
 - in % of normal

COMPOSITION

Coagulation Control AK is prepared from plasma of patients stabilized on long-term oral anticoagulant therapy. Therefore the factors of the prothrombin complex (II, VII, IX and X) have the same coagulation activity as plasma from patients on oral anticoagulant therapy. The PIVKA inhibitors (Proteins induced by Vitamin K Absence) are present in Coagulation Control AK as they are in patient's plasma. Coagulation Control AK contains stabilizers but no bactericide additives.

MATERIAL REQUIRED (not supplied with the kit)

- Pipettes
- Distilled water
- Calibration plasma

REF	5010004	AK-Calibrant	4x1 mL
-----	---------	--------------	--------

WARNING AND PRECAUTIONS

- For Research Use Only
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.
- Each single donor plasma and each lot of Coagulation Control AK are tested and found negative for HbsAg, HIV 1/2 Ab and HCV Ab. However, universal precautions (treating all human source materials as if potentially infectious) should be exercised.

STABILITY AND STORAGE

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8 °C. Stability after reconstitution:

RT*	RT* (Ceveron® alpha) Sample Tray	+2...8 °C	-20 °C
4 hours	4 hours	8 hours	1 month

The vials can be only frozen once. Upon storage, caps should be screwed tightly. * = room temperature

TEST PROCEDURE

PREPARATION AND PERFORMANCE OF THE TEST

- Open the vial carefully and reconstitute the contents in the volume of distilled water specified in the table for the respective test by carefully rotating the vial (avoid frothing).
- Please consider that Coagulation Control AK has to be reconstituted in different volumes of distilled water.
- Allow the reconstituted Coagulation Control AK to stand for 10 minutes at room temperature before use.
- Treat the reconstituted plasma as a citrated patient sample according to the instructions of the respective test.

TABLE OF ASSIGNED VALUES

The stated expected value and the confidence range have been determined by Technoclone.

ANALYSES RESULTS

EVALUATION OF THE RESULTS

The figures in the table are only applicable to the indicated lot number of Coagulation Control AK for Ceveron® alpha.

CEVERON® ALPHA

Before using, [import lot data](#) under Control & Calibrator Lots. You can download easily the lot data under www.ceveron.com or contact your distributor or directly sales@technoclone.com.

STANDARDIZATION

See "Table of assigned values" No international calibrator is established for the INR standardization of control plasmas.

LIMITATION OF THE TEST

The figures in the table are only applicable to the indicated lot number of Coagulation Control AK and the method has to adhere strictly as prescribed for the respective reagent.

LITERATURE

Please contact Technoclone or your local distributor.

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