
# Coagulation Control A f. Ceveron<sup>®</sup> alpha

For Research Use Only

REF 5021035 Coagulation Control A f. Ceveron<sup>®</sup> alpha 5 x 1 mLREF 5021040 Coagulation Control A f. Ceveron<sup>®</sup> alpha 50 x 1 mL

symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Ključova slova / Značenje simbola		
	manufacturer / Hersteller / fabbricante / fabricante / fabricante / fabricant / Tillverkaren / Fabrikanten / Produzent / Κατασκευαστής / Производител / Производител / výrobce / Proizvođač	 expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration / utgångsdatum / udløbsdato / Utløpsdato / Ημερομηνία λήξης / срок на годност / datum expirace / срок годности / datum expirace / Rok trajanja
	storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage / lagringstemperatur / opbevaringstemperatur / Oppbevaringstemperatur / Θερμοκρασία αποθήκευσης / съхранение на / teplota skladování / температура хранения / teplota skladování / Temperatura lagerovanja	 consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation / se användarinstruktioner / følg brugsvejledning / Følg bruksanvisningen / συμβουλευθείτε τις οδηγίες για τη χρήση / прочетете инструкцията за работа / potřeba řídit se instrukcemi / перед использованием читайте инструкцию / sledujte návod k použití / Pročitaj upustvo pre upotrebe
		 determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations / bestämmingar / bestemmelse / Bestemmelse / προσδιορισμοί / брой тестове / stanovení / определний / роčet stanovení / Definicija
<b>AQUA</b>	distilled water / destilliertes Wasser / acqua distillata / agua destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / destilirana voda / destilovaná voda / дистиллированная вода / destilovaná voda / Serija	<b>LOT</b> lot / Charge / lotto / lote / lote / lot / sats / serie / Parti / παρτίδα / партида номер / šarže / лот / šarže / in vitro dijagnostika
<b>BUF</b>	Reaction buffer / Reaktionspuffer / tampone di reazione / tampón de reacció / Tampão de reação / tampon de réaction / Reaktionspuffer / Reaktionsbuffer / Reaktionsbuffer / διάλυμα αντίδρασης / Реакционен буфер / Рабочий буферный раствор / Reakční pufr / Reakcioni pufer	<b>MTP</b> microtiter plate / Mikrotiterplatte / placa microtiter / microplaca / microplaca / microplaques sensibilisées / Mikrotiterplatta / Mikrotiterplade / mikrotiterplate / πλάκα μικροτιτλοδότησης / Микротитрѐна плака / Микропланшет / Mikrotitrační destička / Mikrotitracione ploče
<b>CAL</b>	Calibrator / Kalibrator / Calibratore / calibrador / calibrador / calibreur / Kalibrator / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrator	<b>REF</b> catalogue number / Katalognummer / numero di catalogo / número de catálogo / número de referència / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / καταλογен номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj
<b>CONJ</b>	Conjugate / Konjugat / Coniugato / conjugado / conjugado / conjugaté / Konjugerad / Konjugat / Konjugat / συνδεδετικό / Конюгат / Конъюгат / Konjugát / Konjugat	<b>RTU</b> ready to use / gebrauchsfertig / pronto all'uso / listo para usar / pronto a usar / prêt à l'emploi / færdig att användas / færdig til brug / klar til bruk / έτοιμο προς χρήση / Готов за употреба / готов к использованию / k přímému použití / Razrediti ili rastvoriti
<b>CONT</b>	Control / Kontrolle / controllo / control / control / contrôle / Kontroll / Kontroll / Kontroll / διάλυμα ελέγχου / Контрол / Контрольный образец / Kontrola / Kontrola	<b>STOP</b> stop solution / Stoppløsning / Soluzione di arresto / solución de parada / solução de paragem / solution d'arrêt / Stoppløsning / Stop-opløsning / Stoppløsning / διάλυμα παύσης / Стоп разтвор / Стоп-раствор / Zastavovací roztok / Stop solucija
<b>DIL</b>	dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o dissolver / diluir ou dissolver em / diluer ou dissoudre dans / späd eller upplös i / fortyndes eller opløses i / Fortyndes eller opløses i / αραιωση ή διάλυση σε / разтворете или разрежете с / zředit anebo rozpustit v / разбавить или растворить в / naředte nebo rozpustte v / razrediti ili rastvoriti u	<b>SUB</b> substrate / Substrat / substrato / substrato / substrato / substrat / Substrat / Substrat / Substrat / υπόστρωμα / Субстрат / Субстрат / Substrát / Substrat
<b>INC</b>	incubation buffer / Inkubationspuffer / tampone di incubazione / tampón de incubación / tampão de incubação / tampon d'incubation / Inkubationspuffer / Inkubationsbuffer / Vaskebufferkonsentrat / διάλυμα επώασης / Инкубационен буфер / Буфер для инкубации / Inkubacijski pufr / Inkubacioni pufer	<b>WASH</b> washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado / solución de lavado concentrada / tampão de lavagem concentrado / Tampon de lavage concentré / Vattenlösningskoncentrat / Vaskeopløsningskoncentrat / vaskeløsningskoncentrat / συμπυκνωμένο διάλυμα πλύσης / Концентриран миеш разтвор / Концентрат промывочного раствора / Koncentrát promývacího roztoku / Koncetrat solucije za ispiranje
<b>F I C</b>	For Research Use Only	

**PRODUCT DESCRIPTION**

**INTENDED USE**

For use as accuracy control of screening tests (PT/aPTT/TT/Fibrinogen) and determination of various coagulation factors in the interface between normal and abnormal coagulation activity.

Precision control, i.e. the recording of the day to day variation, is an important procedure in coagulation testing; checking the accuracy of the test system is particularly important.

Coagulation Control A can be used for:

- Accuracy control of all screening tests, i.e. prothrombin time (PT), partial thromboplastin time (aPTT), thrombin time (TT), fibrinogen determination (Clauss), and of various coagulation factors in the interface between normal and abnormal coagulation activity.
- Precision control of all parameters listed in the table.

**COMPOSITION**

Coagulation Control A is prepared from selected citrated plasmas of healthy donors by reducing all clotting factors. Coagulation Control A contains stabilizers but no bactericide additives.

**MATERIAL REQUIRED** (not supplied with the kit)

- Pipettes
- Distilled water

**WARNING AND PRECAUTIONS**

- For Research Use Only
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.
- Each single donor plasma and each lot of Coagulation Control A are tested and found negative for Hb<sub>s</sub>Ag, HIV 1/2 Ab and HCV Ab. However, universal precautions (treating all human source materials as if potentially infectious) should be exercised.

**STABILITY AND STORAGE**

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8 °C.

Stability after reconstitution:

RT*	RT* (Ceveron® alpha) Sample Tray	+2...8 °C	-20 °C
4 hours	4 hours	8 hours	1 month

The vials can be only frozen once.

Upon storage, caps should be screwed tightly.

\*=room temperature <sup>1</sup>

**TEST PROCEDURE**

**PREPARATION AND PERFORMANCE OF THE TEST**

- Open the vial carefully and reconstitute the contents in 1 mL of distilled water by carefully rotating the vial (avoid frothing).
- Allow the reconstituted Coagulation Control A to stand for 10 minutes at room temperature before use.
- Treat the reconstituted plasma as a citrated patient sample according to the instructions of the respective test.

**ANALYSES RESULTS**

**EVALUATION OF THE RESULTS**

The figures in the table are only applicable to the indicated lot number of Coagulation Control A for Ceveron® alpha.

**CEVERON® ALPHA**

Before using, [import lot data](#) under Control & Calibrator Lots. You can download easily the lot data under [www.ceveron.com](http://www.ceveron.com) or contact your distributor or directly [sales@technoclone.com](mailto:sales@technoclone.com).

**STANDARDIZATION AND TRACEABILITY**

The reported values were determined over multiple runs using a specific lot of reagent and against an Internal Reference Standard which is traceable to the current International Standards, identified in the acceptance range table.

For the tests where International Standards are not available, these parameters have been assigned against an Internal Reference Standard which is traceable to a frozen normal plasma pool of 100 donors.

**LIMITATION OF THE TEST**

The figures in the table are only applicable to the indicated lot number of Coagulation Control A and the method has to adhere strictly as prescribed for the respective reagent.

**LITERATURE**

Please contact Technoclone or your local distributor.

<sup>1</sup> Use the reconstituted plasma in single factor determinations within 2 hours stored at room temperature.

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