Factor VIII Inhibitor Plasma HCV neg.

GB

IT

PT

For Research
Use Only

DE

ES

FR

REF 5159010 Factor VIII Inhibitor Plasma HCV neg. 5 x 1 mL

Symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / Explicação dos símbolos / clé des symboles			
	manufacturer / Hersteller / fabbricante / fabricante / fabricant	AQUA	distilled water / destilliertes Wasser / acqua distillata / agua destilada / água destilada / eau distillée
	expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration	DIL	dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o disolver / diluir ou dissolver em / diluer ou dissoudre dans
	storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage	RUO	For Research Use Only
	consult instructions for use / Gebr auchsanweisung beachten / c onsultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation	LOT	lot / Charge / lotto / lote / lote / lot
		REF	catalogue number / Katalognummer / numero di catalogo / número de catálogo / número de referência / réf. de catalogue
	determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations		

PRODUCT DESCRIPTION**INTENDED USE**

For accuracy control of the determination of Factor VIII Inhibitor according to Bethesda. F VIII Inhibitor Plasma HCV neg. may be used as:

- a plasma sample of a hemophilia A-patient with F VIII Inhibitor
- for accuracy control of the determination of F VIII Inhibitor according to Bethesda (using the reagents of the F VIII Inhibitor Reagent Kit and the reagents for F VIII determination of Technoclone).

On the vial label the inhibitor capacity is indicated in Bethesda-Units/mL.

COMPOSITION

F VIII Inhibitor Plasma HCV neg. is an immune-adsorbed lyophilized normal human plasma with a specific inhibitor antibody added directed against F VIII activity (F VIII:C).

MATERIAL REQUIRED (not supplied with the kit)

- Pipette
- Distilled water

WARNING AND PRECAUTIONS

- RUO - For Research Use Only
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.
- This lot of reagents prepared from human blood and each single plasma used for this lot are HB_sAg, HIV 1/2 Ab and HCV Ab negative (see vial label).

STABILITY AND STORAGE

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8 °C.

Stability after reconstitution:

RT* (Ceveron **)	-20°C
2 hours	1 month

Repeated freezing is not recommended.

* = room temperature

** = in the Ceveron[®] alpha in the respective control area in the sample tray.

TEST PROCEDURE**PREPARATION OF REAGENT****CEVERON**

Technoclone provides application sheets for Ceveron[®] alpha. The application sheets contain analyser/assay specific handling and performance information which may differ from that provided in this instruction for use. In this case the information contained in the application sheet supersedes the information in this instruction for use. Please consult the instruction manual of the Ceveron[®] alpha.

MANUAL

- Open the vial carefully and reconstitute the contents in 1 mL of distilled water by carefully rotating the vial (avoid frothing).
- Allow the reconstituted plasma to stand for 10 minutes at room temperature before use. For standardisation the recommended reconstitution is 30 minutes.
- Treat the reconstituted plasma as a citrated patient sample according to the instructions of the respective test.

PERFORMANCE OF THE TEST

Determination of the Bethesda Units is performed with application of the Factor VIII Inhibitor Reagent Kit of Technoclone according to the detailed instructions for the performance of the test.

STANDARDISATION

Standardisation is done by Technoclone, using clotting and chromogenic assays for factor VIII.

LIMITATION OF THE TEST

The values found when testing should be compared to the value given on the data-key for the corresponding lot of Factor VIII Inhibitor Plasma HCV neg. If the results obtained are outside the recommended range, avoid measuring patient samples until the problem is solved.

LITERATURE

1. Carol K. Kaspar et al.: A More Uniform Measurement of Factor VIII Inhibitors; Thrombos. Diathes. haemorrh. 34, (1975), 869
2. Workshop of Inhibitors of Factor VIII and IX, Vienna, January 26. and 27., 1976, Facultas-Verlag

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