Coagulation Reference f. Ceveron®

GB

IT

PT

DE

ES

FR

For Research Use Only

REF 5220130 Coagulation Reference f. Ceveron® 5 x 1 mL

REF 5220135 Coagulation Reference f. Ceveron® 50 x 1 mL

symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Ključova slova / Značenje simbola			
	manufacturer / Hersteller / fabbricante / fabricante / fabricante / fabricant / Tillverkaren / Fabrikanten / Produsent / Κατασκευαστής / Производитель / Производител / Производител / výrobce / Proizvođač		expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration / utgångsdatum / udløbsdato / Utløpsdato / Ημερομηνία λήξης / срок на годност / datum expirace / срок годности / datum expirace / Rok trajanja
	storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage / lagringstemperatur / oppbevaringstemperatur / Oppbevaringstemperatur / Θερμοκρασία αποθήκευσης / съхранение на / teplota skladování / температура хранения / teplota skladování / Temperatura lagerovanja		consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation / se användarinstruktioner / følg brugsvejledning / Følg bruksanvisningen / συμβουλευθείτε τις οδηγίες για τη χρήση / прочетете инструкцията за работа / potřeba řídit se instrukcemi / перед использованием читайте инструкцию / sledujte návod k použití / Pročitaj upustvo pre upotrebe
			determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations / bestämmingar / bestemmelse / Bestemmelse / προσδιορισμοί / брой тестове / stanovení / определний / ročet stanovení / Definicija
AQUA	distilled water / destilliertes Wasser / acqua distillata / agua destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / destilirana voda / destilovaná voda / дистиллированная вода / destilovaná voda / Serija	LOT	lot / Charge / lotto / lote / lote / lot / sats / serie / Parti / παρτίδα / партида номер / šarže / lot / šarže / in vitro diagnostika
BUF	Reaction buffer / Reaktionspuffer / tampone di reazione / tampón de reacció / Tampão de reação / tampon de réaction / Reaktionspuffer / Reaktionspuffer / Reaktionspuffer / διάλυμα αντίδρασης / Реакционен буфер / Рабочий буферный раствор / Reakční pufr / Reakcioni pufer	MTP	microtiter plate / Mikrotiterplatte / placa microtiter / microplaca / microplaca / microplaques sensibilisées / Mikrotiterplatta / Mikrotiterplade / mikrotiterplate / πλάκα μικροτιτλοδότησης / Микротитрѐна плака / Микропланшет / Mikrotitrační destička / Mikrotitracione ploče
CAL	Calibrator / Kalibrator / Calibratore / calibrador / calibrador / calibreur / Kalibrator / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrator	REF	catalogue number / Katalognummer / numero di catalogo / número de catálogo / número de referència / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / каталожен номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj
CONJ	Conjugate / Konjugat / Coniugato / conjugado / conjugado / conjugaté / Konjugerad / Konjugat / Konjugat / συνδεδετικό / Конюгат / Конъюгат / Konjugát / Konjugat	RTU	ready to use / gebrauchsfertig / pronto all'uso / listo para usar / pronto a usar / prêt à l'emploi / færdig att användas / færdig til brug / klar til bruk / έτοιμο προς χρήση / Готов за употреба / готов к использованию / k přímému použití / Razrediti ili rastvoriti
CONT	Control / Kontrolle / controllo / control / control / contrôle / Kontroll / Kontroll / Kontroll / διάλυμα ελέγχου / Контрол / Контрольный образец / Kontrola / Kontrola	STOP	stop solution / Stopplösung / Soluzione di arresto / solución de parada / solução de paragem / solution d'arrêt / Stoppløsning / Stop-opløsning / Stoppløsning / διάλυμα παύσης / Стоп разтвор / Стоп-раствор / Zastavovací roztok / Stop solucija
DIL	dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o dissolver / diluir ou dissolver em / diluer ou dissoudre dans / späd eller upplös i / fortyndes eller opløses i / Fortyndes eller opløses i / αρύωση ή διάλυση σε / разтворете или разреждете с / zředit anebo rozpustit v / разбавить или растворить в / nafedte nebo rozpustte v / razrediti ili rastvoriti u	SUB	substrate / Substrat / substrato / substrato / substrato / substrat / Substrat / Substrat / Substrat / υπόστρωμα / Субстрат / Субстрат / Substrát / Substrat
INC	incubation buffer / Inkubationspuffer / tampone di incubazione / tampón de incubación / tampão de incubação / tampon d'incubation / Inkubationspuffer / Inkubationspuffer / Vaskebufferkonsentrat / διάλυμα επώασης / Инкубационен буфер / Буфер для инкубации / Inkubační pufr / Inkubacioni pufer	WASH	washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado / solución de lavado concentrada / tampão de lavagem concentrado / Tampon de lavage concentré / Vattenlösningsskoncentrat / Vaskeopløsningskoncentrat / vaskeløsningskoncentrat / συμπυκνωμένο διάλυμα πλύσης / Концентриран миеш разтвор / Концентрат промывочного раствора / Koncentrát promývacího roztoku / Koncetrat solucije za ispiranje
RUO	For Research Use Only		

PRODUCT DESCRIPTION**INTENDEND USE**

Calibration plasma for all screening tests and coagulation factors. Coagulation Reference can be used for:

- Establishing calibration curves for screening tests.
- Establishing calibration curves for all coagulation factors given in the table.
- Accuracy control for all screening tests, i.e. prothrombin time (PT), partial thromboplastin time (aPTT), thrombin time (TT) and determination of fibrinogen (Clauss) within the normal range.
- Precision control of all parameters listed in the table.

COMPOSITION

Coagulation Reference is prepared from specially selected citrated plasmas of healthy donors. By establishing this pool the clotting activity of the single donor plasmas is normally distributed. Thereby the "average" presence of all coagulation factors and inhibitors is secured. Coagulation Reference contains stabilizers but no bactericide additives.

MATERIAL REQUIRED (not supplied with the kit)

- Pipettes
- Distilled water

WARNING AND PRECAUTIONS

- RUO For Research Use Only
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.
- Each single donor plasma and each lot of Coagulation Reference are tested and found negative for Hb_sAg, HIV 1/2 Ab and HCV Ab. However, universal precautions (treating all human source materials as if potentially infectious) should be exercised.

STABILITY AND STORAGE

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8 °C.

Stability after reconstitution, without loss of activity:

RT*	RT* (Ceveron®) Sample Tray
4 hours	4 hours

Freezing of the reconstituted reagent is not recommended due to possible loss of activity.

*=room temperature

TEST PROCEDURE**PREPARATION AND PERFORMANCE OF THE TEST**

The lyophilized Coagulation Reference is reconstituted in the volume of distilled water specified for the respective single factor determination in the table. The Coagulation Reference contains 1 Unit/mL of the respective single factor to be determined.

For the indications in g/L the Coagulation Reference is reconstituted with 1 mL of distilled water.

- Open the vial carefully and reconstitute the contents in the volume of distilled water specified in the table for the respective test by carefully rotating the vial (avoid frothing).
- Allow the reconstituted Coagulation Reference to stand for 10 minutes at room temperature before use¹.
- Treat the reconstituted plasma as a citrated sample according to the instructions of the respective test.

LIMITATION OF THE TEST

The figures in the table are only applicable to the indicated lot number of Coagulation Reference and the method has to adhere strictly as prescribed for the respective reagent. The guidelines for evaluation given in this leaflet have been established assuming that a constant sensitivity of different batches is guaranteed by the manufactures of each reagent.

ANALYSES RESULTS**EVALUATION OF THE RESULTS**

The figures in the table are only applicable to the indicated lot number of Coagulation Reference for Ceveron®.

CEVERON

Before using, [import lot data](#) under Control & Calibrator Lots. You can download easily the lot data under www.ceveron.com or contact your distributor or directly sales@technoclone.com.

STANDARDIZATION AND TRACEABILITY

The reported values were determined over multiple runs using a specific lot of reagent and against an Internal Reference Standard which is traceable to the current International Standards, identified in the acceptance range table. For the tests where International Standards are not available, these parameters have been assigned against an Internal Reference Standard which is traceable to a frozen normal plasma pool of 100 donors.

LITERATURE

Please contact Technoclone or your local distributor.

¹ For standardization a reconstitution time of 30 minutes is recommended.

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