COATEST™ APC™ Resistance V - 82 3120 63

Intended use
For determination of resistance to activated protein C (APC), caused by the factor V:Q266 (factor V Leiden) mutation, in plasma from untreated individuals and from patients on oral anticoagulant (OAC) or heparin therapy.

Background and summary
The APC resistance phenotype^{1,2} is in more than 90% of the cases due to a mutation in the factor V gene, resulting in a replacement of Arg⁵⁰⁶ (R) with Gln (Q) in the factor V protein,^{3,4} the selectivity for the factor V:Q266, or other mutations in the factor V gene rendering the protein resistant to inactivation by APC.⁵ It is increased by normalizing the concentrations of other plasma proteins involved in formation and regulation of thrombin. Hence, by performing the APTT-based APC resistance assay in the presence of an excess of factor V deficient plasma, the sensitivity and specificity for the factor V:Q266 mutation is significantly increased. Further, this modification allows for the analysis of plasma from patients who are on OAC therapy.^{1,2}

Measurement principle
Sample plasma is prediluted in V-DEF Plasma and incubated with the APTT reagent for a standard period of time. Coagulation is triggered by the addition of CaCl₂ in the absence and presence of APC and the time for clot formation is recorded.
REAGENTS

- V-DEF Plasma** 4 vials
Stabilized, lyophilized human plasma, with a low level of factor V activity, containing the heparin antagonist Polybrene®. Reconstitute with 4.0 mL of NCCLS type II water ¹². Allow to stand for 30 minutes at 20-25°C. Swirl gently before use.
- CaCl₂** 1 vial
8 mL of calcium chloride, 0.025 mol/L, in Tris buffer containing 0.5% bovine serum albumin.
- APTT reagent** 1 vial
16 mL of purified phospholipids with colloidal silica as contact activator. Contains a preservative. Mix thoroughly on a 30-sec mixer before use.
- APC/CaCl₂** 4 vials
Human activated protein C copolyphlized with CaCl₂. Reconstitute with 2.0 mL of NCCLS type II water ¹². Allow to stand for 30 minutes at 20-25°C. Swirl gently before use.
- Control Plasma Level 1** 1 vial
Lyophilized human plasma. Reconstitute with 1.0 mL of NCCLS type II water ¹². Allow to stand for 30 minutes at 20-25°C. Swirl gently before use.
- Control Plasma Level 2** 1 vial
Lyophilized human plasma. Reconstitute with 1.0 mL of NCCLS type II water ¹². Allow to stand for 30 minutes at 20-25°C. Swirl gently before use.

Reagents 3 and 4 are not interchangeable between lots.
V-DEF Plasma (Art. No 82 3146 63), **Control Plasma Level 1** (Art. No 82 2650 63) and **Control Plasma Level 2** (Art. No 82 2668 63) are also available separately from Chromogenix.
CAUTION: Each donor unit used in the preparation of human source plasma has been tested by FDA approved methods for the presence of Hepatitis B surface antigen and antibodies to HIV 1 and 2 and Hepatitis C and found to be negative. However, since no test can completely rule out the presence of these blood borne diseases, the handling and disposal of human source reagents from this product should be made with care.

Danger
APC/CaCl₂
Hazard class: Resp. sens. 1, H334
Precautionary statements: P261: Avoid breathing dust/fume, P304 + P340: IF INHALED: If breathing is difficult, remove person to fresh air and keep comfortable for breathing, P342 + P311: If experiencing respiratory symptoms: Call a POISON CENTER/doctor. P501: Dispose of contents/container in accordance with local/regional/national/international regulation.

Supplemental hazard information: EUH 208: Contains APC (Activated Protein C). May produce an allergic reaction.
EUH 210: Safety data sheet available on request. Up to 3.97% of the mixture consists of component of unknown acute toxicity (dermal, inhalation) for the human health and unknown hazard to the aquatic environment.
CaCl₂, APTT Reagent, Control Plasma level 1, Control Plasma level 2, V-DEF plasma
Hazard class: none
Hazard statements: none
Precautionary statements: none
Supplemental hazard information:
CaCl₂, EUH 210: Safety data sheet available on request.
APTT Reagent, EUH 208: Contains 1,2-benzisothiazolin-3-one. May produce an allergic reaction. Up to 20.85% of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and unknown hazard to the aquatic environment.
Control Plasma level 1: = 100% of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and unknown hazard to the aquatic environment.
Control Plasma level 2: = 98.75% of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and unknown hazard to the aquatic environment.

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V-DEF plasma: = 98,26% of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and unknown hazard to the aquatic environment.
This product is for *in vitro* diagnostic use.
Materials required but not provided:

- Deionized water, filtered through 0,22 μm or NCCLS type II water.¹³
 - Calibrated pipettes.
 - Automated or semi-automated coagulation instruments which, employ mechanical or optical detection, methods should be used.
- NOTE:** When using automated or semi-automated instruments, always refer to the operator manual from the instrument manufacturer for exact test procedure.
Storage conditions and stability
The sealed reagents are stable at 2-8°C until the expiry date printed on the label.
Avoid contamination of the reagents by microorganisms.
- V-DEF Plasma**
Stability after reconstitution is 8 hours at 15-25°C, 24 hours at 2-8°C or 3 months at -20°C or below when stored in the original vial. *See NOTE.
 - CaCl₂**
Stability reagent in the original vial is stable for 1 week at 15-25°C or 1 month at 2-8°C.
 - APTT reagent**
Opened reagent in the original vial is stable for 1 week at 15-25°C or 1 month at 2-8°C. Do not freeze!
 - APC/CaCl₂**
Stability after reconstitution is 2 hours at 37°C, 8 hours at 15-25°C, 5 days at 2-8°C or 3 months at -20°C or below when stored in the original vial. *See NOTE.
 - Control Plasma Level 1**
Stability after reconstitution is 6 hours at 2-25°C or 3 months at -20°C or below when stored in the original vial. *See NOTE.
 - Control Plasma Level 2**
Stability after reconstitution is 6 hours at 2-25°C or 3 months at -20°C or below when stored in the original vial. *See NOTE.

*NOTE: Frozen reagent should be rapidly thawed at 37°C and gently mixed before use. Do not refreeze.
Quality controls
Control Plasma Level 1 and Level 2 should be used for validation of the assay series. Level 1 shows a normal response to APC whereas Level 2 shows a response consistent with presence of the factor V:Q266 mutation. Ranges of expected APC-V ratios are provided with each batch. If values outside the specified range are obtained, a complete check of reagents and instrument performance should be made and the analysis should be repeated. (See Calibration section for QC use of Control Plasma Level 1 and 2).
Traceability of calibrators and control materials

The reported values were determined over multiple runs on ACL Futura using a specific lot of reagent and against an internal House Standard. As an International Standard is not still available for the APC-V assay, the values have been assigned against a House Standard which is traceable to frozen plasma samples which have been determined to be homozygous or heterozygous respectively for factor V Leiden dependent APC Resistance.

Specimen collection
The patient should be at rest for 10 min. before sampling. Collect blood (9 volumes) in 0.1 mol/L sodium citrate (1 volume) and centrifuge within 24 hours at 2000 x g for 20 min. at room temperature. Take care to avoid contamination from the platelet layer into the plasma when the plasma is separated from the cells. Analyse the plasma within 25 hours from blood sampling.¹³ Alternatively, freeze rapidly at -70°C in aliquots of 1 mL or less and store for not more than 3 years at -70°C. Specimens should not be stored in a self defrosting freezer and not be thawed and refrozen before assay. Treat specimens as potentially infectious. For more information see NCCLS document H21-A3.¹⁴

Procedure

- All reagents must be brought to room temperature before use. Frozen plasma samples should be rapidly thawed at 37°C in a standardized way ensuring negligible loss of activity of labile coagulation factors and absence of cryoprecipitate.
- Pre-warm a sufficient volume of CaCl₂ and APC/CaCl₂ at 37±0.5°C.
- Pre-dilute one volume of sample plasma or Control Plasma with four volumes of V-DEF Plasma. Pre-diluted plasma should be analyzed within 45 minutes.
- Add one volume of plasma to a test tube or cuvette, then add an equal volume of the APTT reagent. Incubate at 37°C for 5 minutes. An instrument with a different, preset, incubation time may be used provided it is at least 3 minutes.
- Add one volume of CaCl₂ and simultaneously begin timing of clot formation. Record the time for clot formation.
- Perform a second analysis on the plasma, exchanging CaCl₂ with APC/CaCl₂ and record the time for clot formation.

Results
Calculate the factor V related APC ratio for the samples and controls:
APC ratio = $\frac{\text{Clot time APC/CaCl}_2}{\text{Clot time CaCl}_2}$

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
ENGLISH - Insert revision 02/2016

APC resistance due to a factor V mutation is indicated when the APC-V ratio is below or equal to the cut-off value (see Calibration). The following results illustrate the applicability of the Coatest APC Resistance V assay on the analysis of plasma from untreated individuals (Untreated), patients on OAC therapy, INR 1.3 - 6.0 (OAC), and patients receiving unfractionated heparin (Heparin) or low molecular weight heparins (shown here are Fragin® and Fraxiparin®).

Sample APC-V ratios are distributed according to the factor V genotype. For a given genotype, similar APC-V ratios are obtained regardless of patient group. Typical median values are 2.6 (normal), 1.7 (heterozygous FV:Q266) and 1.2 (homozygous FV:Q266). Basal APT times for untreated individuals may differ moderately as compared to the original Coatest APC Resistance assay.

Performance Characteristics
PRECISION
APC-V ratios were calculated from 22 single replicate analyses of Control Plasma Level 1 on an ACL instrument using 11 different reagent combinations on 11 different occasions.

Control Plasma Level 1	APC-V ratio 2.4	CV% (Within series)	CV% (Between series)
		≤4 %	≤4 %

Duplicate analyses of plasma from altogether 630 individuals on three different types of instruments (ST-4, Thrombolyzer and MLa/Electra 900) resulted in within series CV ≤ 3% and ≤ 5% for the clotting times obtained in the absence and presence of APC respectively. NCCLS document EP5-T2.¹⁵

Calibration
APC-V ratios obtained with the Coatest APC Resistance 5 method are lower than the APC ratios obtained with the original Coatest APC Resistance method, independent of instrumentation used. It is recommended that each user establishes the performance of his own instrument and determines the factor V related APC resistance cut-off value through the following procedure:

- Perform five independent determinations of the APC-V ratio, using at least triplicates in each series, of a plasma sample with normal APC response. Confirm that the inter and intra assay variation of the APC-V ratio is below 7%. In case a satisfactory performance already has been established with the original Coatest APC Resistance method, this step may be omitted.
- Determine the APC-V ratios for at least 30 plasma samples from healthy individuals in the age range 20 - 65 years. Include Control Plasma Level 1 and Level 2 for assay validation.
- Verify that the APC-V ratios for the Control Plasmas are within their specified ranges.
- Calculate the median APC-V ratio.
- Calculate the factor V related APC resistance cut-off value as 0.8 times the median APC-V ratio when below 2.8 and as 0.75 times the median APC-V ratio when 2.8 or higher.
- The APC-V ratio for Control Plasma Level 1 should be within the normal range. The APC-V ratio for Control Plasma Level 2 should be below the cut-off value.

Sensitivity
Coatest APC Resistance V provides 100% sensitivity for FV:Q266 as determined on Thrombolyzer (n = 447), ACL (n = 295), ST-4 (n = 248) and MLa/Electra (n = 50).
Limitation/interfering factors

No significant differences are obtained between fresh and frozen samples. The sensitivity and specificity for the factor V:Q266 mutation in analysis of plasma from OAC patients is not affected by the INR value. The prescribed assay procedure allows for the analysis of plasma from heparinized patients at heparin levels < 1 IU/mL plasma (unfractionated and low molecular weight heparins). Although the 1+4 predilution strongly decreases interferences, it can not be excluded that analysis of plasma from patients with high inhibitor activity (e.g. phospholipid antibodies) may give an abnormal APT time and thus possibly misleading results. In such cases, increasing the dilution factor (e.g. 1+9 or 1+19) may correct the results. As for any APT based assay, care should be taken to avoid contact activation of samples since this may lead to activation of FVIII and FV.
Reference values

The APC-V ratios obtained from analysis of plasmas from 61 healthy individuals on ACL and ST-4 and from 390 healthy individuals on Thrombolyzer were in the range 2 - 3.5. No difference is found between sexes.

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