# TECHNOVIEW Arixtra®\*

GB

IT

RU

DE

FR

## For Research Use Only

REF	5090010	TECHNOVIEW Arixtra® CAL Set	5 x 1 mL
REF	5090012	TECHNOVIEW Arixtra® CONT Low	6 x 1 mL
REF	5090014	TECHNOVIEW Arixtra® CONT High	6 x 1 mL

symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Klíčova slova / Značenje simbola			
	manufacturer / Hersteller / fabbricante / fabricante / fabricant / Tillverkaren / Fabrikanten / Produzent / Κατασκευαστής / Производитель / Производител / distillerovaná voda / destilovaná voda / Destilovaná voda		expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration / utgångsdatum / udløbsdato / Utløpsdato / Ημερομηνία λήξης / срок на годност / datum expirace/ срок годности / datum expirace / Rok trajanja
	storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservació / temperatura de conservação / température de stockage / lagringstemperatur / Oppbevaringstemperatur / θερμοκρασία αποθήκευσης / съхранение на / teplota skladování / температура хранения / teplota skladování / Temperatura lagerovanja		consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation / se användarinstruktioner / følg bruksvejledning / Følg bruksanvisningen / συμβουλευθείτε τις οδηγίες για τη χρήση / прочетете инструкцията за работа / potfeba říditi se instrukcemi / перед использованием читайте инструкцию / sledujte návod k použití / Pročitaj upustvo pre upotrebe
			determinations / Bestimmungen / determinazioni / determinaciones / determinações / détermination / bestämmingar / bestemmelse / Bestemmelse / προσδιορισμοί / брой тестове / stanovení / определний / počet stanovení / Definicija
<b>AQUA</b>	distilled water / destilliertes Wasser / acqua distillata / agua destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / destilirana voda / destilovaná voda / destilovaná voda / destilovaná voda / Destilovaná voda	<b>LOT</b>	lot / Charge / lotto / lote / lote / lot / sats / serie / Parti / παρτία / партида номер / šarže / lot / šarže / Serija
<b>BUF</b>	Reaction buffer / Reaktionspuffer / tampone di reazione / tampón de reacció / Tampão de reação / tampon de réaction / Reaktionspuffer / Reaktionspuffer / Reaktionspuffer / διάλυμα αντίδρασης / Реакционен буфер / Рабочий буферный раствор / Reakční pufr / Reakcioni pufer	<b>MTP</b>	microtiter plate / Mikrotiterplatte / placa microtiter / microplaca / microplaca / microplaques sensibilisées / Mikrotiterplatta / Mikrotiterplade / mikrotiterplate / πλάκα μικροτιτλοδότησης / Микротитърна пласка / Микропланшет / Mikrotitrační destička / Mikrotitracione ploče
<b>CAL</b>	Calibrator / Kalibrator / Calibratore / calibrador / calibrador / calibreteur / Kalibrator / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrator	<b>REF</b>	catalogue number / Katalognummer / numero di catalogo / número de catálogo / número de referência / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / каталожен номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj
<b>CONJ</b>	Conjugate / Konjugat / Coniugato / conjugado / conjugado / conjugaté / Konjugerad / Konjugat / Konjugat / συνδεδετικό / Конюгат / Конъюгат / Konjugát / Konjugat	<b>RTU</b>	ready to use / gebrauchsfertig / pronto all'uso / listo para usar / pronto a usar / prêt à l'emploi / færdig att användas / færdig til brug / klar til bruk / έτοιμο προς χρήση / Готов за употреба / готов к использованию / k přímému použití / Razrediti ili rastvoriti
<b>CONT</b>	Control / Kontrolle / controllo / control / control / contrôle / Kontroll / Kontroll / Kontroll / διάλυμα ελέγχου / Контрол / Контрольный образец / Kontrola / Kontrola	<b>STOP</b>	stop solution / Stopplösung / Soluzione di arresto / solución de parada / solução de paragem / solution d'arrêt / Stoppløsning / Stop-opløsning / Stoppløsning / διάλυμα παύσης / Стоп разтвор / Стоп-раствор / Zastavovací roztok / Stop solucija
<b>DIL</b>	dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o dissolver / diluir ou dissolver em / diluer ou dissoudre dans / spädd eller upplöst i / fortyndes eller opløses i / Fortyndes eller opløses i / αρωση ή διάλυση σε / разтворете или разрежете с / zředit anebo rozpustit v / разбавить или растворить в / nafedte nebo rozpustit v / razrediti ili rastvoriti u	<b>SUB</b>	substrate / Substrat / substrato / substrato / substrato / substrat / Substrat / Substrat / Substrat / υπόστρωμα / Субстрат / Субстрат / Substrát / Substrat
<b>INC</b>	incubation buffer / Inkubationspuffer / tampone di incubazione / tampón de incubació / Tampão de incubação / tampon d'incubation / Inkubationspuffer/ Inkubationspuffer/ Vaskebufferkonsentrat / διάλυμα επώασης / Инкубационен буфер / Буфер для инкубации / Inkubační pufr / Inkubacioni pufer	<b>WASH</b>	washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado / solución de lavado concentrada / tampão de lavagem concentrado / Tampon de lavage concentré / Vattenlösning / Vaskeopløsningskonsentrat / vaskeløsningskonsentrat / συμπυκνωμένο διάλυμα πλύσης / Концентриран миеш разтвор / Концентрат промывочного раствора / Koncentrat promývacieho roztoku / Koncentrat solucije za ispiranje
<b>RUO</b>	For Research Use Only		

**PRODUCT DESCRIPTION**

**INTENDED USE**

TECHNOVIEW Arixtra Calibrator Set is a set of calibration plasmas and TECHNOVIEW Arixtra Control Plasma is a set of low and high control plasmas for calibration and the quality control of Arixtra (Fondaparinux sodium) measurements, titrated and optimised using anti-Xa assays. Arixtra® can be used as an anticoagulant for curative or preventive indications. Measuring the Arixtra concentration in sample's plasma allows monitoring the therapy and adjusting drug dosage.

**COMPOSITION**

TECHNOVIEW Arixtra Calibrator and Control plasma are prepared from citrated plasmas supplemented with different concentrations of Arixtra®. The plasma contains stabilizers but no bactericide additives.

TECHNOVIEW Arixtra Calibrator Set 5x1 mL [REF] 5090010 contains:

mL	Reagent
1	Calibrator 1 human plasma, freeze dried, without any additive of Arixtra (CAL 1 at 0 µg/mL)
1	Calibrator 2 human plasma, freeze dried, supplemented with Arixtra (CAL 2 at ~ 0.5 µg/mL)
1	Calibrator 3 human plasma, freeze dried, supplemented with Arixtra (CAL 3 at ~ 1 µg/mL)
1	Calibrator 4 human plasma, freeze dried, supplemented with Arixtra (CAL 4 at ~ 1.5 µg/mL)
1	Calibrator 5 human plasma, freeze dried, supplemented with Arixtra (CAL 5 at ~ 2 µg/mL)

TECHNOVIEW Arixtra Control Low 6x1 mL [REF] 5090012 containing 6 vials:

mL	Reagent
1	Control Low human plasma, freeze dried, supplemented with Arixtra (CONT L at ~ 0.4 µg/mL)

TECHNOVIEW Arixtra Control High 6x1 mL [REF] 5090014 containing 6 vials:

mL	Reagent
1	Control High human plasma, freeze dried, supplemented with Arixtra (CONT H at ~ 1.2 µg/mL)

**MATERIAL REQUIRED** (not supplied with the kit)

- Pipettes - Distilled water

**WARNING AND PRECAUTIONS**

- RUO - For Research Use Only
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.
- Each single donor plasma and each lot of Control and Calibration plasma are tested and found negative for Hb<sub>s</sub>Ag, HIV 1/2 Ab and HCV Ab. However, universal precautions (treating all human source materials as if potentially infectious) should be exercised.

**STABILITY AND STORAGE**

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8 °C. Stability after reconstitution:

RT*	+2...8 °C	-20 °C
48 hours	7 days	1 month

The vials can be only frozen once. Upon storage, caps should be screwed tightly.

\*=room temperature <sup>1</sup>

**TEST PROCEDURE**

**PREPARATION AND PERFORMANCE OF THE TEST**

- Open the vial carefully and reconstitute the contents in 1 mL of distilled water by carefully rotating the vial until the product is completely reconstituted (avoid frothing).
- Allow the reconstituted TECHNOVIEW Arixtra Calibrator and Control plasma to stand for 10 minutes at room temperature before use <sup>1</sup>. Invert to mix before use (avoid foaming).
- Treat the reconstituted plasma as a citrated sample according to the instructions of the respective test.

**APPLICATIONS FOR INSTRUMENTS**

Application sheets are available from Technoclone or your local distributor upon request.

**ANALYSES RESULTS**

**EVALUATION OF THE RESULTS**

The figures in the batch table are only applicable to the indicated lot number of TECHNOVIEW Arixtra Calibrator and Control plasma and the method has to adhere strictly as prescribed for the respective reagent.

**STANDARDIZATION**

TECHNOVIEW Arixtra Calibrator and Control plasma are calibrated against an internal Arixtra reference standard. the International Standard preparation 01/608.

**LIMITATION OF THE TEST**

The figures in the batch table are only applicable to the indicated lot number of TECHNOVIEW Arixtra Calibrator and Control plasma and the method has to adhere strictly as prescribed for the respective reagent.

The guidelines for evaluation given in this leaflet have been established assuming that a constant sensitivity of different batches is guaranteed by the manufactures of each reagent.

**PERFORMANCE CHARACTERISTICS**

TECHNOVIEW Arixtra Control plasmas allow validating the calibration curve for the measurements of Arixtra® in plasma, especially with anti-Xa methods. The calibration curve obtained covers the usual concentrations currently observed during Arixtra® therapy. The currently available anti-Xa methods, used for the measurement of heparins and their analogues in plasma, offer a sensitivity threshold of about 0.05 µg/mL.

**LITERATURE**

1. J. Hirsch, "Fondaparinux", BC Decker Inc., Hamilton, 2007.
2. J.M. Walenga, J. Fareed, W.P. Jeske, F.X. Frapaise, R.L. Bick, M.M. Samama, "Development of a Synthetic Heparin Pentasaccharide: Fondaparinux", *Turk J Haematol*, 2002; 19(2):137-150.

<sup>1</sup> For standardization a reconstitution time of 30 minutes is recommended.

\* Arixtra® is a registered trade mark of Glaxo Smith Kline

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