TECHNOVIEW® Edoxaban Controls

3015607 Rev. 001 [02] 03/2018

For Research Use Only

TECHNOVIEW® Edoxaban Controls- English

INTENDED USE

TECHNOVIEW® Edoxaban Controls are plasmas with different concentrations (low, medium and high) of edoxaban to be used for quality control of edoxaban measurements, optimized using TECHNOCHROM® anti-Xa assay on Ceveron® alpha.

SUMMARY

Edoxaban (Savaysa®, Lixiana®) is an oral anticoagulant drug which acts as a direct factor Xa inhibitor.

REAGENTS

The REF 5090251 TECHNOVIEW® Edoxaban CONT Low 5x1mL contains:

mL	Reagent	Description
5 x 1	CONT L	Control low, human plasma, lyophilized, ~30 ng/mL edoxaban

The REF 5090252 TECHNOVIEW® Edoxaban CONT Medium 5x1mL contains:

mL	Reagent	Description
5 x 1	CONT M	Control medium, human plasma, lyophilized, ~125 ng/mL edoxaban

The REF 5090253 TECHNOVIEW® Edoxaban CONT High 5x1mL contains:

mL	Reagent	Description
5 x 1	CONT H	Control high, human plasma, lyophilized, ~400 ng/mL edoxaban

Please consult the batch table provided with the controls for the exact concentrations in ng/mL of each control lot.

Material required (not supplied with the kit)

- Pipettes
- Distilled water

REF 5340250 TECHNOCHROM® anti-Xa Kit 80.T

REF 5090250 TECHNOVIEW® Edoxaban Calibrations Set Low 5x1mL

Warning and precautions

RUO - For Research Use Only

Every single donor plasma and every lot of the calibrators included is tested and found negative for HbSag, HIV 1/2 antibodies and HCV antibodies. However, general precautions should be taken by handling all human source materials as potentially infectious.

All food and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force. Dispose the same way as hospital waste.

A Material Safety Data Sheet for this product is available from www.technoclone.com

Stability and storage

The expiry date printed on the labels is only applicable to storage of the unopened containers at 2...8°C. Stability after reconstitution:

Reagent	18...25°C	2...8°C	<-20°C
CONT L	8 hours	6 days	2 weeks
CONT M	8 hours	6 days	2 weeks
CONT H	8 hours	6 days	2 weeks

Reagent can only be frozen once and in its original vial.

TEST PROCEDURE

Preparation of reagents

Vials with lyophilized controls and distilled water used for reconstitution need to reach room temperature (18...25°C) before reconstitution. Reconstitute each vial (control) with 1 mL of distilled water. Allow the reconstituted material to stand at room temperature for 30 minutes followed by swirling of the vial before use. Avoid foaming.

Performance of the test

The TECHNOVIEW® Edoxaban Calibrator Set is always used in combination with the TECHNOVIEW® Edoxaban Control Low, Medium and High and the TECHNOCHROM® anti-Xa kit.

The edoxaban concentration of the controls may vary from one lot to another but is clearly indicated in the lot specific batch table provided in the control box.

To cover the whole assay range, two applications are established for this product. For standard application, a quality control using control low and control medium is recommended. For the application high, a quality control using control Medium and control high is recommended. Technoclone provides instrument specific application sheets which contain analyser/assay specific handling and performance information.

STANDARDISATION

The TECHNOVIEW® Edoxaban controls are directly traceable to high performance liquid chromatography with mass spectrometry detection (HPLC-MS).

LITERATURE

Please contact Technoclone www.technoclone.com or your local distributor.

Savaysa® and Lixiana® are registered trademarks of Daiichi Sankyo.

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- Pipetten
- Destilliertes Wasser

REF 5340250 TECHNOCHROM® anti-Xa Kit 80.T

REF 5090250 TECHNOVIEW® Edoxaban Calibration Set 5x1mL

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Technoclone Herstellung von Diagnostika und Arzneimitteln GmbH, Brunner Strasse 67, 1230 Vienna, Austria