

arOCell TK 210 Cell Culture Protocol – For research use only. Not for use in diagnostics procedures

Note: Calibration curve and quantification should be carried out as described in the user manual.

Cell Extraction Buffer: 50 mM Tris-HCl, pH = 8.0, 150 mM NaCl, 0.5% NP-40 with 0.1% SDS

1. Spin down the plate 5 min @ 2000 xG

Cell Supernatant

2. Transfer 120 μ L supernatant to an uncoated microtiter plate
3. Add 80 μ L sample dilution buffer ($V_f=200 \mu$ L), briefly mix on orbital shaker
4. Seal and incubate the plate 60 min @ RT
5. Prewash the wells of the coated microtiter plate 4 times with 350 μ L wash buffer.
6. Within 10 minutes of completing the above wash, transfer 150 μ L to prewashed coated microtiter plate
7. Continue with “Joint Detection Protocol”

Cell Pellet/Cell Extract

2. Add 200 μ L prechilled (4°C) wash buffer and spin 5 min @ 2000 xG
3. Remove wash buffer
4. Lyse cells with 120 μ L cell extraction buffer
5. Transfer 80 μ L sample to uncoated microtiter plate
6. Add 80 μ L sample dilution buffer ($V_f=160 \mu$)
7. Seal and incubate the plate 60 min @ RT
8. Prewash the wells of the coated microtiter plate 4 times with 350 μ L wash buffer.
9. Within 10 minutes of completing the above wash, transfer 100 μ L to prewashed coated microtiter plate
10. Continue with “Joint Detection Protocol”

Joint Detection Protocol

- a. Cover and seal. Incubate 2 hours @ RT w/ intermediate shaking
- b. Wash 4 times with 350 μ L wash buffer
- c. Add 100 μ L Biotinylated anti TK1 to each well
- d. Cover and seal. Incubate 1 hour @ RT w/ intermediate shaking
- e. Wash 4 times with 350 μ L wash buffer
- f. Add 100 μ L streptavidin-HRP conjugate
- g. Cover and seal. Incubate 30 minutes @ RT w/ intermediate shaking
- h. Wash 4 times with 350 μ L wash buffer
- i. Add 100 μ L TMB substrate to each well. Incubate 15 minutes @ RT **in the dark** with no rotation
- j. Add 100 μ L Stop solution. Ensure complete mixing
- k. Read within 15 minutes @ 450 nm.