

## Platelet Rich Plasma Protocol

### **PRP preparation:**

1. Divide blood specimen into 2 parts
2. Centrifuge first part of blood specimen at 100 g for 8 minutes.
3. Transfer  $\frac{3}{4}$  of the upper plasma layer to a new tube using a pipette with disposable tips (PRP specimen obtained).
4. Platelets concentration is calculated on hematological analyzer (platelets count)
5. Second part of blood specimen is used for PFP specimen preparation\*
6. PFP is used for PRP dilution to get the desired platelets concentration.

Recommended platelet concentration  $(100-200) \times 10^3 \mu\text{L}^{-1}$

For TD4D measurements (thrombin generation) use PLS kit, but do not add lipids from PLS vial.

In order to increase sensitivity to P2Y<sub>12</sub> inhibitors ADP addition to PRP is used. ADP in concentration 100 uM is added immediately after plasma is mixed with reagent 2. The incubation time was 3 min.

### **\*PFP preparation (as described in User Manual):**

1. Centrifuge blood specimen tube at 1600 g for 15 minutes.
2. Transfer  $\frac{3}{4}$  of the upper plasma layer to a new tube using a pipette with disposable tips.
3. Centrifuge the plasma specimen tube at 10000 g for 5 minutes.
4. Transfer 90% of the specimen liquid volume from the upper PFP layer to a new tube using a pipette with disposable tips.