Intended use

For determination of resistance to activated protein C (APC), caused by the factor V:Q506 (factor V Leiden) mutation, in plasma from untreated individuals and from patients on oral anticoagulant (OAC) or heparin therapy.

Background and summary

The APC resistance phenotype^{1,2} is in more than 90% of the cases due to a mutation in the factor V gene, resulting in a replacement of Arg⁵⁰⁶ (R) with Gln (Q) in the factor V protein.³⁻⁵ The selectivity for the factor V:Q506, or other mutations in the factor V gene rendering the protein resistant to inactivation by APC⁶, is increased by normalizing the concentrations of other plasma proteins involved in formation and regulation of thrombin. Hence, by performing the APTT-based APC resistance assay in the presence of an excess of factor V deficient plasma, the sensitivity and specificity for the factor V:Q506 mutation is significantly increased. Further, this modification allows for the analysis of plasma from patients who are on OAC therapy.⁷⁻¹¹

Measurement principle

Sample plasma is prediluted in V-DEF Plasma and incubated with the APTT reagent for a standard period of time. Coagulation is triggered by the addition of CaCl₂ in the absence and presence of APC and the time for dot formation is recorded.

Reagents

- V-DEF Plasma:** 4 vials - Stabilized, lyophilized human plasma, with a low level of factor V activity, containing the heparin antagonist Polybrene®. Reconstitute with 4.0 mL of NCCLS type II water¹². Allow to stand for 30 minutes at 20-25°C. Swirl gently before use.
- CaCl₂:** 1 vial - 8 mL of calcium chloride, 0.025 mol/L, in Tris buffer containing 0.5% bovine serum albumin.
- APTT reagent:** 1 vial - 16 mL of purified phospholipids with colloidal silica as contact activator. Contains a preservative. Mix thoroughly on a Vortex mixer before use.
- APC/CaCl₂:** 4 vials - Human activated protein C colyophilized with CaCl₂. Reconstitute with 2.0 mL of NCCLS type II water¹². Allow to stand for 30 minutes at 20-25°C. Swirl gently before use.
- Control Plasma Level 1:** 1 vial - Lyophilized human plasma. Reconstitute with 1.0 mL of NCCLS type II water¹². Allow to stand for 30 minutes at 20-25°C. Swirl gently before use.
- Control Plasma Level 2:** 1 vial - Lyophilized human plasma. Reconstitute with 1.0 mL of NCCLS type II water¹². Allow to stand for 30 minutes at 20-25°C. Swirl gently before use.

Reagents 3 and 4 are not interchangeable between lots.

V-DEF Plasma (Cat. No. 82 3146 63), Control Plasma Level 1 (Cat. No. 82 2650 63) and Control Plasma Level 2 (Cat. No. 82 2668 63) are also available separately from Chromogenix.

Caution

APC/Calcium chloride, APC Control Plasma Levels 1 and 2, and Factor V Reagent Plasma contain human source material that tested non-reactive for HIV antibody, Hepatitis B Surface Antigen and Anti-HCV at the donor stage. This product, as with all human based specimens, should be handled with proper laboratory safety procedures to minimize the risk of transmission of infectious disease.

APC/Calcium chloride and Calcium chloride contain bovine material. All donor animals were sourced from BSE-free herds. The cattle received ante- and post-mortem health inspection by a veterinarian, and they were apparently free from infectious and contagious material. However, the material should be treated as potentially infectious.

Calcium Chloride

WARNING
Hazard class: Skin Sens. 1, H317. Aquatic Chronic 3, H412
Hazardous ingredients: Reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1) (CIT/MIT).

Hazard statements: H317: May cause an allergic skin reaction. H412: Harmful to aquatic life with long lasting effects.

Precautionary statements: P272: Contaminated work clothing should not be allowed out of the workplace. P273: Avoid release to the environment. P280: Wear protective clothing, protective gloves. P302+P352: IF ON SKIN: Wash with plenty of water. P333+P313: If skin irritation or rash occurs: Get medical advice/attention. P362+P364: Take off contaminated clothing and wash it before reuse.

APC/CaCl₂

DANGER
Hazard class: Resp. sens. 1, H334
Hazard statements: H334: May cause allergy or asthma symptoms or breathing difficulties if inhaled.

Precautionary statements: P261: Avoid breathing dust/fume. P304 + P340: IF INHALED: IF breathing is difficult, remove person to fresh air and keep comfortable for breathing. P342 + P311: If experiencing respiratory symptoms: Call a POISON CENTER/doctor. P501: Dispose of contents/container in accordance with local/regional/national/international regulation.

Supplemental hazard information: EUH 208: Contains APC (Activated Protein C). May produce an allergic reaction. EUH 210: Safety data sheet available on request. Up to 3.97% of the mixture consists of ingredients of unknown acute toxicity (dermal, inhalation) for the human health and unknown hazard to the aquatic environment.

APTT Reagent, Control Plasma level 1, Control Plasma level 2, V-DEF plasma

Hazard class: none
Hazard statements: none
Precautionary statements: none

Supplemental hazard information:

APTT Reagent: EUH 208: Contains 1,2-benzisothiazolin-3-one. EUH 210: Safety data sheet available on request. May produce an allergic reaction. Up to 20.85% of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and unknown hazard to the aquatic environment.
Control Plasma level 1: Up to 100% of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and unknown hazard to the aquatic environment.
Control Plasma level 2: Up to 98,75% of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and unknown hazard to the aquatic environment.
V-DEF plasma: Up to 98,26% of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and unknown hazard to the aquatic environment.
 This product is for *in vitro* diagnostic use.
Materials required but not provided:
 — Deionized water, filtered through 0.22 µm or NCCLS type II water.¹³
 — Calibrated pipettes.
 — Automated or semi-automated coagulation instruments which, employ mechanical or optical detection, methods should be used.

NOTE: When using automated or semi-automated instruments, always refer to the operator manual from the instrument manufacturer for exact procedures.

Storage conditions and stability

The sealed reagents are stable at 2-8°C until the expiry date printed on the label. Avoid contamination of the reagents by microorganisms.

- V-DEF Plasma:** Stability after reconstitution is 8 hours at 15-25°C, 24 hours at 2-8°C or 3 months at -20°C or below when stored in the original vial. *See NOTE.
- CaCl₂:** Opened reagent in the original vial is stable for 1 week at 15-25°C or 1 month at 2-8°C.
- APTT reagent:** Opened reagent in the original vial is stable for 1 week at 15-25°C or 1 month at 2-8°C. Do not freeze!
- APC/CaCl₂:** Stability after reconstitution is 2 hours at 37°C, 8 hours at 15-25°C, 5 days at 2-8°C or 3 months at -20°C or below when stored in the original vial. *See NOTE.
- Control Plasma Level 1:** Stability after reconstitution is 6 hours at 2-25°C or 3 months at -20°C or below when stored in the original vial. *See NOTE.
- Control Plasma Level 2:** Stability after reconstitution is 6 hours at 2-25°C or 3 months at -20°C or below when stored in the original vial. *See NOTE.

*NOTE: Frozen reagent should be rapidly thawed at 37°C and gently mixed before use. Do not refreeze.

Quality controls

Control Plasma Level 1 and Level 2 should be used for validation of the assay series. Level 1 shows a normal response to APC whereas Level 2 shows a response consistent with presence of the factor V:Q506 mutation. Ranges of expected APC-V ratios are provided with each batch. If values outside the specified range are obtained, a complete check of reagents and instrument performance should be made and the analysis should be repeated. (See Calibration section for QC use of Control Plasma Level 1 and 2).

Traceability of calibrators and control materials

The reported values were determined over multiple runs on ACL Futura using a specific lot of reagent and against an internal House Standard. As an International Standard is not still available for the APC-R V assay, the values have been assigned against a House Standard which is traceable to frozen plasma samples which have been determined to be homozygous or heterozygous respectively for factor V dependent APC Resistance.

Specimen collection


The patient should be at rest for 10 min. before sampling. Collect blood (9 volumes) in 0.1 mol/L sodium citrate (1 volume) and centrifuge within 24 hours at 2000 x g for 20 min. at room temperature. Take care to avoid contamination from the platelet layer into the plasma when the plasma is separated from the cells. Analyse the plasma within 25 hours from blood sampling.¹³ Alternatively, freeze rapidly at -70°C in aliquots of 1 mL or less and store for not more than 3 years at -70°C. Specimens should not be stored in a self defrosting freezer and not be thawed and refrozen before assay. Treat specimens as potentially infectious. For more information see NCCLS document H21-A3.¹⁴

Procedure

- All reagents must be brought to room temperature before use. Frozen plasma samples should be rapidly thawed at 37°C in a standardized way ensuring negligible loss of activity of labile coagulation factors and absence of cryoprecipitate.
- Pre-warm a sufficient volume of CaCl₂ and APC/CaCl₂ at 37±0.5°C.
- Pre-dilute one volume of sample plasma or Control Plasma with four volumes of V-DEF Plasma. Pre-diluted plasma should be analysed within 45 minutes.
- Add one volume of plasma to a test tube or cuvette, then add an equal volume of the APTT reagent. Incubate at 37°C for 5 minutes. An instrument with a different, preset, incubation time may be used provided it is at least 3 minutes.
- Add one volume of CaCl₂ and simultaneously begin timing of clot formation. Record the time for clot formation.
- Perform a second analysis on the plasma, exchanging CaCl₂ with APC/CaCl₂ and record the time for clot formation.

Results

Calculate the factor V related APC ratio for the samples and controls:

Sample APC-V ratios are distributed according to the factor V genotype. For a given genotype, similar APC-V ratios are obtained regardless of patient group. Typical median values are 2.6 (normal), 1.7 (heterozygous FV:Q506) and 1.2 (homozygous FV:Q506). Basal APT times for untreated individuals may differ moderately as compared to the original Coatest APC Resistance assay.

Performance Characteristics

Precision
 APC-V ratios were calculated from 22 single replicate analyses of Control Plasma Level 1 on an ACL instrument using 11 different reagent combinations on 11 different occasions.

Control Plasma Level 1	APC-V ratio 2.4	CV% (Within series)	CV% (Between series)
		≤4%	≤4%

Duplicate analyses of plasma from altogether 630 individuals on three different types of instruments (ST-4, Thrombolyzer and mL/Electra 900) resulted in within series CV ≤ 3% and ≤ 5% for the clotting times obtained in the absence and presence of APC respectively. NCCLS document EP5-T2.¹⁵

Calibration

APC-V ratios obtained with the Coatest APC Resistance 5 method are lower than the APC ratios obtained with the original Coatest APC Resistance method, independent of instrumentation used. It is recommended that each user establishes the performance of his own instrument and determines the factor V related APC resistance cut-off value through the following procedure:

- Perform five independent determinations of the APC-V ratio, using at least triplicates in each series, of a plasma sample with normal APC response. Confirm that the inter and intra assay variation of the APC-V ratio is below 7%. In case a satisfactory performance already has been established with the original Coatest APC Resistance method, this step may be omitted.
- Determine the APC-V ratios for at least 30 plasma samples from healthy individuals in the age range 20-65 years. Include Control Plasma Level 1 and Level 2 for assay validation.
- Verify that the APC-V ratios for the Control Plasmas are within their specified ranges.
- Calculate the median APC-V ratio.
- Calculate the factor V related APC resistance cut-off value as 0.8 times the median APC-V ratio when below 2.8 and as 0.75 times the median APC-V ratio when 2.8 or higher.
- The APC-V ratio for Control Plasma Level 1 should be within the normal range. The APC-V ratio for Control Plasma Level 2 should be below the cut-off value.

Sensitivity

Coatest APC Resistance V provides 100% sensitivity for FV:Q506 as determined on Thrombolyzer (n = 447), ACL (n = 295), ST-4 (n = 248) and mL/Electra (n = 50).

Limitation/interfering factors

No significant differences are obtained between fresh and frozen samples. The sensitivity and specificity for the factor V:Q506 mutation on analysis of plasma from OAC patients is not affected by the INR value. The prescribed assay procedure allows for the analysis of plasma from heparinized patients at heparin levels < 1 IU/mL plasma (unfractionated and low molecular weight heparins). Although the 1+4 predilution strongly decreases interferences, it can not be excluded that analysis of plasma from patients with high inhibitor activity (e.g. phospholipid antibodies) may give an abnormal APT time and thus possibly misleading results. In such cases, increasing the dilution factor (e.g. 1+9 or 1+19) may correct the test result. As for any APTT-based assay, care should be taken to avoid contact activation of samples since this may lead to activation of FVIII and FV.

Reference values

The APC-V ratios obtained from analysis of plasmas from 61 healthy individuals on ACL and ST-4 and from 390 healthy individuals on Thrombolyzer were in the range 2 - 3.5. No difference is found between sexes.

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Bibliography / Literatur / Bibliografia / Bibliographie / Bibliografía

- Dahlbäck B et al. Familial thrombophilia due to a previously unrecognized mechanism characterized by poor anti-coagulant response to activated protein C: Prediction of a cofactor to activated protein C. Proc Natl Acad Sci 90, 1004-1008. (1993)
- Svensson P J and Dahlbäck B. Resistance to activated protein C as a basis for venous thrombosis. N Engl J Med 330, 517-522 (1994).
- Bertina R M et al. Mutation in blood coagulation factor V associated with resistance to activated protein C. Nature 369, 64-67 (1994).
- Voorberg J et al. Association of idiopathic venous thromboembolism with single point mutation at Arg506 of factor V. Lancet 343, 1535-1536 (1994).
- Greengard J S et al. Activated protein C resistance caused by Arg506Gln mutation in factor Va. Lancet 343, 1361-1362 (1994).
- Williamson D et al. Factor V Cambridge: A new mutation (Arg306->Thr) associated with resistance to activated protein C. Blood 91, 1140-1144 (1998).
- Jorquera J I et al. Modified test for activated protein C resistance. Lancet 344, 1162-1163 (1994).
- Trossaert M et al. Modified APC resistance assay for patients on oral anticoagulants. Lancet 344, 1709 (1994).
- Svensson PJ et al. Evaluation of original and modified APC-resistance tests in unselected outpatients with clinically suspected thrombosis and in healthy controls. Thromb Haemost 77, 332-335 (1997).
- Freyburger G et al. Proposal for objective evaluation of the performance of various functional APC-resistance tests in genotyped patients. Thromb Haemost 78, 1360-1365 (1997).
- Hall C et al. Evaluation of a modified APTT-based method for determination of APC resistance in plasma from patients on heparin or oral anticoagulant therapy. Thromb Res 89, 203-209 (1998).
- National Committee for Clinical Laboratory Standards. Specifications for reagent water used in the clinical laboratory. NCCLS Approved Standard: ASC-3.
- Freyburger G et al. Response to activated protein C upon storage of whole blood and plasma. Thromb Res 93, 89-95 (1999).
- National Committee for Clinical Laboratory Standards. Collection, transport and processing of blood specimens for coagulation testing and performance of coagulation assays, NCCLS Document H21-A3; vol. 11 No. 23.
- National Committee for Clinical Laboratory Standards. EP5-T2. User evaluation of clinical chemistry devices. US Patent 5,443,960; EP 0 608 235; EP 0 690 991; Australia 666 484; 690 535; Japan 2562000; Canada 2,119,761; New Zealand 261 190. APC is a trademark of Chromogenix AB ©2022 Instrumentation Laboratory. All rights reserved.

Symbols used / Verwendete Symbole / Símbolos utilizados / Symboles utilisés / Simboli impiegati

IVD In vitro diagnostic medical device In-vitro Diagnostikum De uso diagnóstico in vitro Dispositif médical de diagnostic in vitro Per uso diagnostico in vitro	LOT Batch code Chargen-Bezeichnung Identificación número de lote Designation du lot Numero del lotto	CONTROL Control Kontrollen Control Contrôle Controllo	 Use by Verwendbar bis Caducidad Utilisable jusqu'à Da utilizzare prima del	 Temperature limitation Festgelegte Temperatur Temperatura de Almacenamiento Températures limites de conservation Limiti di temperatura	 Consult instructions for use Beilage beachten Consultar la metódica Lire le mode d'emploi Vedere istruzioni per l'uso	 Change(s) to section since previous revision Der Abschnitt hat sich seit der letzten Überarbeitung geändert La sección ha cambiado desde la revisión anterior La section a changé depuis la révision précédente La sezione è cambiata rispetto alla revisione precedente	 Biological risks Biologisches Risiko Riesgo biológico Risque biologique Rischio biologico	 Manufacturer Hergestellt von Fabricado por Fabricant Prodotto da	EC REP Authorised representative Bevollmächtigter Representante autorizado Mandatata Rappresentanza autorizzata
---	--	---	--	--	---	--	---	--	---

Control Plasma Level 1

Control Plasma Level 2

All Instruments	CaCl ₂ (s)	APC/CaCl ₂ (s)	APC-V ratio
LOT			
			

All Instruments	CaCl ₂ (s)	APC/CaCl ₂ (s)	APC-V ratio
LOT			
			

Expected values **ENGLISH - Insert revision 05/2022**

Different clotting times will be obtained with different types of instruments depending on the clot detection principle. The following ranges cover the expected clotting times when performing analysis of Control Plasma Level 1 with the above batch number using COATEST[™] APC[™] Resistance V / V-S.

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SIZE: 11 x 17" (280 x 432 mm.).

PRINT: Front/Back.

PRINT COLOR: Sheet 1: Front - Top Band Orange PMS 137, Red PMS 179, it is acceptable to print black in lieu of Red PMS 179, all remaining type in black.
Back - Red PMS 179, it is acceptable to print black in lieu of Red PMS 179, all remaining type in black.
Sheet 2: Front - Top Band Orange PMS 137, Red PMS 179, it is acceptable to print black in lieu of Red PMS 179, all remaining type in black.
Back - Red PMS 179, it is acceptable to print black in lieu of Red PMS 179, all remaining type in black.

SPECIAL REQUIREMENTS: Folded

CARTON SIZE: WHITE BOX PCS0003022610 manual pack must have barcodes on the back of pages 2 and 4. Sheet 1 must have "A" at the end of the document part number, i.e., 000304123A. Sheet 2 must have "B" at the end of the document part number, i.e., 000304123B