

Activation kinetics of FIX by specific aPTT reagent explain discrepancies observed in one-stage assay for N9-GP

Rosén P, Bryngelhed P, Rosén S
Rossix AB, Mölndal Sweden



Disclosures for S. Rosén

In compliance with COI policy, ISTH requires the following disclosures to the session audience:

**S. Rosén is Director of Scientific and Medical Affairs
at Rossix AB, Mölndal, Sweden**

This work was financed by NovoNordisk A/S

Introduction

- Multicenter study 2013 directed by NIBSC revealed large discrepancies in assigned FIX potencies among one-stage (OS) clotting methods on assay of new "long-life" rFIX such as N9-GP
- APTT reagents with Silica as activator showed in general much higher assigned FIX potencies than Ellagic acid based reagents with N 9-GP potency ranges:

Ellagic acid (EA) based reagents	2 – 11 IU/mL
Silica (Si) based reagents	3 – 293 IU/mL, median 44 IU/mL

- New sensitive methods for FXIa and FIXa allow evaluation of activation profiles during contact activation and after addition of CaCl_2 in OS methods
- A study was performed on generation of FXIa and FIXa vs time by subsampling in OS APTT based methods.

Materials

FIX preparations

4th IS FIX Conc (07/182)

BeneFIX (Pfizer)

N9-GP (NovoNordisk)

APTT reagents

APTT SP (IL)

SynthAFax (IL)

Actin FS (Siemens)

Chromogenic methods

Rox Factor IX (Rossix)




Rox Factor XIa, modified (Rossix)

Rox FIX-A, modified (Rossix)




Study Design

- Subsampling studies were performed in OS methods using the single sample dilution 1:20 of 1 IU/mL solutions prepared in FIX deficient plasma.
- Each assay series comprised one APTT reagent and all three FIX preparations. Subsampling was performed simultaneously from the three activation mixtures. Two independent runs were made.
- **Contact activation**
Subsampling up to 6 min of 10 μ L into 6.0 mL ice-cold MES buffer pH 5.7.
Assay within about 5 min for determination of FXIa.
- **Coagulation triggering after addition of CaCl₂:**
Subsampling until clot formation (35 - 50 s) of 50 μ L into 50 μ L stop solution with EDTA, Aprotinin and CTI. Mean subsampling time interval about 7 s.
The t=0 sample was withdrawn 10-15 s before addition of CaCl₂.
Assay within about 5 min for determination of FIXa.


Subsampling during contact activation – FXIa generation



Subsampling after addition of CaCl_2 – FIXa generation



Activation kinetics – Rox Factor FIX



Conclusions related to assay of N9-GP

- Differences in generation of FIXa was identified as the cause of discrepancy in FIX potency assignment with OS methods utilizing different aPTT reagents.
- With APTT SP, 50-60% of FIXa was generated already during contact activation, resulting in high FXa generation.
- FIXa generation was impaired with Actin FS.
- The Rox Factor IX kit showed uniform FXa generation vs time for all three FIX preparations.
- PEGylated BeneFIX shows the same FIX activation pattern as N9-GP with all three APTT reagents and thus FIXa formation during contact activation with APTT SP and impaired FIXa formation with Actin FS. Mere addition of PEG had no effect.

Acknowledgements

Dr Mirella Ezban and Dr Egon Persson at Novo Nordisk are acknowledged for valuable suggestions.