

PL CHIP FOR T-TAS[®] 01

A Flow Chamber System to Measure
Primary Hemostatic Function

**T-TAS**[®]
Total Thrombus-formation Analysis System



WHERE CLINICALLY SIGNIFICANT BLEEDING IS POSSIBLE...



Surgery



Interventional Procedures



Trauma



Comprehensive Hemophilia Treatment Center

...it's important to know if the patient's primary hemostasis is normal.

PRIMARY HEMOSTASIS CAN HAVE A MAJOR IMPACT ON:

- Hospital resources and blood product management.^{1,2}
- Pre-surgical wait times and surgical delays.^{3,4}
- Treatment decisions and outcomes.^{1,5}



MEASURING PRIMARY HEMOSTASIS CAN HELP EVALUATE IF:

- The patient might benefit from platelet transfusion.
- Pre-surgical primary hemostasis has been restored.
- Active bleeding is associated with impaired platelet activity.
- Platelet activity is significantly impaired in association with interventional procedures.
- Antiplatelet therapy is showing the desired effect



INTRODUCING THE PL CHIP FOR T-TAS 01

The first ex-vivo flow chamber model of in-vivo primary hemostasis available for clinical use.



Magnified view

THE TECHNOLOGY



Uses physiological arterial shear stress to assess platelet thrombus formation (primary hemostasis) in whole blood.

Contains flow chamber with 26 collagen-coated microcapillaries.



Generates results within 40 minutes of sample collection.

Accommodates two samples per chip.

The test measures primary hemostatic function as the area under the pressure-time curve (AUC). An AUC < 260 suggests abnormal primary hemostatic function.



NORMAL, no primary hemostatic defect identified. Platelet aggregation increases pressure. Flow impeded.




ABNORMAL, abnormal primary hemostatic function. Lack of platelet aggregation. Flow uninterrupted.



CLINICAL PERFORMANCE

The T-TAS 01 PL assay cutoff is sensitive and specific to abnormal primary hemostasis.^{6,9}

AUC RESULTS BY PATIENT TYPE



ASA, aspirin monotherapy, DAPT, dual antiplatelet therapy, VWD, von Willebrand disease, GT, Glanzmann's thrombasthenia

THE T-TAS 01 PL ASSAY SELECTIVELY MEASURES PRIMARY HEMOSTATIC FUNCTION

Not influenced by secondary hemostatic function¹⁰

Significant correlation with intensity of antiplatelet therapy, VWF antigen, VWF activity, and factor VIII activity levels^{10,11}

Dose-response relationship with medications known to affect primary hemostatic function¹⁰

PARAMETER	VALUE	95% CI
Normal*	95.8%	91.1 - 98.0%
ASA**	68.4%	55.5 - 79.0%
Clopidogrel DAPT**	100.0%	81.5 - 100.0%
Prasugrel DAPT**	100.0%	78.2 - 100.0%
Ticagrelor DAPT**	100.0%	76.8 - 100.0%
VWD**	72.0%	50.6 - 87.9%
GT**	100.0%	43.9 - 100.0%

*Negative agreement, **Sensitivity

PERFORMANCE CHARACTERISTICS

Reference Range 270.0 – 447.7 AUC

Cutoff AUC < 260

Assay Imprecision SD ≤ 39 AUC
2.8% CV (high AUC)
19.6% CV (low AUC)

Reportable Range 0.3 – 467.7 AUC


SPECIFICATIONS

Sample Type	BAPA-anticoagulated whole blood
Sample Volume	320 µL
Test Duration	≤ 10 minutes
Sample Stability	Up to 6 hours
Reagent Storage	PL Chip: 4-8 °C BAPA tube: 15-30 °C
Open Pouch Stability	Up to 8 hours
Quality Control	Internal QC External QC (donor blood samples)
Instrument Dimensions (L x W x H)	14.2" x 12.6" x 9.7" (36 x 32 x 24.7 cm)
Instrument Weight	13.2 lbs (6.0 kg)
Operating Conditions	Temperature: 68-86 °F (20-30 °C) Relative Humidity: 20-80%
On-board Storage	Thousands of results

ORDERING INFORMATION

ITEM	CATALOG NUMBER
T-TAS 01 Total Thrombus Formation Analysis System Instrument	18001
PL Chip for T-TAS 01 20 Chips	18002
PL Chip Reservoir Set for T-TAS 01 100 sets	18003
BAPA Tube for T-TAS 01 (3 mL) 50 tubes	18004

**FOR ADDITIONAL INFORMATION
AND TO VIEW PRODUCT VIDEOS,
PLEASE VISIT WWW.T-TAS.INFO**



The T-TAS 01 PL chip is intended for use in the clinical laboratory for the analysis of the platelet thrombus formation process (primary hemostatic function) in patients age 21 and older with a history of conditions associated with impaired primary hemostatic function or use of antiplatelet therapy.

Physicians should use their clinical judgment and experience when deciding how to diagnose and treat patients and in the use of the PL Chip for T-TAS 01 in the treatment of the patient. Please refer to the PL Chip for T-TAS 01 Package Insert and T-TAS 01 User's Manual for full instructions on sample collection and handling, and all other test procedures.

REFERENCES

1. Bobbitt J., et al. Evaluating the effectiveness of pre-operative platelet inhibition testing to reduce costs and pre-operative length of stay. *Applied Nursing Research* 39 2018: 241-243.
2. Orsini S., et al. Bleeding risk of surgery and its prevention in patients with inherited platelet disorders. The Surgery in Platelet disorders And Therapeutic Approach (SPATA) study. *Haematologica*. 2017;102:1192-203.
3. Ferraris V., et al. 2012 update to the Society of Thoracic Surgeons guideline on use of antiplatelet drugs in patients having cardiac and noncardiac operations. *The Annals of thoracic surgery*. 2012: 1761-1781.
4. Khudairy A., et al. Withholding Clopidogrel for 3 to 6 versus 7 Days or More before Surgery in Hip Fracture Patients. *Journal of Orthopaedic Surgery*, 2013, 146-150.
5. Kibler M., et al. Primary Hemostatic Disorders and Late Major Bleeding After Transcatheter Aortic Valve Replacement. *Journal of the American College of Cardiology* 2018: 2139-2148.
6. Oimatsu Y., et al. Total Thrombus-formation Analysis System Predicts Periprocedural Bleeding Events in Patients With Coronary Artery Disease Undergoing Percutaneous Coronary Intervention. *J Am Heart Assoc*. 2017;6(4):e005263.
7. Nogami K., et al. Assessing the clinical severity of type 1 von Willebrand disease patients with a microchip flow-chamber system. *Journal of Thrombosis and Haemostasis*. 2016: 667-674.
8. Ågren A., et al. Monitoring of coagulation factor therapy in patients with von Willebrand disease type 3 using a microchip flow chamber system. *Thromb Haemost*. 2017 Jan 5;117(1):75-85.
9. Ogiwara K., et al. "Comprehensive evaluation of haemostatic function in von Willebrand disease patients using a microchip-based flow chamber system." *Haemophilia*. 2015: 71-80.
10. PL Chip for T-TAS 01 package insert.
11. Arima Y., et al. Assessment of platelet-derived thrombogenicity with the total thrombus-formation analysis system in coronary artery disease patients receiving antiplatelet therapy. *J Thromb Haemost*. 2016 Apr;14(4):850-9.

© 2020 FUJIMORI KOGYO CO., LTD. T-TAS, Zacros, and related logos are registered trademarks of FUJIMORI KOGYO CO., LTD. All rights reserved. PN: TP-SA-004 rev. 4



Total Thrombus-formation Analysis System



www.diapharma.com
info@diapharma.com

FOR USA