

# Ceveron Wash Solution



## English

### Intended Use

Ceveron Wash Solution is a wash solution for automated Ceveron analyzers.

### Summary and principle

The Ceveron Wash Solution is used for cleaning the pipetting system on automated Ceveron analyzers to avoid carry over.

### Reagents

Ceveron Wash solution: buffered aqueous solution containing urea, acids, detergents and preservatives.

REF	CONTENT
9820320-RUO	1 x 30 mL
9820328-RUO	6 x 30 mL

### Precautions and warnings

For research use only.

Universal precautions for the use of chemicals and potentially biohazardous substances must be applied.

Handle waste as potentially biohazardous material and dispose according to accepted laboratory instructions and procedures.

Discard all material in safe and acceptable manner in compliance with relevant local disposal regulations.

Safety data sheet available for professional user on request.

Avoid foam formation in all reagents and sample types (specimens, calibrators and controls).

Hazardous components are classified in accordance with the Regulation (EC) 1272/2008 as follows:



Warning

**Hazardous ingredient:** Acetic acid

### Hazard Statement

H315 Causes skin irritation.

H319 Causes serious eye irritation.

### Precautionary Statement

P264 Wash hands thoroughly after handling.

### Reagent handling

The reagent is ready to use.

Do not use the reagent if you observe any change in appearance of components included in the kit or if you observe any damage in the packaging materials.

Visible turbidity and flocculation are signs of microbial contamination.

Avoid foam formation.

### Reagent storage and stability

Store at 2-25 °C.

Unopened reagents are stable until the stated expiration date indicated on the label.

Store reagents upright in their packaging.

Stability of the opened reagent in original vial:

at 2-25 °C (Closed cap)	8 weeks
at 18-25 °C (Open vial)	5 days

Do not freeze.

Considering the numerous possible combinations of storage conditions (partly at 18-25 °C, partly at 2-8 °C), each laboratory should establish its own stability durations according to its practices if the reagent is stored partly at 2-8 °C. These durations should not exceed the figures mentioned above which have been determined under controlled conditions.

### Procedure

#### Materials provided

See "Reagents" section.

#### Materials required (but not provided)

- General laboratory equipment
- Reagents for performing coagulation assays
- Coagulation analyzer such as Ceveron 100 series

Please note that the application on other analyzers can be validated by the manufacturer under their responsibility as long as the intended purpose and performance are not modified.

### Test procedure

For more details on the use of Ceveron Wash Solution please refer to the user manual of the instrument to be used.

With analyzers of the Ceveron 100 series, after opening of the solution bottle, load it into the instrument according to the recommendations of the instrument manual in positions R4 and R9.

### Symbols

The following symbols and signs are used in addition to those listed in the ISO 15223-1 Standard:

Global Trade Item Number

Content

Research use only

# Ceveron Wash Solution



---

Information and or pictures contained in this document might be protected by copyrights and other intellectual property rights.

Technoclone's logos and product names are registered trademarks.

Major additions, deletions or changes are indicated by lines in the margin.

A point (period/stop) is always used in this instructions for use as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

This document is available in several languages. The translations have been done using the master document in English. In the event of doubts or discrepancies, the wording in the master document in English shall take precedence.

© Technoclone Herstellung von Diagnostika und Arzneimitteln GmbH, 2022. All rights reserved.

RUO



Technoclone Herstellung von Diagnostika und Arzneimitteln GmbH

[www.technoclone.com](http://www.technoclone.com)



+43 1 86 373 0

EN R01\_2023-12