## Special Test Control Level 1 - 0020011000

### Intended use

For the quality control in the low abnormal range of the chromogenic tests (Antithrombin, Plasminogen, Plasmin Inhibitor, Protein C and Factor VIII) and free Protein S assay performed on the IL Coagulation Systems. For the Quality Control of von Willebrand Factor assays (Antigen and Activity) in the normal range on the IL Coagulation Systems.

### Summary and principle

The Special Test Control Level 1 is prepared from human citrated plasma from healthy donors and modified, by means of a dedicated process, to simulate an abnormal coagulation sample.

It is intended for the assessment of precision and accuracy of the chromogenic tests and the free Protein S assay in the range of 50 - 60% activity.

von Willebrand Factor values (Antigen and Activity) are within the normal range.

Use of normal and abnormal controls is recommended for a complete quality control program.

### Composition

The **HemosIL Special Test Control Level 1** consists of:

- Special Test Control Level 1** (Cat. No. 0020011100): 10 x 1.0 mL vials of lyophilized human plasma containing buffer and stabilizers. No preservatives are included.

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### Precautions and warnings:

The material in this product was tested by FDA approved test methods and found nonreactive for Hepatitis B Surface Antigen (HBsAg), Anti-HCV and HIV antibodies. Handle as if potentially infectious.1

**Hazard class:** None

**Hazard statements:** None

**Precautionary statements:** None

**Supplemental hazard information:**

≈ 33.1% of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and unknown hazard to the aquatic environment.

This product is For *in vitro* Diagnostic Use.

### Preparation

Dissolve the contents of each vial with 1.0 mL of CLSI CLRW Type water or equivalent.2 Replace the stopper and swirl gently. Make sure of the complete reconstitution of the product. Keep the control at 15-25°C for 30 minutes and mix gently before use. Do not shake. Avoid foam formation.

### △ Reagent storage and stability

Unopened control is stable until the expiration date shown on the vial when stored at 2-8°C.

Stability after reconstitution: 8 hours at 15-25°C or 24 hours at -20°C in the original vial. Frozen controls may be thawed at 37°C and gently mixed before use. Do not refreeze.

For optimal stability remove control from the system and store it at 2-8°C in the original vials.

-  Low relative humidity is associated with increased evaporation of uncapped reagents and controls, which may decrease on-board stability. For optimum on-board stability, laboratory temperature and humidity should be controlled.

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### ENGLISH - Insert revision 03/2022

### Instrument/test procedure

After reconstitution the Special Test Control Level 1 should be handled in the same manner as fresh citrated plasma.

### Traceability of calibrators and control materials

The reported values were determined over multiple runs on IL Coagulation Systems using a specific lot of reagent and against a Calibration Plasma House Standard which is traceable to the current International Standards, identified in the acceptance range table.

For the tests where International Standards are not available, these parameters (i.e. Plasminogen and Plasmin Inhibitor) have been assigned against a House Standard which is traceable to a frozen normal plasma pool of 100 donors.

### Limitations

This product is designed as an abnormal control for monitoring chromogenic and free Protein S assays and as a normal level for the Quality Control of von Willebrand Factor assay.

The control is subjected to the limitations of the assay system. Deviations in control values may indicate possible problems with one or more components in the test system.

### Performance characteristics/expected values

Refer to specific reagent inserts for performance characteristics. The reported ranges were determined over multiple runs on IL Coagulation Systems using specific lots of IL reagents. The mean of the control range determined in your laboratory may vary due to the lot of reagent used.

Due to differences in reagents and instrumentation, each laboratory should establish its own Target Value and Acceptance Range (mean and standard deviation). However, any properly functioning coagulation system should yield mean values within the Acceptance Range on the package insert.

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- Richmond JY, McKinney RW eds. Biosafety in Microbiological and Biomedical Laboratories, U.S. Dept. of Health and Human Services, Public Health Service, 4<sup>th</sup> Edition, 1999.
- Clinical and Laboratory Standards Institute. Preparation and Testing of Reagent Water in the Clinical Laboratory; Approved Guideline. Fourth, Edition, CLSI Document C3-A4; Vol. 26 No.22.

The Instrumentation Laboratory logo, ACL, ACL AcuStar, ACL Elite, ACL TOP, HemoCell, HemoHub, HemosIL, ReadiplasTin, RecombiPlasTin, SynthAFax and SynthASil are trademarks of Instrumentation Laboratory Company and/or one of its subsidiaries or parent companies and may be registered in the United States Patent and Trademark Office and in other jurisdictions. All other product names, company names, marks, logos and symbols are trademarks of their respective owners.

Patent: <https://www.instrumentationlaboratory.com/us/en/trademarks-and-patents> ©2022 Instrumentation Laboratory. All rights reserved.

**Symbols used / Verwendete Symbole / Símbolos utilizados / Symboles utilisés / Simboli impiegati / Símbolos utilizados**

<b>TARGET VALUE</b>	Target Value / Referenzwert-Zielwert / Valor de referencia / Valeur cible / Valore di riferimento / Valor-Alvo
<b>ACCEPTANCE RANGE</b>	Acceptance range / Referenzbereich / Rango esperado / Valeur usuelles / Intervallo di accettabilità / Intervalo de aceitação
<b>WHO STD</b>	WHO Standard Code / WHO Standard Code / Código del Estándar OMS / Code standard OMS / Codice OMS / Código Padrão WHO
<b>PRODUCT NAME</b>	IL Product name / IL Produktname / Nombre del kit IL / Nom du kit IL / Kit IL / Designação do kit IL
<b>UNITS</b>	Measurement units / Messeinheiten / Unidades de medida / Unités de mesure / Unità di misura / Unidades de Medida

<b>IVD</b> <i>In vitro</i> diagnostic medical device <i>In-vitro</i> Diagnostikum De uso diagnóstico <i>in vitro</i> Dispositif médical de diagnostic <i>in vitro</i> Per uso diagnóstico <i>in vitro</i> Dispositivo médico para utilização em diagnóstico <i>in vitro</i>	<b>LOT</b> Batch code Chargen-Bezeichnung Identificación número de lote Designation du lot Numero del lotto Número de lote	 Use by Verwendbar bis Caducidad Utilizable jusqu'à Da utilizzare prima del Data limite de utilização	 Temperature limitation Festgelegte Temperatur Temperatura de Almacenamiento Températures limites de conservation Limiti di temperatura Limite de temperatura	 Consult instructions for use Beilage beachten Consultar la metodica Lire le mode d'emploi Vedere istruzioni per l'uso Consultar as instruções de utilização	 Change(s) to section since previous revision Der Abschnitt hat sich seit der letzten Überarbeitung geändert La sección ha cambiado desde la revisión anterior La section a changé depuis la révision précédente La sezione è cambiata rispetto alla revisione precedente A seção mudou desde a revisão anterior	 Caution, consult accompanying documents Achtung, Begleitdokumente beachten Atención, veir instrucciones de uso Attention, voir notice d'instructions Attenzione, vedere le istruzioni per l'uso Atenção, consulte a documentação incluída	<b>CONTROL</b>  Control Kontrollen Control Contrôle Controllo Controllo	 Biological risks Biologisches Risiko Riesgo biológico Risque biologique Rischio biologico Risco biológico	 Manufacturer Hergestellt von Fabricado por Fabricant Prodotto da Fabricado por	<b>EC REP</b> Authorised representative Bevollmächtigter Representante autorizado Mandataire Rappresentanza autorizzata Representante autorizado
---	--	---	---	--	--	--	--	--	---	--

**Special Test Control Level 1 0020011000**

---

Printed Insert Sheet: 000303523

Revision: R09

Issued: 03/2022

C.O.: 538568

**LANGUAGES**

---

ENGLISH

GERMAN

SPANISH

FRENCH

ITALIAN

PORTUGUESE

**TECHNICAL SPECS**

---

PAPER: White paper, 50-60 g/m<sup>2</sup> weight, PCS00400073.

SIZE: 11 x 17" (280 x 432 mm.).

PRINT: Front/Back.

PRINT COLOR: Front - Top band Green Pantone 382,  
all remaining type in black.

Back - All type in black.

SPECIAL REQUIREMENTS: Folded