
TECHNOLEIA® D-Dimer

REF	4847200	TECHNOLEIA® D-Dimer Latex Kit	
REF	4847210	TECHNOLEIA® D-Dimer Latex Kit	
REF	4847220	TECHNOLEIA® D-Dimer Latex Reagent	2 x 12 mL
REF	4847222	TECHNOLEIA® D-Dimer Reaction Buffer	2 x 21 mL
REF	4847230	TECHNOLEIA® D-Dimer Control high	5 x 1 mL
REF	4847232	TECHNOLEIA® D-Dimer Control low	5 x 1 mL
REF	4847234	TECHNOLEIA® D-Dimer Calibrator ~ 3000 ng/ml	2 x 1 mL
REF	4847236	TECHNOLEIA® D-Dimer Calibrator 0 ng/ml	2 x 1 mL
REF	4847238	TECHNOLEIA® D-Dimer Latex Reagent	5 x 4 mL
REF	4847240	TECHNOLEIA® D-Dimer Reaction Buffer	5 x 7 mL

symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Klíčova slova / Značenje simbola			
	manufacturer / Hersteller / fabbricante / fabricante / fabricant / Tillverkaren / Fabrikanten / Producent / Κατασκευαστής / Производител / Производитель / výrobce / Proizvodač		expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration / utgångsdatum / udløbsdato / Utløpsdato / Ημερομηνία λήξης / срок на годност / datum expirace/ срок годности / datum expirace / Rok trajanja
	storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage / lagringstemperatur / opbevaringstemperatur / Oppbevaringstemperatur / θερμοκρασία αποθήκευσης / съхранение на / teplota skladování / температура хранения / teplota skladování / Temperatura lagerovanja		consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation / se användarinstruktioner / følg brugsvejledning / Følg bruksanvisningen / συμβουλευθείτε τις οδηγίες για τη χρήση / прочетете инструкцията за работа / potfeba fídit se instrukcemi / перед использованием читайте инструкцию / sledujte návod k použití / Pročitaj upustvo pre upotrebe
			determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations / bestämningar / bestemmelse / Bestemmelse / προδιορισμοί / брой тестове / stanovení / определний / ročet stanovení / Definicija
AQUA	distilled water / destilliertes Wasser / acqua distillata / agua destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / дестилирана вода / destilovaná voda / дестилированная вода / destilovaná voda / Destilisaná Voda	RUO	For research use only
CAL	Calibrator / Kalibrator / Calibratore / calibrador / calibrador / calibreur / Kalibrator / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrator	LOT	lot / Charge / lotto / lote / lote / lot / sats / serie / Parti / παρτίδα / партида номер / šarže / lot / šarže / Serija
DIL	dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o dissolver / diluir ou dissoudre em / diluer ou dissoudre dans / späð eller upplös i / fortyndes eller opløses i / Fortyndes eller oppløses i / αραιωση ή διάλυση σε / разтворете или разрежете с / zredit alebo rozpustit v / разбавить или растворить в / nafedte nebo rozpustte v / razrediti ili rastvoriti u	REF	catalogue number / Katalognummer / numero di catalogo / número de catálogo / número de referência / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / каталожен номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj

PRODUCT DESCRIPTION

INTENDEND USE

Reagents for the quantitative determination of D-Dimer by "latex enhanced immunoassay".

COMPOSITION

TECHNOLEIA® D-Dimer Latex Kit 4847200 (150 T.) / 4847210 (50 T.) contains:

mL	Reagent
1 x 12 / 4	TECHNOLEIA® D-Dimer Latex reagent
1 x 21 / 7	Reaction buffer
1 x 8	NaCl solution 0.9 %
1 x 1	TECHNOLEIA® D-Dimer Calibrator 0 ng/mL
1 x 1	TECHNOLEIA® D-Dimer Calibrator ~3000 ng/mL

TECHNOLEIA® D-Dimer Latex reagent 4847220 and 4847238 contains:

mL	Reagent
2 x 12 / 5 x 4	TECHNOLEIA® D-Dimer Latex reagent

TECHNOLEIA® D-Dimer Reaction buffer 4847222 and 4847240 contains:

mL	Reagent
2 x 21 / 5 x 7	TECHNOLEIA® D-Dimer Reaction buffer

MATERIAL REQUIRED (not supplied with the kit)

- Pipettes - Distilled water - Controls

REF	4847230	TECHNOLEIA® D-Dimer Control high	5 x 1 mL
REF	4847232	TECHNOLEIA® D-Dimer Control low	5 x 1 mL

WARNING AND PRECAUTIONS

- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.
- Controls and calibrators are prepared from human blood and each single plasma used for this lots are HBsAg, HIV 1/2 Ab and HCV Ab negative (see package label and vial label).

STABILITY AND STORAGE

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8 °C.

Stability:

Reagents	Storage	Stability
Sample	+ 2 ... 8 °C	8 hours
	- 20 °C *	1 month
Latex reagent	+ 2 ... 8 °C	4 weeks
	+ 15 °C	1 week
Reaction buffer	+ 2 ... 8 °C	4 weeks
	+ 15 °C	1 week
NaCl Solution	+ 2 ... 8 °C	4 weeks
	+ 15 °C	1 week
D-Dimer Calibrator 0 ng/mL D-Dimer Calibrator ~ 3000 ng/mL D-Dimer Control High D-Dimer Control Low after reconstitution	room temperature + 2 ... 8 °C + 15 °C - 20 °C *	8 hours 1 month

* Frozen sera should be rethawed quickly at 37°C in a water bath.

TEST PROCEDURE

SAMPLES

Preparation of plasma samples

Plasma separation: Mix 9 parts of venous blood and 1 part of Sodium Citrate Solution (0.11 mol/L) and centrifuge for 15 minutes at a RCF of at least 2500 g (corresponding to DIN 58905). The samples should be used as fresh as possible or should be stored before usage at -20°C. (see also STABILITY AND STORAGE).

PREPARATION OF REAGENT

TECHNOLEIA® D-Dimer Latex Reagent and Reaction buffer are ready for use. The latex may sediment during storage

Mix thoroughly before use.
If kept in the analyzer, the reagent should be gently mixed daily.

The TECHNOLEIA® D-Dimer Calibrators and the TECHNOLEIA® D-Dimer Control high and low are lyophilized and have to be reconstituted with 1 mL of the NaCl solution included in the pack at least 10 minutes before use.

PERFORMANCE OF THE TEST

Application sheets for auto analyzers are available on request.

A monoclonal anti D-Dimer antibody coupled to latex particles forms high molecular weight immunocomplexes with D-Dimer in the test sample. The resulting change in turbidity of the sample is detected photometrically by turbidimetric and nephelometric analysers and is directly proportional to the antigen concentration. These absorbance changes are converted into antigen concentration using a reference curve constructed with standards of known concentrations. Depending on the auto analyzer used 10 µL sample, 140 µL Reaction buffer and 80 µL Latex reagent are used per determination.

Measuring wave lengths: 546-800 nm

CEVERON

Technoclon provides application sheets for Ceveron® alpha. The application sheets contain analyser/assay specific handling and performance information which may differ from that provided in this instruction for use. In this case the information contained in the application sheet supersedes the information in this instruction for use. Please consult the instruction manual of the Ceveron® alpha.

ANALYSES RESULTS

REFERENCE RANGE

Elevated levels are found in patients with confirmed deep venous thrombosis (DVT), pulmonary embolism (PE), DIC, and trauma. D-Dimer levels rise during pregnancy and high levels are associated with complications. Samples with a D-Dimer concentration > 250 ng/mL are considered as **pathologic**.

The probability of **absence of thrombotic** events is high if the concentration values are < 135 ng/mL.

Samples **above 3000 ng/mL** (depending on analyzer) should be **diluted** with Saline Solution and re-assayed.

CALCULATION OF THE RESULTS

TECHNOLEIA® D-Dimer Calibrator ~3000 ng/mL is used to construct a reference curve. Dilutions of the standard are made with TECHNOLEIA® D-Dimer Calibrator 0 ng/mL or with NaCl solution (see instrument application sheets).

Example:

	Sample	D-Dimer Calibrator 0 ng/mL	Resulted dilution*
Cal. 1	D-Dimer CAL ~3000 ng/mL	0 µL	3000 ng/mL
Cal. 2	125 µL Cal. 1	125 µL	1500 ng/mL
Cal. 3	125 µL Cal. 2	125 µL	750 ng/mL
Cal. 4	125 µL Cal. 3	125 µL	375 ng/mL
Cal. 5	125 µL Cal. 4	125 µL	187 ng/mL

* The values are lot specific (only example)

Not used Calibrators and Controls can be frozen in portions. The number of calibrator points can vary depending on the instrument and are described in the protocols of the respective auto analyser.

STANDARDISATION

No international calibrator is available for the standardization of D-Dimer. Therefore the D-Dimer calibrator is calibrated by using a commercial ELISA kit.

LIMITATION AND INTERFERENCES OF THE TEST

Accuracy should be checked in each run with TECHNOLEIA® D-Dimer Control high and TECHNOLEIA® D-Dimer Control low. The control plasma is treated in the same way as patient samples and results should be checked against the acceptable confidence limits supplied on the vial label.

The calibration curve can be stored in the instrument until the next batch change (or every 6 months) but must be checked in each run with the control plasma. A new calibration is recommended, if software changes, or following a major service of either instruments or equipment.

Presence of rheumatoid arthritis factor may result in false-positive results (influence not quantified).

No interference is found from CRP (<50mg/L), bilirubin (<0.27 g/L), haemoglobin (<6.7g/L), or triglycerides (<3.3g/L). Highly lipemic samples should be diluted in saline solution and re assayed.

Human anti-mouse antibody testing was not performed with this product. Results from patients with heterophilic antibody should be interpreted with caution since this test kit contains mouse antibodies and interference may occur resulting in falsely elevated or decreased values.

APPLICATIONS FOR INSTRUMENTS

Application Sheets are available from Technoclon or your local distributor upon request.

PERFORMANCE CHARACTERISTICS

Performance data are given below. Results obtained in individual laboratories may differ.

PRECISION

Reproducibility was determined with different samples (in series and day to day).

The following results were obtained:

	Intra assay		Inter assay	
Sample	Sample 1	Sample 2	Sample 1	Sample 2
n	20	50	10	5
MW (µg/L)	3401.2	293.1	2697.2	304.2
SD (%)	65.148	8.389	79.055	4.357
CV (%)	1.9%	2.9%	2.93%	1.43%

COMPARISON OF METHODS OR CORRELATION

Following correlation was obtained in comparing TECHNOLEIA® D-Dimer (Technoclon) with the Auto-Dimer (Trinity) on a:

Thrombolyzer	n=18	y = 0.9794x + 158.81	R = 0.955
Ceveron® alpha	n=24	y = 0.8273x + 19.383	R = 0.928

LINEARITY

D-Dimer method: 101 - 3250 (Activity µg/L)

LITERATURE

Please contact Technoclon or your local distributor.

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