
# TECHNOCLOT® PT PLUS

For Research Use Only

REF 5003220 TECHNOCLOT® PT Plus 5 x 4 mL

REF 5003222 TECHNOCLOT® PT Plus 10 x 10 mL

## symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Klíčova slova / Značenje simbola

	Manufacturer / Hersteller / fabbricante / fabricante / fabricant / Tillverkaren / Fabrikanten / Produsent / Κατασκευαστής / Производител / Производитель / výrobc / Proizvođač		Expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration / utgångsdatum / udløbsdato / Utløpsdato / Ημερομηνία λήξης / срок на годност / datum expirace/ срок годности / datum expirace / Rok trajanja
	Storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage / lagringstemperatur / opbevaringstemperatur / Opbevaringstemperatur / θερμοκρασία αποθήκευσης / съхранение на / teplota skladování / температура хранения / teplota skladování / Temperatura lagerovanja		Consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation / se användarinstruktioner / følg brugsvejledning / Følg bruksanvisningen / συμβουλευθείτε τις οδηγίες για τη χρήση / прочетете инструкцията за работа / potřeba řídit se instrukcemi / перед использованием читайте инструкцию / sledujte návod k použití / Pročitaj upustvo pre upotrebe
			Determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations / déterminations / bestämmingar / bestemmelser / Bestemmelser / προσδιορισμοί / брой тестове / stanovení / определний / роčet stanovení / Definicija
<b>AQUA</b>	Distilled water / destilliertes Wasser / acqua distillata / agua destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / destilirana voda / destilovaná voda / дистиллированная вода / destilovaná voda / Destilovaná Voda	<b>LOT</b>	Lot / Charge / lotto / lote / lote / lot / sats / serie / Parti / ηάρτια / партида номер / šarže / лот / šarže / Serija
<b>BUF</b>	Reaction buffer / Reaktionspuffer / tampone di reazione / tampón de reacció / Tampão de reação / tampon de réaction / Reaktionsbuffert / Reaktionsbuffer / Reaktionsbuffer / διάλυμα αντιδραστή / Реакционен буфер / Рабочий буферный раствор / Reakční pufr / Reakcioni pufer	<b>MTP</b>	Microtiter plate / Mikrotiterplatte / placa microtiter / microplaca / microplaca / microplaques sensibilisées / Mikrotiterplatta / Mikrotiterplade / mikrotiterplate / πλάκα μικροπιπιλοδότησης / Микротительна плака / Микропланшет / Mikrotitrační destička / Mikrotitracione ploče
<b>CAL</b>	Calibrator / Kalibrator / calibratore / calibrador / calibrador / calibreur / Kalibrator / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrátor	<b>REF</b>	Catalogue number / Katalognummer / numero de catalogo / numéro de catálogo / número de referência / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / каталожен номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj
<b>CONJ</b>	Conjugate / Konjugat / Coniugato / conjugado / conjugado / conjugaté / Konjugerad / Konjugat / Konjugat / συνδεδετικό / Конюгат / Конъюгат / Konjugát / Konjugat	<b>RTU</b>	Ready to use / gebrauchsfertig / pronto all'uso / listo para usar / pronto a usar / prêt à l'emploi / færdig att användas / færdig til brug / klar til bruk/ έτοιμο προς χρήση / Готов за употреба / готов к использованию / k přímému použití / Razrediti ili rastvoriti
<b>CONT</b>	Control / Kontrolle / controllo / control / control / contrôle / Kontroll / Kontroll / Kontroll / διάλυμα ελέγχου / Контрол / Контрольный образец / Kontrola / Kontrola	<b>STOP</b>	Stop solution / Stopplösung / Soluzione di arresto / solución de parada / solução de paragem / solution d'arrêt / Stopplösning / Stop-opløsning / Stoppløsning / διάλυμα παύσης / Стоп раствор / Стоп-раствор / Zastavovací roztok / Stop solucija
<b>DIL</b>	Dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o dissolver / diluir ou dissolver em / diluer ou dissoudre dans / späd eller upplös i / fortyndes eller opløses i / Fortyndes eller oppløses i / αραιωση ή διάλυση σε / растворите или разредете с / zředit anebo rozpustit v / разбавить или растворить в / nafedte nebo rozpustite v / razrediti ili rastvoriti u	<b>SUB</b>	Substrate / Substrat / substrato / substrato / substrato / substrat / Substrat / Substrat / Substrat / υπόστρωμα / Субстрат / Субстрат / Substrát / Substrat
<b>INC</b>	Incubation buffer / Inkubationspuffer / tampone di incubazione / tampón de incubación / tampão de incubação / tampon d'incubation / Inkubationsbuffert/ Inkubationsbuffer/ Vaskebufferkonsentrat / διάλυμα επώασης / Инкубационен буфер / Буфер для инкубации / Inkubační pufr / Inkubacioni pufer	<b>WASH</b>	Washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado / solución de lavado concentrada / tampão de lavagem concentrado / Tampon de lavage concentré / Vattenlösningsskoncentrat / Vaskeopløsningskoncentrat / vaskeløsningskoncentrat / συμπυκνωμένο διάλυμα πλύσης / Концентриран миеш разтвор / Концентрат промывочного раствора / Koncentrát promývachio roztoku / Koncentrat solucije za ispiranje
<b>FI C</b>	For Research Use Only		

# TECHNOCLOT® PT PLUS

Ca-thromboplastin for the determination of the prothrombin time (PT, Quick)

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## PRODUCT DESCRIPTION

### INTENDED USE

TECHNOCLOT® PT Plus is used for the determination of the prothrombin time (Quick test). This screening test is used

- for the control of blood coagulation disorders of the extrinsic system,
- for the surveillance of oral anticoagulant therapy,
- for separate determination of individual coagulation factors of the extrinsic system,
- to check the liver's function of synthesizing coagulation factors in hepatopathies.

A heparin neutralizing agent in the reagent permits the assessment of plasma obtained from samples on conventional heparin therapy (1 - 2 U. hep./mL).

### COMPOSITION

**TECHNOCLOT® PT Plus Reagent** is a standardized Thromboplastin reagent obtained from rabbit brain which is characterized by sensitivity to coagulation factors II, V, VII and X. In addition, the reagent contains a heparin neutralizing agent. **TECHNOCLOT® PT Plus Solvent** is an aqueous solvent containing calcium.

### MATERIAL REQUIRED (not supplied with the kit)

- Pipettes
- Distilled water
- Control Plasmas and Calibrators

REF 5020040 Coagulation Control N\*\* 5 x 1 mL

REF 5011050 Coagulation Control AK\*\* 5 x 1 mL

REF 5010004 AK-Calibrant 4 x 1 mL

\*\* or any other package sizes, special Ceveron® alpha or TECHNOCLOT® Control and Calibration reagents of Technoclone.

### WARNING AND PRECAUTIONS

- For Research Use Only
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.

### STABILITY AND STORAGE

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8°C.

Stability after reconstitution:

+37°C 8 hours	RT* 5 days	+2...8°C 5 days	+12°C (Ceveron® alpha) 48 hours
------------------	---------------	--------------------	------------------------------------

Upon storage, caps should be screwed tightly.

\*=room temperature (+18...25°C)

### TEST PROCEDURE

#### PREPARATION OF PLASMA SAMPLES

Plasma separation: Mix 9 parts of freshly drawn venous blood and 1 part of trisodium citrate solution (0,11 mol/L) and centrifuge for 15 min. at an RCF of at least 2.500 g (acc.DIN 58905 or the CLSI document H21-A5) Keep the plasma at room temperature (up to 4 hours).

#### PREPARATION OF REAGENT

Transfer the entire content of one vial of the solvent (white) into one vial of the lyophilized reagent (grey). Allow the reconstituted reagent to stand at room temperature (+18...25°C) for 30 minutes. Then, swirl the reconstituted reagent gently to obtain a homogenous suspension.

#### PERFORMANCE OF THE TEST

The reagent requires a stirred position on automated analyzers, or it is recommended to turn the vial upside down every 4 hours to obtain a homogeneous suspension again.

#### CEVERON

Technoclone provides Application sheets for Ceveron® alpha. The application sheets contain analyser/assay specific handling and performance information which may differ from that provided in these instruction for use. In this case the information contained in the application sheets supersedes the information in this instruction for use. Please consult the instruction manual of the Ceveron® alpha.

#### MANUAL

Pipette the following substances into a test tube preheated to +37°C:

Citrated plasma	0,1 mL
Incube for 2 min at +37°C	
+ TECHNOCLOT® PT Plus +37°C	0,2 mL

Upon addition of TECHNOCLOT® PT Plus, start the stopwatch or the coagulometer resp. and determine the clotting time.

### ANALYSES RESULTS

#### CALCULATION OF THE RESULTS

The prothrombin time is indicated in seconds, in % of normal (Quick), or in INR. Values can be converted by means of the enclosed reference curve in form of a table, considering the indicated method. Reference curves have to be checked for accuracy by means of control plasma prior to use. If the values are beyond the confidence range given for control plasma, a separate reference curve has to be prepared.

For the preparation of the reference curve Technoclone's calibration plasma or the following dilutions of a normal plasma pool (Coagulation Control N) with isotonic saline solution may be used.

% of normal dilution	100% undil.	50% 1+1	25% 1+3	12,5% 1+7	10% 1+9
----------------------	-------------	---------	---------	-----------	---------

These dilutions are to be determined like samples. The mean values of triple determinations are plotted on a reciprocal curve sheet and joined linearly.

The values indicated in % of normal can be transformed to INR according to the enclosed reference curve table. The indication of the prothrombin time in INR, as recommended by the WHO, is based on the following formula:

$$INR = \left( \frac{PT_{\text{sample plasma}} (\text{sec})}{PT_{\text{normal plasma}} (\text{sec})} \right)^{ISI}$$

#### REFERENCE RANGE

Normal time - approximate value: 10 - 14 sec.

Normal range in % of normal: 70 - > 100 %

Therapeutic range of oral anticoagulants: INR of 2.0 - 4.5 of 20 - 42% of normal, resp.

#### STANDARDISATION

The reference batch is calibrated against the WHO Int. Reference Preparation for rabbit brain thromboplastin reagent (plain), RBT/05.

#### LIMITATION OF THE TEST

Inhibitors of the Lupus type anticoagulant can influence prothrombin time and lead to INRs that do not accurately reflect the true level of anticoagulation.

The choice of anticoagulant (i.e. oxalate instead of citrate) and the condition of the specimen (e.g. hemolyzed, lipemic, parenteral feeding, etc.) may affect results. The latter is particularly true of optical instrumentation measurement of the PT.

#### PERFORMANCE CHARACTERISTICS

Performance data are given below. Results obtained in individual laboratories may differ.

#### PRECISION

Reproducibility was determined with different samples (in series and day to day). The following results were obtained:

Sample	Intra assay		Inter assay	
	Sample 1	Sample 2	Sample 1	Sample 2
n	21	21	21	21
MV INR	1.027	2.740	1.018	2.819
SD	0.010	0.037	0.009	0.027
CV (%)	1.01	1.34	0.97	0.96

#### COMPARISON OF METHODS OR CORRELATION

Following correlation (INR) was obtained in comparing TECHNOCLOT® PT PLUS with:

Recombiplastin 2G y = 1.3414x - 0.4592 R² = 0.9656  
Thromborel S y = 1.2007x - 0.2626 R² = 0.9578

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