# AK-Calibrant

GB

IT

PT

DE

ES

FR

For Research Use Only

REF 5010004 AK-Calibrant

4 x 1 mL

## Symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Ključova slova / Značenje simbola

	Manufacturer / Hersteller / fabbricante / fabricante / fabricante / fabricant / Tillverkaren / Fabrikanten / Producent / Κατασκευαστής / Производител / Производитель / výrobc / Proizvođač		Expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration / utgångsdatum / udløbsdato / Utløpsdato / Ημερομηνία λήξης / срок на годност / datum expirace/ срок годности / datum expirace / Rok trajanja
	Storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage / lagringstemperatur / opbevaringstemperatur / Oppbevaringstemperatur / Βερυοκρασία αποθήκευσης / съхранение на / teplota skladování / температура хранения / teplota skladování / Temperatura lagerovanja		Consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation / se användarinstruktioner / følg bruksvejledning / Følg bruksanvisningen / συμβουλευθείτε τις οδηγίες για τη χρήση / прочетете инструкцията за работа / potěba řídit se instrukcemi / перед использованием читайте инструкцию / sledujte návod k použití / Pročitaj uputstvo pre upotrebu
			Determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations / bestämmingar / bestemmelser / Bestemmelser / προδιορισμοί / брой тестове / stanovení / определенний / роčet stanovení / Definicija
<b>AQUA</b>	Distilled water / destilliertes Wasser / acqua distillata / água destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / дестилирана вода / destilovaná voda / дистиллированная вода / destilovaná voda / Destiliana Voda	<b>LOT</b>	Lot / Charge / lotto / lote / lote / lot / sats / serie / Parti / παρτία / партида номер / šarže / lot / šarže / Serija
<b>BUF</b>	Reaction buffer / Reaktionspuffer / tampone di reazione / tampón de reacció / Tampão de reação / tampon de réaction / Reaktionsbuffert / Reaktionsbuffer / Reaktionsbuffer / διάλυμα αντίδρασης / Реакционен буфер / Рабочий буферный раствор / Reakční pufr / Reakcioni pufer	<b>MTP</b>	Microtiter plate / Mikrotiterplatte / placa microtiter / microplaca / microplaca / microplaques sensibilisées / Mikrotiterplatta / Mikrotiterplade / mikrotiterplate / πλάκα μικροτιτλοδότησης / Микротитърна плака / Микропланшет / Mikrotitrační destička / Mikrotitracione ploče
<b>CAL</b>	Calibrator / Kalibrator / Calibratore / calibrador / calibrador / calibreur / Kalibrator / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrator	<b>REF</b>	Catalogue number / Katalognummer / numero di catalogo / número de catálogo / número de referência / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / каталожен номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj
<b>CONJ</b>	Conjugate / Konjugat / Coniugato / conjugado / conjugado / conjugaté / Konjugerad / Konjugat / Konjugat / συνδέτικό / Конюгат / Конъюгат / Konjugát / Konjugat	<b>RTU</b>	Ready to use / gebrauchsfertig / pronto all'uso / listo para usar / pronto a usar / prêt à l'emploi / færdig att användas / færdig til brug / klar til bruk/ έτοιμο προς χρήση / Готов за употреба / готов к использованию / к рíímému použití / Razrediti ili rastvoriti
<b>CONT</b>	Control / Kontrolle / controllo / control / control / contrôle / Kontroll / Kontroll / Kontroll / διάλυμα ελέγχου / Контрол / Контрольный образец / Kontrola / Kontrola	<b>STOP</b>	Stop solution / Stopplösung / Soluzione di arresto / solución de parada / solução de paragem / solution d'arrêt / Stoppløsning / Stop-opløsning / Stoppløsning / διάλυμα παύσης / Стоп разтвор / Стоп-раствор / Zastavovací roztok / Stop solucija
<b>DIL</b>	Dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o disolver / diluir ou dissolver em / diluer ou dissoudre dans / späð eller upplös i / fortyndes eller opløses i / Fortyndes eller oppløses i / αραιωση ή διάλυση σε / растворите или разрежете с / zředí anebo rozpustit v / разбавить или растворить в / nafedite nebo rozpustite v / razrediti ili rastvoriti u	<b>SUB</b>	Substrate / Substrat / substrato / substrato / substrato / substrat / Substrat / Substrat / Substrat / υπόστρωμα / Субстрат / Субстрат / Substrát / Substrat
<b>INC</b>	Incubation buffer / Inkubationspuffer / tampone di incubazione / tampón de incubación / Incubação de incubação / tampon d'incubation / Inkubationsbuffert / Inkubationsbuffer / Vaskebufferkonsentrat / διάλυμα επώασης / Инкубационен буфер / Буфер для инкубации / Inkubační pufr / Inkubacioni pufer	<b>WASH</b>	Washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado / solución de lavado concentrada / tampão de lavagem concentrado / Tampon de lavage concentré / Vattenlösningkoncentrat / Vaskeopløsningskoncentrat / vaskeløsningskonsentrat / συμπυκνωμένο διάλυμα πλύσης / Концентриран миеш разтвор / Концентрат промышленного раствора / Koncentrát promývaciho roztoku / Koncentrat solucije za ispiranje
<b>IVD</b>	In vitro diagnostic use / in vitro Diagnostikum / diagnostico in vitro / diagnóstico en vitro / diagnóstico in vitro / diagnostic in vitro / för in vitro diagnostik / in vitro diagnostik / in vitro diagnostisk bruk / рíση διαγνωστική εντόβ σωλήνα / за ин vitro диагностика / pro in vitro diagnostiku / использовать для диагностики in vitro / diagnosticky prostředek in vitro / in vitro dijagnostika		

**PRODUCT DESCRIPTION**

**INTENDEND USE**

The AK calibration set can be used for the following purposes:

- For direct INR determination of a patient sample (INR reference curve)
- For establishing laboratory, instrument and reagent specific ISI and PT normal plasma values
- For constructing a reference curve in percentage of the normal.

**COMPOSITION**

The calibration set contains 4 calibration plasmas for standardising the PT test. The AK-Calibrant A is a lyophilised normal plasma pool which is produced from the selected citrated plasmas (Austrian source plasma) of at least 100 blood samples of healthy donors. AK-Calibrant B, C, and D are produced exclusively from the plasma of donors on long-term oral anticoagulant treatment. As a result, the prothrombin complex factors (II, VII, IX, X, Protein C and Protein S) are present at the level normally expected in plasma from patients under oral anticoagulant treatment, together with the Proteins induced by vitamin K absence (PIVKA inhibitors). No preservatives are added to the lyophilised plasmas.

mL	Reagent	Description
1 x 1	AK-Calibrant A	lyophilised normal plasma, INR approx.1.0
1 x 1	AK-Calibrant B	lyophilised AK-plasma, INR about 2.0
1 x 1	AK-Calibrant C	lyophilised AK-plasma, INR about. 3.0
1 x 1	AK-Calibrant D	lyophilised AK-plasma, INR about. 4.0

**MATERIAL REQUIRED** (not supplied with the kit)

- Pipettes
- Distilled water
- Control Plasma and reagents\*

REF	5020040	Coagulation Control N	5 x 1 mL
REF	5020050	Coagulation Control N	50 x 1 mL
REF	5011050	Coagulation Control AK	5 x 1 mL
REF	5011060	Coagulation Control AK	50 x 1 mL
REF	5003128	TECHNOCLOT® PT	10 x 4 mL
REF	5003129	TECHNOCLOT® PT	10 x 10 mL
REF	5003220	TECHNOCLOT® PT Plus	5 x 4 mL
REF	5003222	TECHNOCLOT® PT Plus	10 x 10 mL
REF	5005032	TECHNOCLOT® PT Owren	10 x 4 mL
REF	5005037	TECHNOCLOT® PT Owren	10 x 10 mL

\* or any other package sizes, special Ceveron® alpha or TECHNOCLOT® Control and Calibration reagents of Technoclone.

**WARNING AND PRECAUTIONS**

- IVD for *in vitro* diagnostic use
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.
- Each single donor plasma and each lot of AK-Calibrant are tested and found negative for Hb<sub>s</sub>Ag, HIV 1/2 Ab and HCV Ab. However, universal precautions (treating all human source materials as if potentially infectious) should be exercised.

**STABILITY AND STORAGE**

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8 °C.

Stability after reconstitution:

RT**	+15°C	+4°C
4 hours	6 hours	8 hours

The vials can be only frozen once.

Upon storage, caps should be screwed tightly.

\*\*=room temperature

**TEST PROCEDURE**

**PREPARATION AND PERFORMANCE OF THE TEST**

**CEVERON**

Technoclone provides application sheets Ceveron® alpha. The application sheets contain analyser/assay specific handling and performance information which may differ from that provided in this instruction for use. In this case the information contained in the application sheets supersedes the information in this instruction for use. Please consult the instruction manual of the Ceveron® alpha.

**MANUAL**

- Reconstitute the lyophilised plasma in 1 mL distilled water (RT) before performing the test and allow to stand for 10 minutes\*\*\*. For standardising tests a reconstitution time of 30 minutes is recommended.
- The reconstituted plasmas should be treated in the same way as patient samples, using the instructions given for the test.

**ANALYSES RESULTS**

**EVALUATION**

Evaluations in item 1 and 2 are based on the INR definition.

**1. Establishing an INR Reference Curve**

Determine the coagulation times of AK-Calibrant A, B, C and D in duplicate and plot the mean values on double log paper (x-axis: INR values; y-axis: Coagulation time in seconds see Fig. 1). From this reference curve the INR's of patient plasmas can be directly read off.

**2. Establishing Laboratory Specific ISI's**

A linear relationship exists between the log INR value and log PT (sec) as expressed by the equation [log (PT patient plasma) = 1/ISI x log (PT normal plasma)].

A local ISI value can either be calculated by linear regression analysis or read from the graphical plot of the results using the enclosed graph sheet. This sheet contains a pre-set ISI scale printed on the Y2 axis which is derived mathematically from the equation given above.

**3. Graphical Evaluation**

- The INR calibration curve is established as described above in Section 1 using the enclosed graph evaluation sheet and the curve extended to intersect the Y1 and the Y2 axis.
- A second line is then drawn parallel to the calibration curve starting from the co-ordinates INR 1.0/10 sec on the x-axis and continuing through to the Y2 axis. This line is referred to as the standardised straight-line (see Figure 2). The local ISI value for the thromboplastin is read on the printed ISI scale on the Y2 axis where the standardised straight-line intersects this axis (see Figure 2).
- An alternative method to the standardised straight-line is to measure the distance (a) between the calibration curve and x-axis at the co-ordinates INR 1.0/10 sec and then measure the same distance down on the Y2 axis starting from the intersect of the calibration curve on this axis (see Figure 2)

**4. Establishing a Laboratory Specific PT Normal Plasma**

Using the following equation:

$$\text{PT normal plasma (sec)} = \frac{\text{PT AK-Calibrant A (sec)}}{\text{INR AK-Calibrant A}}$$

**5. Setting up the Reference (Calibration) Curve as a Percentage of the Norm**

For thromboplastins, e. g. TECHNOCLOT® PT Plus, the percentages of norm values may be read off from the enclosed table. These values are specific to each thromboplastin and cannot be transferred to other reagents. In order to establish the reference curve, the percentages of the norm are plotted against the mean time values in seconds and the points joined by a straight line (Fig. 3).

**STANDARDISATION**

In order to monitor oral anticoagulant treatment the Prothrombin time (PT) should, according to the WHO, be reported as an International Normalised Ratio (INR). The INR is calculated using the ISI (International Sensitivity Index) by the following formula:

$$\text{INR} = \left[ \frac{\text{TPZ patient plasma}}{\text{TPZ normal plasma}} \right]^{\text{ISI}}$$

The ISI of a thromboplastin is assigned by the manufacturer, and is batch specific. The ISI is also influenced by the type of instrument used and variations can be clearly seen when different instruments are used.

**LIMITATION OF THE TEST**

A new calibration is required for each batch of thromboplastin and for each instrument used. Also, if software changes are introduced or following a major service of either instruments or equipment, recalibration is recommended.

The figures in the table are only applicable to the indicated lot number of AK-Calibrant and the method has to adhere strictly as prescribed for the respective reagent.

**PERFORMANCE CHARACTERISTICS**

Within-run precision was assessed over multiple runs using a specific lot of reagent and AK-Calibrant plasmas. The coefficients of variation obtained in this study were < 3% for PT sec for Level A, and < 6% for Levels B-D.

ISTH and CLSI documents state that the difference between duplicate PT (seconds) as well as the difference between the mean of the run (day- to-day) results must be within 10%.

**TRACEABILITY**

The INR's for AK-Calibrant are determined according to WHO (WHO\_TRS\_889\_A3\_1999) guidelines by Technoclone and/or by external quality assessment.

**LITERATURE**

Please contact Technoclone or your local distributor.

\*\*\* When modifying methods, the reconstitution volumes might change.

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