
DAPTTIN[®] TC

For research
use only

REF	5035060	DAPTTIN [®] TC	5 x 2 mL
REF	5035090	DAPTTIN [®] TC	6 x 10 mL
REF	5035100	DAPTTIN [®] TC	20 x 10 mL

symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Ključova slova / Značenje simbola			
	manufacturer / Hersteller / fabbricante / fabricante / fabricant / fabricant / Tilverkaren / Fabrikanten / Produzent / Κατασκευαστής / Производитель / Производител / Производител / výrobc / Proizvođač		expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration / utgångsdatum / udløbsdato / Utløpsdato / Ημερομηνία λήξης / срок на годност / datum expirace / срок годности / datum expirace / Rok trajanja
	storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage / lagringstemperatur / oppbevaringstemperatur / Oppbevaringstemperatur / θερμοκρασία αποθήκευσης / съхранение на / teplota skladování / температура хранения / teplota skladování / Temperatura lagerovanja		consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation / se användarinstruktioner / følg brugsvejledning / Følg bruksanvisningen / συμβουλευθείτε τις οδηγίες για τη χρήση / прочетете инструкцията за работа / potfeba fídit se instrukcemi / перед использованием читайте инструкцию / sledujte návod k použití / Pročitaj upustvo pre upotrebe
			determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations / bestämmingar / bestemmelser / Bestemmelser / προσδιορισμοί / брой тестове / stanovení / определенний / роčet stanovení / Definicija
AQUA	distilled water / destilliertes Wasser / acqua distillata / agua destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / децилирана вода / destilovaná voda / дистиллированная вода / destilovaná voda / Destiliana Voda	LOT	lot / Charge / lotto / lote / lote / lot / sats / serie / Parti / партиа / партида номер / šarže / lot / šarže / Serija
BUF	Reaction buffer / Reaktionspuffer / tampone di reazione / tampón de reacció / Tampão de reação / tampon de réaction / Reaktionsbuffert / Reaktionsbuffer / Reaktionsbuffer / διάλυμα αντίδρασης / Реакционен буфер / Рабочий буферный раствор / Reakčni pufr / Reakcioni pufer	MTP	microtiter plate / Mikrotiterplatte / placa microtiter / microplaca / microplaca / microplaques sensibilisées / Mikrotiterplatta / Mikrotiterplade / mikrotiterplate / πλάκα μικροτιτλοδότησης / Микротительна плака / Микропланшет / Mikrotitrační destička / Mikrotitracione ploče
CAL	Calibrator / Kalibrator / Calibratore / calibrador / calibrador / calibreur / Kalibrator / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrator	REF	catalogue number / Katalognummer / numero di catalogo / número de catálogo / número de referència / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / каталожен номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj
CONJ	Conjugate / Konjugat / Coniugato / conjugado / conjugado / conjugaté / Konjugerad / Konjugat / Konjugat / συνδεδετικό / Конюгат / Конъюгат / Konjugát / Konjugat	RTU	ready to use / gebrauchsfertig / pronto all'uso / listo para usar / pronto a usar / prêt à l'emploi / færdig att användas / færdig til brug / klar til bruk / έτοιμο προς χρήση / Готов за употреба / готов к использованию / k přímému použití / Razrediti ili rastvoriti
CONT	Control / Kontrolle / controllo / control / control / contrôle / Kontroll / Kontroll / Kontroll / διάλυμα ελέγχου / Контрол / Контрольный образец / Kontrola / Kontrola	STOP	stop solution / Stoppløsning / Soluzione di arresto / solución de parada / solução de paragem / solution d'arrêt / Stoppløsning / Stop-opløsning / Stoppløsning / διάλυμα πάυσης / Стоп разтвор / Стоп-разтвор / Zastavovací roztok / Stop solucija
DIL	dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o disolver / diluir ou dissolver em / diluer ou dissoudre dans / spädd eller upplöst i / fortyndes eller oppløses i / Fortyndes eller oppløses i / αραιωση ή διάλυση σε / разворете или разредете с / zředit anebo rozpustit v / разбавить или растворить в / nafedte nebo rozpustit v / razrediti ili rastvoriti u	SUB	substrate / Substrat / substrato / substrato / substrato / substrat / Substrat / Substrat / Substrat / υπόστρωμα / Субстрат / Субстрат / Substrát / Substrat
INC	incubation buffer / Inkubationspuffer / tampone di incubazione / tampón de incubación / tampão de incubação / tampon d'incubation / Inkubationsbuffert/ Inkubationsbuffer/ Vaskebufferkonsentrat / διάλυμα επώσης / Инкубационен буфер / Буфер для инкубации / Inkubační pufr / Inkubacioni pufer	WASH	washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado / solução de lavado / solución de lavado concentrada / tampão de lavagem concentrado / Tampon de lavage concentré / Vattenlösningskonsentrat / Vaskeopløsningskonsentrat / vaskeløsningskonsentrat / συμπυκνωμένο διάλυμα πλύσης / Концентриран миеш разтвор / Концентрат промывочного раствора / Koncentrát promývacího roztoku / Koncetrat solucije za ispiranje
RUO	For research use only		

PRODUCT DESCRIPTION

INTENDEND USE

The Activated Partial Thromboplastin Time (aPTT) assay is a universally accepted screening procedure used to detect abnormalities in the intrinsic coagulation System. Dapttin® TC can be used to detect deficiencies of Factors II, V, VIII, IX, X, XI, and XII. In addition, Dapttin® TC can be used to detect the Lupus anticoagulant and is recommended when monitoring heparin therapy since it is sensitive to the presence of heparin.

COMPOSITION

Dapttin® TC (double-activated aPTT) is a standardized aPTT-(activated partial thromboplastin time) reagent consisting of silica and sulphatide as surface activators and a blend of highly purified phospholipids. Response to all factors and inhibitors is optimal.

MATERIAL REQUIRED (not supplied with the kit)

-	Pipette: 100 µL	-	Distilled water	
-	Buffers			
REF	5277015	CaCl ₂ 25 mmol/L		100 mL
-	Control Plasma Normal and Abnormal **			
REF	5020040	Coagulation Control N		5 x 1 mL
REF	5021055	Coagulation Control A		5 x 1 mL

** or any other package sizes, special Ceveron® alpha or Technoclot Control and Calibration reagents of Technoclone.

WARNING AND PRECAUTIONS

- Dapttin® TC contains no human material
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.

STABILITY AND STORAGE

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8 °C. Stability after reconstitution:

37 °C	RT*	15 °C	12 °C (Ceveron)	+4 °C	-20 °C
8 hours	3 days	7 days	7 days	15 days	6 months

Upon storage, caps should be screwed tightly. * = room temperature (20...25 °C)
Repeated freezing is not recommended

TEST PROCEDURE

PREPARATION OF PLASMA SAMPLES

Plasma separation: Mix 9 parts of venous blood and 1 part of Sodium Citrate Solution (0.11 mol/L) and centrifuge for 15 minutes at a RCF of at least 2500 g (corresponding to DIN 58905). Store the plasma at room temperature and use within 4 hours. Stability at -20°C: 1 month.

PREPARATION OF REAGENT

The lyophilised reagent should be dissolved in the volume of distilled water (room temperature) indicated. Allow the dissolved reagent to stand for 10 minutes before use.

PERFORMANCE OF THE TEST

CEVERON

Technoclone provides application sheets for Ceveron® alpha. The application sheets contain analyser/assay specific handling and performance information which may differ from that provided in this instruction for use. In this case the information contained in the application sheets supersedes the information in this instruction for use. Please consult the instruction manual of the Ceveron® alpha.

MANUAL

The CaCl₂ 25 mmol/L should be preheated to 37°C.
When using equipment with a mechanical measuring system, it is recommended to preheat Dapttin® TC to 37°C.
Pipetting Scheme:

0.10 mL plasma sample
+ 0.10 mL DAPTIN® TC
Shake briefly and incubate for two minutes at 37°C
+ 0.10 mL CaCl ₂ solution (25 mmol/L; 37°C)
Start timing after the addition of the CaCl ₂ solution and determine the clotting end point

ANALYSES RESULTS

REFERENCE RANGE

29 – 42 seconds
For quality control a normal control plasma (i.e. Coagulation Control N) and abnormal control plasma (i.e. Coagulation Control A) should be tested in each test series.
Dapttin® TC has a low opacity and is suitable for use in all commercially available automated coagulometers. The reagent can be used in conjunction with standard machine aPTT programs and the incubation time varied to suit the machine parameters (range 2 – 5 minutes).

STANDARDISATION

No international calibrator is available for the standardization of aPTT tests. Dapttin® TC documentations are based on studies testing plasma samples from normal subjects and heparinised plasma samples. A house reference batch is established in order to avoid batch-to-batch variation.

LIMITATION OF THE TEST

- Only citrate should be used as anticoagulant (neither EDTA nor heparin)
- The blood samples should be taken in non-activating samples tubes in order to avoid activation of the clotting factors in the plasma sample.
- Plasma samples should be stored at room temperature to avoid cold activation.
- Platelet Poor Plasma samples should be used.
- It is important to use the correct incubation time, when testing both heparinised and normal plasma samples.
- The calcium concentration used for recalcification should be 25 mmol/L when testing normal, abnormal or heparin plasmas.
- The REFERENCE RANGE represents approximate values only. It is recommended that each laboratory should determine its own normal range to take account of laboratory specific factors, such as the type of coagulometer used, reagent temperature, patient population, etc...
- The values found when testing Control Plasma Normal and Control plasma Abnormal should be compared to the value given on the data-key for the corresponding lot of Control plasma.
- If the results obtained are outside the recommended range, avoid measuring patient samples until the problem is solved.
- A new calibration is required for each batch of Dapttin® TC and additional reagents (e.g. deficient plasmas) where a calibration curve is necessary and for each instrument used. Also a new calibration is recommended, if software changes are introduced or following a major service of either instruments or equipment.

PERFORMANCE CHARACTERISTICS

Performance data are given below. Results obtained in individual laboratories may differ.

PRECISION

Reproducibility was determined with different samples (in series and day to day).
The following results were obtained:

Sample	Intra assay		Inter assay	
	Sample 1	Sample 2	Sample 1	Sample 2
n	12	12	10	10
MV sec	34.0	85.7	33.1	80.1
SD (%)	0.394	0.478	0.408	1.44
CV (%)	1.16	0.56	1.23	1.80

COMPARISON OF METHODS OR CORRELATION

Following correlation (%) was obtained in comparing DAPTIN TC (Technoclone) with Actin FS (Siemens) in aPTT method:
aPTT method: n=46 y = 1.10721x + 2.8936 R² = 0.9312

LINEARITY

- FVIII method: 0.8 – 100 (Activity %)
- FIX method: 0.8 – 100 (Activity %)
- FXI method: 1.6 – 200 (Activity %)
- FXII method: 6.25 – 200 (Activity %)

DETECTION LIMIT

- Heparin: UFH: ≤ 1 IU/ml
- Heparin: LMW: ≤ 3 IU/ml
- Triglycride: no detection limit up to 500mg/dl
- Bilirubin: no detection limit up to 0.4mg/dl

LITERATURE

Please contact Technoclone or your local distributor.

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