PRODUCT DESCRIPTION**INTENDED USE**

Siron LS is an activated partial thromboplastin time (aPTT) reagent to be used:

- as a screening test for the intrinsic coagulation system
- in the diagnosis and treatment of haemophilia A and B
- for specific determinations of factors VIII, IX, X, XI and XII when used in combination with the corresponding deficiency plasma, for control of heparin treatment
- for detection of Lupus anticoagulant

Siron LS measures the overall activity of the intrinsic coagulation system. Siron LS is sensitive to the factors XII, XI, X, IX, VIII, V and II, but insensitive to factors VII and XIII.

COMPOSITION

Siron LS is an activated partial thromboplastin time (aPTT) test. This liquid preparation is a stable aqueous suspension of phospholipids. FXII activation of Siron LS is achieved by addition of ellagic acid.

MATERIAL REQUIRED (not supplied with the kit)

- Buffer**

[REF] 5277015 CaCl₂ 25 mmol/L 100 mL

- Control Plasma Normal and Abnormal **

[REF] 5020040 Coagulation Control N 5 x 1 mL

[REF] 5021055 Coagulation Control A 5 x 1 mL

[REF] 5220130 Coagulation Reference 5 x 1 mL

** or any other package sizes, special Ceveron[®] alpha or TECHNOCLOT[®] Control and Calibration reagents of Technoclone.

WARNING AND PRECAUTIONS

- RUO for research use only
- Siron LS contains no human material. FXII activation of Siron LS is achieved by ellagic acid and phospholipids from synthetic origins.
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.

STABILITY AND STORAGE

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8°C.

Stability **opened** vials:

+2...8°C	+12°C (Ceveron [®] alpha)	RT*
4 weeks	10 days	24 hours

Upon storage, caps should be screwed tightly.

*=room temperature (+18...25°C)

TEST PROCEDURE**PREPARATION OF PLASMA SAMPLES**

Plasma separation: Mix 9 parts of venous blood and 1 part of Sodium Citrate Solution (0.11 mol/L or 0.13 mol/L) and centrifuge for 15 minutes at a RCF of at least 2500 g (corresponding to DIN 58905) to obtain platelet poor plasma.

STORAGE OF PLASMA

Samples of plasma are stored in plastic tubes or siliconized vials at room temperature before testing:

- For patients that are not on heparin therapy, within **3 hours**
- For patients on heparin therapy, within **1 hour**
- For longer period - 20°C, 2 weeks
- 70°C, 6 months

PREPARATION OF REAGENT

Each vial contains 4 mL / 10 mL ready to use aPTT reagent. Before each use, let the reagent stabilize for 30 min at room temperature (+18...25°C) while shaking the vial from time to time.

PERFORMANCE OF THE TEST**CEVERON**

Technoclone provides application sheets for Ceveron[®] alpha. The application sheets contain analyser/assay specific handling and performance information which may differ from that provided in this instruction for use. In this case the information contained in the application sheets supersedes the information in this instruction for use. Please consult the instruction manual of the Ceveron[®] alpha.

MANUAL

The CaCl₂ 25 mmol/L should be preheated to +37°C.

Pipetting Scheme:

100 µL plasma sample
+ 100 µL Siron LS in a pre-warmed tube
incubate for 3 minutes
+ 100 µL CaCl ₂ solution (25 mmol/L; +37°C) Start stop-watch simultaneously
Leave the tube for 20 sec. before tilting. Determine the clotting end point when clot is formed. The results are given in seconds.

ANALYSES RESULTS**REFERENCE RANGE**

Normal range: 25 - 38 seconds

The normal range refers to analyzers with optical detection (Ceveron[®] alpha). It is recommended that each laboratory establish and controls its normal range.

Siron LS in control of heparin therapy:

Siron LS seconds of 2 - 4 times the normal has been found to correspond to the accepted therapeutic level of heparin. However, the actual therapeutic level should be decided from the clinical situation.

STANDARDISATION

No international calibrator is available for the standardization of aPTT tests. Siron LS documentation is based on studies testing plasma samples from normal subjects and heparinised plasma samples. An in-house reference batch is established in order to avoid batch-to batch variation.

LIMITATION OF THE TEST

- The blood samples should be taken in non-activating samples tubes in order to avoid activation of the clotting factors in the plasma sample.
- Plasma samples should be stored at room temperature to avoid cold activation.
- Platelet poor plasma samples should be used.
- It is important to use the correct incubation time, which is 3 minutes for Siron LS when testing both heparinised and normal plasma samples.
- The calcium concentration used for re-calcification should be 25 mmol/L when testing normal, abnormal or heparin plasmas.
- The REFERENCE RANGE represents approximate values only. Each laboratory should determine its own normal range to take account of laboratory specific factors, such as the type of coagulometer used, reagent temperature, patient population, etc...
- The values found when testing control plasma normal and control plasma abnormal should be compared to the value given on the data-key for the corresponding lot of Control plasma.
- If the results obtained are outside the recommended range, avoid measuring patient samples until the problem is solved.
- A new calibration is required for each batch of Siron LS and additional reagents (e.g. deficient plasmas) where a calibration curve is necessary and for each instrument used. Also a new calibration is recommended, if software changes are introduced or following a major service of either instruments or equipment.
- Clotting time is not affected in this assay by samples containing ≤ 1000 mg/dL Haemoglobin, ≤46 mg/dL unconjugated bilirubin and ≤ 1000 mg/dL triglycerides. However elevated levels of conjugated bilirubin causes a prolonged clotting time.

PERFORMANCE CHARACTERISTICS

Performance data are given below. Results obtained in individual laboratories may differ.

PRECISION

Reproducibility was determined with different samples (in series and day to day). The following results were obtained:

Sample	Intra assay		Inter assay	
	Sample 1	Sample 2	Sample 1	Sample 2
n	21	21	20	20
Mean (sec)	28.2	73.3	28.2	65.7
CV (%)	0.99	0.92	1.96	2.3

COMPARISON OF METHODS OR CORRELATION

Following correlation was obtained in comparing SIRON LS:

$$\text{Actin}^{\circ} \text{FSL} \quad y = 0.775x + 7.9 \quad R^2 = 0.95$$

LITERATURE

Please contact Technoclone or your local distributor.

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