
Fibrinogen Reagent

For Research Use Only

REF 5138080 Fibrinogen Reagent 5 x 5 mL

REF 5138085 Fibrinogen Reagent 5 x 2 mL

symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Ključova slova / Značenje simbola			
	manufacturer / Hersteller / fabbricante / fabricante / fabricant / Tillverkaren / Fabrikanten / Produzent / Κατασκευαστής / Производитель / Производител / vŕobca / Proizvođač		expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration / utgångsdatum / udløbsdato / Utløpsdato / Ημερομηνία λήξης / срок на годност / datum expirace / срок годности / datum expirace / Rok trajanja
	storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage / lagringstemperatur / opbevaringstemperatur / Oppbevaringstemperatur / Θερμοκρασία αποθήκευσης / съхранение на / teplota skladování / температура хранения / teplota skladování / Temperatura lagerovanja		consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation / se användarinstruktioner / følg brugsvejledning / Følg bruksanvisningen / συμβουλευθείτε τις οδηγίες για τη χρήση / прочетете инструкцията за работа / potfeba fidiť se instrukciami / перед использованием читайте инструкцию / sledujte návod k použití / Pročitaj upustvo pre upotrebe
			determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations / bestämmingar / bestemmelser / Bestemmelser / προσδιορισμοί / брой тестове / stanovení / определний / ročet stanovení / Definicija
AQUA	distilled water / destilliertes Wasser / acqua distillata / agua destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / дестилирана вода / destilovaná voda / дестилированная вода / destilovaná voda / Destilisaná Voda	LOT	lot / Charge / lotto / lote / lote / lot / sats / serie / Parti / παρτία / партида номер / šarže / lot / šarže / Serija
BUF	Reaction buffer / Reaktionspuffer / tampone di reazione / tampón de reacció / Tampão de reação / tampon de réaction / Reaktionsbuffer / Reaktionsbuffer / Reaksjonsbuffer / διάλυμα αντίδρασης / Реакционен буфер / Рабочий буферный раствор / Reakční pufr / Reakcioni pufer	MTP	microtiter plate / Mikrotiterplatte / placa microtiter / microplaca / microplaca / microplaques sensibilisées / Mikrotiterplatta / Mikrotiterplade / mikrotiterplate / πλάκα μικροτιτλοδότησης / Микротитърна плака / Микропланшет / Mikrotitracijska destička / Mikrotitracione ploče
CAL	Calibrator / Kalibrator / Calibratore / calibrador / calibrador / calibreur / Kalibrator / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrator	REF	catalogue number / Katalognummer / numero di catalogo / número de catálogo / número de referência / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / каталожен номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj
CONJ	Conjugate / Konjugat / Coniugato / conjugado / conjugado / conjugaté / Konjugerad / Konjugat / Konjugat / συνδεδετικό / Конюгат / Конъюгат / Konjugát / Konjugat	RTU	ready to use / gebrauchsfertig / pronto all'uso / listo para usar / pronto a usar / prêt à l'emploi / færdig att användas / færdig til brug / klar til bruk / έτοιμο προς χρήση / Готов за употреба / готов к использованию / k přímému použití / Razrediti ili rastvoriti
CONT	Control / Kontrolle / controllo / control / control / contrôle / Kontroll / Kontroll / Kontroll / διάλυμα ελέγχου / Контрол / Контрольный образец / Kontrola / Kontrola	STOP	stop solution / Stoppløsning / Soluzione di arresto / solución de parada / soluçào de paragem / solution d'arrêt / Stoppløsning / Stop-opløsning / Stoppløsning / διάλυμα παύσης / Стоп разтвор / Стоп-раствор / Zastavovací roztok / Stop solucija
DIL	dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o disolver / diluir ou dissolver em / diluer ou dissoudre dans / spädd eller upplöst i / fortyndes eller opløses i / Fortyndes eller oppløses i / αραιωση ή διάλυση σε / разтворете или разрежете с / zředit anebo rozpustit v / разбавить или растворить в / nafedte nebo rozpustit v / razrediti ili rastvoriti u	SUB	substrate / Substrat / substrato / substrato / substrato / substrat / Substrat / Substrat / Substrat / υπόστρωμα / Субстрат / Субстрат / Substrát / Substrat
INC	incubation buffer / Inkubationspuffer / tampone di incubazione / tampón de incubación / tampão de incubação / tampon d'incubation / Inkubationsbuffer / Inkubationsbuffer / Vaskebufferkonsentrat / διάλυμα επώασης / Инкубационен буфер / Буфер для инкубации / Inkubační pufr / Inkubacioni pufer	WASH	washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado / solución de lavado concentrada / tampão de lavagem concentrado / Tampon de lavage concentré / Vattenlösningskoncentrat / Vaskeopløsningskoncentrat / vaskeløsningskoncentrat / συμπυκνωμένο διάλυμα πλύσης / Концентриран миеш разтвор / Концентрат промывочного раствора / Koncentrat promývacieho roztoku / Koncentrat solucije za ispiranje
RUO	for research use only		

Fibrinogen Reagent

PRODUCT DESCRIPTION

INTENDED USE

For the determination of Fibrinogen based on a modified Clauss method. The method for Fibrinogen determination is based on a modification of the method as described by Clauss. Citrated plasma is diluted, mixed with an excess of Thrombin (approx. 80 I.U./mL) and the coagulation time is determined. There is a linear correlation between the logarithm of the coagulation time and the logarithm of the fibrinogen concentration.

COMPOSITION

Fibrinogen Reagent Kit contains:

mL	Reagent	Other data
5 x 2	Fibrinogen Reagent	~ 80 I.U. bovine Thrombin/mL
5 x 5	Fibrinogen Reagent	~ 80 I.U. bovine Thrombin/mL

MATERIALS REQUIRED (not supplied with the kit)

- Solutions / Buffer: Distilled water
- REF 5410012 Imidazole Buffer 90 mL
- Control Plasmas and Calibrators**
- REF 5020040 Coagulation Control N 5 x 1 mL
- REF 5021055 Coagulation Control A 5 x 1 mL
- REF 5220110 Coagulation Reference 5 x 1 mL

** or any other package sizes, special Ceveron® alpha or TECHNOCLOT® Control and Calibration reagents of Technoclone.

WARNING AND PRECAUTIONS

- RUO for research use only
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.

STABILITY AND STORAGE

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8 °C. Stability after reconstitution:

RT*	+2...8 °C / +12 °C (Ceveron)
5 days	12 days

Upon storage, caps should be screwed tightly. Avoid contamination with microorganisms.

*room temperature

TEST PROCEDURE

PREPARATION OF PLASMA SAMPLES

Mix 9 parts of venous blood and 1 part of sodium citrate solution (0.11 mol/L) and centrifuge for 15 min at the RCF of at least 2500g (corresponding to DIN 58905). Store the plasma at room temperature (up to 1 day). Keep the plasma in a plastic tube to avoid contact activation. For quality control a normal control plasma (e.g. Coagulation Control N) and abnormal control plasma (e.g. Coagulation Control A) should be tested in each test series.

PREPARATION OF REAGENT

The lyophilized reagents should be reconstituted in the volume of distilled water (room temperature) indicated. Allow the dissolved reagents to stand for 10 minutes at room temperature. For standardized tests a reconstitution time of 30 minutes is recommended.

PERFORMANCE OF THE TEST

CEVERON

Technoclone provides Application sheets for Ceveron® alpha. The Application sheets contain analyser/assay specific handling and performance information which may differ from that provided in this instruction for use. In this case the information contained in the Application sheets supersedes the information in this instruction for use. Please consult the instruction manual of the Ceveron® alpha.

MANUAL

Patient plasma and Control Plasmas are diluted 1:10 with Imidazole Buffer (1 volume plasma and 9 volumes of buffer). The optimum range of determination is between 9 and 40 seconds. If the time value is below 9 seconds the determination should be repeated using a higher dilution, e.g. 1:20, and if the time value is above 40 seconds with a lower dilution, e.g. 1:5

0.20 mL	diluted citrated plasma
	incubate for 1 minute at 37°C
+ 0.20 mL	Fibrinogen reagent (room temperature)
	determine the point of coagulation


ANALYSES RESULTS

TABLE FOR DETERMINATION OF FIBRINOGEN

Each package contains a lot specific table (calibration curve) for determination of fibrinogen. The values for the above mentioned controls are read off this table and have to be in range for the test to be valid. Alternatively, an own calibration curve can be established.

CALIBRATION CURVE

An own calibration curve for a 1:10 dilution of patient's plasma is established using the Coagulation Reference, not included in the package: Reconstitute the Coagulation Reference with 1 mL distilled water, prepare a geometric dilution of 1:5 (0.5 mL Coagulation Reference and 2.0 mL buffer) and from this dilution a 1:10, 1:20, and 1:40 dilution (see below) using imidazole Buffer.


- The coagulation times of these dilutions are determined.
- Calculate the fibrinogen concentration of the dilutions from the given value (mg Fibrinogen/100 mL) of the Coagulation Reference by multiplying by 10 and dividing by the appropriate dilution factor.
- Plot the values on double log paper x-axis: concentration in mg / 100 mL, y-axis: coagulation time in seconds

Example:

Dilutions	sec	mg Fibrinogen / 100 mL*
1:5	7.2	490
1:10	14.4	245
1:20	28.3	123
1:40	57.0	61.3

* mg Fibrinogen / 100 mL of the Coagulation Reference (e.g. 245) multiplied by 10 and divided by the dilution factors 5, 10, 20 and 40.

Example of a calibration curve for a 1:10 dilution of patient's plasma:

EVALUATION

Time values of the 1:10 diluted plasma samples can be read off the calibration curve directly. Using other dilutions the Fibrinogen values found have to be calculated as follows:

$$\frac{\text{fibrinogen value (table or calibration curve)}}{10} \times \text{dilution factor} = \frac{\text{mg fibrinogen}}{100 \text{ mL plasma}}$$

REFERENCE RANGE

180-450 mg Fibrinogen / 100mL

STANDARDISATION

The standardization of the lot specific table for determination was done by Technoclone, by using Coagulation Reference.

LIMITATION OF THE TEST

The tilting method is not suited for determination of Fibrinogen based on the Clauss method. We, therefore, recommend the hooking method or a coagulometer to determine the coagulation times. Some auto analyzers, based on optical clot detection, need special applications.

PERFORMANCE CHARACTERISTICS

Performance data are given below. Results obtained in individual laboratories may differ.

PRECISION

Reproducibility was determined with different samples (in series and day to day). The following results were obtained:

Sample	Intra assay		Inter assay	
	Sample 1	Sample 2	Sample 1	Sample 2
n	12	12	12	12
MV %	284.2	101.7	267.4	122.7
SD (%)	13.79	7.177	12.27	4.437
CV (%)	4.85	7.06	4.59	3.62

COMPARISON OF METHODS OR CORRELATION

Following correlation (%) was obtained in comparing Fibrinogen determination based on a modified Clauss method with Fibrinogen (Technoclone) and Fibrinogen (Diagnostica STAGO):
Fibrinogen method: n=192 $y = 0.986821x + 4.61536$ $R^2 = 0.911212$

LINEARITY

Fibrinogen method: 60 - 700 mg/dl

INTERFERENCES

No interference for: Heparin: UFH: ≤ 2 IU/ml
CBN Fibrinogen Fragment: ≤ 10 mg/ml
Triglyceride: ≤ 500 mg/dl
Bilirubin: ≤ 0.4 mg/dl

LITERATURE

Please contact Technoclone or your local distributor

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