Factor VIII Deficient Plasma immunads.

For Research Use Only

REF 5154002 Factor VIII Deficient Plasma, immunads.

5 x 1 mL

REF 5154004 Factor VIII Deficient Plasma, immunads.

50 x 1 mL

symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / Explicação dos símbolos / explication des symboles / Symbolförklaring / symbolforklaring / Tegnforklaring / символы			
	manufactured by / Hergestellt von / prodotto da / fabricado por / fabricado por / fabriqué par / Producerat av / fremstillet af / Produsert av / производитель	AQUA	distilled water / destilliertes Wasser / acqua distillata / agua destilada / Água destilada / eau distillée / Destilleret vand / destilleret vand / Destillert vann / дистиллированная вода
	expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / Data de validade / date d'expiration / Forfallsdag / udløbsdato / Utløpsdato / срок годности	DIL	dilute or dissolve in / verdünnen oder lösen in / a diluire o a sciogliere in / diluir o dissolver / diluir ou dissolver em / à diluer ou à dissoudre / Fortyndes eller opløses i / fortyndes eller opløses i / Fortynnes eller oppløses i / разбавить или растворить в
	storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / Temperatura de conservação / température de conservation / Förvaringstemperatur / opbevaringstemperatur / Oppbevaringstemperatur / температура хранения	RUO	For Research Use Only
	consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / Consultar o manual de instruções / consulter la notice d'utilisation / Følg bruksanvisningen / følg brugsvejledning / Følg bruksanvisningen / перед использованием прочитайте инструкцию	LOT	lot / Charge / lotto / lote / Lote / lot / Parti / serie / Parti / лот
		REF	catalogue number / Katalognummer / numero di catalogo / número de catálogo / Número de referência / référence / Artikelnummer / katalognummer / Katalognummer / каталожный номер
	determinations / Bestimmungen / determinazioni / determinaciones / Determinações / déterminations / Bestämningar / bestemmelse / Bestemmelser / определений		

Factor VIII Deficient Plasma, immunads.

GB

PRODUCT DESCRIPTION

INTENDEND USE

Factor VIII deficient plasma immunads. is used in the determination of Coagulation Factor VIII by one-stage method based on the Activated Partial Thromboplastin Time (aPTT).

COMPOSITION

The Factor VIII deficient plasma immunads. is an immune-adsorbed lyophilised, stabilised human plasma with a Factor VIII content of <1%, prepared from HIV 1/2 Ab negative plasmas.

MATERIAL REQUIRED (not supplied with the kit)

- Pipettes - Distilled water - Solutions/buffers:

REF	5410010	Imidazole buffer	50 mL
REF	5277015	CaCl ₂ 25 mmol/L solution	100 mL

- Reagents**

REF	5035060	Daptin® TC	5 x 2 mL
REF	5035105	Siron (aPTT liquid)	2 x 4 mL
- Control Plasmas and Calibrators***			
REF	5020040	Coagulation Control N	5 x 1 mL
REF	5021055	Coagulation Control A	5 x 1 mL
REF	5220110	Coagulation Reference	5 x 1 mL

** or any other package sizes of Technoclone.

*** or any other package sizes, special Ceveron® alpha or Technodot Control and Calibration reagents of Technoclone.

WARNING AND PRECAUTIONS

- RUO For Research Use Only

- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.
- This lot of reagents prepared from human blood and each single plasma used for this lot are HBsAg, HIV 1/2 Ab and HCV Ab negative (see package label and vial label).

STABILITY AND STORAGE

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8 °C. Stability after reconstitution:

RT* (Ceveron****) 3 hours	-20 °C 1 month
------------------------------	-------------------

Upon storage, caps should be screwed tightly. * = room temperature

**** = in the Ceveron® alpha in the respective control area in the sample tray.

When storing, the tubes must be securely capped.

Deep frozen reagent must be thawed for at least 10 minutes at 37°C and mixed thoroughly before use. Repeated freezing is not recommended

TEST PROCEDURE

PREPARATION OF PLASMA SAMPLES

Mix 9 parts of venous blood and 1 part of Sodium Citrate solution (0.11 mol/L) and centrifuge for 15 min at a RCF of at least 2500 (corresponding to DIN 58905). Store the plasma at room temperature (up to 4 hours). If quickly frozen, samples can be stored at -20°C up to 15 days or at -80°C up to 1 month. Thaw the sample at 37°C. Before carrying out the test, the sample plasma is diluted 1:5 (0.1 mL + 0.4 mL) using the imidazole buffer. For very low and very high Factor VIII levels, however, other dilutions should be used. If the presence of inhibitors is suspected, various dilutions of the plasma sample should be tested.

PREPARATION OF REAGENT

Reconstitute the lyophilised reagents in the prescribed quantities of distilled water and allow them to stand for 10 min at room temperature.¹ Prewarm the CaCl₂ 25 mmol/L solution to 37°C.

PERFORMANCE OF THE TEST

CEVERON

Technoclone provides application sheets for Ceveron® alpha. The application sheets contain analyser/assay specific handling and performance information which may differ from that provided in this instruction for use. In this case the information contained in the application sheet supersedes the information in this instruction for use. Please consult the instruction manual of the Ceveron® alpha.

MANUAL

Pipetting scheme:

0.1 mL	Factor VIII Deficient Plasma, immune-adsorbed
+ 0.1 mL	diluted plasma samples
+ 0.1 mL	PTT reagent (shake briefly and incubate for 3 min. at 37 °C)
+ 0.1 mL	CaCl ₂ 25 mmol/L solution 37°C
determine the point of coagulation	

ANALYSES RESULTS

REFERENCE RANGE

70-150 % of normal

CALCULATION OF THE RESULTS

The predilution 1:5 of the plasma sample is not to be considered in the evaluation.

The plasma samples prediluted at a ratio of 1:5 may be read off directly from the calibration curve. If dilution ratios other than 1:5 are used, the % F VIII values of the calibration curve may be converted by using the following formula:

$$\frac{\% \text{ FVIII content}}{5} \times \text{dilution} = \% \text{ FVIII of the sample}$$

CALIBRATION CURVE

Reconstitute the Coagulation Reference as indicated in the table. A predilution of 1:5 is attained by using an imidazole buffer in the ratio 1:5 (1 part Coagulation Reference plus 4 parts buffer). Prepare a geometric series of dilutions (1:1 to 1:128) of the predilution (1:5). The 1:1 ratio corresponds to the predilution 1:5.

Predilution	Dilutions of calibration curve							
1:5 =	1:1	1:2	1:4	1:8	1:16	1:32	1:64	1:128
% FVIII	100	50	25	12.5	6.3	3.13	1.56	0.78

Determine the coagulation times of the geometric series of dilutions and plot them on semilog paper (x-axis: log activity %; y-axis: coagulation time in seconds).

QUALITY CONTROL

In order to verify the accuracy of the results, abnormal control plasma (i.e. Coagulation Control A) and normal control plasma (i.e. Coagulation Control N) should always be tested in the same way as the sample plasma.

LIMITATION OF THE TEST

Incorrect sample handling can lead to partial activation of the coagulation factors and to falsely elevated single factor determinations.

Lupus anticoagulant can affect the apparent factor activity in single factor determination.

A new calibration is required for each batch of reagents where a calibration curve is necessary and for each instrument used. Also a new calibration is recommended, if software changes are introduced or following a major service of either instruments or equipment.

APPLICATIONS FOR INSTRUMENTS

Application sheets are available from Technoclone or your local distributor upon request

PERFORMANCE CHARACTERISTICS

Performance data are given below. Results obtained in individual laboratories may differ.

PRECISION

Reproducibility was determined with different samples (in series and day to day).

The following results were obtained:

	Intra assay		Inter assay	
Sample	Sample 1	Sample 1	Sample 2	Sample 2
n	56	6	6	
MV (sec)	36.5	47.3	75.6	
SD (%)	0.261	0.728	3.08	
CV (%)	0.7%	1.54%	4.08%	

COMPARISON OF METHODS OR CORRELATION

Following correlation (%) was obtained in comparing TECHNOCHROM® FVIII:C (Technoclone) with the FVIII deficient Plasma immunads in the Imidazole buffer method and deficient plasma method :

Imidazole buffer method: n=26 y = 1.1137x – 6.9503 R² = 0.9754
Deficient plasma method: n=20 y = 1.1475x – 14.32 R² = 0.9809

LINEARITY

Imidazole buffer method: 0.8 – 100 (Activity %)

Deficient plasma method: 0.4 – 100 (Activity %)

DETECTION LIMIT

Imidazole buffer method: 0.8 % (Activity %)

Deficient plasma method: 0.4 % (Activity %)

LITERATURE

Please contact Technoclone or your local distributor.

¹ For standardisation a reconstitution time of 30 min is recommended.

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