
Fletcher Trait Plasma

For Research Use Only

REF 5205006 Fletcher Trait Plasma

2 x 1 mL

Symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Ključova slova / Značenje simbola			
	Manufacturer / Hersteller / fabbricante / fabricante / fabricante / fabricant / Tillverkaren / Fabrikanten / Produzent / Κατασκευαστής / Производитель / Производител / výroba / Proizvođač		Expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration / utgångsdatum / udløbsdato / Utløpsdato / Ημερομηνία λήξης / срок на годност / datum expirace/ срок годности / datum expirace / Rok trajanja
	Storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage / lagringstemperatur / opbevaringstemperatur / Oppbevaringstemperatur / θερμοκρασία αποθήκευσης / съхранение на / teplota skladování / температура хранения / teplota skladování / Temperatura lagerovanja		Consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation / se användarinstruktioner / følg brugsvejledning / Følg bruksanvisningen / συμβουλευθείτε τις οδηγίες για τη χρήση / прочетете инструкцията за работа / potfeba řidit se instrukcemi / перед использованием читайте инструкцию / sledujte návod k použití / Pročitaj upustvo pre upotrebe
			Determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations / bestämningar / bestemmelser / Bestemmelser / προσδιορισμοί / брой тестове / stanovení / определний / počet stanovení / Definicija
AQUA	Distilled water / destilliertes Wasser / acqua distillata / agua destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / децилирана вода / destilovaná voda / дистиллированная вода / destilovaná voda / Destilisaná Voda	LOT	Lot / Charge / lotto / lote / lote / lot / sats / serie / Parti / партида номер / šarže / лот / šarže / Serija
BUF	Reaction buffer / Reaktionspuffer / tampone di reazione / tampón de reacció / Tampão de reação / tampon de réaction / Reaktionsbuffert / Reaktionsbuffer / Reaktionsbuffer / διάλυμα αντίδρασης / Реакционен буфер / Рабочий буферный раствор / Reakční pufr / Reakcioni pufer	MTP	Microtiter plate / Mikrotiterplatte / placa microtiter / microplaca / microplaca / microplaques sensibilisées / Mikrotiterplatta / Mikrotiterplade / mikrotiterplate / πλάκα μικροτιτλοδότησης / Микротитр-на плака / Микропланшет / Mikrotitrační destička / Mikrotitracione ploče
CAL	Calibrator / Kalibrator / Calibratore / calibrador / calibrador / calibreur / Kalibrator / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrator	REF	Catalogue number / Katalognummer / numero di catalogo / número de catálogo / número de referência / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / каталожен номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj
CONJ	Conjugate / Konjugat / Coniugato / conjugado / conjugado / conjugaté / Konjugerad / Konjugat / Konjugat / συνδεδετικό / Конюгат / Конъюгат / Konjugát / Konjugat	RTU	Ready to use / gebrauchsfertig / pronto all'uso / listo para usar / pronto a usar / prêt à l'emploi / færdig att användas / færdig til brug / klar til bruk / έτοιμο προς χρήση / Готов за употреба / готов к использованию / k přímému použití / Razrediti ili rastvoriti
CONT	Control / Kontrolle / controllo / control / control / contrôle / Kontroll / Kontroll / Kontroll / διάλυμα ελέγχου / Контрол / Контрольный образец / Kontrola / Kontrola	STOP	Stop solution / Stopplösung / Soluzione di arresto / solución de parada / solução de paragem / solution d'arrêt / Stopplösning / Stop-opløsning / Stop-opløsning / διάλυμα πάυσης / Стоп раствор / Стоп-раствор / Zastavovací roztok / Stop solucija
DIL	Dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o dissolver / diluir ou dissolver em / diluer ou dissoudre dans / spädd eller upplös i / fortyndes eller opløses i / Fortyndes eller opløses i / αραιωση ή διάλυση σε / растворите или разредете с / zředit anebo rozpustit v / разбавить или растворить в / nafedte nebo rozpustite v / razrediti ili rastvoriti u	SUB	Substrate / Substrat / substrato / substrato / substrato / substrat / Substrat / Substrat / Substrat / υπόστρωμα / Субстрат / Субстрат / Substrát / Substrat
INC	Incubation buffer / Inkubationspuffer / tampone di incubazione / tampón de incubación / tampão de incubação / tampon d'incubation / Inkubationsbuffert/ Inkubationsbuffer/ Vaskebufferkonsentrat / διάλυμα επώασης / Инкубационен буфер / Буфер для инкубации / Inkubační pufr / Inkubacioni pufer	WASH	Washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado / solución de lavado concentrada / tampão de lavagem concentrado / Tampon de lavage concentré / Vattenlösningsskoncentrat / Vaskeopløsningskoncentrat / vaskeløsningskoncentrat / συμπυκνωμένο διάλυμα πλύσης / Концентриран миеш разтвор / Концентрат промывочного раствора / Koncentrát promývacieho roztoku / Koncentrat solucije za ispiranje
RUO	For Research Use Only		

PRODUCT DESCRIPTION

INTENDENT USE

Fletcher Trait Plasma is used for measuring pre-kallikrein deficiency. A deficiency in pre-kallikrein will prolong the activated partial thromboplastin time (aPTT) considerably but does not usually result in the expected haemorrhagic tendency. This abnormality in the activating phase of the aPTT can be corrected in vitro by adding Factor XII Deficient Plasma or Fitzgerald Trait Plasma (HMW-Kininogen Deficient Plasma). In addition a prekallikrein deficiency will also affect fibrinolysis, kinin formation and the complement system.

COMPOSITION

Fletcher Trait Plasma is a lyophilized and stabilized human pre-kallikrein deficiency plasma.

MATERIAL REQUIRED (not supplied with the kit)

-	Pipettes	- Distilled water	- Solutions/buffers:
REF	5279025	CaCl ₂ 50 mmol/L Solution	100 mL
REF	5400045	Citrate-Natriumchloride-Buffer	60 mL
-	Reagents		
REF	5035060	Daptin® TC	5 x 2 mL
REF	5035090	Daptin® TC	6 x 10 mL
-	Control and Calibration Plasma		
REF	5020040	Coagulation Control N	5 x 1 mL
REF	5021055	Coagulation Control A	5 x 1 mL
REF	5220110	Coagulation Reference	5 x 1 mL

WARNING AND PRECAUTIONS

- For Research Use Only
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.
- This lot of reagents prepared from human blood and each single plasma used for this lot are HB_eAg, HIV 1/2 Ab and HCV Ab negative (see package label and vial label).

STABILITY AND STORAGE

The expiry date printed on the labels applies to storage of the unopened bottles at +2...+8 °C. Stability after reconstitution:

RT*	RT* Ceveron® alpha	-20 °C
2 hours	2 hours	1 month

Upon storage, caps should be screwed tightly.
* = room temperature: +18...+25 °C

TEST PROCEDURE

PREPARATION OF PLASMA SAMPLES

Mix 9 parts of venous blood and 1 part of Sodium Citrate solution (0.11 mol/L) and centrifuge for 15 min at a RCF of at least 2500 (corresponding to DIN 58905). Store the plasma at room temperature. Before carrying out the test, the sample plasma is diluted 1:100 using the Citrate-Natriumchloride-Buffer; dilution is to be carried out in two stages on a 1:10 scale (for instance 0.05 mL and 0.45 mL each). For very low and very high pre-kallikrein levels, however, other dilutions should be used.

PREPARATION OF REAGENT

The lyophilized reagent is reconstituted with 1 mL of distilled water and allowed to stand for 10 min at room temperature¹. Prewarm the CaCl₂ solution to +37 °C.

SUGGESTED PERFORMANCE OF THE TEST

The results obtained with the test arrangement as suggested depend to a large extent on the respective lot of the PTT reagent used. Should the values not seem plausible, the pre-dilution of the plasma sample and/or the incubation period should be varied accordingly.

CEVERON

Technoclone provides application sheets for Ceveron® alpha. The application sheets contain analyser/assay specific handling and performance information which may differ from that provided in this instruction for use. In this case the information contained in the application sheets supersedes the information in this instruction for use. Please consult the instruction manual of the Ceveron® alpha.

MANUAL Pipetting scheme:

0,10 mL	Fletcher Trait Plasma
+ 0,10 mL	diluted plasma sample
+ 0,10 mL	PTT reagent (f.i. Daptin TC) (shake briefly and incubate for 1 min. at +37 °C)
+ 0,10 mL CaCl ₂ 50 mmol/L solution + 37 °C	
determine the point of coagulation	

ANALYSES RESULTS

REFERENCE RANGE

70-120 % of normal

CALCULATION OF THE RESULTS

The predilution 1:100 of the plasma sample is not to be considered in the evaluation. The plasma sample prediluted at a ratio of 1:100 may be read off directly from the calibration curve. If dilution ratios other than 1:100 are used, the % pre-kallikrein values of the calibration curve may be converted by using the following formula:

$$\frac{\% \text{ pre-kallikrein}}{100} \times \text{dilution} = \% \text{ pre-kallikrein of the sample}$$

CALIBRATION CURVE

The Coagulation Reference is a lyophilised normal plasma pool to be dissolved in 1 mL distilled water or accordance with data indicated on the table. With Citrate-Natriumchloride-Buffer predilution is effected in line with the following instructions: 1:100 (twice one part Coagulation Reference plus 9 parts Citrate-Natriumchloride-Buffer). Prepare a geometric series of dilutions (1:1 to 1:32) of the predilution (1:100). The 1:1 ratio corresponds to the predilution 1:100.

Predilution	Dilution of calibration curve					
1:100	1:1	1:2	1:4	1:8	1:16	1:32
% pre-kallikrein	100	50	25	12,5	6,3	3,1

Determine the coagulation time of the different dilutions and plot on log-log paper (x-axis: log activity in %; y-axis: log coagulation time in sec).

QUALITY CONTROL

In order to verify the accuracy of the results, Coagulation Control A and N should always be tested in the same way as the sample plasma.

LIMITATION OF THE TEST

- Incorrect sample handling can lead to partial activation of the coagulation factors and to falsely elevated single factor determinations.
- A new calibration is required for each batch of reagents where a calibration curve is necessary and for each instrument used. Also a new calibration is recommended, if software changes are introduced or following a major service of either instruments or equipment.

APPLICATIONS FOR INSTRUMENTS

Application Sheets are available from Technoclone or your local distributor upon request.

LITERATURE

Please contact Technoclone or your local distributor.

¹ For standardisation a reconstitution time of 30 min is recommended.

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