## REAAADS COAGULATION CONTROL 2

**REAAADS Coagulation Control 2 (CC-2)** is an assayed control plasma for *in vitro* diagnostic use in coagulation studies.

### PRINCIPLE OF THE PROCEDURE

Coagulation Control 2 is tested in the same manner as citrated patient plasma samples to assess the performance of each assay run, for the parameters listed on the reverse side. Testing variables in each laboratory, including equipment, reagents, and technique may influence control recovery. Although an expected range is provided for each parameter, laboratories should establish their own expected range for their particular instrument-reagent system.

### REAGENT

**Contents:** lyophilized pool of citrated human plasma from healthy donors, with stabilizers, to be reconstituted with 1.0 mL deionized water.

### Precautions:

- CC-2 has been tested and found negative for HBsAg, HIV-I/II, and HCV antibodies. However, all human blood derivatives, including patient samples, should be treated as potentially infectious.

Biological Risk  Warning .

- Signs of deterioration: the unreconstituted material should appear as a light yellow, dry plug or straw-colored pieces. Notify manufacturer before use if appearance of product has changed.

### Instructions for use:

**Reagent preparation:** reconstitute with 1.0 mL of deionized water. Swirl gently to mix. Allow to stand 10 minutes before use for complete dissolution.

**Storage and stability:** lyophilized product is stable until expiration date printed on the vial, when stored at 2 - 8°C; reconstituted product is stable for 4 hours when maintained at 2 - 8°C.

### PERFORMANCE CHARACTERISTICS:

Assayed values for the various coagulation parameters should be used as a guideline. Reagents from Corgenix, Inc. or Helena Laboratories were used to assign values traceable to World Health Organization standards, where available. For best results, each laboratory should establish mean values and expected ranges for their own particular coagulation system.

### WARRANTY

This product is warranted to perform as described in this package insert. Corgenix, Inc. disclaims any implied warranty of merchantability or fitness for a particular use, and in no event shall Corgenix, Inc. be liable for consequential damage.

**For Technical or Customer Service in the United States, phone 1-800-729-5661. Outside the United States, phone (303) 457-4345, fax (303) 457-4519, or contact a Corgenix authorized distributor. [techsupport@corgenix.com](mailto:techsupport@corgenix.com) or [customersupport@corgenix.com](mailto:customersupport@corgenix.com) [www.corgenixonline.com](http://www.corgenixonline.com)**

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© 2016, Corgenix, Inc.
MT Promedy Consulting GmbH