## COATEST™ APC™ Resistance Control Plasma Level 1 - 82 2650 63 ENGLISH - Insert revision 03/2013

### Intended use

For control of accuracy when using the Coatest APC Resistance and Coatest APC Resistance V kits.

### Reagent

**Control Plasma Level 1** 5 vials

Lyophilized preparation of citrated human plasma.

**CAUTION:** Each donor unit used in the preparation of human source reagent has been tested by FDA approved methods for the presence of Hepatitis B surface antigen and antibodies to HIV 1 and 2 and Hepatitis C and found to be negative. However, since no test can completely rule out the presence of these blood borne diseases, the handling and disposal of human source reagents from this product should be made with care.

Hazard class: **None**

Risk phrases: **None**

Safety phrases: **None**

This product is for *in vitro* diagnostic use.

### Procedure

Reconstitute each vial with 1.0 mL deionized water, filtered through 0.22 µm or NCCLS type II water. Allow to stand at room temperature for 30 minutes. Swirl gently before use.

### Storage conditions and stability

Stability after reconstitution is 6 hours at 2-25°C or 3 months at -20°C or below. The frozen reagent should be thawed at 37°C and gently mixed before use. Do not refreeze.

### Expected values

Different clotting times will be obtained with different types of instruments depending on the clot detection principle. The following ranges cover the expected clotting times when performing analysis of Control Plasma Level 1 with the above batch number using the Coatest APC Resistance and Coatest APC Resistance V kits.

### Traceability of calibrators and control materials

The reported values were determined over multiple runs on ACL Futura using a specific lot of reagent and against an internal House Standard. As an International Standard is not still available for the APCR-V assay, the values have been assigned against a House Standard which is traceable to frozen plasma pool.

## COATEST™ APC™ Resistance Control Plasma Level 1 - 82 2650 63 DEUTSCH - Beipackzettel Version 03/2013

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ACCEPTANCE RANGE			
COATEST™ APC™ Resistance			
TP (e.g. ACL, MLA)	CaCl2 (s)	APC/CaCl2 (s)	APC ratio
EM (e.g. ST 4, KC, STA)	CaCl2 (s)	APC/CaCl2 (s)	APC ratio
COATEST™ APC™ Resistance V			
ALL Instruments	CaCl2 (s)	APC/CaCl2 (s)	APC ratio

LOT

---

***Coatest APC Resistance Control Plasma Level 1      82 2650 63***

---

Printed Insert Sheet: 302198  
Revision: R3  
Issued: March 2013  
C.O.: 433694

---

***LANGUAGES***

---

ENGLISH  
DEUTSCH  
ESPAÑOL  
FRANÇAIS  
ITALIANO  
DANSK  
SWENSK  
GREEK  
PORTUGUÊS

---

***TECHNICAL SPEC'S***

---

PAPER: White paper, 50-60 g/m<sup>2</sup> weight.  
SIZE: 11 x 8.5" (280 x 216 mm.).  
PRINT: Front/Back.  
PRINT COLOR: Front - Top rule Orange Pantone 137,  
all remaining type in black.  
Back - All type in black.