COATEST™ APC™ Resistance Control Plasma Level 2 - 82 2668 63 ENGLISH - Insert revision 12/2016

Intended use

For control of accuracy when using the Coatest APC Resistance and Coatest APC Resistance V kits.

Reagent

Control Plasma Level 2 5 vials Lyophilized preparation of citrated human plasma from selected donors.

CAUTION: Each donor unit used in the preparation of human source reagent has been tested by FDA approved methods for the presence of Hepatitis B surface antigen and antibodies to HIV 1 and 2 and Hepatitis C and found to be negative. However, since no test can completely rule out the presence of these blood borne diseases, the handling and disposal of human source reagents from this product should be made with care.

Hazard class: None

Hazard statements: None

Precautionary statements: None

Supplemental hazard information: ≈ 98,7% of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and unknown hazard to the aquatic environment.

This product is for *in vitro* diagnostic use.

Procedure

Reconstitute each vial with 1.0 mL deionized water, filtered through 0.22 µm or NCCLS type II water. Allow to stand at room temperature for 30 minutes. Swirl gently before use.

Storage conditions and stability

Stability after reconstitution is 6 hours at 2-25°C or 3 months at -20°C or below. The frozen reagent should be thawed at 37°C and gently mixed before use. Do not refreeze.

Expected values

Different clotting times will be obtained with different types of instruments depending on the clot detection principle. The following ranges cover the expected clotting times when performing analysis of Control Plasma Level 2 with the above batch number using the Coatest APC Resistance and Coatest APC Resistance V kits.

Traceability of calibrators and control materials

The reported values were determined over multiple runs on ACL Futura using a specific lot of reagent and against an internal House Standard. As an International Standard is not still available for the APCR-V assay, the values have been assigned against a House Standard which is traceable to frozen plasma samples which have been determined to be homozygous or heterozygous respectively for factor V dependent APC Resistance.

COATEST™ APC™ Resistance Control Plasma Level 2 - 82 2668 63 DEUTSCH - Beipackzettel Version 12/2016

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ACCEPTANCE RANGE			
COATEST™ APC™ Resistance			
TP (e.g. ACL, MLA)	CaCl2 (s)	APC/CaCl2 (s)	APC ratio
EM (e.g. ST 4, KC, STA)	CaCl2 (s)	APC/CaCl2 (s)	APC ratio
COATEST™ APC™ Resistance V			
ALL Instruments	CaCl2 (s)	APC/CaCl2 (s)	APC ratio

LOT

Coatest APC Resistance Control Plasma Level 2 82 2668 63

Printed Insert Sheet: 302204
Revision: R4
Issued: December 2016
C.O.: 473324

LANGUAGES

ENGLISH
DEUTSCH
ESPAÑOL
FRANÇAIS
ITALIANO

TECHNICAL SPEC'S

PAPER: White paper, 50-60 g/m² weight.
SIZE: 11 x 8.5" (280 x 216 mm.).
PRINT: Front/Back.
PRINT COLOR: Front - Top rule Orange Pantone 137,
all remaining type in black.
Back - All type in black.