
PRODUCT INFORMATION SHEET
Monoclonal antibodies detecting human antigens

Anti-Hemoglobin S (HbS)

FITC ASR REF IQP-574F CONC 50 µg/mL

ASR Analyte Specific Reagent, analytical and performance characteristics are not established.
The product complies with the ASR definitions of the U.S. Food & Drug Administration (FDA).



Description

Clone 57-8

Isotype IgG1

Specificity The HbS antibody recognizes the β -chain of the human hemoglobin. Reactivity is seen with human erythrocytes and nucleated red cell precursors containing HbS. No binding is observed with HbA, HbA2, HbC, HbE and HbF.

Applications Flow cytometry

Representative data

57-8 was analyzed by flow cytometry on a blood sample using adult red blood cells (RBCs). RBCs were fixed, permeabilized and stained with anti-HbS FITC.




Figure: Flow cytometry analysis of whole blood stained with anti-Hemoglobin S. Histograms: uninterrupted line represents HbS^{high} sample, interrupted line represents blanc sample.



Handling and Storage

Antibodies are supplied in 0.01 M sodium phosphate, 0.15 M NaCl; pH 7.3, 0.2% BSA, 0.09% sodium azide (NaN₃). Store the vials at 2-8 °C. Monoclonal antibodies should be protected from prolonged exposure to light. Reagents are stable for the period shown on the vial label when stored properly.

Warranty Products sold hereunder are warranted only to conform to the quantity and contents stated on the label at the time of delivery to the customer. There are no warranties, expressed or implied, which extend beyond the description on the label of the product. IQ Products is not liable for property damage, personal injury, or economic loss caused by the product.

Characterization To ensure consistently high-quality reagents, each batch of monoclonal antibody is tested for conformance with characteristics of a standard reagent. Representative flow cytometric data is included in this data sheet.

Warning All products contain sodium azide. This chemical is poisonous and hazardous. Handling should be done by trained staff only.

Explanation of used symbols

	Description
	Catalogue number
	Caution, consult accompanying document
	Keep away from (sun)light
	Biological risks
	Temperature limitation (°C)
	Analyte Specific Reagent
	Batch number
	Use by yyyy-mm-dd
	Manufacturer

	Conjugates	Ex -max (nm)	Em -max (nm)
P	PURE Unconjugated antibody		
F	FITC Fluorescein Isothiocyanate	488	519
R	R-PE R-Phycoerythrin	488, 532	578
C	CyQ Tandem conjugate of R-PE-and Cy5.18	488, 532	667
A	APC Allophycocyanin	595, 633, 635, 647	660
D	Dy-410 Violet Dye 410	405	460
PC	PerCP Peridinin-chlorophyll-protein	488, 532	678
PCC	PerCP-Cy5.5 Tandem conjugate of PerCP-and Cy5.5	488, 532	695



IQ Products BV
Rozenburglaan 13a
9727 DL Groningen, The Netherlands

+31 (0)50 57 57 000
 +31 (0)50 57 57 002
 Technical marketing@iqproducts.nl
 Orders orders@iqproducts.nl
 www.iqproducts.nl