COAMATIC® Factor VIII - 82 2585 63

Intended use of the kit

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Bibliography / Literatur / Bibliografia / Bibliographie / Bibliografía

- ROSEN S et al: Clinical application of a chromogenic substrate method for determination of factor VIII activity. *Thromb Haemost* 54, 818-823 (1985).
- MERTENS K et al: The role of factor VIII in the activation of blood coagulation factor X by activated factor IX. *Thromb Haemost* 54, 654-660 (1985).
- LEHAGEN S et al: Clinical application of the chromogenic assay of factor VIII in hemophilia A and different variants of von Willebrand's disease. *Scand J Haematol* 37, 448-453 (1986).

Symbols used / Verwendete Symbole / Símbolos utilizados / Symboles utilisés / Simboli impiegati

IVD In Vitro Diagnostic Medical Device In Vitro Diagnostikum Dispositif médical de diagnostic in vitro Producto sanitario para diagnóstico in vitro Dispositivo medicodiagnóstico in vitro	LOT Batch code Chargenbezeichnung Codigo de lote Code du lot Codice del lotto	 Use by Verwendbar bis Fecha de caducidad Utiliser jusque Utilizzare entro	 Temperature limitation Zulässiger Temperaturbereich Limite de temperatura Limites de température Limite di temperatura	 Consult instructions for use Gebrauchsanweisung beachten Consulte las instrucciones de uso Consulte las instrucciones d'utilisation Consultare le istruzioni per l'uso	 Manufacturer Hersteller Fabricant Fabricante Fabbriicante	EC REP Authorised representative Bevollmächtigter Representante autorizado Mandatario Rappresentanza autorizzata
--	---	---	--	--	---	--

Additional reagents

- Deionized water, filtered through 0.22 µm or NCCLS type II water.
- Calibration plasma, calibrated against an International Standard
- Control Plasma Abnormal and Normal, calibrated against an International Standard for Factor VIII
- Saline (0.9% NaCl).
- Acetic acid 20% or citric acid 2% (end-point methods).

Material required but not provided

- Spectrophotometer, 405 nm (and 490 nm for microplate procedure)
- Incubator 37°C ±0.2
- Microplate or semi-micro cuvettes
- Centrifuge, 2000 x g
- Plastic test tubes
- Stopwatch
- Vortex mixer
- Calibrated pipettes
- Linear graph paper

Quality control

Appropriate controls for plasma or concentrates calibrated against an International Standard for Factor VIII should be used. Periodically within each run a control should be analyzed. The control material should be treated in the same way as a test sample. A range of allowable variation should be established for controls in each laboratory. If a value outside the established control range is obtained, a complete check of calibration, reagents and instrument performance should be made.

Results

Faktor VIII results are reported in IU/mL. 1.0 IU/mL of FVIII is equivalent to 100% FVIII.

Expected values

The Faktor VIII levels measured in 61 healthy individuals, 28 males and 33 females, aged between 21 and 55, were in the range of 0.5-2 IU/mL. Median 1.13 IU/mL, SD 0.39 IU/mL. (Microplate method).

Procedures

All conditions included in this package insert are referred to Microplate method. Detailed instrument settings including instructions for preparation of the reagents for a variety of automated instruments are available on request from Chromogenix.

Calibration

A standard curve should be performed with every run, prepared by different dilution of Calibration Plasma which should be traceable to the International Standard. A normal plasma can be also used for preparation of standard dilutions (it should be calibrated against an International Standard for plasma factor VIII). In case the normal plasma does not contain exactly 1 IU/mL (100%) factor VIII, the values of the standards must be recalculated accordingly, in order to obtain a correct factor VIII potency assignment.

This plasma is diluted according to a two-step procedure (predilution followed by final dilution) described in the below table. For factor VIII levels below 0.05 IU/mL (hemophilia A patients), the low range should be used.

	Plasma	Pre-dilution	Buffer working solution µL	Plasma predilution µL	Final dilution	Buffer working solution µL
Faktor VIII IU/mL	µL	µL	µL	µL	µL	µL
Normal range						
1.42	undiluted	—	25	25	1400	
1.00	undiluted	—	25	25	2000	
0.50	100	100	25	25	2000	
0.25	50	150	25	25	2000	
0	—	—	—	—	2000	
Low range						
0.050	100	1900	25	25	2000	
0.024	50	2000	25	25	2000	
0.012	25	2000	25	25	2000	
0.006	25	4000	25	25	2000	
0	—	—	—	—	2000	

Assay condition for microplate and test tube techniques

The test should be performed within 30 minutes after sampling or thawing of the plasma samples.

Dilution of samples and controls:

Samples/controls	25 µL
Buffer working solution	2000 µL
Mix well	

For extremely low levels of factor VIII, 0.005-0.05 IU/mL, a special range with increased reaction times is used in order to secure optimal resolution.

	Assay range
Diluted samples/controls/standards incubate at 37°C	0-1.5 IU/mL 50 µL
Faktor reagent (pre-warm at 37°C) incubate at 37°C	0-0.05 IU/mL 50 µL
S-2765 + I-2581 (pre-warm at 37°C)	3-4 min 50 µL
A. Kinetic method: read ΔA/min	2 min 4 min 50 µL
B. End-point method: proceed as described below incubate at 37°C	2 min 10 min 50 µL
Acetic acid 20% or 2% citric acid	50 µL 50 µL
A. Kinetic method: read the absorbance change at 405 nm for 30-120 S.	
B. End-point method: read the absorbance against buffer at 405 nm.	

2. For manual method in plastic tubes, use 200 µL instead of 50 µL for all pipetting steps in the above scheme

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