## COATEST™ APC™ Resistance V - S 0082313863

CE

**Intended use**

For determination of resistance to activated protein C (APC), caused by the factor V:Q<sup>506</sup> (factor V Leiden) mutation, in plasma from untreated individuals and from patients on oral anticoagulant (OAC) or heparin therapy.

**Background and summary**

The APC resistance phenotype<sup>1,2</sup> is in more than 90% of the cases due to a mutation in the factor V gene, resulting in a replacement of Arg<sup>506</sup> (R) with Gln (Q) in the factor V protein.<sup>3-5</sup> The selectivity for the factor V:Q<sup>506</sup>, or other mutations in the factor V gene rendering the protein resistant to inactivation by APC<sup>6</sup>, is increased by normalizing the concentration of other plasma proteins involved in formation and regulation of thrombin. Hence, by performing the APTT-based APC resistance assay in the presence of an excess of factor V deficient plasma, the sensitivity and specificity for the factor V:Q<sup>506</sup> mutation is significantly increased. Further, this modification allows for the analysis of plasma from patients who are on OAC therapy.<sup>7-11</sup>

**Measurement principle**

Sample plasma is prediluted in V-DEF Plasma and incubated with the APTT reagent for a standard period of time. Coagulation is triggered by the addition of CaCl<sub>2</sub> in the absence and presence of APC and the time for clot formation is recorded.

**Reagents**

**1. V-DEF Plasma:** 2 vials - Stabilized, lyophilized human plasma, with a low level of factor V activity, containing the heparin antagonist Polysorbate®. Reconstitute with 4.0 mL of NCCLS type II water.<sup>12</sup> Allow to stand for 30 minutes at 20-25°C. Swirl gently before use.

**2. CaCl<sub>2</sub>:** 2 vials - 2 mL of calcium chloride, 0.025 mol/L, in Tris buffer containing 0.5% bovine serum albumin.

**3. APTT reagent:** 2 vials - 4 mL of purified phospholipids with colloidal silica as contact activator. Contains a preservative. Mix thoroughly on a Vortex mixer before use.

**4. APC/CaCl<sub>2</sub>:** 2 vials - Human activated protein C colyophilized with CaCl<sub>2</sub>. Reconstitute with 2.0 mL of NCCLS type II water.<sup>12</sup> Allow to stand for 30 minutes at 20-25°C. Swirl gently before use.

**5. Control Plasma Level 1:** 1 vial - Lyophilized human plasma. Reconstitute with 1.0 mL of NCCLS type II water.<sup>12</sup> Allow to stand for 30 minutes at 20-25°C. Swirl gently before use.

**6. Control Plasma Level 2:** 1 vial - Lyophilized human plasma. Reconstitute with 1.0 mL of NCCLS type II water.<sup>12</sup> Allow to stand for 30 minutes at 20-25°C. Swirl gently before use.

**Reagents 3 and 4 are not interchangeable between lots.**

V-DEF Plasma (Cat. No. 00082314663), Control Plasma Level 1 (Cat. No. 00082265063) and Control Plasma Level 2 (Cat. No. 00082266863) are also available separately from Chromogenix.

**Caution**

APC/Calcium chloride, APC Control Plasma Levels 1 and 2, and Factor V Reagent Plasma contain human source material that tested non-reactive for HIV antibody, Hepatitis B Surface Antigen and Anti-HCV at the donor stage. This product, as with all human based specimens, should be handled with proper laboratory safety procedures to minimize the risk of transmission of infectious disease.

APC/Calcium chloride and Calcium chloride contain bovine material. All donor animals were sourced from BSE-free herds. The cattle received ante- and post-mortem health inspection by a veterinarian, and they were apparently free from infectious and contagious material. However, the material should be treated as potentially infectious.

**Calcium Chloride****WARNING**

Hazard class: Skin Sens. 1, H317. Aquatic Chronic 3, H412

Hazardous ingredients: Reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1) (CIT/MIT).

**Hazard statements:** H317: May cause an allergic skin reaction. H412: Harmful to aquatic life with long lasting effects.

**Precautionary statements:** P272: Contaminated work clothing should not be allowed out of the workplace. P273: Avoid release to the environment. P280: Wear protective clothing, protective gloves. P302+P352: IF ON SKIN: Wash with plenty of water. P333+P313: If skin irritation or rash occurs: Get medical advice/attention. P362+P364: Take off contaminated clothing and wash it before reuse.

**APC/CaCl<sub>2</sub>****DANGER**

Hazard class: Resp. sens. 1, H334

**Hazard statements:** H334: May cause allergy or asthma symptoms or breathing difficulties if inhaled.

**Precautionary statements:** P261: Avoid breathing dust/fume. P304 + P340: IF INHALED: If breathing is difficult, remove person to fresh air and keep comfortable for breathing. P342 + P311: If experiencing respiratory symptoms: Call a POISON CENTER/doctor. P501: Dispose of contents/container in accordance with local/regional/national/international regulation.

**Supplemental hazard information:** EUH208: Contains APC (Activated Protein C). May produces an allergic reaction. EUH210: Safety data sheet available on request. Up to 3.97% of the mixture consists of ingredients of unknown acute toxicity (dermal, inhalation) for the human health and unknown hazard to the aquatic environment.

**APTT Reagent, Control Plasma level 1, Control Plasma level 2, V-DEF plasma**

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**APTT Reagenz, Control Plasma level 1, Control Plasma level 2, V-DEF plasma**

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**Supplemental hazard information:**

**APTT Reagenz:** EUH 208: Contains 1,2-benzisothiazolin-3-one. EUH 210: Safety data sheet available on request. May produce an allergic reaction. Up to 20.85% of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and unknown hazard to the aquatic environment.

**Control Plasma level 1:** Up to 100% of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and unknown hazard to the aquatic environment.

**Control Plasma level 2:** Up to 98.75% of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and unknown hazard to the aquatic environment.

**V-DEF plasma:** Up to 98.26% of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and unknown hazard to the aquatic environment.

This product is for *in vitro* diagnostic use.

**Materials required but not provided:**

— Deionized water, filtered through 0.22 µm or NCCLS type II water.<sup>12</sup>

— Calibrated pipettes.

— Automated or semi-automated coagulation instruments which, employ mechanical or optical detection, methods should be used.

**NOTE:** When using automated or semi-automated instruments, always refer to the operator manual from the instrument manufacturer for exact procedures.

**Storage conditions and stability**

The sealed reagents are stable at 2-8°C until the expiry date printed on the label. Avoid contamination of the reagents by microorganisms.

**1. V-DEF Plasma:** Stability after reconstitution is 8 hours at 15-25°C, 24 hours at 2-8°C or 3 months at -20°C or below when stored in the original vial. \*See NOTE.

**2. CaCl<sub>2</sub>:** Opened reagent in the original vial is stable for 1 week at 15-25°C or 1 month at 2-8°C.

**3. APTT reagent:** Opened reagent in the original vial is stable for 1 week at 15-25°C or 1 month at 2-8°C. Do not freeze!

**4. APC/CaCl<sub>2</sub>:** Stability after reconstitution is 2 hours at 37°C, 8 hours at 15-25°C, 5 days at 2-8°C or 3 months at -20°C or below when stored in the original vial. \*See NOTE.

**5. Control Plasma Level 1:** Stability after reconstitution is 6 hours at 2-25°C or 3 months at -20°C or below when stored in the original vial. \*See NOTE.

**6. Control Plasma Level 2:** Stability after reconstitution is 6 hours at 2-25°C or 3 months at -20°C or below when stored in the original vial. \*See NOTE.

\*NOTE: Frozen reagent should be rapidly thawed at 37°C and gently mixed before use. Do not refreeze.

**Quality controls**

Control Plasma Level 1 and Level 2 should be used for validation of the assay series. Level 1 shows a normal response to APC whereas Level 2 shows a response consistent with presence of the factor V:Q<sup>506</sup> mutation. Ranges of expected APC-V ratios are provided with each batch. If values outside the specified range are obtained, a complete check of reagents and instrument performance should be made and the analysis should be repeated. (See Calibration section for QC use of Control Plasma Level 1 and 2).

**Traceability of calibrators and control materials**

The reported values were determined over multiple runs on ACL Futura using a specific lot of reagent and against an internal House Standard. As an International Standard is not still available for the APC-V assay, the values have been assigned against a House Standard which is traceable to frozen plasma samples which have been determined to be homozygous or heterozygous respectively for factor V dependent APC Resistance.

**Specimen collection**

The patient should be at rest for 10 min. before sampling. Collect blood (9 volumes) in 0.1 mol/L sodium citrate (1 volume) and centrifuge within 24 hours at 2000 x g for 20 min. at room temperature. Take care to avoid contamination from the platelet layer into the plasma when the plasma is separated from the cells. Analyse the plasma within 25 hours from blood sampling.<sup>13</sup> Alternatively, freeze rapidly at -70°C in aliquots of 1 mL or less and store for not more than 3 years at -70°C. Specimens should not be stored in a self defrosting freezer and not be thawed and refrozen before assay. Treat specimens as potentially infectious. For more information see NCCLS document H21-A3.<sup>14</sup>

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
Calculate the factor V related APC ratio for the samples and controls:

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## Control Plasmas Level 1 and 2

## Control Plasma Level 1

All Instruments	CaCl <sub>2</sub> (s)	APC/CaCl <sub>2</sub> (s)	APC-V ratio
-----------------	-----------------------	---------------------------	-------------

 LOT

## Control Plasma Level 2

All Instruments	CaCl <sub>2</sub> (s)	APC/CaCl <sub>2</sub> (s)	APC-V ratio
-----------------	-----------------------	---------------------------	-------------

 LOT

## Expected values

## ENGLISH - Insert revision 05/2022

Different clotting times will be obtained with different types of instruments depending on the clot detection principle. The following ranges cover the expected clotting times when performing analysis of Control Plasma Level 1 with the above batch number using COATEST™ APC™ Resistance V / V-S.

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SIZE: 11 x 17" (280 x 432 mm.).

PRINT: Front/Back.

PRINT COLOR: Sheet 1: Front - Top Band Orange PMS 137, Red PMS 179, it is acceptable to print black in lieu of Red PMS 179, all remaining type in black.

Back - Red PMS 179, it is acceptable to print black in lieu of Red PMS 179, all remaining type in black.

Sheet 2: Front - Top Band Orange PMS 137, Red PMS 179, it is acceptable to print black in lieu of Red PMS 179, all remaining type in black.

Back - Red PMS 179, it is acceptable to print black in lieu of Red PMS 179, all remaining type in black.

SPECIAL REQUIREMENTS: Folded

CARTON SIZE: WHITE BOX PCS0003022610 manual pack must have barcodes on the back of pages 2 and 4. Sheet 1 must have "A" at the end of the document part number, i.e., 000304123A. Sheet 2 must have "B" at the end of the document part number, i.e., 000304123B