

# S-2251™

X

# For Laboratory Use Only

For General Laboratory Use

S-2251 is a chromogenic substrate for plasmin and streptokinase-activated plasminogen.

## COMPOSITION

Each vial contains chromogenic substrate S-2251 25 mg and mannitol 60 mg as a bulking agent.

PRECAUTIONS AND WARNINGS:

Hazard class: None

Precautionary statements: None

Supplemental Hazard Information: ≈ 2% of the mixture consists of component of unknown acute toxicity (dermal, inhalation) for the human health and unknown hazard to the aquatic environment

## CHEMISTRY

Formula:

Mol. wt:

Chemical name: H-D-Valyl-L-leucyl-L-lysine-

p-Nitroaniline dihydrochloride H-D-Val-Leu-Lys-pNA · 2 HCl

551.6

ε<sub>316 nm</sub>: 1.27 · 10<sup>4</sup> mol<sup>-1</sup> · L · cm<sup>-1</sup>

Solubility: > 40 mmol/L in H<sub>2</sub>O

Stability: Substance: Stable until expiry date

if stored at 2-8°C. Avoid exposure to light. The substance is hygroscopic and should be stored dry. Solution: 3 mmol/L in H<sub>2</sub>O is

Solution: 3 mmol/L in H<sub>2</sub>O is stable for at least 6 months at 2-8°C. Contamination by microorganisms may cause

hvdrolvsis.

Suitable stock solution:

3-4 mmol/L in H<sub>2</sub>O.

## PRINCIPLE

### Enzyme

### KINETIC DATA

Plasmin (human): K<sub>m</sub> = 3 · 10<sup>-4</sup> mol/L

 $V = 0.5 \cdot 10^{-6} \text{ mol/min} \cdot CU$ 

Plasminogen  $^{-}$  SK:  $K_m = 2 \cdot 10^{-4} \text{ mol/L},$  $V = 1 \cdot 10^{-6} \text{ mol/min} \cdot \text{mL plasma}.$ 

Determined at 37°C in 2.5 mL 0.05 mol/L Tris buffer pH 7.4, I 0.5.



# **CHROMOGENIX**



Instrumentation Laboratory Company - Bedford, MA 01730-2443 (USA)

Instrumentation Laboratory SpA-V.le Monza 338-20128 Milano (Italy)

Tel. +1 800 678 0710

## STANDARDIZATION

An activity of  $\Delta A/min=0.05$  (37°C) is obtained by using a substrate concentration of 2 . km and:

- 1. 0.010 CU/mL of human plasmin from Chromogenix AB. 2. 0.0011 U/mL of the plasmin standard from NIBSC.
- 0.0011 U/mL of the plasmin standard from NIBSC, Potters Bar, London.
- 0.0078 CU/mL of SK-activated human plasminogen from Chromogenix AB.

The substrate is insensitive to kallikrein (glandular and plasma) and urokinase.

## APPLICATIONS

The substrate has been used for the determination of:

- 1. Antiplasmin in plasma (1,2,3,5)
- 2. Plasminogen in plasma (4,5,6,7)



- EDY J et al.: Inhibition of plasmin by normal and antiplasmindepleted human plasma. Thromb Res 8, 513-518 (1976)
- TEGER-NILSSON A-C et al.: Determination of a new rapid plasmin inhibitor in human blood by means of a plasmin specific tripeptide substrate. Scand J Clin Lab Invest 37, 403-409, (1977).
- Chromogenix AB. Determination of antiplasmin in plasma with S-2251.
  Laboratory Instruction.
- SORIA J et al.: Dosage du plasminogène a l'aide d'un substrat chromogène tripeptidique. Pathologie Biologie 24, 725-729 (1976).
- FRIBERGER P.: Methods for the determination of plasmin, antiplasmin and plasminogen by means of the substrate S-2251. Haemostasis 7, 138-145, (1975)
- PHILO R D and GAFFNEY P J.: Some observation on the assay of plasminogen using streptokinase and S-2251. Haemostasis 11 suppl. 1. 66 (1982).
- Chromogenix AB. Determination of plasminogen in plasma with S-2251. Laboratory Instruction.



# S-2251

Printed Insert Sheet: 301933

Revision: R3

Issued: 12/2016 C.O.: 473479

# **LANGUAGES**

**ENGLISH** 

# TECHNICAL SPECS

PAPER: White paper,

50-60 g/m<sup>2</sup> weight.

SIZE: 4.1 x 5.9" (104 x 150 mm.).

PRINT: Front/Back.

PRINT COLOR: All type in black.