TECHNOZYM® PAP Complex ELISA Kit

For research use only

REF	TC12060	TECHNOZYM® PAP Complex ELISA Kit	
REF	TC12062	TECHNOZYM® PAP Calibrator Set	5 x 0.5 mL
REF	TC12064	TECHNOZYM® PAP Control Set	2 x 0.5 mL

Symbols key / Symbolschlüssel / clé des symboles

	Manufactured by / Hergestellt von / fabriqué par	INC	Incubation buffer / Inkubationspuffer / tampon d'incubation
	Expiry date / Verfallsdatum / date d'expiration	LOT	Lot / Charge / lot
	Storage temperature / Lagertemperatur / température de conservation	MTP	Microtiter plate / Mikrotiterplatte / microplaques sensibilisées
	Consult instructions for use / Gebrauchsanweisung beachten / consulter la notice d'utilisation	REF	Catalogue number / Katalognummer / référence
	Determinations / Bestimmungen / déterminations	RTU	Ready to use / gebrauchsfertig / prêt à l'emploi
AQUA	Distilled water / destilliertes Wasser / eau distillée	RUO	For research use only / nur für Forschungszwecke / Pour la recherche uniquement
CONJ	Conjugate / Konjugat / conjugate	SUB	Substrate / Substrat / substrat
CONT	Control / Kontrolle / contrôle	STOP	Stop solution / Stopplösung / solution d'arrêt
CAL	Calibrator / Kalibrator / calibrateur	WASH	Washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado
DIL	Dilute or dissolve in / verdünnen oder lösen in / à diluer ou à dissoudre	STAB	Stable peroxide solution / Stabile Peroxidlösung / solution stable de peroxyde

PRODUCT DESCRIPTION

INTENDED USE

The TECHNOZYM® Plasmin- α -2-antiplasmin (PAP) complex ELISA can be used to detect elevated levels of this complex which may occur during thrombotic events, in cases of endogenous hyperfibrinolysis and during thrombolytic therapy.

COMPOSITION

- ELISA test strips (12): with 8 wells each, precoated with a monoclonal antibody directed against the neoantigen in the PAP complex; the drying agent is supplied in an aluminium bag
- Washing buffer concentrate (PBS; pH 7.3); containing detergent; 2 vials, 20 mL each.
- Dilution buffer concentrate (PBS; pH 7.3); contains stabiliser protein and protease inhibitors; 1 vial; 20 mL
- Calibrators (Standards) numbered 1-5; lyophilised; 1 vial each; 500 μ L.
Concentrations are lot-dependent; consult the label on the vial.
Control plasmas "low level" and "high level", lyophilised; 1 vial each.
Concentrations are lot-dependent; consult the label on the vial.
- Conjugate: anti-Plasminogen POX; dyed blue; 1 vial, 0.3 mL
- Chromogenic substrate TMB (tetramethylbenzidine); 1 vial; 12 mL; ready to use
- Stopping solution: sulphuric acid 0.45 mol/L; 1 vial; 12 mL; ready to use.
- Adhesive film: for ELISA test strips; 2 pieces

MATERIAL REQUIRED (but not supplied with the kit)

- Distilled water
- Test tubes for diluting calibrator and samples
- Measuring cylinder (1000 mL)
- Precision pipettes (10, 100 and 1000 μ L)
- Variable Pipette (1000 μ L)
- Multichannel and/or dispensing Pipettes (100 and 200 μ L)
- ELISA washer or multichannel pipette
- ELISA reader with 450 nm filter, with a 620 nm reference filter if available.
- Incubator (+37 °C)
- Vortex mixer

WARNING AND PRECAUTIONS

- For research use only.
- All human blood or plasma products as well as samples must be considered as potentially infectious. They have to be handled with appropriate care and in strict observance of safety regulations. The rules pertaining to disposal are the same as applied to disposing hospital waste.
- The calibrator made from human blood and any individual plasma involved in the procedure is HBsAg, HIV 1/2 Ab and HCV-Ab-negative (see labels on kit and/or bottles).
- Stopping solution (sulphuric acid) may irritate the skin. Should acid get into your eyes, wash out immediately with water and consult a doctor.
- The reagents sometimes contain preserving agents. Beware of swallowing! Avoid contact with skin or mucous membranes.

STABILITY AND STORAGE

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8°C.

Stability after reconstitution/opening:

Material/Reagent	State	Storage	Stability
Dilution buffer	After opening*	+2...8 °C	3 months
Wash buffer	After opening*	+2...8 °C	3 months
POX AB	After opening*	+2...8 °C	6 months
	Working solution	RT (+18...25 °C)	60 min
Calibrators, control plasmas	After reconstitution*	-20°C	6 months
		RT (+18...25 °C)	4 hours

* only if any contamination is avoided

TEST PROCEDURE

PREPARATION OF TEST SAMPLES

In general citrated plasma can be used. For samples undergoing thrombolytic therapy "inhibitor cocktail" plasma samples (citrae plasma containing 2000 KIU/mL aprotinin and 20 mM benzamidine) should be used. For this purpose collect the blood of samples to be tested in pre-coated plastic or siliconized tubes already containing 3.2 % buffered citrate in 1:10 ratio to blood as anticoagulant as well as aprotinin and benzamidine in a ratio yielding the final concentrations of 2000 KIU/mL and 20 mM, respectively. As alternative, commercially available PPACK tubes can be used. After filling the tubes, samples should be gently mixed by inverting the tube 5 times, and, then are placed in a crushed ice-water mixture. Centrifuge the blood within 30 min. after the puncture at 2000 xg for 30 min. at 4 °C (preferably in a cold centrifuge with swing out rotor). Immediately after centrifugation, plasma should be carefully pipetted off and collected in a small plastic tube. After carefully mixing of the pooled plasma, precooled and prelabelled plastic tubes for storage are filled. The freezing operation is very important: snap-freezing, that means almost instantaneous freezing, is highly recommended. Storage temperature should be kept constantly below a temperature of -30 °C and preferably at -70 °C. The total time between blood collection and plasma freezing should not exceed 90 minutes. Thawing and refreezing of plasma aliquots is not recommended. Thawing for assay is achieved rapidly using a water bath at 37 °C. After thawing, plasma samples are placed in a crushed ice-water mixture until analysis. Haemolytic and lipaemic plasmas may be used if no other samples are available.

PREPARATION OF REAGENTS

- Before starting the test, all required components are to be brought to room temperature.
- Preparation of washing buffer: Mix 20 mL of washing buffer concentrate with 230 mL of distilled water (1:12.5). Mix well!! There may be crystalline precipitations which will dissolve at 37°C within 10 minutes.
- Preparation of dilution buffer: Mix 20 mL of dilution buffer concentrate with 30 mL of distilled water (1:2.5). Mix well!!
- Reconstitution of calibrators and control plasmas: Calibrators and control plasmas are reconstituted with **500 μ L** distilled water and mixed for 10 seconds after a reconstitution time of 15 min. (vortex mixer).
Reconstituted components are clear to slightly turbid.
Calibrators and Control plasmas are used undiluted.
- Sample dilution:
 - For determination of LOW PAP LEVELS
Sample dilution 1:10
25 μ L plasma + 225 μ L dilution buffer
 - For determination of HIGH PAP LEVELS
Sample dilution 1:100 or 1:1000
Predilution 1:10:20 μ L plasma + 180 μ L dilution buffer
1:100: 25 μ L predilution + 225 μ L dilution buffer
1:1000: 10 μ L predilution + 990 μ L dilution buffer

- Preparation of conjugate working solution (1+50): Dilute 1 part by volume conjugate with 50 parts by volume incubation buffer.

For 8 test wells: Mix 20 μ L conjugate with 1000 μ L dilution buffer

PERFORMANCE OF THE TEST

SAMPLE INCUBATION (reference 1,2)	Pipette calibrators, control plasmas, diluted samples into test wells; cover test strips with film	100 μ L
	Incubate at +2...8 °C	Over night
WASHING (reference 1,3,4)	Washing buffer	4 x 200 μ L
CONJUGATE REACTION (reference 1,2)	Pipette conjugate working solution into wells, cover test strip with film	100 μ L
	Incubate at +37 °C	2 hours
WASHING (reference 1,3,4)	Washing buffer	4 x 200 μ L
SUBSTRATE REACTION (reference 1,2)	Pipette substrate solution into test wells, cover test strip with film	100 μ L
	Incubate at room temperature	20 minutes
STOPPING (reference 1,2)	Pipette stopping solution into wells	100 μ L
MEASURING (reference 5)	ELISA-Reader, 450 nm	shake 10 sec., measure within 10 min.

References:

- Reagents of different lots must not be combined
- Precision and performance, among others, essentially depend on the following factors:
 - Thorough mixing of all substances used for dilution, 10 sec. with Vortex Mixer
 - Test calibrators and samples in duplicates
 - Incubate at indicated temperature (Room temperature: 18...25 °C)
 - Strict observance of the order of pipetting and of the time element as indicated.
 - The time for sample incubation, conjugate and substrate reaction as indicated starts after pipetting the last sample. Incubation times shall not vary by more than \pm 5 %.
 - During sample incubation and conjugate reaction, the time for pipetting calibrators, samples and conjugate solutions must not exceed 60 seconds per ELISA test strip (8 wells).
 - During substrate reaction and at stopping, the time needed for pipetting the substrate and/or the stopping solution must not exceed 10 seconds per ELISA test strip. Short pipetting times may be secured by using multichannel- and dispensing pipettes.
- Label/number strips with a water resistant pen in case the strips accidentally fall out of the frame during testing.
- After the last washing, wells must be aspirated thoroughly, turned upside down and positioned on a blotting paper; by gentle tapping, the last remnants must be removed.
- By measuring the difference in wave length at 450 and 620 nm the precision of the test is increased.

LIMITATION OF THE TEST

Samples collected in citrated plasma without inhibitors may be used only if samples are not undergoing thrombolysis therapy. Otherwise PAP complexes can form ex vivo and lead to false high levels of PAP.

ANALYSES RESULTS

CALCULATION OF THE RESULTS

Setting up a reference curve:


X axis: concentration PAP complexes (ng/mL)
Y axis: extinction at 450 nm

Graph plot is linear – linear with a best fit.

Assessment of reference curve

- The validity of the test may be checked on the basis of the calculated control values.

Example of standard curve:

Measuring concentration of samples

- Read off the concentration from the reference curve.
Do not forget to multiply the results by a factor 10, 100 or 1000 for diluted samples.

REFERENCE RANGE

Normal range of PAP levels are 0 – 514 ng/mL. PAP-levels above 120 ng/mL during lytic therapy indicate complete consumption of α -2-antiplasmin which may lead to generalised plasminemia.

STANDARDISATION

The standard was calibrated against PAP complexes made in vitro from purified components.

PERFORMANCE CHARACTERISTICS

Performance data are given below. Results obtained in individual laboratories may differ.

PRECISION

Reproducibility was determined with different samples (in series and day to day). The following results were obtained.

Sample	Intra assay variation		Inter assay variation	
	Sample 1	Sample 2	Sample 3	Sample 4
N	10	10	3	3
Mean (ng/mL)	132.9	20.4	131.9	21.6
SD (ng/mL)	1.71	0.83	2.47	1.44
CV (%)	1.29	4.06	1.87	6.65

ASSAY RANGE

0.6 ng/mL – 225 ng/mL

(or up to the actual value of calibrator 1 and calibrator 5)

DETECTION LIMIT

2.7 ng/mL

LITERATURE

Please contact Technoclone or your local distributor.

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