# CMV Brite™ Turbo

*CMV antigenemia detection test*

INTERNATIONAL PACKAGE INSERT  
ENGLISH • DEUTSCH • FRANÇAIS

## CMV Brite™ Turbo Kit

Rapid CMV pp65 Antigenemia for the detection of active CMV infection

**Product code: VIR-CMV 110**

CE0344 In Vitro Diagnostic medical device

**Product code: VIR-CMV110 BDC**

510(k) #991650 In Vitro Diagnostic medical device

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# CMV Brite™ Turbo Kit

## Intended Use

The *CMV Brite™ turbo kit* Kit is intended for the rapid qualitative detection of Cytomegalovirus (CMV) lower matrix protein pp65 by indirect immunofluorescence using microscopy in isolated peripheral blood leukocytes obtained from ethylenediaminetetraacetic acid (EDTA) or heparin anti-coagulated human peripheral blood. The detection of CMV pp65 in human peripheral blood cells aids in the diagnosis of acute or reactivated CMV infection. This product is not FDA cleared (approved) for use in testing (i.e. screening) of blood or plasma donors.

## Summary and Explanation

Human cytomegalovirus (CMV) is a DNA virus, approximately 200 nm in diameter, belonging to the herpes virus family. It has a double-stranded DNA core, a capsid composed of 162 capsomeres, and a surrounding envelope (1). Infection with CMV is common and is usually subclinical in the healthy host. However, CMV infection in the immuno-compromised host and the developing fetus may result in either localized or disseminated disease. Its clinical manifestations include pneumonia, retinitis, hepatitis, enteritis, and neurologic disease. Despite improved treatment modalities, CMV infection may result in significant morbidity and mortality in transplant patients, AIDS patients and cancer patients, particularly those with leukemia or lymphoma. Patients are at risk from both primary CMV infection and reactivation of latent infection (2-4).

Early and rapid diagnosis of CMV infection is of great importance in avoiding over treatment with immunosuppressive drugs and in guiding antiviral therapy.

However, distinguishing asymptomatic viral shedding from significant disease that requires treatment can be difficult. Isolation of CMV from blood specimens (CMV viremia) has been shown to correlate best with symptomatic disease, whereas isolation from saliva and urine are commonly found without apparent illness. Conventional methods for isolation of CMV from blood are sensitive, but may not yield results for several weeks. Shell vial cultures provide a result within 1 to 2 days, but are not sensitive for detection of CMV in blood specimens. In contrast, detection of CMV antigen in peripheral blood polymorphonuclear (PMN) cells (CMV antigenemia) is both sensitive and rapid (5, 6). This technique uses monoclonal antibodies to detect the CMV lower matrix phosphoprotein (pp65), an early antigen in virus replication, that is abundantly present in antigen-positive PMNs (6-9). The *CMV Brite™ Turbo* assay can be completed within 2 hours of sample collection. The CMV antigenemia assay is valuable in the diagnosis and monitoring of active CMV infection in solid organ and bone marrow transplant patients (1, 2, 10-13, 19). It has also been reported that CMV antigenemia shows correlation with CMV disease in AIDS patients (4, 18).

## Principles of the Procedure

The CMV antigenemia assay has been developed using a cocktail of two monoclonal antibodies (C10/C11) directed against CMV lower matrix protein pp65 (6). The assay uses the C10/C11 cocktail in an indirect immunofluorescence staining of cytospin preparations of peripheral blood leukocytes. The **CMV Brite™ Turbo** antigenemia assay is completed within two hours of blood collection, which saves time and means a rapid answer for the clinician.

### *The CMV Brite™ Turbo method consists of:*

- a. Direct lysis of peripheral blood erythrocytes
- b. Preparation of cytospin slides
- c. Fixation and permeabilization
- d. Indirect immunofluorescence staining using monoclonal antibodies directed against CMV pp65 protein
- e. Reading and evaluation of results

The first step in the **CMV Brite™ Turbo** method involves direct lysis of the peripheral blood erythrocytes (22). Following lysis the leukocytes are cytocentrifuged onto a slide, fixed and permeabilized to allow subsequent detection of CMV pp65 antigen. The presence of the CMV pp65 antigen is detected by the C10/C11 antibody cocktail and visualized by means of a specific secondary FITC-labeled antibody. CMV antigen-positive leukocytes exhibit homogeneous yellow-green polylobate nuclear staining when observed using a fluorescence microscope. The number of CMV antigen-positive cells is counted per duplicate stain.

The whole procedure can be performed in approximately 2 hours.

## Reagents

*Materials and reagents supplied in the CMV Brite™ Turbo Kit*

*Reagents for 110 tests including 5 control slides.*

<b>CMV Brite™ Turbo Reagent A</b>	Erythrocyte Lysing Solution - Ammonium chloride solution containing < 0.1% sodium azide. (10x). (Dilute before use in demineralized water).  <b>WARNING</b>	200 ml
<b>CMV Brite™ Turbo Reagent B</b>	Fixative Solution - Formaldehyde, in PBS, containing < 0.1% sodium azide. (Dilute before use in PBS). (5x)  <b>DANGER</b>	290 ml
<b>CMV Brite™ Turbo Reagent C</b>	Permeabilization Solution - Igepal Ca-630 and Newborn Calf Serum in PBS containing < 0.1% sodium azide.(5x) (Dilute before use in PBS).	290 ml
<b>CMV Brite™ Turbo Reagent D</b>	Monoclonal antibody C10/C11 (IgG1/IgG1) against CMV lower matrix protein pp65 containing < 0.1% sodium azide. (Ready to use)	4.0 ml
<b>CMV Brite™ Turbo Reagent E</b>	FITC-conjugated sheep anti-mouse Immunoglobulins - with Evans Blue containing < 0.1% sodium azide. (Ready to use)	4.0 ml
<b>CMV Brite™ Turbo Control Slides</b>	CMV Antigenemia Control Microscope Slides - Control slide in a sealed pouch with desiccant. (Ready to use).	5 slides

Each control slide contains a negative control spot (fixed CMV antigen-negative leukocytes) and a positive control spot (fixed CMV lower matrix protein pp65 antigen-positive cells mixed with antigen-negative leukocytes). Prior to opening, allow control slides to equilibrate to room temperature.

Upon receipt, store reagents at 2-8 °C. Avoid direct sunlight. Reagents stored according to

stated storage instructions are stable until the expiration date indicated on the label.

The **CMV Brite™ Turbo Kit** is available as a 110 test kit and contains 5 control slides. The reagents are sufficient for a minimal of 110 tests, or, 100 patient tests and 5 control slides. The recommended procedure is to run the antigenemia assay in duplicate for each patient blood sample.

## Warnings and Precautions

### FOR IN VITRO DIAGNOSTIC USE

#### CMV Brite™ Turbo Control Slides

In the preparation of the **CMV Brite™ Turbo Control Slides**, leukocytes have been used obtained from a healthy human blood donor. Each donor sample has been tested by a FDA approved method and found non-reactive for the presence of antibodies to HIV-1, HIV-2, HCV and CMV as well as for HBsAg. While these methods are highly accurate, they do not guarantee that all infected samples will be detected. This product may also contain other human source material for which there is no approved test. Because no known test method can offer complete assurance that HIV-1, HIV-2, HCV, hepatitis B virus, CMV or other infectious agents are absent, all these reagents should be handled in accordance with good laboratory practices using appropriate precautions as described in "Biosafety in Microbial and Biomedical Laboratories," 2nd ed., 1988. HHS Publication No (CDC) 88-8395, Centers for Disease Control. The control slides supplied with the **CMV Brite™ Turbo Kit** have already been fixed and permeabilized.

#### Caution:

- Reagents containing sodium azide may react with lead or copper plumbing to form explosive metal azides. On disposal, flush with large amounts of water to prevent azide build-up.
- All reagents should be handled in accordance with good laboratory practices using appropriate precautions. In addition, handle all patient samples with appropriate precautions. Do not pipette by mouth.
- Reagent A (Erythrocyte Lysing Solution) contains Ammonium chloride.
  - ⚠ **WARNING**
  - H302 Harmful if swallowed; H319 Causes serious eye irritation. P305 + P351 + P338
  - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

- Reagent B (Fixative Solution) contains less than 9,5 % formaldehyde. Formaldehyde is a highly toxic allergenic and potentially carcinogenic reagent. Avoid skin or eye contact.

#### **DANGER**

- H301 + H311 + H331 Toxic if swallowed, in contact with skin or if inhaled; H314 Causes severe skin burns and eye damage; H317 May cause an allergic skin reaction; H335 May cause respiratory irritation; H351 Suspected of causing cancer; H370 Causes damage to organs. P260 Do not breathe dust/ fume/ gas/ mist/ vapours/ spray; P280 Wear protective gloves/ protective clothing/ eye protection/Face protection; P301 + P310 IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician; P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing; P310 Immediately call a POISON CENTER or doctor/ physician.
- Reagent E contains Evans Blue. Avoid contact with eyes, skin, clothing. Keep container tightly closed. Wash thoroughly after handling.

## Specimen Collection and Preparation

### *Processing of the blood sample*

Collect between 3 to 5 ml venous blood into an EDTA-treated tube, using aseptic venipuncture. Send the sample to the laboratory without delay. The blood sample should be kept at room temperature (20-25 °C) until processing. Processing should be performed within 6 to 8 hours of sample collection since it has been shown that a decrease of antigen-positive cells (collected in heparin) can be found after storage (15-17). Therefore, immediate preparation of the slides should always be attempted. In patients with severe neutropenia (absolute neutrophil count less than 200/mm<sup>3</sup>) at least 10 ml of blood may be required.

*Materials and reagents supplied*

• CMV Brite™ Turbo	Reagent A	200 ml
• CMV Brite™ Turbo	Reagent B	290 ml
• CMV Brite™ Turbo	Reagent C	290 ml
• CMV Brite™ Turbo	Reagent D	4.0 ml
• CMV Brite™ Turbo	Reagent E	4.0 ml
• CMV Brite™ Turbo	Control Slides	5 slides

*Laboratory material and equipment required**(but not included in the CMV Brite™ Turbo Kit)*

- Laboratory centrifuge
- 50 ml conical bottomed centrifuge tubes; sterile
- Pipettes
- Micropipette and tips
- Phosphate buffered saline (PBS); pH 7.4; Ca, Mg free
- Demineralized (distilled) water
- Hemocytometer or automated cell counter
- Cytocentrifuge slides
- Cytocentrifuge (such as Shandon Southern Products, Ltd., model Cytospin 3)
- Laboratory marker pen (Dako #52002)
- Coplin jars or histology staining jars
- Humid chamber
- 37 °C incubator
- Fluorescence microscope capable of 250x to 1000x magnification
- Mounting medium (non-fluorescent, such as CITI FLUOR, Glycerol PBS solution, UKC Chem Lab; or Immunoconcept Mounting Medium, Catalogue # 111)
- Fume Hood
- Micro cover glass
- Stop watch/Timer

## TEST PROCEDURE

*Note: Unless stated otherwise, reagents (including PBS) should be at room temperature when used in the test procedure. Room temperature is defined as 20-25 °C.*

*I. Preparation of leukocyte suspension*

*Note: Appropriate controls should be included for the specimen preparation process.*

- Dilute Reagent A 1:10 in demineralized (distilled) water and allow to cool to 4°C.
- Mix 2 ml blood with 30 ml of cold (+4°C) diluted erythrocyte lysing solution in a 50 ml conical bottomed tube and incubate for 5 minutes at +4°C.
- Centrifuge for 2 minutes at 1000xg (2500 rpm). Discard the supernatant.
- Repeat the cold lysing step if lysing of erythrocytes is not sufficient after the first time.
- Resuspend the cell pellet in 30 ml PBS
- Centrifuge for 2 minutes at 1000xg (2500 rpm). Discard the supernatant
- Resuspend the pellet in 1 ml PBS. (In patients with severe neutropenia re-suspension in as low as 0.2 ml may be required).

*II. Cell Counting*

- Count cells using a Hemocytometer or automated cell counter.
- Adjust concentration to  $2.0 \times 10^6$  cells/ml by diluting in PBS.

*III. Preparation of cytocentrifuge slides*

- Centrifuge 100 µl of the  $2.0 \times 10^6$  cells/ml suspension at approx. 600 rpm (54x g) for 4 minutes onto glass slides by cytocentrifuge.
- Prepare at least 3 slides per patient specimen (2 for testing, plus an additional slide for backup).

- Let slides dry for approximately 5 minutes.
- Circle cell area on the slide using a laboratory marker pen.
- Slides can be kept at room temperature overnight before fixation.

#### *IV. Fixation and permeabilization*

- Dilute fixative (reagent B) 1:5 in PBS in a fume hood prior to use.
- Dilute permeabilization solution (reagent C) 1: 5 in PBS in a fume hood prior to use. (Do not reuse).
- Immerse 2 slides in diluted reagent B for 5 minutes at room temperature in the fume hood.
- Dip slides 3 times in PBS (washing solution) and leave in the washing solution for 3 minutes.
- Immerse slides in diluted reagent C for 1 minute at room temperature.
- Dip slides 3 times in washing solution and place in fresh washing solution for 5 minutes (or any time up to 60 minutes). If staining with monoclonal antibody is to follow directly then proceed to step V.

To store slides:

If slides are to be stored then rinse in demineralized (distilled) water for 15 seconds. Let slides dry under the ventilator for approximately 20 minutes.

Once dry, slides should be packed in aluminum foil and stored at 4°C for 24 hours, or frozen at -80 °C.

#### *V. Immunofluorescence staining*

- From this point on, do not allow the cell preparations to dry out during the remainder of the staining procedure.
- **CMV Brite™ Turbo** Control slide. Rehydrate control slide in PBS for 1 or 2 minutes
- Remove one slide at a time from the washing solution, carefully dry the area surrounding the cell spot.
- Apply 35 µl of C10/C11 MoAb solution (reagent D) incubate for 20 minutes at 37 °C in a humid chamber.
- Dip slides 3 times in washing solution (PBS) and put in fresh washing solution for 3 minutes.
- Remove one slide at a time from the washing solution, carefully dry the area surrounding the cell spot.
- Apply 35 µl conjugate (reagent E), incubate for 20 minutes at 37 °C in a humid chamber.
- Wash twice in fresh PBS and carefully rinse with tap water (3 times) and mount with mounting medium and a micro cover glass.

#### *VI. Reading*

Perform reading as soon as possible. Slides may be stored for up to 8 hours at 4°C covered tightly in order to minimize fading. Perform microscopic evaluation using an immunofluorescence microscope at 400 x magnification. A higher magnification of 1000x may be used to increase resolution. (See note on recommendations for microscope and fluorescence cube under “Laboratory material and equipment”).

Scan the whole surface of the spot. Two spots should be scanned per patient.

Positive cells show homogenous yellow-green polylobate nuclear staining.

Negative cells show no yellow-green nuclear staining.

*Equivocal readings due to artifacts may occur rarely (less than 5%):*

- The most common artifact is due to eosinophils. They are recognized by a nonspecific cytoplasmatic staining, with the nucleus appearing as black holes, often spectacle-shaped. The cytoplasm may appear dull and more yellow.
- If all cells appear greenish, this represents an artifact, which may be associated with some types of patients and is very rare.
- A greenish staining of cells at the periphery of the spot only, is not considered positive. This may occur if the spot has started to dry out during incubation with either monoclonal antibody solution or conjugate solution. Be sure to incubate in a humid chamber and check that the whole cell spot is covered with the reagent.

If slides are not interpretable due to equivocal readings of artifacts, stain backup slide(s) or if unresolved obtain and test an additional specimen.

## Quality Control

The positive and negative control slides are provided with the kit. The slides are only used to check the staining procedure and do not influence the diagnostic value of the kit. The positive control must demonstrate appropriate

staining before patient specimens can be evaluated. The positive control should exhibit homogeneous yellow-green staining of positive cells with round morphology (nucleus is not visible). The negative control should show no yellow-green staining. Despite the fact that the quality procedures are followed strictly, single positive cells could be observed in the negative control. In such cases the staining run should not be considered invalid. The single positive cells on the negative control do not mean that the slide with patient material can not be interpreted.

## Results

### 1. Reading and Interpretation

Results of evaluation of patient specimen slides are qualitative. The clinical significance of an increase or decrease in the number of antigen-positive cells during time in one patient has not been determined. Therefore results should be reported only as Positive or Negative, according to the following definitions:

**Positive** - one or more CMV antigen-positive cells present per duplicate stain.

**Negative** - no CMV antigen-positive cells present per duplicate stain.

A minimum of approximately 50,000 specimen cells should be present in order to determine that a result is negative. If equivocal readings due to artifacts occur on both duplicate stains, the results should not be interpreted. Stain back-up slides or obtain and test an additional patient specimen.

## 2. Choice of Microscope

Selecting the correct filter cube configuration, which matches the fluorochrome in use, is important. This will vary for each microscope manufacturer. Select the correct combination for FITC use, such as the Olympus BX-40 or BX-60 microscope with:

<i>Fluorescence cube:</i>	<i>U-MNB</i>
<i>Excitation wavelength:</i>	<i>470-490 nm</i>
<i>Emission wavelength:</i>	<i>&gt;515 nm</i>

## Limitations of the Procedure

1. Detection of CMV pp65 lower matrix, early structural protein should be performed by laboratories experienced in immunocytochemical techniques.
2. Leukocyte preparation should be performed by personal experienced in aseptic techniques.
3. The efficacy of the **CMV Brite™ Turbo Kit** with samples other than human blood leukocytes has not been established.
4. Test performance characteristics have been established using the **CMV Brite™ Kit** (EDTA and heparin specimens) and validated in internal studies using the **CMV Brite™ Turbo Kit** (EDTA specimens only).
5. The **CMV Brite™ Turbo Kit** is intended for use only with immunofluorescent microscopy and not for use with flow cytometry.
6. The type and condition of the instrumentation used will influence the visual appearance of the image obtained. The reaction may vary due to the type of microscope employed, the light source, age of the bulb, filter assembly and filter thickness, differences in sensitivity of the antigen substrate, or the assay procedure used. Each laboratory should establish its own criteria for the reading of patient specimens using appropriate controls.
7. The detection and confirmation of CMV pp65 lower matrix, early structural protein in peripheral blood leukocytes is not diagnostic of symptomatic illness, since CMV pp65 may be present for a significant period following acute infection and patients with CMV antigenemia (especially low levels) may have asymptomatic illness. There are many published reports of patients who are antigenemia-negative and viremia positive and some of these patients may have CMV disease. Thus, a negative antigenemia test does not absolutely exclude CMV infection or disease.
8. The **CMV Brite™ Turbo Kit** is not intended for antiviral drug monitoring.
9. Test performance characteristics have not been established using neonatal specimens.
10. A decrease can be noted of antigen-positive cells (collected in heparin) per slide after storage (15-17). Therefore, immediate preparation of the slides should always be attempted. The maximal specimen storage time at room temperature has not been determined.
11. In patients with severe neutropenia (absolute neutrophil count less than 200/mm<sup>3</sup>) at least 10 ml of blood may be required.
12. Since the monoclonal antibodies have been prepared using a prototype strain, they may not detect all antigenemic or new strains of CMV. For example, monoclonal antibodies may fail to detect strains of CMV which have undergone minor amino acid changes in the target epitope region.
13. Test performance characteristics have not been established using specimens collected in heparin.

## Expected Values

The CMV Brite™ Kit was evaluated at two geographically distinct clinical virology laboratories in bone marrow or solid organ transplant or AIDS treatment centers (20). A total of 300 venous blood samples were collected from transplant recipients (159 allogeneic bone marrow and 47 solid organ), 85 HIV positive/AIDS patients, and 9 immunocompetent patients. 92 samples (30.7%) were determined to be CMV positive by the CMV Brite™ Kit. In random blood donor populations, the number of specimens found repeatedly reactive for CMV antigenemia by the CMV Brite™ Kit has typically been 0%.

## Specific Performance Characteristics

*Note: The performance characteristics have been established using the CMV Brite™ Kit and validated in internal studies using the CMV Brite™ Turbo Kit.*

The CMV Brite™ Kit was compared to CMV virus detection in conventional culture (CC) and shell vial culture (SV) using human peripheral blood leukocytes. 300 clinical samples were evaluated in comparison to conventional culture (CC) and shell vial culture (SV).

*Comparison of CMV Brite™ Kit with CC and SV (Clinical Site # 1)*

### Test CMV Brite™

CC/SV	+	-	Total numbers
+	3	0	3
-	45	103	148
Total numbers	48	103	151

*Comparison of CMV Brite™ Kit with CC and SV (Clinical Site # 2)*

### Test CMV Brite™

CC/SV	+	-	Total numbers
+	28	3	31
-	16	102	118
Total numbers	44	105	149

*Comparison of CMV Brite™ Kit with CC and SV (Two Sites Combined)*

### Test CMV Brite™

CC/SV	+	-	Total numbers
+	31	3	34
-	61	205	266
Total numbers	92	208	300

- Sensitivity  $31/34 = 91.2\%$   
(95% CI = 76.3 to 98.1%)
- Specificity  $205/266 = 77.1\%$   
(95% CI = 70.2 to 81.1%)

The imprecision of CMV culture (conventional and shell vial) when compared to CMV antigenemia has been documented in the literature (20).

Discrepant results appear to be the consequence of antiviral therapy or sample variability (12). Further, detection of CMV antigenemia has been accepted as clear evidence of active CMV infection in addition to shell vial culture and conventional culture (4th International CMV Workshop, Paris, April 19-21, 1993).

Discrepant results were resolved by looking for evidence of CMV infection at other sites in the patient, e.g. positive cultures of urine, throat, broncho alveolar lavage (BAL) samples and reviewing patient history for other related clinical findings evident of CMV infection. Additionally, evidence of CMV infection was obtained from histology and cultures of biopsies, and/or CMV IgG or IgM serology (not for bone marrow transplant patients or AIDS patients), and/or reference (in house validated, not commercially available) CMV antigenemia assay results.

*Comparison of CMV Brite™ Kit with other methods for detection of CMV infection (Clinical Site # 1)*

#### Test CMV Brite™

CMV expected	+	-	Total numbers
+	47	6	53
-	1	97	98
Total numbers	48	103	151

*Comparison of CMV Brite™ Kit with other methods for detection of CMV infection (Clinical Site # 2)*

#### Test CMV Brite™

CMV expected	+	-	Total numbers
+	42	3	45
-	2	102	104
Total numbers	44	105	149

*Comparison of CMV Brite™ Kit with other methods for detection of CMV infection (Two Sites Combined)*

19 of the 64 discrepant results were resolved on the basis of the reference CMV antigenemia assay result in conjunction with CMV serology.

#### Test CMV Brite™

CMV expected	+	-	Total numbers
+	89	9	98
-	3	199	202
Total numbers	92	208	300

- Resolved sensitivity  $89/98 = 90.8\%$   
(95% CI = 83.3 to 95.7%)
- Resolved specificity  $199/202 = 98.5\%$   
(95% CI = 95.7 to 99.7%)
- Positive Predictive Value  $89/92 = 96.7\%$   
(95% CI = 90.8 to 99.3%)
- Negative Predictive Value  $199/208 = 95.7\%$   
(95% CI = 91.6 to 97.9%)

The performance of the **CMV Brite™ Turbo Kit** was validated in a study of 183 patient samples. The patient samples included 173 organ transplant patients, 9 immunocompetent patients and 1 HIV positive patient sample. Each patient was tested in parallel with the **CMV Brite™** and **CMV Brite™ Turbo Kits**. The results of this validation study are provided in the table below.

Test CMV Brite™

CMV Brite™ Turbo Kit	+	-	Total numbers
+	43	7	50
-	6	127	133
Total numbers	49	134	183

## Cross Reactivity

Cross reactivity studies were performed using the **CMV Brite™** monoclonal antibody cocktail and conjugate, and single isolates of viruses (except for HSV-1 which was tested with 7 isolates) related to CMV, inducing similar clinical symptoms and/or isolated from peripheral blood. Unless otherwise noted, clinical isolates of the following viruses were tested.

- Herpes virus type 1
- Herpes virus type 2
- Varicella-zoster virus
- Adenovirus 2
- Adenovirus 4
- Adenovirus 5
- Parainfluenza virus 1
- Parainfluenza virus 2
- Parainfluenza virus 3
- Respiratory syncytial virus
- Poliovirus 3 (wild type 3)

- Echovirus 11
- Echovirus 30
- Epstein-Barr virus (laboratory strain P3HR1)\*
- Human herpes virus type 6 (laboratory strain GS) \*
- Human immunodeficiency virus (type 1 commercially available IF slides) \*\*

\* National Institute of Public Health and Environmental Protection, The Netherlands

\*\* Virion, Switzerland

No cross reactivity was observed except for a weak positive reaction with six of the seven HSV-1 isolates in the shell vial assay. The weak staining observed with the six HSV-1 isolates was cytoplasmic (i.e., seen only outside the nucleus, in a limited number of small foci), similar to the IF pattern usually found with HSV-1 monoclonal antibodies. The weak staining observed may have been due to Fc receptors being expressed by the HSV-1 infected cells in the shell vial assay. Subsequent analysis concluded that there is no evidence of cross reactivity between HSV and the CMV antigenemia assay using the C10/C11 monoclonal antibody cocktail.

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Product code: VIR-CMV 110  
CE0344 In Vitro Diagnostikum

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# Bibliography

1. Van Son, W.J. and The, T.H. (1989). Cytomegalovirus infection after organ transplantation: an update with special emphasis on renal transplantation. *Transpl. Int.* 2, 147-164.
2. The, T.H., Van der Bij, W., Van den Berg, A.P., Van der Giessen, M., Weits, J., Sprenger, H. and Van Son, W.J. (1990). Cytomegalovirus antigenemia. *Rev. Inf. Dis.* 12, S737-744.
3. Boeckh, M., Bowden, R.A., Goodrich, J.M., Pettinger, M., Meyers, J.D. (1992). Cytomegalovirus antigen detection in peripheral blood leukocytes after allogeneic marrow transplantation. *Blood* 80, 1358 - 1364.
4. Mazzulli, T., Rubin, R.H., Ferraro, M.J., D'Aquila, R.T., Doveikis, S.A., Smith, B.R., The, T.H. and Hirsch, M.S. (1993). Cytomegalovirus antigenemia: clinical correlations in transplant recipients and in persons with AIDS. *J.Clin. Microbiology* 31, 2824-2827.
5. Van der Bij, W., Schirm, J., Torensma, R., Van Son, W.J., Tegzess, A.M., The, T.H. (1988) Comparison between viremia and antigenemia for detection of cytomegalovirus in blood. *J. Clin. Microbiology* 26, 2531 - 2535.
6. Grefte, J.M.M., Van der Gun, B.T.F., Schmolke, S., Van der Giessen, M., Van Son, W.J. and The, T.H. (1992). The lower matrix protein pp65 is the principal viral antigen present in peripheral blood leukocytes during an active cytomegalovirus infection. *J.Gen.Virol.* 73, 2923-2932.
7. Gerna, G., Revello, M.G., Percivalle, E., Zavattoni, M., Parea, M., Battaglia, M. (1990). Quantification of human cytomegalovirus viremia by using monoclonal antibodies to different viral proteins. *J. Clin. Microbiology* 28, 2681 - 2688.
8. Revello, M.G., Percivalle, E., Di Matteo, A., Morini, F., Gerna, G. (1992). Nuclear expression of the lower matrix protein of human cytomegalovirus in peripheral blood leukocytes of immunocompromised viremic patients. *J. Gen Virol* 73, 437 - 442.
9. Gerna, G., Revello, M.G., Percivalle, E., Morini, F. (1992). Comparison of different immunostaining techniques and monoclonal antibodies to the lower matrix phosphoprotein (pp65) for optimal quantitation of human cytomegalovirus antigenemia. *J. Clin. Microbiology* 30, 1232 - 1237.
10. Van den Berg, A.P., Van der Bij, W. Van Son, W.J., Anema, J., Van der Giessen, M., Schirm, J., Tegzess, A.M., The, T.H. (1989). Cytomegalovirus antigenemia as a useful marker of symptomatic cytomegalovirus infection after renal transplantation: a report of 130 consecutive patients. *Transplantation* 48, 991 - 995.
11. Van den Berg, A.P., Klompmaker, I.J., Haagsma, E.B., Scholten-Sampson, A., Bijleveld, C.M.A., Schirm, J., Van der Giessen, M., Slooff, M.J.H., The, T.H. (1991). Antigenemia in the diagnosis and monitoring of active cytomegalovirus infection after liver transplantation. *J. Infect. Dis.* 164, 265 - 270.

12. Gerna, G., Zipeto, D., Parea, M., Revello, M.G., Sillini, E., Percivalle, E., Zavattoni, M., Milanese, G. (1991). Monitoring of human cytomegalovirus infections and ganciclovir treatment in heart transplant recipients by determination of viremia antigenemia and DNAemia. *J. Infect. Dis.* 164, 488 - 498.
13. The, T.H., van der Ploeg, M., Van der Berg, A.P., Vlieger, A.M., Van der Giessen, M., Van Son, W.J. (1992). Direct detection of cytomegalovirus in peripheral blood leukocytes - a review of the antigenemia assay and the polymerase chain reaction. *Transplantation* 54, 193 - 198.
14. Landry, M.L. and Ferguson D. (1993). Comparison of quantitative cytomegalovirus antigenemia assay with culture methods and correlation with clinical disease. *J.Clin.Microbiology* 31, 2851-2856.
15. Boeckh, M., Woogerd, P.M., Stevens-Ayers, T., Ray, C.G., and Bowden, R.A. (1994). Factors influencing detection of quantitative Cytomegalovirus antigenemia. *J.Clin.Microbiology* 32, 832-834.
16. Landry, M.L., Ferguson D., Cohen, S., Huber, K., and Wetherill, P. (1995). Effect of delayed specimen processing on Cytomegalovirus antigenemia test results. *J.Clin.Microbiology* 33, 257-259.
17. Niubò, J., Pérez, J.L., Carvajal, A., Ardanuy, C., and Martin, R. (1994). Effect of delayed processing of blood samples on performance of cytomegalovirus antigenemia assay. *J.Clin.Microbiology* 32, 1119-1120.
18. Gerna, G., Parea, M., Percivalle, E., Zipeto, D., Sillini, E., Barbarini, G., and Milanese, G. (1990). Human cytomegalovirus viraemia in HIV-1 seropositive patients at various clinical stages of infection. *AIDS* 4, 1027 - 1031.
19. Koskinen, P.K., Nieminen, M.S., Mattila, S.P., Häyry, P.J., and Lautenschlager, I.T. (1993). The correlation between symptomatic CMV infection and CMV antigenemia in heart allograft recipients. *Transplantation*, 55, 547-551.
20. Landry, M.L., Ferguson, D., Stevens-Ayers, T., De Jonge, M.W.A., and Boeckh, M. (1996). Evaluation of CMV Brite™ Kit for detection of Cytomegalovirus pp65 antigenemia in peripheral blood leukocytes by immunofluorescence. *J.Clin. Microbiology* 34, 1337-1339.
21. Erice, A., Holm, M.A., Gill, P.C., Henry, S.A. et al (1992). Cytomegalovirus (CMV) antigenemia assay is more sensitive than shell vial cultures rapid detection of CMV in polymorphonuclear blood leukocytes. *J. Clin. Microbiology* 30, 2822-2825.
22. Ho, S.K.N., Lo, C-Y., Cheng, J.K.P. and Chan, T-M., (1998) Rapid cytomegalovirus pp65 antigenemia assay by direct Erythrocyte Lysis and Immunofluorescence Staining. *J. Clin.Microbiology.* 36, 638-640.


# Transplantation Products

PRODUCT	APPLICATION	NO. OF TESTS	PRODUCT CODE	REGULATORY STATUS
CMV Brite™ Turbo Kit	Complete kit for the diagnosis of an active CMV infection (2 hours)	110	VIR-CMV 110 VIR-CMV110 BDC	IVD C €0344 510(k) #k991650 USA
CMV Brite™ Kit	Complete kit for the diagnosis of an active CMV infection (4 hours)	100	VIR-CMV 100 VIR-BDC100	IVD C €0344 510(k) #k951550 USA
CMV C10/C11 cocktail	Mix of 2 CMV pp65 monoclonal antibodies	200	VIR-CMV C10/C11	IVD C €0344
FITC Conjugate	Sheep anti-mouse-FITC	200	VIR-FITC	IVD C €0344
CMV Control slides	CMV pp65 positive and negative cytopots	set of 5	VIR-CMV CS05	IVD C €0344


Figure 1

Human peripheral blood leukocytes from a patient with an active CMV infection, stained with the **cmv Brite™ Turbo** antigenemia Kit. Immunofluorescence staining of CMV pp65 antigen positive polymorphonuclear cells. Negative cells are counterstained (red) with Evans Blue.

Figure 2

Control slide with the **cmv Brite™ Turbo kit** showing green immunofluorescence staining of a CMV pp65 antigen positive cell. The negative cells are normal human leukocytes counterstained with Evans Blue.

IQ Products

Rozenburglaan 13a

9727 DL Groningen

Netherlands

Tel: +31 (0) 50 5757 000

Fax: +31 (0) 50 5757 002

Email: [marketing@iqproducts.nl](mailto:marketing@iqproducts.nl)

[www.iqproducts.nl](http://www.iqproducts.nl)