
	ANALYSENZERTIFIKAT CERTIFICATE OF ANALYSIS	DQ6MBB0
---	---	---------


TECHNOCLOT® PT Owren manual 10x10mL

REF 5005037
LOT 6M45BB0.02
 31-Jan-2017

Analytische Prozedur Analytical Procedure	Anforderung Requirement	Ergebnis Result		
pH	7.4 – 7.8	7.7	OK	
Verwendbarkeit Usability KC 10	CCN LOT 1P37B00	90.2 – 122.0 % 0.84 – 1.14 INR	114.5 1.00	OK
	CCAK LOT 2B35B00	13.5 – 18.3 % 2.69 – 3.65 INR	16.2 3.11	OK
Etiketten und Barcode Überprüfung Label and bar code Verification		OK		

**Durchführung von
Prepared by**


Lucas ROSENMAYR


Datum / Date:

**Freigabe von
Approved by**


Fikret NASUFI


Datum / Date: