

# Factor IXa Calibrator

Art No 9599

Lot# **GSL1516**

Reconstitute in 2.0 mL purified water

---

**For In Vitro Research Use Only – Not for Diagnostic Use**

---

## PROPERTIES

Lyophilized preparation of human Factor IXa. Factor IXa activity after reconstitution in 2.0 mL purified water:

**6.8 mIU/mL**

The above Factor IXa activity apply only to **lot# GSL1516**

The calibration has been made against the 1<sup>st</sup> IS FIXa 97/562 using the Rox FIX-A kit (Art No 95 00 30)



Keep at -20°C or below



Use by 2019 04

## PACKAGE

10 x 2.0 mL

## STABILITY

A reconstituted vial is stable for 8 hours at 2-8°C.