

English

PL Chip Reservoir set for T-TAS[®] 01 Package Insert

Contains accessories for use with the PL Chip for T-TAS 01

INTENDED USE:

The PL Chip Reservoir set for T-TAS 01 contains accessories for use with the PL Chip for T-TAS 01.

SUMMARY:

The PL Chip Reservoir set for T-TAS 01 contains disposable accessories that are used with the PL Chip for T-TAS 01 for measuring primary hemostatic ability of whole blood samples.

For additional details about the use of the PL Chip Reservoir set for T-TAS 01, please consult the PL Chip for T-TAS instructions for use and the T-TAS 01 User's Manual.

MATERIALS PROVIDED:

The PL Chip Reservoir set for T-TAS 01 contains 100 sets of reservoirs, caps, and over-caps.

Item	Contents	Catalog Number
PL Chip Reservoir set for T-TAS 01	100 sets	18003

WARNINGS AND PRECAUTIONS:

- For in vitro diagnostic use only.
- For professional use only.
- Used reservoirs, caps and over-caps are potentially infectious. Proper handling and disposal methods should be followed in accordance with local, state and federal regulations.
- Carefully follow the instructions and procedures described in this package insert.
- Do not use materials that are bent, deformed, broken, or contain foreign matter.

STORAGE AND HANDLING REQUIREMENTS:

Store unused PL Chip Reservoir set contents at room temperature. Do not remove the reservoir, caps, or over-caps from their bag until ready for use.

ASSISTANCE:

For assistance, please contact your local distributor.


DEFINITION OF SYMBOLS:

Symbol	Definition
	Do Not Re-use
	CE Mark
	Consult Instructions for use

LOT	Batch Code
REF	Catalogue number
IVD	<i>In Vitro</i> diagnostic medical device
	Manufacturer
EC REP	Authorized representative In the European Community
	Contains sufficient for <n> tests
Rx only	This product is restricted to sale by or on the order of a licensed healthcare practitioner
	Importer



FUJIMORI KOGYO CO., LTD.
1-1-1 Koishikawa, Bunkyo-ku,
Tokyo 112-0002 Japan



Medical Device Safety Service GmbH
Schiffgraben 41 30175 Hannover, Germany



MedEnvoy
Prinses Margrietplantsoen 33 - Suite 123
2595 AM The Hague
The Netherlands

T-TAS® is a registered trademark of FUJIMORI KOGYO CO., LTD.

©2022 FUJIMORI KOGYO CO., LTD.

All rights reserved

