




Count on DiaPharma for your von Willebrand Factor (vWF) Research Products

Reagent-complete kits for vWF Activity, Antigen, Collagen Binding, Propeptide, and more.

Von Willebrand Disease (vWD) is the most frequent bleeding disorder. **Von Willebrand factor (vWF)** has two important functions in hemostasis: it promotes adhesion of platelets to the injured vessel wall, and functions as a protective carrier protein for factor VIII. vWF is a glycosylated protein organized in multimers, its subunits contain binding sites for platelet glycoprotein receptors and collagen.



Explore our full portfolio of von Willebrand Disease research products.



Ordering Information

Product	Description	REF
REAADS vWF Activity ELISA	vWF Activity	K10826
REAADS vWF Antigen ELISA	vWF Antigen	K034-001
TECHNOZYM® vWF:CBA ELISA*	vWF:CB	5450301
TECHNOZYM® vWF:CBA Calibrator Set*	Calibrator	5450310
TECHNOZYM® vWF:CBA Control Set*	Control	5450312
TECHNOZYM® vWF:Ag ELISA Kit*	vWF Antigen	5450201
TECHNOZYM® vWF:Ag Calibrator Set*	Calibrator	5450210
TECHNOZYM® vWF:Ag Control Set*	Control	5450212
Inter Array vWFpp ELISA*	vWF propeptide (vWFpp) Antigen	13.2.095.0096
T-TAS® 01	Microfluidic Clot Occlusion Instrument	18001
Coamatic® Factor VIII	Chromogenic Factor VIII	K822585
Coatest® SP Factor VIII	Chromogenic Factor VIII	K824086
Coatest® SP4 Factor VIII	Chromogenic Factor VIII	K824094
Rox Factor VIII*	Chromogenic Factor VIII	800070
Technochrom® Factor VIII*	Chromogenic Factor VIII	5344101

*For research use only. Not for use in diagnostic procedures.



VWF ACTIVITY MEASUREMENT

REAADS vWF Activity ELISA


- Mouse monoclonal anti-human vWF antibody recognizes the functional A1 domain of vWF

vWF propeptide ELISA*

- For research related to Type 1C vWD

Technozym® vWF:CBA ELISA*


- For research related to Type 2 vWD



VWF ANTIGEN MEASUREMENT

REAADS vWF Antigen ELISA

- Easily quantify vWF Antigen



FACTOR VIII ACTIVITY

Coamatic® FVIII

Coatest® SP/SP4 FVIII

ROX FVIII*


Technochrom® FVIII*

- Chromogenic Factor VIII kits featuring bovine Factor reagents

PLATELET FUNCTION TESTING

T-TAS® 01

- Total platelet function, simplified
- Global evaluation of overall hemostatic ability using physiologic, flow conditions



It's not just support. It's solutions.

Explore our wide portfolio of hemostasis assays at diapharma.com.

800.526.5224 | info@diapharma.com | diapharma.com



diapharma